



LTV Seat®
Owner's Manual



Package Contents

- LTV® Cushion
- LTV Cover
- Owner's manual
- Hand pump
- Repair kit

LTV Seat® Owner's Manual

To Adjust Your LTV Seat:

1. Turn the valve counterclockwise to open. Inflate with the hand pump until all of the air cells are visibly filled. Turn the valve clockwise to close, and then remove the hand pump.
2. Place the LTV Seat on your sitting surface: The non-skid side of the LTV Cover should face down, and the zippered portion of the cover should face the back. The valve should be on the front left corner.
3. Sit on the LTV Seat in your normal sitting position.
4. Open the valve to release air. As air is released, you will slowly sink into the LTV Seat. Release air until you get close to your desired comfort level, and then turn the valve clockwise to close. Release air or add more air until you reach your desired comfort level.

To clean the LTV Cover:

Remove the cover from the cushion.

If the cover is Quilted Fabric: Machine wash cold, gentle cycle, with like colors. Air dry.

If the cover is Ultraleather™: Hand wash with mild detergent. Wipe or air dry.



To clean the LTV Cushion:

1. Close the valve and place the cushion in a large basin. Use a liquid dishwashing detergent, laundry detergent, or another multipurpose detergent that is used for general cleaning. Follow the surface-cleaning directions on the product label.
 2. Use a soft plastic bristle brush, sponge, or cloth to gently scrub all surfaces.
 3. Rinse the cushion with clean water and air dry or towel dry.
- DO NOT expose to direct sunlight.
 - DO NOT use automatic dishwasher detergents, which are too caustic for personal safety.

Caution:

- Protect the LTV Seat from sharp objects and hot ashes.
- DO NOT use the LTV Seat as a water-flotation device. It will NOT support you in water.

LTV® Cushion made in the U.S.A.

© 2005, 2015 ROHO, Inc. LTV®, LTV Seat®, ROHO®, and shape fitting technology® are registered trademarks of ROHO, Inc. Ultraleather™ is a registered trademark of Ultrafabrics.

Limited Warranty

What is covered? ROHO, Inc. covers any defects in materials or workmanship of the product. Any substitutions of supplied components will void the warranty.

For how long? 12 months from the date the product was originally purchased, with the exception of the cover, which has a 6-Month Limited Warranty.

What is not covered? A product that has been misused, damaged by an accident, or damaged from "an act of Nature," e.g., flood, tornado, earthquake, fire. Any substitutions of supplied product accessories or modification to the product will void the warranty.

What we will do: Within a reasonable amount of time after you return the product to us, we will repair or replace at our option any defect in material and workmanship and ship it back to you, all free of charge.

What the customer must do: Within the U.S.A., the customer must contact the Warranty Service Department at warranty@roho.com or call toll-free 1-800-851-3449. Outside the U.S.A., contact your country's ROHO distributor. Note: Any items sent to ROHO, Inc. without prior authorization will not be returned to the customer.

Questions? See Contact Information.

Disclaimer: Any implied warranties, including the MERCHANTABILITY and FITNESS FOR A PARTICULAR PURPOSE are limited to the same time period of the product types referenced above in the "For how long?" section, and any action for a breach of such warranties or any express warranties herein must be commenced within this time. ROHO, Inc. shall not be responsible for any incidental or consequential damages with respect to economic loss to property, whether as a result of a breach of express or implied warranties, or otherwise.

Your Rights and the Effect of this Warranty upon Them: This warranty restricts the duration of all implied warranties as noted above and limits or excludes recovery for incidental or consequential damages. Some jurisdictions do not allow limitations on how long implied warranties will last, and some jurisdictions do not allow the exclusion or limitation of consequential or incidental damages, so those provisions may not apply to you. While this warranty gives you specific legal rights, you may also have other rights which vary from jurisdiction to jurisdiction.

Contact Information	
In the U.S.A.	Phone: 800-851-3449 • Fax: 888-551-3449 Warranty Service: warranty@roho.com Customer Care: cc@roho.com
Outside the U.S.A.	intl@roho.com Phone: 618-277-9150 • Fax: 618-277-6518 Or visit roho.com for a list of ROHO International distributors.
If directed to return any items, please send to:	ROHO, Inc. 1501 South 74th Street Belleville, IL 62223
ROHO, Inc. has a policy of continual product improvement and reserves the right to amend this document. A current version of this document is available at roho.com .	
This product may be covered by one or more U.S. and foreign patents and trademarks. For further details, please refer to roho.com .	



Manufactured by ROHO, Inc.

100 N. Florida Avenue, Belleville, IL 62221-5429 USA

U.S.: 800-851-3449 • Fax 888-551-3449

Outside the U.S.: 618-277-9150 • Fax 618-277-6518

Mailing/Shipping Address:

1501 South 74th Street, Belleville, IL 62223-5900 USA

TheROHOStore.com



T36000

Revision: 08/2015

Qty:

Print Date: