

permobil
POWER WHEELCHAIRS


REFINED DESIGN



C350
CORPUS 3G



Smaller turning radius than other rear wheel chairs allowing for greater maneuverability.



CORPUS SEATING

Power Elevate | Power Tilt | Power Recline | Power Legs



SMOOTH RIDE

4-wheel adjustable independent suspension for the ultimate ride quality.

COLOR OPTIONS



Vintage Wine



Extreme Yellow



Crystal Blue



Platinum Silver



Onyx Black

FABRIC OPTIONS



Leatherette



Mesh

SPECIFICATIONS

Maximum User Weight	300 lbs
Maximum Speed	5.5 mph, 6.5 mph
Range	13 mi*
Base Width	24.5"
Base Length	36"
Minimum Turning Radius	26"
Ground Clearance	3"
Recommended Battery Type	Group 34
Weight Including Batteries	367 lbs
Drive Electronics	90A /120A R-Net
Maximum Slope	10° Up, 10° Down†
Seat-to-Floor Height	17.5" (elevator) 17.5", 18.25", 19" (seat tube)

Optional Seat Elevator	8"
Tilt Options	0° – 50° (power) -5° – 45° (power)
Recline Options	85° – 175° (power) 85° – 120° (manual)
Legrest Elevation	85° – 170° (power) 85° – 170° (manual)
Backrest Height	20", 23" – 28"
Seat Depth	14" – 22"
Seat Width	17", 19", 21", 23"
Armrest Pad Lengths	10", 13", 16", 18"

*Actual driving range will vary based on driving and battery conditions.
†Refer to Owner Manual for safe driving practices.