

permobil

POWER WHEELCHAIRS

REFINED DESIGN




C350

CORPUS 3G



Smaller turning radius than other rear wheel chairs allowing for greater maneuverability.



CORPUS SEATING

Power Elevate | Power Tilt | Power Recline | Power Legs



SMOOTH RIDE

4-wheel adjustable independent suspension for the ultimate ride quality.

COLOR OPTIONS



Vintage Wine



Extreme Yellow



Crystal Blue



Platinum Silver



Onyx Black

FABRIC OPTIONS



Leatherette



Mesh

SPECIFICATIONS

| | |
|----------------------------|--|
| Maximum User Weight | 300 lbs |
| Maximum Speed | 5.5 mph, 6.5 mph |
| Range | 13 mi* |
| Base Width | 24.5" |
| Base Length | 36" |
| Minimum Turning Radius | 26" |
| Ground Clearance | 3" |
| Recommended Battery Type | Group 34 |
| Weight Including Batteries | 367 lbs |
| Drive Electronics | 90A /120A R-Net |
| Maximum Slope | 10° Up, 10° Down† |
| Seat-to-Floor Height | 17.5" (elevator) 17.5", 18.25", 19" (seat tube) |

| | |
|------------------------|---|
| Optional Seat Elevator | 8" |
| Tilt Options | 0° – 50° (power) -5° – 45° (power) |
| Recline Options | 85° – 175° (power) 85° – 120° (manual) |
| Legrest Elevation | 85° – 170° (power) 85° – 170° (manual) |
| Backrest Height | 20", 23" – 28" |
| Seat Depth | 14" – 22" |
| Seat Width | 17", 19", 21", 23" |
| Armrest Pad Lengths | 10", 13", 16", 18" |

*Actual driving range will vary based on driving and battery conditions.
†Refer to Owner Manual for safe driving practices.