
## M300 PS Jr







# \* Indicates a Required Field Dealer Information



* Contact:
Dealer Code:
* Dealer Name:
Address:
* City:
* State:
* Phone#:
* Fax#:
PO#:
Email Address:
Client Information
Permobil recommends that the client is evaluated by a
certified rehab specialist.
* First Name:
* Last Name:
Diagnosis:
Funding Source:
Client Age:
Client Measurements
* Weight:
*Height:
A) Top of Shoulders:
B) Chest Depth:
C) Chest Width:
D) Seat Depth:
E) Top of Head:
F) Elbow to Hand:
G) Seat Pan to Elbow:
H) Hip Width:
*I) Knee to Foot:

Permobil Australia Pty Ltd. Unit 3A / 6 Boundary Rd. Northmead, NSW 2152 Australia www.permobil.com



## **Chair Model and Colors**

CHOOSE BASE: Part Number    110857	<b>Description</b> M300 Base PS Jr	
	Price Includes: M300 Power Base, Shock Absorbing Suspension System, R-net Controller, Tie-Down Hardware For Strap Systems, Charger, Flat Free Tires, and Standard Positioning Belt. Must Choose Flat Free or Pneumatic Tires Below. 10kmh, Weight Capacity = 75kg, STFH = 17".	
CHOOSE WHEEL PA	ACKAGE:	
Part Number	Description Croy Tire Poekers MM/D	Price AUD
☐ I11054	Grey Tire Package - MWD  Includes grey tires on silver casters and drive wheel rims. Hub cap color is selected below.	No Charge
	Please select caster size and drive wheel type: 6" Casters: 8" Casters:	
□ 1440EE	Flat free drive tires: Pneumatic drive tires:	No Charas
☐ I11055	Black Tire Package - MWD Includes black tires on silver casters and drive wheel rims. Hub cap color is selected below.	No Charge
	Please select caster size and drive wheel type: 6" Casters: 8" Casters:	
□ 1440E6	Flat free drive tires: Pneumatic drive tires:	
☐ I11056	Blackout Wheel Package - MWD  Custom Option. Includes black tires and custom matte black casters, drive wheels and hub	
	caps. Please select caster size and drive wheel type:	
	6" Casters: 8" Casters: Flat free drive tires: Pneumatic drive tires:	
OPTIONAL DRIVE TI Part Number	RES: Description	Price AUD
110782	Tri-Spoke Split Rim w/ Snow Tire Pneumatic 3x8"	FIICE AUD
	Snow tires are designed for use in cold weather conditions. The nylon compound in these tires remains flexible at low temperatures, but wears more easily. Therefore, the snow tires (sold as wheel assemblies) are not recommended for year-long use and must be ordered along with	
	standard drive wheels. Not available with black tires or black rims.	
CHOOSE SHROUD	COLOR:	
Part Number	Description	Price AUD
☐ I105239-99-0	M300 Shroud - Onyx Black	No Charge
☐ I105243-99-0 ☐ I105242-99-0	M300 Shroud - Platinum Silver M300 Shroud - Sublime Green	No Charge No Charge
☐ I106172-99-0	M300 Shroud - Cobalt Blue	No Charge
☐ I106172-99-0	M300 Shroud - Phantom Grey	No Charge
☐ I106161-99-0	M300 Shroud - Pop Star Pink	No Charge
☐ I106173-99-0	M300 Shroud - Purple Rain	No Charge
I106162-99-0	M300 Shroud - Volcano Red	No Charge
_		J
CHOOSE HUBCAP C	Description	
T dit Humber	Match Hubcap Color to Shroud:	
	Note: Hub caps will be matched to shroud color if no selection is made.	
☐ I10871	Solid Color Hub Cap - Pair	No Charge
	Platinum Silver: Onyx Black: Cobalt Blue: Phantom Gray: Pop Star Pink: Volcano Red: Purple Rain: Titan Blue:	
	Sublime Green: Crystal Blue: Vintage Wine: Extreme Yellow:	
	Chili Pepper Red:	
CHOOSE ONE OF TI	HE FOLLOWING:	
Part Number	Description	Price AUD
<b>√</b> 110901	Australian Freight & Handling Fee	



## **Base Options**

MUST CHOOSE EITH Part Number	HER SEAT ELEVATOR OR FIXED SEAT TUBE: Description	Price AUD
☐ I101994-99-0	Seat Tube Fixed	No Charge
1106958-99-0	Power Adjustable Seat Height - 8" Travel	140 Onargo
MUST CHOOSE BAT	TERIES:	
Part Number	Description	Price AUD
☐ I107031-99-0	Batteries, Grp 24 (72Ah) Sealed Gel, x2, Installed	
☐ I107032-99-0	Batteries, Grp 34 (60Ah) Sealed Gel, x2, Installed	
OPTIONAL ITEMS:		
Part Number	Description	Price AUD
☐ I105428-99-0	Enhanced Steering Perform Unit (ESP), R-net	
—	Recommended item with switched drive inputs for improved tracking.	
I105297-99-0	LED Light Kit for M300 / M400 w/ R-net Electronics	
☐ ID51113	Programming Key for R-net - Permobil Version	
☐ I106223-99-0	Stealth BPAP Vent Tray - MWD Base	
	Designed for BPAP systems, the tray has a $12"W \times 7"D$ surface that can support up to 40 lbs. A $1"$ lip and adjustable strap help secure contents. The tray mounts to the rear tie-down loops, is tapered to fit between the rear swing arms and extends $1"$ beyond the trailing casters.	
Joystick Options		
	PLEASE NOTE:	
	For additional input devices, refer to Alterative Drive Systems section.	
☐ I105246-99-0	R-net LED Remote Joystick Module	No Charge
	This Joystick does NOT have a Display Screen and does NOT contain mono jacks. Limited to	
☐ I107492-99-0	five profiles with this Joystick Module. R-net Remote Color JS w/ Mono Jacks, No Lights	
☐ I107491-99-0	R-net Remote Color JS w/ Toggle Switches, No Lights	
☐ I107491-99-0 ☐ I104923-99-0	R-net Remote Color JS w/ Mono Jacks, Lights	
1107493-99-0	R-net Remote Color JS w/ North Jacks, Lights R-net Remote Color JS w/ Toggle Switches, Lights	
—		
☐ I105700-99-0	R-net Remote Color JS w/ Lights, Heavy Duty Gimbal	
	Heavy-Duty modification enables gimbal to hold up to extreme operating forces. Reinforced housing is attached to a standard Color LCD R-net Joystick. 3" long joystick knob has 3/4"	
	diameter. Not compatible with alternative joystick handles.	
☐ I103814-99-0	Compact Advanced Joystick for R-net	
_	Compact sized joystick that does not contain an LCD Display. Contains two mono jacks for	
	mode & on/off. Does NOT contain a wrist support pad. If you want a wrist support pad, please	
☐ I10621	select I10621 below. Add Wrist Pad to Compact Joystick Advanced	
	Add Whist Fad to Compact Joystick Advanced	
CHOOSE JOYSTICK		
Part Number	Description One of the stick Keeks Orandeed	Price AUD
☐ I10045	Cone Shaped Joystick Knob - Standard	No Charge
1100960-99-0	Chin Cup For Joystick Knob	
I100961-99-0	"T" Handle For Joystick Knob	
I100962-99-0	"Large Ball" For Joystick Knob	
I100963-99-0	"Softball" For Joystick Knob	
I100964-99-0	"Mushroom" For Joystick Knob	
☐ I100965-99-0	"Stick" For Joystick Knob	
☐ I105442-99-0	1.75" Black Foam Ball For Joystick Knob	
☐ I105846-99-0	Bodypoint J/S Handle 3" U-Shaped w/ Flex-Shaft	
☐ I105847-99-0	Bodypoint J/S Handle 4" U-Shaped w/ Flex-Shaft	
1105848-99-0	Bodypoint J/S Handle Dome Shaped, Rubber	



## **Joystick Options**

CHOOSE JOYSTICK I	MOUNT:	
Part Number	Description	Price AUD
☐ I10412	Joystick Mount Right - Fixed, R-net Only	No Charge
☐ I10413	Joystick Mount Left - Fixed, R-net Only	No Charge
☐ I102537-99-0	Joystick Mount Right - Swing Away	
☐ I102538-99-0	Joystick Mount Left - Swing Away	
☐ I103820-99-0	Retractable Joystick Mount Right	
☐ I103819-99-0	Retractable Joystick Mount Left	
☐ I103472-99-0	Mount Standard Joystick on Backrest Shell	
	Standard R-net Joystick is mounted to backrest. Can be chosen instead of PN ID50882.  Mount on Left: Mount on Right:	
☐ ID50882	Attendant Control for R-net Control System Only	
☐ I107350-99-0	R-net Bluetooth iDevice Module	
☐ I107367-99-0	Allows input device commands (such as joystick nudges or switch clicks) to be used to operate an iPhone, iPad or iPod. By assigning each input to an iOS function, users can make phone calls, send texts and emails and access applications. Kit includes a Bluetooth iDevice Module, an R-net Bus Cable, and two mono jack cables for direct switch inputs (switches not included). Must be ordered with a Color JSM or Omni for mode display.  R-net Bluetooth Mouse Module	
_	Enables mouse control through any R-net input device. It can operate the mouse on a computer and on select Android devices (with R-net for Android app). Kit includes a Bluetooth Mouse Module, an R-net Bus Cable, two mono jack cables for L/R click (switches not included) and a USB BT Adapter for user's computer. Note: BT Adapter may not work with all computers. Must be ordered with a Color JSM or Omni for mode display.	
☐ I107368-99-0	Add a second Bluetooth Mouse Module	
	Allows R-net to simultaneously pair to a second compatible computer or Android device. Each module is accessible through a separate R-net mode. Kit includes a unique Bluetooth Mouse Module, an R-net Bus Cable, two mono jack cables for L/R click and a USB BT Adapter.	
☐ I10582	Output Module (IOM) for R-net Control System	
	This allows a user to control an external device (ex: Communication Device, ECU, Door Opener, etc.) through the drive input (JSM, Head Array, etc.). Each module can provide up to 4 switch closures. May require specialty cable or wiring. Up to three output modules can be installed on one wheelchair.	
Seating System		
CHOOSE ONE OF TH	IF FOLLOWING:	
Part Number	Description	Price AUD
☐ I103712-99-0	PS Junior Seat Frame - Only Weight limit 75kg.	No Charge
☐ I103715-99-0	PS Jr Seat w/ 45° Manual Tilt  Must be ordered with Stroller Handles (I103811-99-0). Weight limit 75kg.	
☐ I103716-99-0	PS Junior Seat Frame with 45° Power Tilt  Weight limit 75kg.	
OPTIONAL ITEMS:		
Part Number	Description	Price AUD
☐ I103811-99-0	Stroller Handle for PS Jr Seating System	
	Angle adjustable. This item must be selected when ordering manual tilt (1103715-99-0).	

## **Seating System Options**

☑ I10585 Seat Function Control Kit for R-net

Please indicate below if you want push buttons or toggles on the ICS Switchbox.



## **Seating System Options**

CHOOSE PUSH BUTT Part Number  ☐ 1102798-99-0 ☐ 1102799-99-0	TONS OR TOGGLES:  Description  Push Buttons on ICS Switchbox  Toggle Switches on ICS Switcbox	Price AUD No Charge
1105409-99-0	Do Not Install ICS Switchbox, Send as Part	No Charge
	POSITION OF ICS SWITCHBOX:  Description  Mount ICS Switchbox onto R-net Joystick Bracket  ICS Switchbox will be mounted to the same bracket as the R-net Joystick Module. If this item is selected, the ICS Switchbox will swing away with the Joystick Module.	<b>Price AUD</b> No Charge
☐ I104651-99-0 ☐ I106256-99-0	Mount ICS Switchbox onto Armrest  Mount on Left: Mount on Right:  Mount ICS Switchbox on Backrest	No Charge
	Mount on Left: Mount on Right:  Description  ICS Alternative Switchbox for B. not Electronics	Price AUD
☐ I104196-99-0 ☐ I106743-99-0	ICS Alternative Switchbox for R-net Electronics  Seat function interface box that accommodates up to 8 momentary switches. Provides an additional method for accessing the seat functions. 1/8" mini-mono switches not included.  Remote Stop System - All Electronics  Consists of a Radio Frequency (RF) handheld transmitter and an RF receiver mounted on the wheelchair. Custom charges may be required when ordering with a PS/PS Jr seat, a low backrest height or a pediatric wheelchair. If ordering on a K450, push handles must also be selected.	
<b>Backrest Options</b>		
CHOOSE ONE OF THE Part Number ☑ I103717-99-0	E FOLLOWING:  Description  PS Jr Back Canes, Height & Angle Adjustable  Must also select PS Jr Style Armrests if ordering this backrest. Choose PS Jr Style Backrest  Type (Metal Back or Omit Back) below. Backrest width should match seat width. Back height is adjustable from 16.5"T - 22.5"T, measured from the seat pan to the top of the back.	<b>Price AUD</b> No Charge
CHOOSE ONE OF THE Part Number  110448	Description PS Jr Curved Metal Backrest Backrest width will match seat width. Back height is adjustable from 16.5"T - 22.5"T, measured from the seat pan to the top of the backrest. Please note: As back height is increased, a proportional gap is created between the seat pan and the bottom of the backrest. PS Jr - Omit Backrest	Price AUD No Charge
☐ I104616-99-0	Backrest Adjustment Quick Lock, PS / PS Jr	
Seat Sizing		
CHOOSE ONE OF THE Part Number ☐ I10432 ☐ I10433 ☐ I10434	E FOLLOWING:  Description PS Jr Seat Width 11" PS Jr Seat Width 12" PS Jr Seat Width 13"	Price AUD No Charge No Charge No Charge

No Charge

PS Jr Seat Width 16"

T 110437



#### Seat Sizing CHOOSE ONE OF THE FOLLOWING: **Part Number** Description **Price AUD** SET SEAT DEPTH AT: \_ □ I10985 No Charge Seat depth is adjustable from 5" to 18". Please round desired seat depth to nearest 0.5". CHOOSE ONE OF THE FOLLOWING: **Part Number** Description **Price AUD** ☐ I10276 Miniflex Seat Cushion Leatherette Approximately 2" thick. Please select cushion size: 10"W x 10"D: \_\_\_ \_\_\_ 12"W x 12"D: \_\_\_ 14"W x 14"D: Custom sizes are available upon request. □ I10999 Ergonomic Seat Cushion Leatherette Please select cushion size: 16.5"W x 15"D: 17.5"W x 16.5"D: \_\_\_\_\_ 17.5"W x 19.5"D: \_\_\_ Custom sizes available upon request. **Armrest Options** CHOOSE ONE OF THE FOLLOWING: **Price AUD Part Number** Description 110824 Flip Back, Height Adjustable Armrests Mounted off the backrest canes. Includes: 45° of angle adjustability, rotational adjustment, 3" offset adjustment both medially and laterally, and and height adjustment of 6" - 14" from the seat pan. In one inch increments, please indicate where you would like the armrests set. SET ARMRESTS AT (6" - 14"): **110825** Flip Back, Height Adjustable Armrests - NC No Charge No charge when ordering power tilt. In one inch increments, please indicate where you would like the armrests set. SFT ARMRESTS AT (6" - 14"): \_ CHOOSE ONE OF THE FOLLOWING: **Price AUD Part Number** Description 1105100-99-0 4"x10" Corpus Armrest Pad Right - Leatherette No Charge 1105104-99-0 4"x13" Corpus Armrest Pad Right - Leatherette No Charge 1107212-99-0 4"x10" Corpus Armrest Pad Right - Black Fabric No Charge 1106047-99-0 4"x13" Corpus Armrest Pad Right - Black Fabric No Charge 1106486-99-0 4"x10" Corpus Arm w/ Stealth Gel Pad - Right 4"x13" Corpus Arm w/ Stealth Gel Pad - Right 1106488-99-0 3.5"x8" Unitrack Armrest Pad Right (C) 1104951-99-0 No Charge 1104950-99-0 3.5"x12.5" Unitrack Armrest Pad Right (G) No Charge CHOOSE ONE OF THE FOLLOWING: **Price AUD Part Number** Description 1105098-99-0 4"x10" Corpus Armrest Pad Left - Leatherette No Charge 1105102-99-0 4"x13" Corpus Armrest Pad Left - Leatherette No Charge 1107213-99-0 4"x10" Corpus Armrest Pad Left - Black Fabric No Charge 1106046-99-0 4"x13" Corpus Armrest Pad Left - Black Fabric No Charge 4"x10" Corpus Arm w/ Stealth Gel Pad - Left 1106487-99-0 1106489-99-0 4"x13" Corpus Arm w/ Stealth Gel Pad - Left 1104952-99-0 3.5"x8" Unitrack Armrest Pad Left (C) No Charge 1104949-99-0 3.5"x12.5" Unitrack Armrest Pad Left (G) No Charge **OPTIONAL ITEMS:** Description Price AUD **Part Number** 1107339-99-0 Armrest Pouch - Right

Page 6 of 13

Armrest Pouch - Left

1107338-99-0



**Price AUD** 

No Charge

No Charge

#### **Armrest Options**

OPTIONAL ITEMS:		
Part Number	Description	Price AUD
☐ I107840-99-0	UT Armrest Side Support - Right Outer	
_	Includes 2"T x 6.5"W tapered pad.	
☐ I107863-99-0	UT Armrest Side Support - Left Outer	
	Includes 2"T x 6.5"W tapered pad.	
☐ I107864-99-0	UT Armrest Side Support - Right Inner	
	Includes 2"T x 6.5"W tapered pad. May require an increase to the width between the armrests	
—	to ensure armrests will flip-up.	
☐ I107865-99-0	UT Armrest Side Support - Left Inner	
	Includes 2"T x 6.5"W tapered pad. May require an increase to the width between the armrests	
	to ensure armrests will flip-up.	
☐ I107839-99-0	UT Elbow Support - Right	
	Must also select PN I10881.	
	UT Elbow Support - Left	
	Must also select PN I10882.	
☐ I10881	Adjustable Elbow Support Hardware - Right	
☐ I10882	Adjustable Elbow Support Hardware - Left	
<del></del>		

#### **Legrest Options**

**Part Number** 

CHOOSE ONE C	OF THE F	FOLLOWING:
--------------	----------	------------

Description

□ 14.0E96	70° Swing Away Lagrage w Multi Avia Adi Fastalata
☐ I10586	70° Swing Away Legrests w Multi Axis Adj Footplate
	Includes 5"W x 8"D Multi Axis Footplate. Range of adjustment from the seat pan to the top of each footplate is: 6 - 17".
☐ I10823	Aeron Legrest Assembly for PS Jr Seat - Narrow
	90° Legrest Hanger with Manual Adjustment Legrest Elevation. Price includes 6"W x 9"D Corpus Style, Angle Adjustable, Two Piece Footplates. Range of adjustment from the seat pan to the top of the footplate is: 12 - 17" (with calf pads) and 6.5 - 17" (without calf pads). Aeron Legrest Assembly for PS Jr Seat - Narrow
	90° Legrest Hanger with Manual Adjustment Legrest Elevation. Price includes 12.5"W x 10"D, Angle Adjustable, One Piece Footplate. Range of adjustment from the seat pan to the top of the footplate is: 12 - 17" (with calf pads) and 6.5 - 17" (without calf pads).

#### **OPTIONAL ITEMS:**

1104686-99-0

Description	Price AUD
UT Calf Support Kit (B) w Mesh Cover & Adj Hardware	No Charge
3.5"W x 4.5"T Pad. Recommended pad to be used with lower leg lengths of 12" or less.	
	UT Calf Support Kit (B) w Mesh Cover & Adj Hardware

Bracket designed for use with Aeron Center Mount Legrest Assembly. Cannot be ordered with Swing Away Legs (I10586 & I10587).

UT Calf Support Kit (J) w Mesh Cover & Adj Hardware

5"W x 6.5"T Pad. Bracket designed for use with Aeron Center Mount Legrest Assembly.

Cannot be ordered with Swing Away Legs (I10586 & I10587).

I107503-99-0 BP Ankle Huggers, Med Side Release, Adj (Pair)

This item will be sent as a part with the wheelchair. It will not be installed.

#### **Accessories**

#### MUST CHOOSE POSITIONING BELT:

Part Number	Description
	NOTE:
	See supplemental order form for more Bodypoint options.
<b>√</b> 1104518-99-0	Bodypoint Non-Padded Hip Belt - Medium

odypoint Non-Padded Hip Belt - Medium No Charge

Bodypoint Padded Hip Belt - Extra Small

#### M300 PS Jr



## Accessories

CHOOSE ONE OF		
Part Number	Description (DC) to the state of the state o	Price AUD
☐ I100749-99-0	MF Headrest for Koala/PS Jr - Leatherette	
	Must also select PN I10552.	
☐ I10552	Adj/Removable Headrest Hardware	
OPTIONAL ITEMS:		
Part Number	Description	Price AUD
1103077-99-0	UT Hip Suppt (B) w/ Mesh Cover - Low (Pair)	
_	4.5"W x 3.5"T Pad w/ 2" Tall Mounting Bracket. Must also select PN I10550.	
1103079-99-0	UT Hip Suppt (A) w/ Mesh Cover - Low (Pair)	
_	6"W x 5"T Pad w/ 2" Tall Mounting Bracket. Must also select PN I10550.	
☐ I103095-99-0	UT Thigh/Hip Suppt (C) w/ Mesh Cover - High (Pair)	
_	8"W x 3.5"T Pad w/ 4" Tall Mounting Bracket. Must also select PN I10550.	
☐ I10550	Adj Removable - Thigh Support Hardware (Pair)	
I107499-99-0	Universal Headrest Adapter	
	Please Note: Recline is limited to 160° or less (depending on backrest height) when ordering a	
	Universal Headrest Adapter on a 23" or 24" tall Corpus backrest.	
☐ I106210-99-0	RAM Self-Leveling Cup Holder	
☐ I106806-99-0	RAM X-Grip Phone Holder	
1107100-99-0	RAM X-Grip 10" Tablet Holder w/Swing Away Mount	
_	Compatible with most 10" screen tablets with or without a protective case. Includes adjustable	
	swing away mount (I103476-99-0). Note: Customization may be required if ordered along with an Omni.	
1107222-99-0	USB Charger	
_	Provides up to 1.5A (at 5V) to charge most smartphones and tablets. The charger features an isolated circuit design for improved safety; and by pulling power from the ICS, prevents battery drainage when the chair is powered off. Includes mounting hardware and ICS connector.  Please select armrest for mounting, keeping in mind that charger may need to be relocated.  Mount on left armrest: Mount on right armrest:	
☐ I105805-99-0	Medical Necessities Bag Hooks	
☐ I102059-99-0	Please Note: Recline is limited to 165° or less (depending on backrest height) when ordering with 23" or 24" tall Corpus backrest.  Medical Necessities Bag	
☐ ISH212	Includes hardware to attach bag to backrest. Please Note: Recline is limited to 165° or less (depending on backrest height) when ordering with a 23" or 24" tall Corpus backrest.  BP Shoulder Harness Medium, Front Pull	
ш -	This item will be sent as a part with the wheelchair. It will not be installed.	
☐ I100578-99-0	Tray, Small, Angle Adjustable	
	Must also select PN I10556.	
☐ I105221-99-0	Tray, Angle Adjustable, 14" W	
_	Must also select PN I10556.	
☐ I105222-99-0	Tray, Angle Adjustable, 16" W	
_	Must also select PN I10556.	
T106863-99-0	Tray, UniTrack, Angle Adjustable	
	Side mounted tray that attaches to the UniTrack Accessory Rail. Dimensions are: 19.5"W x 13"D. Can only be ordered with chairs that offer the UniTrack System. Must also select PN I10556.	
☐ I10556	Adj Removable - Tray Hardware (1 Each)	
☐ I103798-99-0	After Market Seating System Adapter	
	Used to connect after market seating systems to the UniTrack Accessory Rail on the PS Jr Seat Frame. Includes Mounting Hardware.	

## **Stealth Accessories**

Rev 2017-1



## **Stealth Accessories**

OPTIONAL ITEMS: Part Number  1106222-99-0	Description Stealth 3" x 5" Lateral Supports - PS Jr Must also select I10723. Includes THIN Stealth hardware and specialized mounting bracket for Flip Back Armrest assembly. Only compatible with Flip Back Armrest hardware. Adj Removable Stealth Lateral Hardware (Pair)	Price AUD
Alternative Drive S	Systems	
	ORTIONAL JOYSTICK SYSTEMS:	
Part Number  [] 110751	Description Mini-Joystick for Omni / IOM (9-pin connector) Operating Force = 13 g Neutral to Full Deflection = 5mm	Price AUD
☐ I103893-99-0	Compact Joystick for Omni / IOM (9-pin connector)  Operating Force = 220 g Neutral to Full Deflection = 28mm Joystick includes two 3.5mm mono jacks for Mode and On/Off. Includes Standard Joystick Knob and 1.75" (42mm) Foam Ball Knob.	
☐ I103894-99-0	Compact Joystick LITE for Omni / IOM (9-pin conn)  Operating Force = 50 g Neutral to Full Deflection = 17mm Joystick includes two 3.5mm mono jacks for Mode and On/Off. Includes Standard Joystick Knob and 1.75" (42mm) Foam Ball Knob.	
ALT PROPORTIONAL Part Number	JOYSTICK MOUNTING HARDWARE:  Description  NOTE:	
☐ I104690-99-0	Please see the PERMOfix Mounting Solutions, within the Alt Drive Accessories section, for additional mounting options.  Alt. JS Swing Away Mount Kit - Armrest	
☐ I107135-99-0	Mount in front of Left Arm: Mount in front of Right Arm: Alt JS/Omni Combo Swing Away Mount Kit - Armrest Used to mount Omni and Alternative Proportional Joystick to the same armrest. Includes	
	Swing Away Omni/Tablet Mount, with a longer mounting rod, two PERMOfix modules, one 3" rod and one 4" rod.  Mount in front of Left Arm: Mount in front of Right Arm:	
☐ I104695-99-0	Alt. JS Adjustable Mount Kit - Midline  Mount on Left Arm: Mount on Right Arm:	
☐ I104143-99-0	Alt. JS Adjustable Mount Kit - Chin Boom	
☐ I10752	Swing-Away / Removable Chin Boom mounted on right of backrest shell. Custom charge may be applied if ordering lower backrest heights.  Alt. JS Adjustable Mount Kit - Haims Harness, Lrg	
☐ I10753	For neck width up to 19cm. Alt. JS Adjustable Mount Kit - Haims Harness, Sml	
	For neck width up to 14cm.	
ALT PROPORTIONAL Part Number	JOYSTICK ACCESSORIES: Description	Price AUD
☐ I104904-99-0	Dual Switches for Compact Joystick - Fixed Mount Includes two pushbutton switches mounted on rigid Switches can be used for Mode & On/Off.	
☐ I104641-99-0	Switch operation force = 130 g.  Mount at front of Joystick: Mount at rear of Joystick:  Mount at left of Joystick: Mount at right of Joystick:  Dual Switches for Compact Joystick - Gooseneck Mnt  Includes two pushbutton switches mounted on adjustable gooseneck arms. Switches are	
01D & DUEE 0.40TE: 1	positioned on Left and Right sides of Compact Joystick. Switches can be used for Mode & On/Off. Switch operation force = 130 g	

SIP & PUFF SYSTEMS:

Part Number Description

Price AUD



## **Alternative Drive Systems**

SIP & PUFF SYSTE		Daire AUD
Part Number	Description Sip & Puff System for Omni	Price AUD
	Requires Omni Control Module. Cannot use this part with an Input Module. Includes breath tube kit.	
☐ I10574	Sip & Puff System with Fiber Optic Opti-Stop	
	Requires Omni Control Module. Cannot use this part with an Input Module. Includes breath tube kit.	
☐ I106807-99-0	Sip & Puff Kit for Total Control Head Array	
	For use with Total Control Head Array. Requires Omni Control Module. Cannot be used with an Input Module. Includes breath tube kit.	
T 110762	Adjustable Mounting Hrdw - Sip & Puff	
_	Required when selecting PN I10573, I10574 or I106807-99-0.	
PROPORTIONAL FO	OOT CONTROL:	
Part Number	Description	Price AUD
☐ I10767	Proportional Foot Control	
<b>—</b>	Mount on Left: Mount on Right:	
☐ I10766	Adjustable Mounting Hrdw - Foot Control	
	Required when selecting PN I10767.	
Total Control Hea	ad Array	
CONFIGURE HEAD	ARRAY:	
Part Number	Description	Price AUD
☐ I106805-99-0	Permobil Head Array - 3-Switch System	
☐ I11072	Features an occipital pad sensor and two stainless steel arms configured with proximity sensors. The arms are mounted to the upper attachment points on the head array. Longer arms are also included to extend the sensor adjustment range. ESP recommended. Choose Omni Control Module (ID51154IR) separately.  Program Head Array as a 5-Switch System	No Charge
	Select additional arms/sensors below and indicate the intended use of each sensor in the Notes & Comments section. Choose Omni Control Module (ID51154IR) or Input Module (I10581) separately.	ŭ
	CHOOSE OCCIPITAL PAD SIZE (REQUIRED):	
☐ I105920-99-0	Small Occipital Pad for Head Array - 4.25"W x 3"H	No Charge
1105921-99-0	Large Occipital Pad for Head Array - 6"W x 4"H	No Charge
	CHOOSE OCCIPITAL PAD MOUNT (REQUIRED):	
☐ I10847	Fixed Mount for Occipital Pad - Head Array	No Charge
1105936-99-0	Adj. Swivel Mount for Head Array Occipital Pad	_
_	Enables Occipital Pad to be angled up/down and left/right independent of the lateral arms.	
	CHOOSE BACKREST MOUNTING HARDWARE (REQUIRED):	
T105927-99-0	Stealth TWB480 Backrest Mounting Hdwr with Links	
_	Please select the length of the vertical stem below:	
<b>—</b>	Standard length (10"): Short length (6"):	
1105941-99-0	Stealth TWB480 Backrest Mtng Hdwr with Flipback	
1105925-99-0	Stealth HMO475 Backr Mtng Hdwr w 17mm Ball	
☐ I106319-99-0	Stealth TWB480 Backrest Mnt Hdwr w SubOcptl Mount	
☐ I106282-99-0	Stealth 600 SubOccipital Pad - 8.5"W x 2.75"H	
	Must select when PN I106319-99-0 is ordered.	
	SELECT ADDITIONAL LEFT ARM PACKAGE (OPTIONAL):	



## **Total Control Head Array**

CONFIGURE HEAD A	RRAY: Description	Price AUD
☐ I10976	Add Left Lateral Arm w Proximity Sensor	
	Choose Lateral Arm Adjustment Range: Short (5.5" - 10"): Long (5.5" - 14"):	
☐ I10977	Add Left Lateral Arm w Egg Switch	
_	Choose Lateral Arm Adjustment Range:	
☐ I10978	Short (5.5" - 10"): Long (5.5" - 14"): Add Left Lateral Arm w Pad (No Sensor or Switch)	
	Choose Lateral Arm Adjustment Range:	
	Short (5.5" - 10"): Long (5.5" - 14"):	
	SELECT ADDITIONAL RIGHT ARM PACKAGE (OPTIONAL):	
☐ I10979	Add Right Lateral Arm w Proximity Sensor	
	Choose Lateral Arm Adjustment Range: Short (5.5" - 10"): Long (5.5" - 14"):	
☐ I10980	Add Right Lateral Arm w Egg Switch	
_	Choose Lateral Arm Adjustment Range:	
<b>—</b> 14.0004	Short (5.5" - 10"): Long (5.5" - 14"):	
☐ I10981	Add Right Lateral Arm w Pad (No Sensor or Switch)  Choose Lateral Arm Adjustment Range:	
	Short (5.5" - 10"): Long (5.5" - 14"):	
Alt Drive Systems	Satur	
All Drive Systems	Setup	
CHOOSE OMNI OR II	NPUT MODULE (REQUIRED): Description	Price AUD
☐ ID51154IR	Omni Control Module for R-net w/ Infrared	
_	Must also select mounting hardware, either I103476-99-0 or I107135-99-0. Please specify which type of alternative drive input you will be using if not ordering it from Permobil. Input Device:	
	Adjustable Swing-Away Hardware - Omni/Tablet	
_	Mount in front of Left Arm: Mount in front of Right Arm:	
☐ I10581	Input Module (IOM) for R-net Control System	
	This device allows a proportional or switched drive input to be connected to an R-net Control System without the use of an Omni. It must be ordered in conjunction with an R-net Remote	
	Joystick. It does NOT support any of the following input devices: Sip-n-Puff, 3-Switch Head	
	Array or RIM control.	
	POSITION OF SEAT CONTROL SWITCHBOX:	
Part Number	Description  Mount Seat Control Switchbox on Backrest	Price AUD No Charge
□ 110000	Would Seat Control Switchbox on Backlest	No Charge
Alt Drive Accesso	ries	
PERMOFIX MOUNTIN	IC SOLUTIONS:	
Part Number	Description	Price AUD
I106655-99-0	PERMOfix UniTrack Fixed Mount Kit	
_	Includes two base joints, two 4" rods and one 8" rod.	
☐ I106656-99-0	PERMOfix UniTrack Removable Mount Kit	
	Includes quick release mount and thigh support hardware, two base joints and three different length rods.	
☐ I106657-99-0	PERMOfix UniTrack Swing Away Mount Kit	
	Includes PERMOfix Swing Away module, two base joints, one 2" rod and one 4" rod.	



### **Alt Drive Accessories**

PERMOFIX MOUNTI	NG SOLUTIONS: Description	Price AUD
☐ I105907-99-0	Proximity Sensor + Swivel Assm - Not Mounted	
	Only available with the Total Control Head Array. Swivel assembly allows sensor to be mounted with PERMOfix or directly to lateral arm of Total Control Head Array. Since Proximity Sensors require power, they CANNOT be used as On/Off switches.	
☐ I105919-99-0	Egg Switch + Swivel Mount Assm - Not Mounted Swivel assembly allows switch to be mounted with PERMOfix or directly to lateral arm of Total Control Head Array.	
☐ I105912-99-0	Solid Lateral Pad+Swivel Joint - Not Mounted	
	Swivel assembly allows pad to be mounted with PERMOfix or directly to lateral arm of Total Control Head Array.	
☐ I106675-99-0	Oval Plate + Swivel Assm - No Switch - Not Mounted	
	Blank plate provides mounting option for Micro Lite, Buddy Button or other switch options. Swivel assembly allows pad to be mounted with PERMOfix or directly to lateral arm of Total Control Head Array.	
ODTIONAL ITEMS:		
OPTIONAL ITEMS:		
Part Number	Description	Price AUD
<b>Part Number</b> ☐ I101769-99-0	Micro Lite Switch Qty:	Price AUD
Part Number ☐ 1101769-99-0 ☐ 1101767-99-0	Micro Lite Switch Qty: Buddy Button Qty:	Price AUD
Part Number ☐ I101769-99-0 ☐ I101767-99-0 ☐ I106293-99-0	Micro Lite Switch Qty: Buddy Button Qty: Egg Switch Black Qty:	Price AUD
Part Number  ☐ I101769-99-0 ☐ I101767-99-0 ☐ I106293-99-0 ☐ I101771-99-0	Micro Lite Switch Qty: Buddy Button Qty: Egg Switch Black Qty: Mini Cup Switch Qty:	Price AUD
Part Number  ☐ I101769-99-0 ☐ I101767-99-0 ☐ I106293-99-0 ☐ I101771-99-0 ☐ I107500-99-0	Micro Lite Switch Qty:  Buddy Button Qty:  Egg Switch Black Qty:  Mini Cup Switch Qty:  Proximity Switch Qty:	Price AUD
Part Number  ☐ I101769-99-0 ☐ I101767-99-0 ☐ I106293-99-0 ☐ I107771-99-0 ☐ I107500-99-0 ☐ I107495-99-0	Micro Lite Switch Qty: Buddy Button Qty: Egg Switch Black Qty: Mini Cup Switch Qty:	Price AUD
Part Number	Micro Lite Switch Qty:  Buddy Button Qty:  Egg Switch Black Qty:  Mini Cup Switch Qty:  Proximity Switch Qty:  Mini Joystick 5 Switch OC for Omni/IOM (9-pin)  Touch Contact with CR2032 Battery Qty:	Price AUD
Part Number	Micro Lite Switch Qty:  Buddy Button Qty:  Egg Switch Black Qty:  Mini Cup Switch Qty:  Proximity Switch Qty:  Mini Joystick 5 Switch OC for Omni/IOM (9-pin)  Touch Contact with CR2032 Battery Qty:  3.5mm Ext. Cable for Sensor/Switch, 0.3m Long	Price AUD
Part Number	Micro Lite Switch Qty:  Buddy Button Qty:  Egg Switch Black Qty:  Mini Cup Switch Qty:  Proximity Switch Qty:  Mini Joystick 5 Switch OC for Omni/IOM (9-pin)  Touch Contact with CR2032 Battery Qty:  3.5mm Ext. Cable for Sensor/Switch, 0.3m Long  3.5mm Ext. Cable for Sensor/Switch, 1m Long	Price AUD
Part Number	Micro Lite Switch Qty: Buddy Button Qty: Egg Switch Black Qty: Mini Cup Switch Qty: Proximity Switch Qty: Mini Joystick 5 Switch OC for Omni/IOM (9-pin) Touch Contact with CR2032 Battery Qty: 3.5mm Ext. Cable for Sensor/Switch, 0.3m Long 3.5mm Ext. Cable for Sensor/Switch, 1m Long 3.5mm Ext. Cable for Sensor/Switch, 2m Long	Price AUD
Part Number	Micro Lite Switch Qty: Buddy Button Qty: Egg Switch Black Qty: Mini Cup Switch Qty: Proximity Switch Qty: Mini Joystick 5 Switch OC for Omni/IOM (9-pin) Touch Contact with CR2032 Battery Qty: 3.5mm Ext. Cable for Sensor/Switch, 0.3m Long 3.5mm Ext. Cable for Sensor/Switch, 1m Long 3.5mm Ext. Cable for Sensor/Switch, 2m Long 3.5mm Ext. Cable for Sensor/Switch, 3m Long	Price AUD
Part Number	Micro Lite Switch Qty: Buddy Button Qty: Egg Switch Black Qty: Mini Cup Switch Qty: Proximity Switch Qty: Mini Joystick 5 Switch OC for Omni/IOM (9-pin) Touch Contact with CR2032 Battery Qty: 3.5mm Ext. Cable for Sensor/Switch, 0.3m Long 3.5mm Ext. Cable for Sensor/Switch, 1m Long 3.5mm Ext. Cable for Sensor/Switch, 2m Long	Price AUD

Additional Order Instructions:	
Chair Order Policy:	with ANY quantions regarding proper completion of our order forms. The order connect he processed
without the following:	with ANY questions regarding proper completion of our order forms. The order cannot be processed
	authorized by your company to acknowledge the items selected. Permobil is not
responsible for configuration or size discrepanci	es resulting from customer errors on the order form. I summarizing the items selected. Please review this carefully.
-Once the order is shipped, any changes are subj	
Order Acknowledgement:	
I,, am an agent of the	e medical equipment provider named on this order form and I have the authority to contract for the
accurate to the best of my knowledge. I further ackn	on behalf of said provider. I acknowledge that I have reviewed this order and that it is complete and owledge that any changes to this order after submission of this order form are subject to the returns policy
and may carry additional charges.	Signed:
Printe	ed Name:
	Title:
	Date:
NOTE: All specifications and prices are subject to cha parts ordered for service or future alterations will carry	nge without notice. Please note that prices displayed are only valid if ordered with the wheelchair. Any different pricing and warranty.