




AERO T

PURE
GEOMETRY.



Take in the curves and lines of TiLite's most advanced aluminum chair. The Aero T's dual-tube frame minimizes weight while maximizing performance. Add to that full adjustability and TiFit, and you've got something special. It's all in the ride.



AERO T

TECHNICAL SPECIFICATIONS.

Frame Features	Aluminum Dual-Tube Rigid TiFit Full Adjustability
Frame Options	Titanium Frame Upgrade SuperLite Edition Heavy Duty Frame Power Adaptable Frame Reinforced Frame
Transport Weight*	12 lb.
Weight Limit	265 lb. Standard 350 lb. Heavy Duty
Total Chair Width*	23"
Seat Width	12" - 20"
Front Seat Width (Seat Taper)	3 Options
Seat Depth	12" - 20"
Front Seat Height	15" - 21"
Rear Seat Height	13" - 21"
Camber	0°, 2°, 4°, 6°, 8°, or 12°
Front Frame Angle	70°, 80°, 85°, or 90°
Footrest Width Options	Standard Taper or V-Front End
Center of Gravity Adjustment	5"
Front Wheels	23 Options
Rear Wheels	67 Options
Handrims	15 Options
Convenience	Folding Angle Adjustable Seat Back
Aesthetics	5 Aluminum Frame Finishes 16 Frame Colors 6 Color Anodize Packages 7 Ultrasuede Colors 8 Spinergy Wheel Colors
HCPCS Code	K0005
Warranty	Lifetime on Frame One Year on TiLite Components

*Weights and widths based on 16" x 16" chair with 0° camber, minimum rear wheel spacing and all standard features. Rear wheels are not included in transport weight. For reference, a pair of standard 24" rear wheels weighs 7.2 lb.



TIFIT. PERSONAL ERGONOMICS.

The difference between a TiFit® chair and all other chairs is that a TiFit chair is a wheeled prosthetic: Every frame is made to measure. The optimized fit of a TiFit chair brings unparalleled performance and function. You will feel the difference in everything you do.

CUSTOM CONFIGURATION.

Custom Configuration allows you to add options and components to your chair to meet your needs and to optimize your riding experience. This includes choices of wheels, handrims, tires, wheel locks, color, tattoos, and more. TiLite offers more options than any other manufacturer!



TiLite | The Ultimate Ride™ | © 2015 TiLite, LLC.
800.545.2266 | tilite.com | MS0020_RevB_AeroT
TiLite wheelchairs are shown with available options.