



AERO Z

OPEN FRAME.
OPEN ROAD.



The open frame of the aluminum Aero Z calls out to the open road. The Aero Z is the only TiFit, mono-tube aluminum frame on the market. In a class by itself, the Aero Z is lightweight and fully adjustable and it won't break the bank!

AERO Z

TECHNICAL SPECIFICATIONS.

Frame Features	Aluminum Mono-Tube Rigid TiFit Full Adjustability
Frame Options	Titanium Frame Upgrade SuperLite Edition Power Adaptable Frame Reinforced Frame
Transport Weight*	12.4 lb.
Weight Limit	265 lb.
Total Chair Width*	23"
Seat Width	10" - 20"
Front Seat Width (Seat Taper)	3 Options
Seat Depth	10" - 20"
Front Seat Height	16" - 21"
Rear Seat Height	12.5" - 20"
Camber	0°, 2°, 4°, 6°, 8°, or 12°
Front Frame Angle	70°, 75°, 80°, 85°, or 90°
Footrest Width Options	Standard Taper or V-Front End
Center of Gravity Adjustment	6"
Front Wheels	23 Options
Rear Wheels	67 Options
Handrims	15 Options
Convenience	Transit Tie-Down Option Folding Angle Adjustable Seat Back Depth Adjustable Back/Frame
Aesthetics	5 Aluminum Frame Finishes 16 Frame Colors 6 Color Anodize Packages 7 Ultrasuede Colors 8 Spinergy Wheel Colors
HCPCS Code	K0005
Warranty	Lifetime on Frame One Year on TiLite Components

*Weights and widths based on 16" x 16" chair with 0° camber, minimum rear wheel spacing and all standard features. Rear wheels are not included in transport weight. For reference, a pair of standard 24" rear wheels weighs 7.2 lb.



TIFIT. PERSONAL ERGONOMICS.

The difference between a TiFit® chair and all other chairs is that a TiFit chair is a wheeled prosthetic: Every frame is made to measure. The optimized fit of a TiFit chair brings unparalleled performance and function. You will feel the difference in everything you do.

CUSTOM CONFIGURATION.

Custom Configuration allows you to add options and components to your chair to meet your needs and to optimize your riding experience. This includes choices of wheels, handrims, tires, wheel locks, color, tattoos, and more. TiLite offers more options than any other manufacturer!



TiLite | The Ultimate Ride™ | © 2015 TiLite, LLC.
800.545.2266 | tilite.com | MS0026_RevC_AeroZ
TiLite wheelchairs are shown with available options.