


C350 Corpus 3G

permobil



Dealer Information

* Contact: _____

Dealer Code: _____

* Dealer Name: _____

Address: _____

* City: _____

* Province: _____

* Phone#: _____

* Fax#: _____

PO#: _____

Email Address: _____

Client Information

Permobil recommends that the client is evaluated by a certified rehab specialist.

* First Name: _____

* Last Name: _____

Diagnosis: _____

Funding Source: _____

Client Age: _____

Client Measurements

* Weight: _____

* Height: _____

A) Top of Shoulders: _____

B) Chest Depth: _____

C) Chest Width: _____

D) Seat Depth: _____

E) Top of Head: _____

F) Elbow to Hand: _____

G) Seat Pan to Elbow: _____

H) Hip Width: _____

*I) Knee to Foot: _____

Display ADP Codes on order/quote:

Three column pricing on order/quote:

Please send order/quote to fax#: (800) 231-3256

Email to: sales@permobil.com

Permobil Ltd.
12 - 75 Mary St.
Aurora, Ontario
L4G 1G3
Tel: (800) 736-0925
Fax: (800) 231-3256
www.permobil.com

Prices effective September 1, 2017 and payable in Canadian dollars.

Chair Model and Colors

CHOOSE BASE:

<input type="checkbox"/>	I105206-99-0	C350 Base Corpus 3G R-net - C300 PS2 <i>Price Includes: RWD C350 Power Base, Shock Absorbing Suspension System, R-net Basic Controller Set Up with 5mph Motor Programming, Tie-Down, Hardware for Strap Systems, Anti-Tippers, Charger, Flat Free Tires, and Standard Bodypoint Positioning Belt. Must Choose Flat Free or Pneumatic Tires Below. 5mph, Weight Capacity = 300lbs, STFH = 17.5".</i>	8,069.00	WPA31P030
<input type="checkbox"/>	I10604	6.5mph High Speed Upgrade Package - C350 <i>Upgrade option for Group 3 C350 Power Base.</i>	1,663.00	
<input type="checkbox"/>	I10742	C350 Corpus 3G 6.5mph - C350 HS PS2 <i>Price Includes: RWD C350 Power Base, Shock Absorbing Suspension System, 6.5mph Motor Programming, Tie-Down, Hardware for Strap Systems, Anti-Tippers, Charger, Flat Free Tires, and Standard Bodypoint Positioning Belt. Must Choose Flat Free or Pneumatic Tires Below. 6.5mph, Weight Capacity = 300lbs, STFH = 17.5".</i>	9,613.00	

CHOOSE DRIVE TIRES:

Part Number	Description	MSRP - CA	ADP Code
<input type="checkbox"/>	I103499-99-0 Tri-Spoke Split Rim w/ Flat Free Tire 3x8" Grey	No Charge	
<input type="checkbox"/>	I103498-99-0 Tri-Spoke Split Rim w/ Pneumatic Tire 3x8" Grey	No Charge	

OPTIONAL DRIVE TIRES:

Part Number	Description	MSRP - CA	ADP Code
<input type="checkbox"/>	I10782 Tri-Spoke Split Rim w/ Snow Tire Pneumatic 3x8" <i>Snow tires are designed for use in cold weather conditions. The nylon compound in these tires remains flexible at low temperatures, but wears more easily. Therefore, the snow tires (sold as wheel assemblies) are not recommended for year-long use and must be ordered along with standard drive wheels. Not available with black tires or black rims.</i>	524.00	

CHOOSE SHROUD COLOR:

Part Number	Description	MSRP - CA	ADP Code
<input type="checkbox"/>	I103532-99-0 C350 Shroud Vintage Wine	No Charge	
<input type="checkbox"/>	I103533-99-0 C350 Shroud Crystal Blue	No Charge	
<input type="checkbox"/>	I103529-99-0 C350 Shroud Onyx Black	No Charge	
<input type="checkbox"/>	I103530-99-0 C350 Shroud Platinum Silver	No Charge	
<input type="checkbox"/>	I103531-99-0 C350 Shroud Extreme Yellow	No Charge	

CHOOSE HUBCAP COLOR:

Part Number	Description	MSRP - CA	ADP Code
	Match Hubcap Color to Shroud: _____ <i>Note: Hub caps will be matched to shroud color if no selection is made.</i>		
<input type="checkbox"/>	I10871 Solid Color Hub Cap - Pair <i>Platinum Silver: _____ Onyx Black: _____ Cobalt Blue: _____ Phantom Gray: _____ Pop Star Pink: _____ Volcano Red: _____ Purple Rain: _____ Titan Blue: _____ Sublime Green: _____ Crystal Blue: _____ Vintage Wine: _____ Extreme Yellow: _____ Chili Pepper Red: _____</i>	No Charge	

Base Options

MUST CHOOSE EITHER SEAT ELEVATOR OR FIXED SEAT TUBE:

Part Number	Description	MSRP - CA	ADP Code
<input type="checkbox"/>	I106753-99-0 Power Adjustable Seat Height - 8" Travel	3,008.00	
<input type="checkbox"/>	I105421-99-0 Seat Tube Fixed	No Charge	

MUST CHOOSE BATTERIES:

Part Number	Description	MSRP - CA	ADP Code
<input checked="" type="checkbox"/>	I107032-99-0 Batteries, Grp 34 (60Ah) Sealed Gel, x2, Installed	953.00	

OPTIONAL ITEMS:

Part Number	Description	MSRP - CA	ADP Code
<input type="checkbox"/>	I103496-99-0 Light Kit for C350	1,152.00	

Base Options

OPTIONAL ITEMS:

Part Number	Description	MSRP - CA	ADP Code
<input type="checkbox"/> I10789	Seat Stabilizer Option	73.00	
<input type="checkbox"/> I10715	Compact Laptop Ventilator Tray - Corpus Seat <i>Mounts to rear of Corpus backrest hinge. Accommodates Pulmonetic Systems / CareFusion LTV Series and similar sized Laptop Ventilators. The tray is 330mm Tall, adjustable Width between 200mm - 330mm, and adjustable Depth between 50mm - 165mm. Overall backrest angle limited to 150° when mounting this ventilator tray. Custom Option.</i>	1,682.00	
<input type="checkbox"/> I10787	Battery Tray for Laptop Ventilator Tray <i>Mounts to rear of Laptop Ventilator Tray. Measures 290mm W x 170mm H x 42mm D. Accommodates Pulmonetic Systems Lithium-Ion Battery Pack, CareFusion SprintPack Battery System, or Richardson PowerTech Vent Power Center. Custom Option.</i>	269.00	
<input type="checkbox"/> I10788	O2 Holder for Laptop Ventilator Tray <i>Mounts to Left or Right Side of Laptop Vent Tray. O2 Holder is 175mm Tall. Width infinitely adjustable 90mm - 170mm. Depth infinitely adjustable 95mm - 125mm. Custom Option.</i>	380.00	WEPX
<input type="checkbox"/> ID51113	Programming Key for R-net - Permobil Version	524.00	

Joystick Options

CHOOSE JOYSTICK TYPE:

Part Number	Description	MSRP - CA	ADP Code
	PLEASE NOTE: <i>For all remote joystick modules shown below, you must also select an expandable controller. For additional input devices, refer to Alternative Drive Systems Section.</i>		
<input type="checkbox"/> I105246-99-0	R-net LED Remote Joystick Module <i>This Joystick does NOT have a Display Screen and does NOT contain mono jacks. Limited to five profiles with this Joystick Module.</i>	No Charge	
<input type="checkbox"/> I107492-99-0	R-net Remote Color JS w/ Mono Jacks, No Lights	1,803.00	
<input type="checkbox"/> I107491-99-0	R-net Remote Color JS w/ Toggle Switches, No Lghts	2,042.00	WEX1
<input type="checkbox"/> I104923-99-0	R-net Remote Color JS w/ Mono Jacks, Lights	1,803.00	
<input type="checkbox"/> I107493-99-0	R-net Remote Color JS w/ Toggle Switches, Lights	1,990.00	WEX1
<input type="checkbox"/> I105700-99-0	R-net Remote Colour JS w/ Lights, Heavy Duty Gimbal <i>Heavy-Duty modification enables gimbal to hold up to extreme operating forces. Reinforced housing is attached to a standard Colour LCD R-net Joystick. 3" long joystick knob has 3/4" diameter. Not compatible with alternative joystick handles.</i>	2,042.00	WEX1
<input type="checkbox"/> I103814-99-0	Compact Advanced Joystick for R-net <i>Compact sized joystick that does not contain an LCD Display. Contains two mono jacks for mode & on/off. Does NOT contain a wrist support pad. If you want a wrist support pad, please select I10621 below.</i>	1,990.00	WEX1
<input type="checkbox"/> I10621	Add Wrist Pad to Compact Joystick Advanced	253.00	

CHOOSE JOYSTICK HANDLE TYPE:

Part Number	Description	MSRP - CA	ADP Code
<input type="checkbox"/> I10045	Cone Shaped Joystick Knob - Standard	No Charge	
<input type="checkbox"/> I100960-99-0	Chin Cup For Joystick Knob	119.00	
<input type="checkbox"/> I100961-99-0	"T" Handle For Joystick Knob	176.00	
<input type="checkbox"/> I100962-99-0	"Large Ball" For Joystick Knob	119.00	
<input type="checkbox"/> I100963-99-0	"Softball" For Joystick Knob	160.00	
<input type="checkbox"/> I100964-99-0	"Mushroom" For Joystick Knob	160.00	
<input type="checkbox"/> I100965-99-0	"Stick" For Joystick Knob	160.00	
<input type="checkbox"/> I105442-99-0	1.75" Black Foam Ball For Joystick Knob	119.00	
<input type="checkbox"/> I105846-99-0	Bodypoint J/S Handle 3" U-Shaped w/ Flex-Shaft	153.00	
<input type="checkbox"/> I105847-99-0	Bodypoint J/S Handle 4" U-Shaped w/ Flex-Shaft	153.00	
<input type="checkbox"/> I105848-99-0	Bodypoint J/S Handle Dome Shaped, Rubber	153.00	

Part Number	Description	MSRP - CA	ADP Code
<input type="checkbox"/> I10664	Joystick Mount Right - Fixed	No Charge	

Joystick Options

Part Number	Description	MSRP - CA	ADP Code
<input type="checkbox"/> I10665	Joystick Mount Left - Fixed	No Charge	
<input type="checkbox"/> I103820-99-0	Retractable Joystick Mount Right	452.00	WEPM
<input type="checkbox"/> I103819-99-0	Retractable Joystick Mount Left	452.00	WEPM
<input type="checkbox"/> I105210-99-0	Mount Standard Joystick on Corpus Back Hinge <i>Standard R-net Joystick is mounted to backrest hinge. Can be chosen instead of PN I105408-99-0. Mount on Left: _____ Mount on Right: _____</i>	314.00	
<input type="checkbox"/> I101260-99-0	Midline Control Mount <i>Mount Midline on Left Armrest: _____ Mount Midline on Right Armrest: _____</i>	392.00	WE3D
<input type="checkbox"/> I106266-99-0	Stealth Gatlin Midline Mount - 8 Position <i>Features eight locking positions each 45° apart. Adjustments are made by pulling out on mounting arm and rotating to desired angle. Joystick bracket sold separately (select from the brackets listed below). Please indicate if you want it mounted on left or right armrest: Mount Midline on Left Armrest: _____ Mount Midline on Right Armrest: _____</i>	690.00	WE3D
<input type="checkbox"/> I106265-99-0	Stealth Gatlin Midline Mount - 2 Position <i>Designed for easy swing away. Mount can be set to swing from horizontal to 95° or from horizontal to 180°. Joystick bracket sold separately (select from the brackets listed below). Please indicate if you want it mounted on left or right armrest: Mount Midline on Left Armrest: _____ Mount Midline on Right Armrest: _____</i>	690.00	WE3D
<input type="checkbox"/> I106267-99-0	Stealth Gatlin Bracket for R-net/Compact JS Adv	283.00	
<input type="checkbox"/> I106295-99-0	Stealth Gatlin Bracket for Compact Joystick	283.00	
<input type="checkbox"/> I106296-99-0	Stealth Gatlin Bracket for Mini Joystick	283.00	
<input type="checkbox"/> I106297-99-0	Stealth Gatlin Swivel Bracket for Omni	322.00	
<input type="checkbox"/> I10842	PermoFix Joint for Gatlin (0027-7027a) <i>Includes a shim that allows the PERMOFix Joint to be mounted on the 1/2" Stealth Gatlin arm. Provides an adjustable mount for a Mini Joystick or Compact Joystick.</i>	283.00	
<input type="checkbox"/> I106268-99-0	Stealth Height Adjustable Mount <i>Replaces standard UT panel bracket. Provides 1.75" of height adjustability to all Permobil joystick mounts and Stealth Gatlin Midline mounts.</i>	252.00	
Part Number	Description	MSRP - CA	ADP Code
<input checked="" type="checkbox"/> I10411	Expandable Controller - R-net <i>Must be ordered with all expandable/remote joysticks.</i>	No Charge	
<input checked="" type="checkbox"/> I10513	Harness for Expandable Controller <i>Must be ordered with any Expandable Controller.</i>	No Charge	
<input type="checkbox"/> I105225-99-0	Co-Pilot Attendant Control R-net - Corpus	5,034.00	
<input type="checkbox"/> I105408-99-0	Attendant Control R-net	1,307.00	
<input type="checkbox"/> I107350-99-0	R-net Bluetooth iDevice Module <i>Allows input device commands (such as joystick nudges or switch clicks) to be used to operate an iPhone, iPad or iPod. By assigning each input to an iOS function, users can make phone calls, send texts and emails and access applications. Kit includes a Bluetooth iDevice Module, an R-net Bus Cable, and two mono jack cables for direct switch inputs (switches not included). Must be ordered with a Color JSM or Omni for mode display.</i>	1,055.00	
<input type="checkbox"/> I107367-99-0	R-net Bluetooth Mouse Module <i>Enables mouse control through any R-net input device. It can operate the mouse on a computer and on select Android devices (with R-net for Android app). Kit includes a Bluetooth Mouse Module, an R-net Bus Cable, two mono jack cables for L/R click (switches not included) and a USB BT Adapter for user's computer. Note: BT Adapter may not work with all computers. Must be ordered with a Color JSM or Omni for mode display.</i>	1,055.00	
<input type="checkbox"/> I107368-99-0	Add a second Bluetooth Mouse Module <i>Allows R-net to simultaneously pair to a second compatible computer or Android device. Each module is accessible through a separate R-net mode. Kit includes a unique Bluetooth Mouse Module, an R-net Bus Cable, two mono jack cables for L/R click and a USB BT Adapter.</i>	1,055.00	

Joystick Options

Part Number	Description	MSRP - CA	ADP Code
<input type="checkbox"/> I10582	Output Module (IOM) for R-net Control System <i>This allows a user to control an external device (ex: Communication Device, ECU, Door Opener, etc.) through the drive input (JSM, Head Array, etc.). Each module can provide up to 4 switch closures. May require specialty cable or wiring. Up to three output modules can be installed on one wheelchair.</i>	1,749.00	

Seating System

CHOOSE SEAT MOUNTING PLATE TYPE:

Part Number	Description	MSRP - CA	ADP Code
<input type="checkbox"/> I105055-99-0	Standard (0°) Tilt Mounting Plate <i>Range of Tilt will be 0° to +50°.</i>	No Charge	
<input type="checkbox"/> I105068-99-0	5° Anterior Tilt Mounting Plate <i>Range of Tilt will be -5° to +45°. This option will raise the minimum seat to floor height by about 1" when the seat is at 0° of tilt.</i>	344.00	

CHOOSE SEATING SYSTEM:

Part Number	Description	MSRP - CA	ADP Code
PLEASE NOTE: <i>Weight Limit 300lbs. For the Corpus 3G Seat on the C350, there is a half speed limitation when the Total Backrest Angle (TBA) > 120°. TBA = Tilt Angle + Backrest Recline Angle. Example 1: 15° Tilt + 100° Recline = 115° TBA (Full Speed) Example 2: 35° Tilt + 110° Recline = 145° TBA (Half Speed)</i>			
<input type="checkbox"/> I11116	C3G Seat w/ Manual Recline <i>Manual Recline angle limited to 120°.</i>	1,293.00	
<input type="checkbox"/> I10744	C3G Seat w/ No Tilt & 150° Power Recline	8,724.00	
<input type="checkbox"/> I10649	C3G Seat w/ 50° DuraCore Tilt & Manual Recline <i>Manual Recline angle limited to 120°.</i>	7,824.00	
<input type="checkbox"/> I10745	C3G Seat w/ 50° DuraCore Tilt & 150° Power Recline	11,932.00	

Seating System Options

CHOOSE SEAT FUNCTION CONTROL TYPE:

Part Number	Description	MSRP - CA	ADP Code
<input type="checkbox"/> I10418	Single Seat Function Control Switch for R-net <i>Please choose either Push Buttons or Toggle Switches below.</i>	No Charge	
<input type="checkbox"/> I10419	Single Seat Function Control Kit for R-net <i>For operation of a single seat function through the wheelchair electronics via the drive control and/or ICS Switchbox. Order with R-net Controller or remote joystick options. Please indicate below if you want push buttons or toggles on the ICS Switchbox.</i>	1,661.00	
<input type="checkbox"/> I10585	Multiple Seat Function Control Kit for R-net <i>MUST be selected if ordering two or more seat functions. This device is used to operate two or more seat functions through the wheelchair electronics via the drive control and/or ICS switch box. Order with R-net electronics or remote joystick options. Please indicate below if you want push buttons or toggles on the ICS Switchbox.</i>	2,768.00	SEC1P1005

CHOOSE PUSH BUTTONS OR TOGGLES:

Part Number	Description	MSRP - CA	ADP Code
<input type="checkbox"/> I102798-99-0	Push Buttons on Seat Function Control Kit	No Charge	
<input type="checkbox"/> I102799-99-0	Toggle Switches on Seat Function Control Kit	296.00	
<input type="checkbox"/> I105409-99-0	Do Not Install ICS Switchbox, Send as Part	No Charge	

CHOOSE MOUNTING POSITION OF ICS SWITCHBOX:

Part Number	Description	MSRP - CA	ADP Code
<input type="checkbox"/> I10626	Mount ICS Switchbox onto R-net Joystick Bracket <i>ICS Switchbox will be mounted to the same bracket as the R-net Joystick Module. If this item is selected, the ICS Switchbox will swing away with the Joystick Module.</i>	No Charge	
<input type="checkbox"/> I104651-99-0	Mount ICS Switchbox onto Armrest <i>Mount on Left: _____ Mount on Right: _____</i>	No Charge	

Seating System Options

CHOOSE MOUNTING POSITION OF ICS SWITCHBOX:

Part Number	Description	MSRP - CA	ADP Code
<input type="checkbox"/> I106256-99-0	Mount ICS Switchbox on Backrest Mount on Left: _____ Mount on Right: _____	166.00	

OPTIONAL ITEMS:

Part Number	Description	MSRP - CA	ADP Code
<input type="checkbox"/> I107705-99-0	Virtual Seating Coach Package <i>The Virtual Seating Coach (VSC) Package includes a Bluetooth dongle, a USB charger, a RAM X-grip Phone Holder and access to VSC coaching and regimen management. The VSC app can be downloaded for free from the App Store or Google Play. Requires an iPhone (iOS 9 or higher) or an Android phone (OS 4. 1. 2 or higher).</i> Mount VSC to Left Armrest: _____ <i>To provide access, USB Charger and Phone Holder may require relocation.</i> Mount VSC to Right Armrest: _____ <i>To provide access, USB Charger and Phone Holder may require relocation.</i>	2,261.00	
<input type="checkbox"/> I104196-99-0	ICS Alternative Switchbox for R-net Electronics <i>Seat function interface box that accommodates up to 8 momentary switches. Provides an additional method for accessing the seat functions. 1/8" mini-mono switches not included.</i>	1,307.00	
<input type="checkbox"/> I106743-99-0	Remote Stop System - All Electronics <i>Consists of a Radio Frequency (RF) handheld transmitter and an RF receiver mounted on the wheelchair. Custom charges may be required when ordering with a PS/PS Jr seat, a low backrest height or a pediatric wheelchair. If ordering on a K450, push handles must also be selected.</i>	1,247.00	
<input type="checkbox"/> I10750	Memory Seat Program, Corpus w/ R-net	121.00	
<input type="checkbox"/> I106209-99-0	Legrest Movement Linked to Recline Electronically <i>Power Legrest will raise when Power Recline moves backward and will lower when back moves forward. If Power Legrest is operated, Power Recline will NOT automatically reposition. Feature works through ICS Switchbox and Driver Control.</i>	49.00	
<input type="checkbox"/> I10898	Enable Independent Repositioning Mode <i>IRM allows for sequential movement of power seat functions into a posterior position (35° tilt, 115° recline, 135° legrest) and returns to upright position (5° tilt, 100° recline, 110° legrest). This can be accessed using a standard drive control or a single switch. On the Permobil Joystick Module, IRM is accessible through the mode/profile jack (long switch hold). Please see the Alt Drive Accessories section for switch options.</i>	No Charge	

Backrest Options

CHOOSE FIXED OR SLIDING BACKREST KIT:

Part Number	Description	MSRP - CA	ADP Code
<input type="checkbox"/> I105067-99-0	Sliding Backrest Kit - 5" Mechanical Shear	No Charge	
<input type="checkbox"/> I105066-99-0	Fixed Backrest Kit - No Mechanical Shear	53.00	

CHOOSE BACKREST TYPE:

Part Number	Description	MSRP - CA	ADP Code
BACKREST HEIGHT NOTE: <i>Backrest heights are measured from the seat pan to the top of the backrest shell with the backrest set at 90°. Ergo Backs include lateral wedges and lumbar pads.</i>			
<input type="checkbox"/> I105113-99-0	Corpus 3G Ergo Backrest - 14"W x 20"T (C3GERGBCA) <i>Power Recline will be limited to 150° when ordering this backrest.</i>	1,072.00	
<input type="checkbox"/> I105114-99-0	Corpus 3G Ergo Backrest - 16"W x 20"T (C3GERGBCA) <i>Power Recline will be limited to 150° when ordering this backrest.</i>	1,072.00	
<input type="checkbox"/> I105109-99-0	Corpus 3G Ergo Backrest - 18"W x 20"T (C3GERGBCA) <i>Power Recline will be limited to 150° when ordering this backrest.</i>	1,072.00	
<input type="checkbox"/> I105852-99-0	Corpus 3G Ergo Backrest - 20"W x 20"T (C3GERGBCA) <i>Power Recline will be limited to 150° when ordering this backrest.</i>	1,072.00	

PLEASE NOTE:

Backrest Options

CHOOSE BACKREST TYPE:

Part Number	Description	MSRP - CA	ADP Code
PLEASE NOTE:			
<i>All taller backrests listed below include: lateral wedges, lumbar pads, a UT rail mounted at the top for lateral adjustment of the headrest, and split backrest feature with height adjustment in 1" increments. Please indicate below where you want the height set at so that a properly sized backrest cushion can be provided.</i>			
<input type="checkbox"/> I105110-99-0	Corpus 3G Ergo Back -14"W x 23 - 28"T (C3GERGBCA)	1,072.00	
<input type="checkbox"/> I105111-99-0	Corpus 3G Ergo Back -16"W x 23 - 28"T (C3GERGBCA)	1,072.00	
<input type="checkbox"/> I105112-99-0	Corpus 3G Ergo Back -18"W x 23 - 28"T (C3GERGBCA)	1,072.00	
<input type="checkbox"/> I105853-99-0	Corpus 3G Ergo Back -20"W x 23 - 28"T (C3GERGBCA)	1,072.00	
	Set Backrest at 23" _____		
	Set Backrest at 24" _____		
	Set Backrest at 25" _____		
	Set Backrest at 26" _____		
	Set Backrest at 27" _____		
	Set Backrest at 28" _____		

CHOOSE BACKREST CUSHION:

Part Number	Description	MSRP - CA	ADP Code
<input type="checkbox"/> I10659	Leatherette Cushion - Corpus Ergo Back	No Charge	
<input type="checkbox"/> I10660	Black Mesh Cushion - Corpus Ergo Back	55.00	
<input type="checkbox"/> I105117-99-0	Omit Backrest Cushion	No Charge	

OPTIONAL ITEMS:

Part Number	Description	MSRP - CA	ADP Code
<input type="checkbox"/> I105782-99-0	Tubular Backrest Frame - Corpus (1" Canes) <i>Provides 1" diameter, vertical canes for mounting aftermarket backrests. Frame can be adjusted to support backrest widths from 14" to 20". Canes are attached with 2" long UniTrack linkage disks that can rotate to provide space for backrest hardware. The tops of the canes can be adjusted up to a height of 27" above the seat pan. Depending on the aftermarket backrest, several inches of seat depth may be lost when mounting this system.</i>	477.00	
<input type="checkbox"/> I106072-99-0	Direct Backrest Frame - Corpus <i>Used to mount ROHO Agility, Jay 3, Matrx Elite, Varilite Icon and Ride Java backrests ranging in width from 14" to 20". Slotted brackets and 1" UniTrack linkage disks provide a close, adjustable mount. Backrest pitch can be tuned by rotating the linkage disks. Some seat depth may be lost due to the position and orientation of the backrest.</i>	477.00	
<input type="checkbox"/> I10725	Split Backrest Cushion - Corpus <i>Custom Option.</i>	355.00	
<input type="checkbox"/> I10786	Tailored Backrest Cushion - Corpus <i>Sides of cushion are trimmed to provide additional clearance for lateral trunk supports. Must also select a backrest cushion when ordering this custom option.</i>	111.00	

Seat Sizing

CHOOSE SEAT WIDTH:

Part Number	Description	MSRP - CA	ADP Code
<input type="checkbox"/> I107023-99-0	Set Seat Width at 17": _____	No Charge	
<input type="checkbox"/> I107024-99-0	Set Seat Width at 19": _____	No Charge	
<input type="checkbox"/> I107025-99-0	Set Seat Width at 21": _____	No Charge	
<input type="checkbox"/> I107026-99-0	Set Seat Width at 23": _____	No Charge	

CHOOSE SEAT DEPTH:

Part Number	Description	MSRP - CA	ADP Code
PLEASE NOTE:			
<i>There are three seat depth kits, each providing three seat depth settings: (1) 14" - 16", (2) 17" - 19" and (3) 20" - 22". Please select the required seat depth and the appropriate kit will be provided.</i>			
<input type="checkbox"/> I106232-99-0	Seat Depth 14"	No Charge	

Seat Sizing

CHOOSE SEAT DEPTH:

Part Number	Description	MSRP - CA	ADP Code
<input type="checkbox"/> I106233-99-0	Seat Depth 15"	No Charge	
<input type="checkbox"/> I106239-99-0	Seat Depth 16"	No Charge	
<input type="checkbox"/> I106235-99-0	Seat Depth 17"	No Charge	
<input type="checkbox"/> I106237-99-0	Seat Depth 18"	No Charge	
<input type="checkbox"/> I106241-99-0	Seat Depth 19"	No Charge	
<input type="checkbox"/> I106242-99-0	Seat Depth 20"	No Charge	
<input type="checkbox"/> I106243-99-0	Seat Depth 21"	No Charge	
<input type="checkbox"/> I106244-99-0	Seat Depth 22"	No Charge	

CHOOSE SEAT CUSHION:

Part Number	Description	MSRP - CA	ADP Code
NOTE: A Corpus II Ergonomic Seat Cushion will be provided for select seat sizes.			
<input type="checkbox"/> I10651	Corpus Ergo Cushion - Leatherette (C3GERGSCA)	605.00	SES1P0005
<input type="checkbox"/> I10652	Corpus Ergo Cushion - Black Mesh (C3GERGSCA)	605.00	SES1P0005
<input type="checkbox"/> I10965	Corpus Ergo Air Cushion (C3GRHOSCA)	642.00	SE0001267
<i>Features ROHO Dry Floatation technology. The removable ROHO cell pad provides adjustable skin protection while the contoured foam and thigh support wedges help maintain proper positioning and stability. Includes removable stretch fabric cover. Available in widths of 17", 19" and 21" and depths of 17", 19" and 21", as well as 17" x 23".</i>			
<input type="checkbox"/> I11113	Incontinence Cover for Ergo Cushion	192.00	SEMCF106L
<i>Made from a stretch Dartex fabric and suitable for all Ergo cushions including the Corpus Ergo Air. The cover is cut to the contours of the Ergo cushion and is designed to be used as an outer layer of protection. Cover size will be matched to cushion size. Sent as a part. Part Number will be provided on the quote / order confirmation based on the Seat Width and Seat Depth selected.</i>			
<input type="checkbox"/> I105118-99-0	Omit Seat Cushion	No Charge	

Armrest Options

CHOOSE ARMBAR KIT AND MOUNTING POSITION:

Part Number	Description	MSRP - CA	ADP Code
PLEASE NOTE: Narrow Armbar Kit can be adjusted from 15"W to 19"W between the armrest pads without any additional parts. Wide Armbar Kit can be adjusted from 19"W to 23"W between the armrest pads without any additional parts. Please select both Armbar Kit and Width Setting below.			
<input type="checkbox"/> I106543-99-0	Narrow Armbar Kit, Inf Adjustable 15"W - 19"W	No Charge	
<input type="checkbox"/> I106544-99-0	Wide Armrest Kit, Inf Adjustable 19"W - 23"W	No Charge	
<i>Required when selecting a 20" wide backrest. Set the Width Between the Armrests to: _____ Must be wider than backrest width to ensure proper fit. Please round to nearest 1".</i>			

CHOOSE ARMREST HEIGHT MOUNTING POSITION:

Part Number	Description	MSRP - CA	ADP Code
<input type="checkbox"/> I105069-99-0	Height Adjustable Armrests - NC	No Charge	
SET ARMREST HEIGHT AT: _____ <i>Corpus Armrests are infinitely adjustable from 7" - 15", measured from the seat pan to the top of the armrest pad. No additional parts are needed to get this full range of adjustment. In one inch increments (ex: 8", 9", etc.), please indicate where you would like the armrests set.</i>			

CHOOSE ARMREST SWIVEL LOCKED OR ENABLED:

Part Number	Description	MSRP - CA	ADP Code
<input type="checkbox"/> I10645	Left Armrest Bracket - Fixed	No Charge	
<i>Medial Swivel not available with this option. Required when ordering arm trough adapters.</i>			
<input type="checkbox"/> I10646	Right Armrest Bracket - Fixed	No Charge	
<i>Medial Swivel not available with this option. Required when ordering arm trough adapters.</i>			

Armrest Options

CHOOSE ARMREST SWIVEL LOCKED OR ENABLED:

Part Number	Description	MSRP - CA	ADP Code
<input type="checkbox"/> I105073-99-0	Left Armrest Bracket - Medial Swivel Enabled	No Charge	
<input type="checkbox"/> I105072-99-0	Right Armrest Bracket - Medial Swivel Enabled	No Charge	

CHOOSE ARMREST PAD RIGHT:

Part Number	Description	MSRP - CA	ADP Code
<input type="checkbox"/> I105100-99-0	4"x10" Corpus Armrest Pad Right - Leatherette	No Charge	
<input type="checkbox"/> I105104-99-0	4"x13" Corpus Armrest Pad Right - Leatherette	No Charge	
<input type="checkbox"/> I105105-99-0	4"x16" Corpus Armrest Pad Right - Leatherette	No Charge	
<input type="checkbox"/> I105999-99-0	4"x18" Corpus HD Armrest Pad Right - Leatherette	No Charge	
<input type="checkbox"/> I107212-99-0	4"x10" Corpus Armrest Pad Right - Black Fabric	No Charge	
<input type="checkbox"/> I106047-99-0	4"x13" Corpus Armrest Pad Right - Black Fabric	No Charge	
<input type="checkbox"/> I107214-99-0	4"x16" Corpus Armrest Pad Right - Black Fabric	No Charge	
<input type="checkbox"/> I105995-99-0	4"x18" Corpus HD Armrest Pad Right - Black Fabric	No Charge	
<input type="checkbox"/> I106486-99-0	4"x10" Corpus Arm w/ Stealth Gel Pad - Right	191.00	
<input type="checkbox"/> I106488-99-0	4"x13" Corpus Arm w/ Stealth Gel Pad - Right	197.00	
<input type="checkbox"/> I106490-99-0	4"x16" Corpus Arm w/ Stealth Gel Pad - Right	202.00	
<input type="checkbox"/> I106281-99-0	4"x18" Corpus HD Arm w/ Stealth Gel Pad - Right	210.00	

CHOOSE ARMREST PAD LEFT:

Part Number	Description	MSRP - CA	ADP Code
<input type="checkbox"/> I105098-99-0	4"x10" Corpus Armrest Pad Left - Leatherette	No Charge	
<input type="checkbox"/> I105102-99-0	4"x13" Corpus Armrest Pad Left - Leatherette	No Charge	
<input type="checkbox"/> I105107-99-0	4"x16" Corpus Armrest Pad Left - Leatherette	No Charge	
<input type="checkbox"/> I105998-99-0	4"x18" Corpus HD Armrest Pad Left - Leatherette	No Charge	
<input type="checkbox"/> I107213-99-0	4"x10" Corpus Armrest Pad Left - Black Fabric	No Charge	
<input type="checkbox"/> I106046-99-0	4"x13" Corpus Armrest Pad Left - Black Fabric	No Charge	
<input type="checkbox"/> I107215-99-0	4"x16" Corpus Armrest Pad Left - Black Fabric	No Charge	
<input type="checkbox"/> I105994-99-0	4"x18" Corpus HD Armrest Pad Left - Black Fabric	No Charge	
<input type="checkbox"/> I106487-99-0	4"x10" Corpus Arm w/ Stealth Gel Pad - Left	191.00	
<input type="checkbox"/> I106489-99-0	4"x13" Corpus Arm w/ Stealth Gel Pad - Left	197.00	
<input type="checkbox"/> I106491-99-0	4"x16" Corpus Arm w/ Stealth Gel Pad - Left	202.00	
<input type="checkbox"/> I106303-99-0	4"x18" Corpus HD Arm w/ Stealth Gel Pad - Left	210.00	

OPTIONAL ITEMS:

Part Number	Description	MSRP - CA	ADP Code
<input type="checkbox"/> I105463-99-0	Corpus Arm Trough Adapter 1" - Pair <i>For aftermarket armrest pads. Please indicate below if they should be installed or on chair sent as parts. Please Note: Medial Swivel must be must be locked when selecting this option.</i>	347.00	
<input type="checkbox"/> I105462-99-0	Corpus Arm Trough Adpater 7/8" - Pair <i>For aftermarket armrest pads. Please indicate below if they should be installed or on chair sent as parts. Please Note: Medial Swivel must be must be locked when selecting this option.</i>	347.00	
	Install Arm Trough Adapters: _____		
	Send Arm Trough Adapters as Parts: _____		
<input type="checkbox"/> I107339-99-0	Armrest Pouch - Right	252.00	
<input type="checkbox"/> I107338-99-0	Armrest Pouch - Left	252.00	
<input type="checkbox"/> I107840-99-0	UT Armrest Side Support - Right Outer <i>Includes 2"T x 6.5"W tapered pad.</i>	361.00	
<input type="checkbox"/> I107863-99-0	UT Armrest Side Support - Left Outer <i>Includes 2"T x 6.5"W tapered pad.</i>	361.00	
<input type="checkbox"/> I107864-99-0	UT Armrest Side Support - Right Inner <i>Includes 2"T x 6.5"W tapered pad. May need to increase width between the armrests to ensure armrests will flip-up.</i>	361.00	

Armrest Options

OPTIONAL ITEMS:

Part Number	Description	MSRP - CA	ADP Code
<input type="checkbox"/> I107865-99-0	UT Armrest Side Support - Left Inner <i>Includes 2"T x 6.5"W tapered pad. May need to increase width between the armrests to ensure armrests will flip-up.</i>	361.00	
<input type="checkbox"/> I107839-99-0	UT Elbow Support - Right <i>Must also select PN I10881.</i>	194.00	
<input type="checkbox"/> I107838-99-0	UT Elbow Support - Left <i>Must also select PN I10882.</i>	194.00	
<input type="checkbox"/> I10881	Adjustable Elbow Support Hardware - Right	283.00	
<input type="checkbox"/> I10882	Adjustable Elbow Support Hardware - Left	283.00	

Legrest Options

CHOOSE LEGREST TYPE:

Part Number	Description	MSRP - CA	ADP Code
<input type="checkbox"/> I107764-99-0	Corpus Manual Elevating Legrest <i>No charge when ordering power tilt and/or power recline.</i>	No Charge	
<input type="checkbox"/> I10948	Corpus Manual Elevating Legrest	797.00	
<input type="checkbox"/> I107762-99-0	Corpus Power Elevating Legrest - 170° <i>Includes mechanical articulation.</i>	3,549.00	

PLEASE NOTE

If choosing the Power Extending Legrest or the Power Transfer Legrest, footplate height can be adjusted from 12"-16". For the Power Transfer Legrest, there is no pre-set position within the range. Footplate height is determined solely by user adjustment.

<input type="checkbox"/> I10889	Power Extending Legrest <i>Provides 4" of independent legrest extension. Restricts ability to adjust each footplate independently. Compatible with Manual and Power Elevating Legrests.</i>	3,856.00	
<input type="checkbox"/> I105078-99-0	Power Transfer Legrest <i>A single input command extends the footplate to the ground and then raises the seat for easier transfer. Footplate can be extended up to 27" from seat pan to footplate. Compatible with Manual and Power Elevating Legrests. Seat elevator is required for full functionality. Please Note: Calf support pads are only available as a custom option. Customization may also be required when ordered with a 5° anterior tilt mounting plate.</i>	4,242.00	

CHOOSE FOOTPLATE STYLE:

Part Number	Description	MSRP - CA	ADP Code
<input type="checkbox"/> I10711	2pc Wide Corpus Footplates - 9"D x 7.5"W <i>Not compatible with Power Transfer Legrest. Note: Depending on legrest angle, the footplates can come in contact with the caster wheels.</i>	379.00	
<input type="checkbox"/> I10712	2pc Narrow Corpus Footplates - 9"D x 6"W <i>Not compatible with Power Transfer Legrest.</i>	No Charge	
<input type="checkbox"/> I106760-99-0	1pc Narrow Corpus Footpl - 9.75"D x 12.5"W <i>Footplate is 12.5"W in front, 8"W in back.</i>	379.00	
<input type="checkbox"/> I105793-99-0	1pc Narrow Power Transfer Footplate - 10"D x 13"W <i>Footplate is 13"W in front, 8"W in back.</i>	No Charge	

FOOTPLATE HEIGHT POSITION:

Part Number	Description	MSRP - CA	ADP Code
	SET FOOTPLATE HEIGHT AT: _____ <i>This is measured from the SEAT PAN to the top of the FOOTPLATE, when the legrest assembly is at 90°. DO NOT INCLUDE CUSHION THICKNESS IN YOUR MEASUREMENT. In one inch increments, please indicate desired footplate height. NOTE: Corpus footplate height range is 6"-17" and Corpus VS height range is 9"-17". Height may need to be reduced for shipping. Calf Pads are not available at heights less than 9".</i>		

OPTIONAL ITEMS:

Part Number	Description	MSRP - CA	ADP Code
-------------	-------------	-----------	----------

Legrest Options

OPTIONAL ITEMS:

Part Number	Description	MSRP - CA	ADP Code
<input type="checkbox"/> I105252-99-0	UT Calf Support Kit (B) w Mesh Cover & Adj Hardware 3.5"W x 4.5"T Pad.	No Charge	
<input type="checkbox"/> I105258-99-0	UT Calf Support Kit (J) w Mesh Cover & Adj Hardware 5"W x 6.5"T Pad.	No Charge	
<input type="checkbox"/> I105250-99-0	UT Calf Support Kit (A) w Mesh Cover & Adj Hardware 5.5"W x 6"T Pad.	No Charge	
<input type="checkbox"/> I105256-99-0	UT Calf Support Kit (H) w Mesh Cover & Adj Hardware 7"W x 7"T Curved Pad.	392.00	
<input type="checkbox"/> I105254-99-0	UT Calf Support Kit (D) w Mesh Cover & Adj Hardware 8"W x 8"T Pad.	392.00	
<input type="checkbox"/> I106530-99-0	Stealth Gel Pads - Narrow Footplates (Pair)	399.00	
<input type="checkbox"/> I106225-99-0	Stealth Gel Pads - Wide Footplates (Pair)	411.00	

Accessories

<input type="checkbox"/> I104518-99-0	Bodypoint Non-Padded Hip Belt - Medium	No Charge	
<input type="checkbox"/> I104520-99-0	Bodypoint Padded Hip Belt - Medium	144.00	SERND0030
<input type="checkbox"/> I104514-99-0	Bodypoint Non-Padded Hip Belt - Large	116.00	SERND0030
<input type="checkbox"/> I104516-99-0	Bodypoint Padded Hip Belt - Large	197.00	SERND0030
<input type="checkbox"/> I100519-99-0	Positioning Belt Retractable	298.00	SERND0030
<input type="checkbox"/> I100600-99-0	Positioning Belt Retractable Long	323.00	SERND0030

CHOOSE ONE OF THE FOLLOWING:

Part Number	Description	MSRP - CA	ADP Code
<input type="checkbox"/> I107470-99-0	Permobil 7"W Headrest w/ Link Hardware <i>Must also select I107480-99-0 or I107481-99-0.</i>	439.00	SEH1P0005
<input type="checkbox"/> I107469-99-0	Permobil 10"W Headrest w/ Link Hardware <i>Must also select I107480-99-0 or I107481-99-0.</i>	439.00	SEH1P0005
<input type="checkbox"/> I107468-99-0	Permobil 14"W Headrest w/ Link Hardware <i>Must also select I107480-99-0 or I107481-99-0.</i>	439.00	SEH1P0005
<input type="checkbox"/> I107480-99-0	Adj/Removable Headrest Link Hardware - Standard	305.00	
<input type="checkbox"/> I107481-99-0	Adj/Removable Headrest Link Hardware - Extended <i>Includes an additional 3" link to extend headrest pad range.</i>	329.00	
<input type="checkbox"/> I106331-99-0	Permobil 10"W Headrest, Leatherette <i>Must also select PN I10552.</i>	621.00	SEH1P0005
<input type="checkbox"/> I106333-99-0	Permobil Ergo Shape Headrest, Black Fabric <i>Must also select PN I10552.</i>	513.00	SEH1P0005
<input type="checkbox"/> I10552	Adj/Removable Headrest Hardware	133.00	SEHND0020
<input type="checkbox"/> I107499-99-0	Universal Headrest Adapter <i>Please Note: Recline is limited to 160° or less (depending on backrest height) when ordering a Universal Headrest Adapter on a 23" or 24" tall Corpus backrest.</i>	250.00	

CHOOSE ONE OF THE FOLLOWING:

Part Number	Description	MSRP - CA	ADP Code
<input type="checkbox"/> I107488-99-0	Permobil 3" x 5" Lateral Supports - Pair <i>Must also select I107490-99-0. Only compatible with Corpus backrest shells.</i>	459.00	SELND0010
<input type="checkbox"/> I107489-99-0	Permobil 5" x 7" Lateral Supports - Pair <i>Must also select I107490-99-0. Only compatible with Corpus backrest shells.</i>	459.00	SELND0010
<input type="checkbox"/> I107490-99-0	Adj/Removable Lateral Link Hardware - Pair	629.00	SEICF203L

OPTIONAL ITEMS:

Part Number	Description	MSRP - CA	ADP Code
<input type="checkbox"/> I103070-99-0	UT Amp. Support Pad (H) w/ Mesh Cover - (1 Each) 7"W x 7"T Curved Pad. <i>Must also select PN I10553.</i>	373.00	SEFND0005

Accessories

OPTIONAL ITEMS:

Part Number	Description	MSRP - CA	ADP Code
<input type="checkbox"/> I103071-99-0	UT Amp. Support Pad (D) w/ Mesh Cover - (1 Each) <i>8"W x 8"T Flat Pad. Must also select PN I10553.</i>	373.00	SEFN0005
<input type="checkbox"/> I10553	Adj Removable - Amp Support Hardware (1 Each) <i>Please indicate if you would like the amputation support mounted on the right or left side of the seat frame: Mount Amp Suppt on Left: _____ Mount Amp Suppt on Right: _____</i>	133.00	
<input type="checkbox"/> I106104-99-0	UT Hip Suppt (B) w/ Mesh Cover - High (Pair) <i>4.5"W x 3.5"T Pad w/ 4" Tall Mounting Bracket. Must also select PN I10550.</i>	304.00	
<input type="checkbox"/> I106103-99-0	UT Hip Suppt (A) w/ Mesh Cover - High (Pair) <i>6"W x 5.5"T Pad w/ 4" Tall Mounting Bracket. Must also select PN I10550.</i>	304.00	
<input type="checkbox"/> I103094-99-0	UT Thigh/Hip Suppt (J) w/ Mesh Cover - High (Pair) <i>6.5"W x 5"T Pad w/ 4" Tall Mounting Bracket. Must also select PN I10550.</i>	493.00	SELND0015
<input type="checkbox"/> I103095-99-0	UT Thigh/Hip Suppt (C) w/ Mesh Cover - High (Pair) <i>8"W x 3.5"T Pad w/ 4" Tall Mounting Bracket. Must also select PN I10550.</i>	493.00	SELND0015
<input type="checkbox"/> I102977-99-0	UT Thigh/Hip Suppt (G) w/ Mesh Cover - High (Pair) <i>12"W x 3.5"T Pad w/ 4" Tall Mounting Bracket. Must also select PN I10550.</i>	493.00	SELND0015
<input type="checkbox"/> I10550	Adj Removable - Thigh Support Hardware (Pair)	133.00	
<input type="checkbox"/> I106762-99-0	Abductor Kit Corpus - Leatherette <i>Must also select PN I10555.</i>	369.00	
<input type="checkbox"/> I10555	Adj Removable - Abductor Hardware (1 Each)	133.00	
<input type="checkbox"/> I105260-99-0	Transfer Handles, Corpus Seat	513.00	
<input type="checkbox"/> I105292-99-0	Push Handles, Corpus Seat	452.00	
<input type="checkbox"/> I106210-99-0	RAM Self-Leveling Cup Holder	88.00	
<input type="checkbox"/> I106806-99-0	RAM X-Grip Phone Holder	88.00	
<input type="checkbox"/> I107100-99-0	RAM X-Grip 10" Tablet Holder w/Swing Away Mount <i>Compatible with most 10" screen tablets with or without a protective case. Includes adjustable swing away mount. Note: Customization may be required if ordered along with an Omni.</i>	489.00	
<input type="checkbox"/> I107222-99-0	USB Charger <i>Provides up to 1.5A (at 5V) to charge most smartphones and tablets. The charger features an isolated circuit design for improved safety; and by pulling power from the ICS, prevents battery drainage when the chair is powered off. Includes mounting hardware and ICS connector. Please select armrest for mounting, keeping in mind that charger may need to be relocated. Mount on left armrest: _____ Mount on right armrest: _____</i>	209.00	
<input type="checkbox"/> I105805-99-0	Medical Necessities Bag Hooks <i>Please Note: Recline is limited to 165° or less (depending on backrest height) when ordering with 23" or 24" tall Corpus backrest.</i>	96.00	
<input type="checkbox"/> I102059-99-0	Medical Necessities Bag <i>Includes hardware to attach bag to backrest. Please Note: Recline is limited to 165° or less (depending on backrest height) when ordering with a 23" or 24" tall Corpus backrest.</i>	272.00	
<input type="checkbox"/> I105840-99-0	Bodypoint Monoflex Medium Chest Support Belt Kit <i>18.25"L x 4"W Pad Size. Item will be shipped as a part.</i>	187.00	SERND0001
<input type="checkbox"/> I105839-99-0	Bodypoint Monoflex Large Chest Support Belt Kit <i>21.5"L x 4.5"W Pad Size. Item will be shipped as a part.</i>	187.00	SERND0001
<input type="checkbox"/> I105850-99-0	Bodypoint Monoflex X-Lrge Chest Support Belt Kit <i>24.5"L x 5"W Pad Size. Item will be shipped as a part.</i>	187.00	SERND0001
<input type="checkbox"/> I10843	Upper Extremity Support Small, Angle Adjustable <i>Must also select PN I10556. Mount on Left: _____ Mount on Right: _____</i>	455.00	SETND0015
<input type="checkbox"/> I105221-99-0	Upper Extremity Support Angle Adjustable. 14" W <i>Must also select PN I10556.</i>	455.00	SETND0010
<input type="checkbox"/> I105222-99-0	Upper Extremity Support Angle Adjustable. 16" W <i>Must also select PN I10556.</i>	455.00	SETND0010
<input type="checkbox"/> I105220-99-0	Upper Extremity Support Angle Adjustable. 18" W <i>Must also select PN I10556.</i>	455.00	SETND0010

Accessories

OPTIONAL ITEMS:

Part Number	Description	MSRP - CA	ADP Code
<input type="checkbox"/> I105262-99-0	Upper Extremity Support Joystick Cutout 16" W <i>Must also select PN I10556.</i>	455.00	SETND0010
<input type="checkbox"/> I105263-99-0	Upper Extremity Support Joystick Cutout 18" W <i>Must also select PN I10556.</i>	455.00	SETND0010
<input type="checkbox"/> I106863-99-0	Upper Extremity Support Angle Adjust - UniTrack <i>Must also select PN I10556.</i>	455.00	
<input type="checkbox"/> I10556	Adj Removable - Tray Hardware (1 Each)	133.00	
<input type="checkbox"/> I105224-99-0	Oxygen holder	331.00	WEPX
<input type="checkbox"/> I105223-99-0	Crutch holder	344.00	

Alternative Drive Systems

ALTERNATIVE PROPORTIONAL JOYSTICK SYSTEMS:

Part Number	Description	MSRP - CA	ADP Code
<input type="checkbox"/> I10751	Mini-Joystick for Omni / IOM (9-pin connector) <i>Operating Force = 13 g Neutral to Full Deflection = 5mm</i>	3,031.00	
<input type="checkbox"/> I103893-99-0	Compact Joystick for Omni / IOM (9-pin connector) <i>Operating Force = 220 g Neutral to Full Deflection = 28mm Joystick includes two 3.5mm mono jacks for Mode and On/Off. Includes Standard Joystick Knob and 1.75" (42mm) Foam Ball Knob.</i>	1,479.00	
<input type="checkbox"/> I103894-99-0	Compact Joystick LITE for Omni / IOM (9-pin conn) <i>Operating Force = 50 g Neutral to Full Deflection = 17mm Joystick includes two 3.5mm mono jacks for Mode and On/Off. Includes Standard Joystick Knob and 1.75" (42mm) Foam Ball Knob.</i>	1,692.00	

ALT PROPORTIONAL JOYSTICK MOUNTING HARDWARE:

Part Number	Description	MSRP - CA	ADP Code
	NOTE: <i>Please see the PERMOfix Mounting Solutions, within the Alt Drive Accessories section, for additional mounting options.</i>		
<input type="checkbox"/> I104690-99-0	Alt. JS Swing Away Mount Kit - Armrest <i>Mount in front of Left Arm: _____ Mount in front of Right Arm: _____</i>	287.00	
<input type="checkbox"/> I107135-99-0	Alt JS/Omni Combo Swing Away Mount Kit - Armrest <i>Used to mount Omni and Alternative Proportional Joystick to the same armrest. Includes Swing Away Omni/Tablet Mount, with a longer mounting rod, two PERMOfix modules, one 3" rod and one 4" rod.</i> <i>Mount in front of Left Arm: _____ Mount in front of Right Arm: _____</i>	582.00	
<input type="checkbox"/> I104695-99-0	Alt. JS Adjustable Mount Kit - Midline <i>Mount on Left Arm: _____ Mount on Right Arm: _____</i>	873.00	
<input type="checkbox"/> I104143-99-0	Alt. JS Adjustable Mount Kit - Chin Boom <i>Swing-Away / Removable Chin Boom mounted on right of backrest shell. Custom charge may be applied if ordering lower backrest heights.</i>	931.00	
<input type="checkbox"/> I10752	Alt. JS Adjustable Mount Kit - Haims Harness, Lrg <i>For neck width up to 19cm.</i>	665.00	
<input type="checkbox"/> I10753	Alt. JS Adjustable Mount Kit - Haims Harness, Sml <i>For neck width up to 14cm.</i>	665.00	

ALT PROPORTIONAL JOYSTICK ACCESSORIES:

Part Number	Description	MSRP - CA	ADP Code
<input type="checkbox"/> I104904-99-0	Dual Switches for Compact Joystick - Fixed Mount <i>Includes two pushbutton switches mounted on rigid Switches can be used for Mode & On/Off.</i> <i>Switch operation force = 130 g</i> <i>Mount at front of Joystick: _____ Mount at rear of Joystick: _____</i> <i>Mount at left of Joystick: _____ Mount at right of Joystick: _____</i>	276.00	

Alternative Drive Systems

ALT PROPORTIONAL JOYSTICK ACCESSORIES:

Part Number	Description	MSRP - CA	ADP Code
<input type="checkbox"/> I104641-99-0	Dual Switches for Compact Joystick - Gooseneck Mnt <i>Includes two pushbutton switches mounted on adjustable gooseneck arms. Switches are positioned on Left and Right sides of Compact Joystick. Switches can be used for Mode & On/Off. Switch operation force = 130 g</i>	502.00	

SIP & PUFF SYSTEMS:

Part Number	Description	MSRP - CA	ADP Code
<input type="checkbox"/> I10573	Sip & Puff System for Omni <i>Requires Omni Control Module. Cannot use this part with an Input Module. Includes breath tube kit.</i>	1,778.00	WEX5
<input type="checkbox"/> I10574	Sip & Puff System with Fiber Optic Opti-Stop <i>Requires Omni Control Module. Cannot use this part with an Input Module. Includes breath tube kit.</i>	2,366.00	WEX5
<input type="checkbox"/> I106807-99-0	Sip & Puff Kit for Total Control Head Array <i>For use with Total Control Head Array. Requires Omni Control Module. Cannot be used with an Input Module. Includes breath tube kit.</i>	982.00	
<input type="checkbox"/> I10762	Adjustable Mounting Hrdw - Sip & Puff <i>Required when selecting PN I10573, I10574 or I106807-99-0.</i>	287.00	

ALTERNATIVE PROPORTIONAL CONTROL SYSTEMS:

Part Number	Description	MSRP - CA	ADP Code
<input type="checkbox"/> I10767	Proportional Foot Control <i>Mount on Left: _____ Mount on Right: _____</i>	3,121.00	
<input type="checkbox"/> I10766	Adjustable Mounting Hrdw - Foot Control <i>Required when selecting PN I10767.</i>	287.00	

Total Control Head Array

CONFIGURE HEAD ARRAY:

Part Number	Description	MSRP - CA	ADP Code
<input type="checkbox"/> I106805-99-0	Permobil Head Array - 3-Switch System <i>Features an occipital pad sensor and two stainless steel arms configured with proximity sensors. The arms are mounted to the upper attachment points on the head array. Longer arms are also included to extend the sensor adjustment range. ESP recommended. Choose Omni Control Module (ID51154IR) separately.</i>	5,485.00	WEX4
<input type="checkbox"/> I11072	Program Head Array as a 5-Switch System <i>Select additional arms/sensors below and indicate the intended use of each sensor in the Notes & Comments section. Choose Omni Control Module (ID51154IR) or Input Module (I10581) separately.</i>	No Charge	
CHOOSE OCCIPITAL PAD SIZE (REQUIRED)			
<input type="checkbox"/> I105920-99-0	Small Occipital Pad for Head Array - 4.25"W x 3"H	No Charge	
<input type="checkbox"/> I105921-99-0	Large Occipital Pad for Head Array - 6"W x 4"H	No Charge	
CHOOSE OCCIPITAL PAD MOUNT (REQUIRED):			
<input type="checkbox"/> I10847	Fixed Mount for Occipital Pad - Head Array	No Charge	
<input type="checkbox"/> I105936-99-0	Adj. Swivel Mount for Head Array Occipital Pad <i>Enables Occipital Pad to be angled up/down and left/right independent of the lateral arms.</i>	491.00	
CHOOSE BACKREST MOUNTING HARDWARE (REQUIRED):			
<input type="checkbox"/> I105927-99-0	Stealth TWB480 Backrest Mounting Hdwr with Links <i>Please select the length of the vertical stem below: Standard length (10"): _____ Short length (6"): _____</i>	305.00	
<input type="checkbox"/> I105941-99-0	Stealth TWB480 Backrest Mtng Hdwr with Flipback	504.00	
<input type="checkbox"/> I105925-99-0	Stealth HMO475 Backr Mtng Hdwr w 17mm Ball	258.00	
<input type="checkbox"/> I106319-99-0	Stealth TWB480 Backrest Mnt Hdwr w SubOcpl Mount	548.00	

Total Control Head Array

CONFIGURE HEAD ARRAY:

Part Number	Description	MSRP - CA	ADP Code
<input type="checkbox"/> I106282-99-0	Stealth 600 SubOccipital Pad - 8.5"W x 2.75"H <i>Must select when PN I106319-99-0 is ordered.</i>	111.00	
SELECT LEFT ARM PACKAGE (OPTIONAL):			
<input type="checkbox"/> I10976	Add Left Lateral Arm w Proximity Sensor <i>Choose Lateral Arm Adjustment Range: Short (5.5" - 10"): _____ Long (5.5" - 14"): _____</i>	416.00	
<input type="checkbox"/> I10977	Add Left Lateral Arm w Egg Switch <i>Choose Lateral Arm Adjustment Range: Short (5.5" - 10"): _____ Long (5.5" - 14"): _____</i>	277.00	
<input type="checkbox"/> I10978	Add Left Lateral Arm w Pad (No Sensor or Switch) <i>Choose Lateral Arm Adjustment Range: Short (5.5" - 10"): _____ Long (5.5" - 14"): _____</i>	169.00	
SELECT RIGHT ARM PACKAGE (OPTIONAL):			
<input type="checkbox"/> I10979	Add Right Lateral Arm w Proximity Sensor <i>Choose Lateral Arm Adjustment Range: Short (5.5" - 10"): _____ Long (5.5" - 14"): _____</i>	416.00	
<input type="checkbox"/> I10980	Add Right Lateral Arm w Egg Switch <i>Choose Lateral Arm Adjustment Range: Short (5.5" - 10"): _____ Long (5.5" - 14"): _____</i>	277.00	
<input type="checkbox"/> I10981	Add Right Lateral Arm w Pad (No Sensor or Switch) <i>Choose Lateral Arm Adjustment Range: Short (5.5" - 10"): _____ Long (5.5" - 14"): _____</i>	169.00	

Alt Drive Systems Setup

CHOOSE OMNI OR INPUT MODULE (REQUIRED):

Part Number	Description	MSRP - CA	ADP Code
<input type="checkbox"/> ID51154IR	Omni Control Module for R-net w/ Infrared <i>Must also select mounting hardware, either I103476-99-0 or I107135-99-0. Please specify which type of alternative drive input you will be using if not ordering it from Permobil. Input Device: _____</i>	2,478.00	
<input type="checkbox"/> I103476-99-0	Adjustable Swing-Away Hardware - Omni/Tablet <i>Mount in front of Left Arm: _____ Mount in front of Right Arm: _____</i>	282.00	
<input type="checkbox"/> I10581	Input Module (IOM) for R-net Control System <i>This device allows a proportional or switched drive input to be connected to an R-net Control System without the use of an Omni. It must be ordered in conjunction with an R-net Remote Joystick. It does NOT support any of the following input devices: Sip-n-Puff, 3-Switch Head Array or RIM control.</i>	1,749.00	

CHOOSE MOUNTING POSITION OF SEAT CONTROL SWITCHBOX:

Part Number	Description	MSRP - CA	ADP Code
<input type="checkbox"/> I10080	Mount Seat Control Switchbox on Backrest	No Charge	

Alt Drive Accessories

PERMOFIX MOUNTING SOLUTIONS:

Part Number	Description	MSRP - CA	ADP Code
<input type="checkbox"/> I106655-99-0	PERMOfix UniTrack Fixed Mount Kit <i>Includes two base joints, two 4" rods and one 8" rod.</i>	309.00	
<input type="checkbox"/> I106656-99-0	PERMOfix UniTrack Removable Mount Kit <i>Includes quick release mount and thigh support hardware, two base joints and three different length rods.</i>	543.00	
<input type="checkbox"/> I106657-99-0	PERMOfix UniTrack Swing Away Mount Kit <i>Includes PERMOfix Swing Away module, two base joints, one 2" rod and one 4" rod.</i>	783.00	

Alt Drive Accessories

PERMOFIX MOUNTING SOLUTIONS:

Part Number	Description	MSRP - CA	ADP Code
<input type="checkbox"/> I105907-99-0	Proximity Sensor + Swivel Assm - Not Mounted <i>Only available with the Total Control Head Array. Swivel assembly allows sensor to be mounted with PERMOfix or directly to lateral arm of Total Control Head Array. Since Proximity Sensors require power, they CANNOT be used as On/Off switches.</i>	404.00	
<input type="checkbox"/> I105919-99-0	Egg Switch + Swivel Mount Assm - Not Mounted <i>Swivel assembly allows switch to be mounted with PERMOfix or directly to lateral arm of Total Control Head Array.</i>	269.00	
<input type="checkbox"/> I105912-99-0	Solid Lateral Pad+Swivel Joint - Not Mounted <i>Swivel assembly allows pad to be mounted with PERMOfix or directly to lateral arm of Total Control Head Array.</i>	164.00	
<input type="checkbox"/> I106675-99-0	Oval Plate + Swivel Assm - No Switch - Not Mounted <i>Blank plate provides mounting option for Micro Lite, Buddy Button or other switch options. Swivel assembly allows pad to be mounted with PERMOfix or directly to lateral arm of Total Control Head Array.</i>	142.00	

OPTIONAL ITEMS:

Part Number	Description	MSRP - CA	ADP Code
<input type="checkbox"/> I101769-99-0	Micro Lite Switch Qty: _____	178.00	
<input type="checkbox"/> I101767-99-0	Buddy Button Qty: _____	134.00	
<input type="checkbox"/> I106293-99-0	Egg Switch Black Qty: _____	148.00	
<input type="checkbox"/> I101771-99-0	Mini Cup Switch Qty: _____	131.00	
<input type="checkbox"/> I107500-99-0	Proximity Switch Qty: _____	621.00	
<input type="checkbox"/> I107495-99-0	Mini Joystick 5 Switch OC for Omni/IOM (9-pin)	1,055.00	
<input type="checkbox"/> I107498-99-0	Touch Contact with CR2032 Battery Qty: _____	299.00	
<input type="checkbox"/> I10772	Dome Switch, Backrest Mnt w Gooseneck Qty: _____ <i>Mount on Left: _____ Mount on Right: _____ Mounts at top of backrest shell. Custom charge may be applied if ordering with short backrest heights. Kit Includes: One PERMOfix joint with Flat Plate, two PERMOfix joints with 45° swivel adjustment, one 150mm rod, one 200mm rod and one 225mm gooseneck switch with 1.2m cable and mounting clamp.</i>	553.00	
<input type="checkbox"/> I10846	Dome Switch, Backrest Swg Awy Mnt w Gsneck Qty: _____ <i>Mount on Left: _____ Mount on Right: _____ w 90° Swing-Away Function. Mounts at top of back. Custom charge applies if ordered with short back. Kit Includes: One PERMOfix joint with Flat Plate, one PERMOfix joint with 45° swivel adjustment, two 6" rods, one 90° Swing-Away joint, one 9" gooseneck switch w 50" cable and mounting clamp.</i>	830.00	
<input type="checkbox"/> I105935-99-0	3.5mm Ext. Cable for Sensor/Switch, 0.3m Long	38.00	
<input type="checkbox"/> I105931-99-0	3.5mm Ext. Cable for Sensor/Switch, 1m Long	42.00	
<input type="checkbox"/> I105932-99-0	3.5mm Ext. Cable for Sensor/Switch, 2m Long	46.00	
<input type="checkbox"/> I105933-99-0	3.5mm Ext. Cable for Sensor/Switch, 3m Long	50.00	
<input type="checkbox"/> I105934-99-0	3.5mm Y-splitter Cable for Sensor/Switch <i>Allows two Proximity Sensors/Mechanical Switches to connect to the same direction jack.</i>	42.00	

Notes & Additional Instructions:

Chair Order Policy:

You are highly encouraged to contact Permobil Sales with ANY questions regarding proper completion of our order forms at 1-800-736-0925. The order cannot be processed without the following:

- A client name or code, the client's measurements, the client's weight
- The order form must also be signed by a person authorized by your company to acknowledge the items selected. Permobil, LTD. is not responsible for configuration or size discrepancies resulting from customer errors on the order form.
- Order confirmations will be provided by Permobil summarizing the items selected. Please review this carefully.
- Once the order is shipped, any changes are subject to the returns policy and may be prohibited.

Order Acknowledgement:

I, _____, am an agent of the medical equipment provider named on this order form and I have the authority to contract for the purchase of powered wheelchairs and related parts on behalf of said provider. I acknowledge that I have reviewed this order and that it is complete and accurate to the best of my knowledge. I further acknowledge that any changes to this order after submission of this order form are subject to the returns policy and may carry additional charges.

Signed: _____

Printed Name: _____

Title: _____

Date: _____

NOTE: All specifications and prices are subject to change without notice. Please note that prices displayed are only valid if ordered with the wheelchair. Any parts ordered for service or future alterations will carry different pricing and warranty. *The ADP codes provided are not intended to be billing or legal advice, rather our interpretation of the code definitions. Use of the codes does not ensure coverage or payment for the item. For coverage information, verify the policy of the appropriate payer.*