

permobil

POWER WHEELCHAIRS

TOUGH ENOUGH



M300
CORPUS **HD**




RUGGED STYLE, REFINED DETAILS

M300 HD takes all the things you love about the Permobil's Mid Wheel Drive and amplifies them. Outfitted with stronger actuators, more powerful motors, and a reinforced seating system; then completed with elegant features like fully adjustable leg rests and up to four power seat functions, the M300 HD is tough enough.

COLOR OPTIONS



FABRIC OPTIONS



SPECIFICATIONS

Maximum User Weight	450 lbs
Maximum Speed	5 mph
Range	17.8 mi*
Base Width	26"
Base Length	38"
Minimum Turning Radius	22"
Ground Clearance	3"
Recommended Battery Type	Group 24 & Group 34
Weight Including Batteries	445 lbs with Group 24 425 lbs with Group 34
Drive Electronics	R-Net 120A
Maximum Slope	7.5° Up, 7.5° Down†

Seat-to-Floor Height	18", 19", 20"
Optional Seat Elevator	8"
Power Tilt	0° - 45° (power)
Recline Options	85° - 120° (manual) 85° - 150° (power)
Legrest Elevation	85° - 175°
Backrest Height	21", 24" - 29" (by 1" incr.)
Seat Depth	19" - 25" (by 1" incr.)
Seat Width	19", 21", 23", 25"
Armrest Pad Lengths	13", 16", 18"

*Actual driving range will vary based on driving and battery conditions.
†Refer to Owner Manual for safe driving practices.