

permobil
POWER WHEELCHAIRS

ADVENTURE
AWAITS...

K450
MX



EXPLORE MORE

Uniquely engineered seating system that brings the world to your fingertips.

KID-TOUGH TIRES

Large flat tires improve traction and reduce maintenance.

ALL-WHEEL SUSPENSION

Suspension system designed for increased balance, improved ride and safety.

High or low, explore more.



SPECIFICATIONS

Maximum User Weight	125 lbs
Maximum Speed	5 mph
Range	17 mi*
Base Width	25"
Base Length	35"
Minimum Turning Radius	28"
Ground Clearance	2.5"
Recommended Battery Type	Group 34
Weight Including Batteries	337 lbs
Drive Electronics	90A R-Net
Maximum Slope	10° Up, 10° Down†
Seat-to-Floor Height	3" – 26" (elevator)
Seat Elevator	8"
Power Tilt	0° – 45°
Manual Recline	85° – 105°
Manual Legrest Elevation	90° – 125°
Backrest Height	14" – 22.5"
Seat Depth	10" – 18"
Seat Width	10", 12", 14", 16"
Armrest Pad Lengths	8", 12.5"

*Actual driving range will vary based on driving and battery conditions.
†Refer to Owner Manual for safe driving practices.

COLOR OPTIONS



Chili
Pepper
Red



Pop Star
Pink



Extreme
Yellow



Titan
Blue



Platinum
Silver



Onyx
Black