

permobil

READY OR NOT...



**Koala**  
MINIFLEX





#### R-NET ELECTRONICS

Programmable R-net electronics come standard



#### MINIFLEX POWER SEATING

Designed to provide valuable repositioning



#### SAFETY FIRST

Optional bumper available— the ideal safety accommodation for beginner drivers



#### REMOTE STOP

Optional remote stop control gives any parent that extra layer of protection

### AVAILABLE SHROUD COLORS



Volcano Red



Pop Star Pink



Extreme Yellow



Green Apple



Cobalt Blue



Purple Rain



Onyx Black



Platinum Silver

### SPECIFICATIONS

Maximum User Weight	80 lbs
Maximum Speed	4 mph
Range	15 mi*
Base Width	23"
Base Length	32"
Minimum Turning Radius	21"
Ground Clearance	2.5"
Recommended Battery Type	M40
Drive Electronics	R-net
Maximum Slope	10° Up, 6° Down†
Seat-to-Floor Height (STFH)	15"

Optional Seat Elevator	8"
Tilt Options	No Manual Option 0° – 25° (power)
Recline Options	85° – 120° (manual) No Power Option
Legrest Elevation	90° – 150° (manual) No Power Option
Backrest Height	14", 16", 20"
Seat Depth	10", 12", 14"
Seat Width	10", 12", 14"
Armrest Pad Lengths	8", 10"

\*Actual driving range will vary based on driving and battery conditions.  
†Refer to Owner Manual for safe driving practices.