

**permobil**

POWER WHEELCHAIRS

HOW SWEET IT IS

**M300**  
PS JR.







Highly adjustable seat that grows as fast as you do.



**CATLINK SUSPENSION**

Provides a smoother ride, even over rough terrain, to help you handle the adventures ahead.



**HIGHLY CONFIGURABLE**


Allows you to maintain proper positioning as you grow.



**MULTI-MOTION ARMREST**

Perfectly positioned armrest with unique rotation and flip back movements that allow for easier transfers.

**COLOR OPTIONS**



**SPECIFICATIONS**

Maximum User Weight	165 lbs	Optional Seat Elevator	8"
Maximum Speed	6 mph	Tilt Options	0° – 45° (power) 0° – 45° (manual)
Range	13 mi*	Manual Recline	85° – 120°
Base Width	24"	Legrest Elevation	70° (swing-away) 90° – 150° (center mount elevating)
Base Length	36.5"	Backrest Height	16.5" – 22.5"
Minimum Turning Radius	20"	Seat Depth	10" – 18"
Ground Clearance	3"	Seat Width	11" – 16"
Recommended Battery Type	Group 24, Group 34	Armrest Pad Lengths	10", 13"
Weight Including Batteries	313 lbs G34 / 333 lbs G24	*Actual driving range will vary based on driving and battery conditions. †Refer to Owner Manual for safe driving practices.	
Drive Electronics	120A R-Net		
Maximum Slope	9° Up, 9° Down†		
Seat-to-Floor Height	17" (elevator) 17", 17.75", 18.5" (seat tube)		