

**permobil**

TILITE MANUAL WHEELCHAIRS

**TR**

## THE BEST FROM EVERY ANGLE

The sweeping curves of the titanium dual-tube frame are the perfect fusion of striking design and driven performance. Lighter, more efficient, and more elegant than any chair anywhere. As our flagship rigid chair, the TR is simply the best of the best.





### TIFIT. PERSONAL ERGONOMICS.

The difference between a TiFit® chair and all other chairs is that a TiFit chair is a wheeled prosthetic: Every frame is made to measure. A TiFit chair brings unparalleled performance and function. You will feel the difference in everything you do. TiFit...not standard.

### CUSTOM CONFIGURATION.

Custom Configuration allows you to add options and components to your chair to meet your needs and to optimize your riding experience. This includes choices of wheels, handrims, tires, wheel locks, color, tattoos, and more. TiLite offers more options than any other manufacturer!

### TECHNICAL SPECIFICATIONS

Frame Features	Titanium Dual-Tube Rigid TiFit Limited Adjustability	Camber	0°, 2°, 4°, 6°, 8° or 12°
Total Chair Width*	22"	Center of Gravity Adjustment	5"
Seat Width	12" - 20"	Front Seat Width (Seat Taper)	0", 1", 2", 3"
Seat Depth	12" - 20"	Ergonomic Seat	4" - 10"
Front Seat Height	15" - 21"	Transit Options	Q-Strait Lock Down Approved
Rear Seat Height	12.5" - 21"	Transport Weight*	9.3 lb.
Front Frame Angle	60°, 65°, 70°, 75°, 80°, 85°, or 90°	Weight Limit	265 lb.
Seat Back Type	<i>Folding Seat Back</i> Adj. Height 8.5" - 21" Adj. Angle 80° - 101°  <i>Fixed Seat Back</i> Fixed Height 5" - 20" Fixed Angle 80° - 101°	Warranty	Five Years on Frame One Year on TiLite Components

\*Weights and widths based on 16" x 16" chair with 0° camber, minimum rear wheel spacing and all standard features. Rear wheels are not included in transport weight. For reference, a pair of standard 24" rear wheels weighs 7.2 lb.


### COLOR OPTIONS

#### HIGH GLOSS PAINT



TiLite Red    Burgundy    Neon Pink    Super Wet Black

#### MATTE PAINT



Deep Blue Ice    Titanium Grey    Black Lava

#### CANDY PAINT



Ruby Red    Azure Blue    Electric Plum

#### NATURAL



Satin Titanium    Polished Titanium

#### METALLIC PAINT



Salsa Red Metallic    Tangerine Metallic    Acid Green Metallic    Ocean Blue Metallic

#### ANODIZED COLORS



Red Pepper    Burnt Orange    Sapphire    Ultraviolet    Fuchsia    Black

#### PEARLESCENT PAINT



Midnight Blue Pearl    White Pearl