



permobil

TILITE MANUAL WHEELCHAIRS

TRA

**ROLL BETTER.
LIVE BETTER.**

With the power and style of its titanium dual-tube frame, TiFit precision, and full adjustability, the TRA just rolls better. The result is better performance and comfort built around you.



TIFIT. PERSONAL ERGONOMICS.

The difference between a TiFit® chair and all other chairs is that a TiFit chair is a wheeled prosthetic: Every frame is made to measure. A TiFit chair brings unparalleled performance and function. You will feel the difference in everything you do. TiFit...not standard.

CUSTOM CONFIGURATION.

Custom Configuration allows you to add options and components to your chair to meet your needs and to optimize your riding experience. This includes choices of wheels, handrims, tires, wheel locks, color, tattoos, and more. TiLite offers more options than any other manufacturer!

TECHNICAL SPECIFICATIONS

Frame Features	Titanium Dual-Tube Rigid TiFit Limited Adjustability	Center of Gravity Adjustment	5"
Total Chair Width*	23"	Front Seat Width (Seat Taper)	0", 1", 2", 3"
Seat Width	12" – 20"	Transit Options	Q-Straint Lock Down Approved
Seat Depth	12" – 20"	Transport Weight*	11.1 lb.
Front Seat Height	15" – 21"	Weight Limit	265 lb. Standard 350 lb. Heavy Duty
Rear Seat Height	13" – 21"	HCPCS Code	K0005
Front Frame Angle	70°, 75°, 80°, 85°, or 90°	Warranty	Five Years on Frame One Year on TiLite Components
Seat Back Type	Folding Seat Back Adj. Height 8.5" – 21" Adj. Angle 80° – 101°	* Weights and widths based on 16" x 16" chair with 0° camber, minimum rear wheel spacing and all standard features. Rear wheels are not included in transport weight. For reference, a pair of standard 24" rear wheels weighs 7.2 lb.	
Camber	0°, 2°, 4°, 6°, 8° or 12°		


COLOR OPTIONS

HIGH GLOSS PAINT



TiLite Red Burgundy Neon Pink Super Wet Black

MATTE PAINT




Deep Blue Ice Titanium Grey Black Lava

CANDY PAINT




Ruby Red Azure Blue Electric Plum

NATURAL



Satin Titanium Polished Titanium

METALLIC PAINT



Salsa Red Metallic Tangerine Metallic Acid Green Metallic Ocean Blue Metallic

ANODIZED COLORS



Red Pepper Burnt Orange Sapphire Ultraviolet Fuchsia Black

PEARLESCENT PAINT



Midnight Blue Pearl White Pearl