



TR

FORM & FUNCTION
AS ONE.



The sweeping curves of the titanium dual-tube frame are the perfect fusion of striking design and driven performance. Lighter, more efficient, and more elegant than any chair anywhere. As TiLite's flagship rigid chair, the TR is simply the best of the best.



TECHNICAL SPECIFICATIONS.

| | |
|-------------------------------|----------------------------------------------------------------------------------------------------------------------------|
| Frame Features | Titanium Dual-Tube Rigid TiFit Limited Adjustability |
| Frame Options | Black Edition SuperLite Edition Power Adaptable Frame Reinforced Frame |
| Transport Weight* | 9.3 lb. |
| Weight Limit | 265 lb. |
| Total Chair Width* | 22" |
| Seat Width | 12" - 20" |
| Front Seat Width (Seat Taper) | 4 Options |
| Seat Depth | 12" - 20" |
| Ergonomic Seat | 7 Options |
| Front Seat Height | 15" - 21" |
| Rear Seat Height | 12.5" - 21" |
| Camber | 0°, 2°, 4°, 6°, 8°, or 12° |
| Front Frame Angle | 60°, 65°, 70°, 75°, 80°, 85°, or 90° |
| Footrest Width Options | Standard Taper or V-Front End |
| Center of Gravity Adjustment | 5" |
| Front Wheels | 23 Options |
| Rear Wheels | 67 Options |
| Handrims | 15 Options |
| Convenience | Folding Angle Adjustable Seat Back |
| Aesthetics | 6 Titanium Frame Finishes 16 Frame Colors 6 Color Anodize Packages 7 Ultrasuede Colors 8 Spinergy Wheel Colors |
| Warranty | Lifetime on Frame One Year on TiLite Components |

*Weights and widths based on 16" x 16" chair with 0° camber, minimum rear wheel spacing and all standard features. Rear wheels are not included in transport weight. For reference, a pair of standard 24" rear wheels weighs 7.2 lb.



TIFIT. PERSONAL ERGONOMICS.

The difference between a TiFit® chair and all other chairs is that a TiFit chair is a wheeled prosthetic: Every frame is made to measure. The optimized fit of a TiFit chair brings unparalleled performance and function. You will feel the difference in everything you do.

CUSTOM CONFIGURATION.

Custom Configuration allows you to add options and components to your chair to meet your needs and to optimize your riding experience. This includes choices of wheels, handrims, tires, wheel locks, color, tattoos, and more. TiLite offers more options than any other manufacturer!

To learn more about TiLite options that can take your chair's performance to the next level, visit tilite.com/optimize.



TiLite | The Ultimate Ride™ | © 2015 TiLite, LLC.
800.545.2266 | tilite.com | MS0016_RevB_TR3
TiLite wheelchairs are shown with available options.