


# C350 Corpus 3G

# permobil



Three column pricing on order/quote:

Display HCPCS Codes on order/quote:

Please send order/quote to fax#: **(800) 231-3256**

Email to: [sales@permobil.com](mailto:sales@permobil.com)

**\* Indicates a Required Field**



## Dealer Information

\* ATP: \_\_\_\_\_

Dealer Code: \_\_\_\_\_

\* Dealer Name: \_\_\_\_\_

Address: \_\_\_\_\_

\* City: \_\_\_\_\_

\* State/Zip: \_\_\_\_\_

\* Phone#: \_\_\_\_\_

\* Fax#: \_\_\_\_\_

PO#: \_\_\_\_\_

Email Address: \_\_\_\_\_

## Client Information

*Permobil recommends that the client is evaluated by a certified rehab specialist.*

\* First Name: \_\_\_\_\_

\* Last Name: \_\_\_\_\_

Diagnosis: \_\_\_\_\_

Funding Source: \_\_\_\_\_

Client Age: \_\_\_\_\_

## Client Measurements

\* Weight: \_\_\_\_\_

\* Height: \_\_\_\_\_

A) Top of Shoulders: \_\_\_\_\_

B) Chest Depth: \_\_\_\_\_

C) Chest Width: \_\_\_\_\_

D) Seat Depth: \_\_\_\_\_

E) Top of Head: \_\_\_\_\_

F) Elbow to Hand: \_\_\_\_\_

G) Seat Pan to Elbow: \_\_\_\_\_

H) Hip Width: \_\_\_\_\_

\* I) Knee to Foot: \_\_\_\_\_

**Permobil Inc.**

300 Duke Dr.

Lebanon, TN 37090

Tel: (800) 736-0925

Fax: (800) 231-3256

[www.permobil.com](http://www.permobil.com)

Prices effective **March 1, 2017**

# C350 Corpus 3G

## Chair Model and Colors

### CHOOSE BASE:

Part Number	Description	Price	HCPCS Code
<input type="checkbox"/> I10740	C350 Corpus 3G 5mph - C350 PS0 <i>Price Includes: RWD C350 Power Base Set Up No Power Seat Functions, Shock Absorbing Suspension System, R-net Controller, Tie-Down Hardware for Strap Systems, Anti-Tippers, Charger, Flat Free Tires, and Standard Positioning Belt. Must Choose Flat Free or Pneumatic Tires Below. 5.0mph, Weight Capacity = 300lbs, STFH = 17.5".</i>	7,885.00	K0848
<input type="checkbox"/> I10741	C350 Corpus 3G 5mph - C350 PS1 <i>Price Includes: RWD C350 Power Base Set Up For One Power Seat Function, Shock Absorbing Suspension System, R-net Controller, Tie Down Hardware for Strap Systems, Anti-Tippers, Charger, Flat Free Tires, and Standard Positioning Belt. Must Choose Flat Free or Pneumatic Tires Below. 5.0mph, Weight Capacity = 300lbs, STFH = 17.5".</i>	7,885.00	K0856
<input type="checkbox"/> I105206-99-0	C350 Corpus 3G 5mph - C350 PS2 <i>Price Includes: RWD C350 Power Base Set Up For Multiple Power Seat Functions, Shock Absorbing Suspension System, Tie-Down Hardware for Strap Systems, Anti-Tippers, Charger, Flat Free Tires, and Standard Positioning Belt. Must Choose Flat Free or Pneumatic Tires Below. 5.0mph, Weight Capacity = 300lbs, STFH = 17.5".</i>	7,885.00	K0861
<input type="checkbox"/> I10604	6.5mph High Speed Upgrade Package - C350 <i>Upgrade option for Group 3 C350 Power Base.</i>	3,010.00	
<input type="checkbox"/> I10742	C350 Corpus 3G 6.5mph - C350 HS PS2 <i>Price Includes: RWD C350 Power Base Set Up For Multiple Power Seat Functions, Shock Absorbing Suspension System, Tie-Down Hardware for Strap Systems, R-net Controller, Anti-Tippers, Charger, Flat Free Tires, and Standard Positioning Belt. Must Choose Flat Free or Pneumatic Tires Below. 6.5mph, Weight Capacity = 300lbs, STFH = 17.5".</i>	10,114.00	K0884

### CHOOSE DRIVE TIRES:

Part Number	Description	Price	HCPCS Code
<input type="checkbox"/> I103499-99-0	Tri-Spoke Split Rim w/ Flat Free Tire 3x8" Grey	No Charge	
<input type="checkbox"/> I103498-99-0	Tri-Spoke Split Rim w/ Pneumatic Tire 3x8" Grey	No Charge	

### OPTIONAL DRIVE TIRES:

Part Number	Description	Price	HCPCS Code
<input type="checkbox"/> I10782	Tri-Spoke Split Rim w/ Snow Tire Pneumatic 3x8" <i>Snow tires are designed for use in cold weather conditions. The nylon compound in these tires remains flexible at low temperatures, but wears more easily. Therefore, the snow tires (sold as wheel assemblies) are not recommended for year-long use and must be ordered along with standard drive wheels. Not available with black tires or black rims.</i>	456.00	

### CHOOSE SHROUD COLOR:

Part Number	Description	Price	HCPCS Code
<input type="checkbox"/> I103532-99-0	C350 Shroud Vintage Wine	No Charge	
<input type="checkbox"/> I103533-99-0	C350 Shroud Crystal Blue	No Charge	
<input type="checkbox"/> I103529-99-0	C350 Shroud Onyx Black	No Charge	
<input type="checkbox"/> I103530-99-0	C350 Shroud Platinum Silver	No Charge	
<input type="checkbox"/> I103531-99-0	C350 Shroud Extreme Yellow	No Charge	

### CHOOSE HUBCAP COLOR:

Part Number	Description	Price	HCPCS Code
<input type="checkbox"/> I10871	Match Hub Cap Color to Shroud: _____ <i>Note: Hub caps will be matched to shroud color if no selection is made.</i> Solid Color Hub Caps - Pair <i>Platinum Silver: _____ Onyx Black: _____ Cobalt Blue: _____ Phantom Grey: _____ Pop Star Pink: _____ Volcano Red: _____ Purple Rain: _____ Titan Blue: _____ Sublime Green: _____ Crystal Blue: _____ Vintage Wine: _____ Extreme Yellow: _____ Chili Pepper Red: _____</i>	No Charge	

## Base Options

### MUST CHOOSE EITHER SEAT ELEVATOR OR FIXED SEAT TUBE:

Part Number	Description	Price	HCPCS Code
<input type="checkbox"/> I106753-99-0	Power Adjustable Seat Height - 8" Travel	3,145.00	E2300
<input type="checkbox"/> I105421-99-0	Seat Tube Fixed	No Charge	

# C350 Corpus 3G

## Base Options

### MUST CHOOSE BATTERIES:

Part Number	Description	Price	HCPCS Code
<input checked="" type="checkbox"/> I107032-99-0	Batteries, Grp 34 (60Ah) Sealed Gel, x2, Installed	932.00	E2359

### OPTIONAL ITEMS:

Part Number	Description	Price	HCPCS Code
<input type="checkbox"/> I103496-99-0	Light Kit for C350	1,178.00	
<input type="checkbox"/> I10789	Seat Stabilizer Option	62.00	
<input type="checkbox"/> I10715	Compact Laptop Ventilator Tray - Corpus Seat <i>Mounts to rear of Corpus backrest hinge. Accommodates Pulmonetic Systems / CareFusion LTV Series and similar sized Laptop Ventilators. The tray is 13"T, adjustable in width from 8" - 13", and adjustable in depth from 2" - 6.5". Overall backrest angle limited to 150° when mounting this ventilator tray. Custom Option.</i>	1,560.00	E1029
<input type="checkbox"/> I10787	Battery Tray for Laptop Ventilator Tray <i>Measures 11 3/8" W x 6 1/2" H x 1 5/8" D. Accommodates Pulmonetic Systems Lithium-Ion Battery Pack, CareFusion SprintPack Battery System, or Richardson PowerTech Vent Power Center. Custom Option.</i>	261.00	K0108
<input type="checkbox"/> I10788	O2 Holder for Laptop Ventilator Tray <i>Mounts to Left or Right Side of Laptop Vent Tray. O2 Holder is 7" Tall. Width infinitely adjustable 3.5" - 7". Depth infinitely adjustable 3.75" - 5". Custom Option.</i>	331.00	K0108
<input type="checkbox"/> ID51113	Programming Key for R-net - Permobil Version	511.00	

## Joystick Options

### CHOOSE JOYSTICK TYPE:

Part Number	Description	Price	HCPCS Code
PLEASE NOTE: <i>For additional input devices, please refer to the Alternative Drive Systems Section.</i>			
<input type="checkbox"/> I105246-99-0	R-net LED Remote Joystick <i>This Joystick does NOT have a Display Screen and does NOT contain mono jacks. Limited to five profiles with this Joystick Module.</i>	No Charge	
<input type="checkbox"/> I10746	R-net Remote Color JS w/ Mono Jacks, No Lights, NC	No Charge	
<input type="checkbox"/> I107492-99-0	R-net Remote Color JS w/ Mono Jacks, No Lights	860.00	
<input type="checkbox"/> I107491-99-0	R-net Remote Color JS w/ Toggle Switches, No Lghts	1,082.00	
<input type="checkbox"/> I10747	R-net Remote Color JS w/ Mono Jacks, Lights, NC	No Charge	
<input type="checkbox"/> I104923-99-0	R-net Remote Color JS w/ Mono Jacks, Lights	860.00	
<input type="checkbox"/> I107493-99-0	R-net Remote Color JS w/ Toggle Switches, Lights	1,082.00	
<input type="checkbox"/> I105700-99-0	R-net Remote Color JS w/ Lights, Heavy Duty Gimbal <i>Heavy-Duty modification enables gimbal to hold up to extreme operating forces. Reinforced housing is attached to a standard Color LCD R-net Joystick. 3" long joystick knob has 3/4" diameter. Not compatible with alternative joystick handles.</i>	1,196.00	
<input type="checkbox"/> I103814-99-0	Compact Advanced Joystick for R-net <i>Compact sized joystick that does not contain an LCD Display. Must order Expandable Controller with this part. Contains two mono jacks for mode &amp; on/off. Does NOT contain a wrist support pad. If you want a wrist support pad, please select I10621 below.</i>	823.00	E2373
<input type="checkbox"/> I10621	Add Wrist Pad to Compact Joystick Advanced	258.00	

### CHOOSE JOYSTICK HANDLE TYPE:

Part Number	Description	Price	HCPCS Code
<input type="checkbox"/> I10045	Cone Shaped Joystick Knob - Standard	No Charge	
<input type="checkbox"/> I100960-99-0	Chin Cup For Joystick Knob	123.00	E2324
<input type="checkbox"/> I100961-99-0	"T" Handle For Joystick Knob	174.00	E2323
<input type="checkbox"/> I100962-99-0	"Large Ball" For Joystick Knob	123.00	E2323
<input type="checkbox"/> I100963-99-0	"Softball" For Joystick Knob	162.00	E2323
<input type="checkbox"/> I100964-99-0	"Mushroom" For Joystick Knob	162.00	E2323
<input type="checkbox"/> I100965-99-0	"Stick" For Joystick Knob	162.00	E2323
<input type="checkbox"/> I105442-99-0	1.75" Black Foam Ball For Joystick Knob	123.00	E2323
<input type="checkbox"/> I105846-99-0	Bodypoint J/S Handle 3" U-Shaped w/ Flex-Shaft	162.00	E2323
<input type="checkbox"/> I105847-99-0	Bodypoint J/S Handle 4" U-Shaped w/ Flex-Shaft	162.00	E2323

# C350 Corpus 3G

## Joystick Options

### CHOOSE JOYSTICK HANDLE TYPE:

Part Number	Description	Price	HCPCS Code
<input type="checkbox"/> I105848-99-0	Bodypoint J/S Handle Dome Shaped, Rubber	162.00	E2323

### CHOOSE JOYSTICK MOUNT:

Part Number	Description	Price	HCPCS Code
<input type="checkbox"/> I10664	Joystick Mount Right - Fixed	No Charge	
<input type="checkbox"/> I10665	Joystick Mount Left - Fixed	No Charge	
<input type="checkbox"/> I102537-99-0	Joystick Mount Right - Swing Away	363.00	E1028
<input type="checkbox"/> I102538-99-0	Joystick Mount Left - Swing Away	363.00	E1028
<input type="checkbox"/> I103820-99-0	Retractable Joystick Mount Right	563.00	E1028
<input type="checkbox"/> I103819-99-0	Retractable Joystick Mount Left	563.00	E1028
<input type="checkbox"/> I105210-99-0	Mount Standard Joystick on Corpus Back Hinge <i>Standard R-net Joystick is mounted to backrest hinge. Can be chosen instead of PN I105408-99-0. Mount on Left: _____ Mount on Right: _____</i>	334.00	
<input type="checkbox"/> I101260-99-0	Midline Control Mount <i>Mount Midline on Left Armrest: _____ Mount Midline on Right Armrest: _____</i>	332.00	E1028
<input type="checkbox"/> I106266-99-0	Stealth Gatlin Midline Mount - 8 Position <i>Features eight locking positions each 45° apart. Adjustments are made by pulling out on mounting arm and rotating to desired angle. Joystick bracket sold separately (select from the brackets listed below). Please indicate if you want it mounted on left or right armrest: Mount Midline on Left Armrest: _____ Mount Midline on Right Armrest: _____</i>	623.00	E1028
<input type="checkbox"/> I106265-99-0	Stealth Gatlin Midline Mount - 2 Position <i>Designed for easy swing away. Mount can be set to swing from horizontal to 95° or from horizontal to 180°. Joystick bracket sold separately (select from the brackets listed below). Please indicate if you want it mounted on left or right armrest: Mount Midline on Left Armrest: _____ Mount Midline on Right Armrest: _____</i>	623.00	E1028
<input type="checkbox"/> I106267-99-0	Stealth Gatlin Bracket for R-net/Compact JS Adv	259.00	E1028
<input type="checkbox"/> I106295-99-0	Stealth Gatlin Bracket for Compact Joystick	259.00	E1028
<input type="checkbox"/> I106296-99-0	Stealth Gatlin Bracket for Mini Joystick	259.00	E1028
<input type="checkbox"/> I106297-99-0	Stealth Gatlin Swivel Bracket for Omni	293.00	E1028
<input type="checkbox"/> I10842	PERMOFix Joint for Gatlin (0027-7027a) <i>Includes a shim washer that allows the PermoFix joint to be mounted on the 1/2" Stealth Gatlin arm. Provides an adjustable mount for a Mini Joystick or Compact Joystick.</i>	259.00	E1028
<input type="checkbox"/> I106268-99-0	Stealth Height Adjustable Mount <i>Replaces standard UT panel bracket. Provides 1.75" of height adjustability to all Permobil joystick mounts and Stealth Gatlin Midline mounts.</i>	232.00	E1028

### OPTIONAL ITEMS:

Part Number	Description	Price	HCPCS Code
<input type="checkbox"/> I10411	Expandable Controller - R-net <i>This item is required when ordering power tilt, power recline, and power elevating legrests together OR when ordering an Alternative Drive Control.</i>	679.00	E2377
<input type="checkbox"/> I10513	Harness for Expandable Controller <i>Must be ordered with any Expandable Controller (I10411 or I11068).</i>	602.00	E2313
<input type="checkbox"/> I105225-99-0	Co-Pilot Attendant Control R-net - Corpus	4,966.00	E2331
<input type="checkbox"/> I105408-99-0	Attendant Control R-net	1,332.00	E2331
<input type="checkbox"/> I107350-99-0	R-net Bluetooth iDevice Module <i>Allows input device commands (such as joystick nudges or switch clicks) to be used to operate an iPhone, iPad or iPod. By assigning each input to an iOS function, users can make phone calls, send texts and emails and access applications. Kit includes a Bluetooth iDevice Module, an R-net Bus Cable, and two mono jack cables for direct switch inputs (switches not included). Must be ordered with a Color JSM or Omni for mode display.</i>	934.00	K0108
<input type="checkbox"/> I107367-99-0	R-net Bluetooth Mouse Module <i>Enables mouse control through any R-net input device. It can operate the mouse on a computer and on select Android devices (with R-net for Android app). Kit includes a Bluetooth Mouse Module, an R-net Bus Cable, two mono jack cables for L/R click (switches not included) and a USB BT Adapter for user's computer. Note: BT Adapter may not work with all computers. Must be ordered with a Color JSM or Omni for mode display.</i>	934.00	K0108



# C350 Corpus 3G

## Joystick Options

### OPTIONAL ITEMS:

Part Number	Description	Price	HCPCS Code
<input type="checkbox"/> I107368-99-0	Add a second Bluetooth Mouse Module <i>Allows R-net to simultaneously pair to a second compatible computer or Android device. Each module is accessible through a separate R-net mode. Kit includes a unique Bluetooth Mouse Module, an R-net Bus Cable, two mono jack cables for L/R click and a USB BT Adapter.</i>	934.00	K0108
<input type="checkbox"/> I10582	Output Module (IOM) for R-net Control System <i>Enables a user to control an external device (ex: Communication Device, ECU, Door Opener, etc.) through the drive input (JSM, Head Array, etc.). Each module can provide up to 4 switch closures. May require specialty cable or wiring. Up to three output modules can be installed on one wheelchair.</i>	1,254.00	K0108

## Seating System

### CHOOSE SEAT MOUNTING PLATE TYPE:

Part Number	Description	Price	HCPCS Code
<input type="checkbox"/> I105055-99-0	Standard (0°) Tilt Mounting Plate <i>Range of Tilt will be 0° to +50°.</i>	No Charge	
<input type="checkbox"/> I105068-99-0	5° Anterior Tilt Mounting Plate <i>Range of Tilt will be -5° to +45°. This option will raise the minimum seat fo floor height by about 1" when the seat is at 0° of tilt.</i>	327.00	K0108

### CHOOSE SEATING SYSTEM:

Part Number	Description	Price	HCPCS Code
PLEASE NOTE: <i>Weight Limit 300lbs. For the Corpus 3G Seat on the C350, there is a half speed limitation when the Total Backrest Angle (TBA) &gt; 120°. TBA = Tilt Angle + Backrest Recline Angle. Example 1: 15° Tilt + 100° Recline = 115° TBA (Full Speed) Example 2: 35° Tilt + 110° Recline = 145° TBA (Half Speed)</i>			
<input type="checkbox"/> I10647	Corpus 3G Seat w/ No Tilt <i>Backrest angle can be manually adjusted up to 120°</i>	7,338.00	K0108
<input type="checkbox"/> I10744	C3G Seat w/ No Tilt & 150° Power Recline	7,795.00	E1004
<input type="checkbox"/> I10649	Corpus 3G Seat w/ 50° DuraCore Tilt <i>Backrest angle can be manually adjusted up to 120°</i>	7,575.00	E1002
<input type="checkbox"/> I10745	C3G Seat w/ 50° DuraCore Tilt & 150° Power Recline	11,995.00	E1007

## Seating System Options

### CHOOSE SEAT FUNCTION CONTROL TYPE:

Part Number	Description	Price	HCPCS Code
<input type="checkbox"/> I10418	Single Seat Function Control Switch for R-net <i>Please choose either Push Buttons or Toggle Switches below.</i>	No Charge	E2310
<input type="checkbox"/> I10419	Single Seat Function Control Kit for R-net <i>For operation of a single seat function through the wheelchair electronics via the drive control and/or ICS Switchbox. Order with R-net Controller or remote joystick options. Please indicate below if you want push buttons or toggles on the ICS Switchbox.</i>	1,867.00	E2310
<input type="checkbox"/> I10585	Multiple Seat Function Control Kit for R-net <i>MUST be selected if ordering two or more seat functions. This device is used to operate two or more seat functions through the wheelchair electronics via the drive control and/or ICS switch box. Order with R-net electronics or remote joystick options. Please indicate below if you want push buttons or toggles on the ICS Switchbox.</i>	3,533.00	E2311

### CHOOSE PUSH BUTTONS OR TOGGLES:

Part Number	Description	Price	HCPCS Code
<input type="checkbox"/> I102798-99-0	Push Buttons on ICS Switchbox (1.4m)	No Charge	
<input type="checkbox"/> I102799-99-0	Toggle Switches on ICS Switchbox (1.4m)	316.00	
<input type="checkbox"/> I105409-99-0	Do Not Install ICS Switchbox, Send as Part	No Charge	

### CHOOSE MOUNTING POSITION OF ICS SWITCHBOX:

Part Number	Description	Price	HCPCS Code
<input type="checkbox"/> I10626	Mount ICS Switchbox onto R-net Joystick Bracket <i>ICS Switchbox will be mounted to the same bracket as the R-net Joystick Module. If this item is selected, the ICS Switchbox will swing away with the Joystick Module.</i>	No Charge	

# C350 Corpus 3G

## Seating System Options

### CHOOSE MOUNTING POSITION OF ICS SWITCHBOX:

Part Number	Description	Price	HCPCS Code
<input type="checkbox"/> I104651-99-0	Mount ICS Switchbox onto Armrest <i>Mount on Left: _____ Mount on Right: _____</i>	No Charge	
<input type="checkbox"/> I106256-99-0	Mount ICS Switchbox on Backrest <i>Mount on Left: _____ Mount on Right: _____</i>	152.00	

### OPTIONAL ITEMS:

Part Number	Description	Price	HCPCS Code
<input type="checkbox"/> I107705-99-0	Virtual Seating Coach Package w/ Phone Holder <i>The Virtual Seating Coach (VSC) Package includes a Bluetooth dongle, a USB charger, a RAM X-grip Phone Holder and access to VSC coaching and regimen management. The VSC app can be downloaded for free from the App Store or Google Play. Requires an iPhone (iOS 9 or higher) or an Android phone (OS 4. 1. 2 or higher).</i> Mount VSC to Left Armrest: _____ <i>To provide access, USB Charger and Phone Holder may require relocation.</i> Mount VSC to Right Armrest: _____ <i>To provide access, USB Charger and Phone Holder may require relocation.</i>	2,055.00	K0108
<input type="checkbox"/> I104196-99-0	ICS Alternative Switchbox for R-net Electronics <i>Seat function interface box that accommodates up to 8 momentary switches. Provides an additional method for accessing the seat functions. 1/8" mini-mono switches not included.</i>	1,235.00	K0108
<input type="checkbox"/> I106743-99-0	Remote Stop System - All Electronics <i>Consists of a Radio Frequency (RF) handheld transmitter and an RF receiver mounted on the wheelchair. Custom charges may be required when ordering with a PS/PS Jr seat, a low backrest height or a pediatric wheelchair. If ordering on a K450, push handles must also be selected.</i>	1,221.00	E2321
<input type="checkbox"/> I10750	Memory Seat Program, Corpus w/ R-net	128.00	
<input type="checkbox"/> I106209-99-0	Legrest Movement Linked to Recline Electronically <i>Power Legrest will raise when Power Recline moves backward and will lower when back moves forward. If Power Legrest is operated, Power Recline will NOT automatically reposition. Feature works through ICS Switchbox and Driver Control.</i>	48.00	
<input type="checkbox"/> I10898	Enable Independent Repositioning Mode <i>IRM allows for sequential movement of power seat functions into a posterior position (35° tilt, 115° recline, 135° legrest) and returns to upright position (5° tilt, 100° recline, 110° legrest). This can be accessed using a standard drive control or a single switch. On the Permobil Joystick Module, IRM is accessible through the mode/profile jack (long switch hold). Please see the Alt Drive Accessories section for switch options.</i>	No Charge	

## Backrest Options

### CHOOSE FIXED OR SLIDING BACKREST KIT:

Part Number	Description	Price	HCPCS Code
<input type="checkbox"/> I105067-99-0	Sliding Backrest Kit - 5" Mechanical Shear	No Charge	
<input type="checkbox"/> I105066-99-0	Fixed Backrest Kit - No Mechanical Shear	38.00	

### CHOOSE BACKREST TYPE:

Part Number	Description	Price	HCPCS Code
<b>BACKREST HEIGHT NOTE:</b> <i>Backrest heights are measured from the seat pan to the top of the backrest shell with the backrest set at 90°. Ergo Backs include lateral wedges and lumbar pads.</i>			
<input type="checkbox"/> I105113-99-0	Corpus 3G Ergo Backrest - 14"W x 20"T (C3GERGBCA) <i>Power Recline will be limited to 150° when ordering this backrest.</i>	1,054.00	E2620
<input type="checkbox"/> I105114-99-0	Corpus 3G Ergo Backrest - 16"W x 20"T (C3GERGBCA) <i>Power Recline will be limited to 150° when ordering this backrest.</i>	1,054.00	E2620
<input type="checkbox"/> I105109-99-0	Corpus 3G Ergo Backrest - 18"W x 20"T (C3GERGBCA) <i>Power Recline will be limited to 150° when ordering this backrest.</i>	1,054.00	E2620
<input type="checkbox"/> I105852-99-0	Corpus 3G Ergo Backrest - 20"W x 20"T (C3GERGBCA) <i>Power Recline will be limited to 150° when ordering this backrest.</i>	1,054.00	E2620
<input type="checkbox"/> I10667	Solid Back 16"W X 20"T - Corpus 3G Seat <i>Power Recline will be limited to 150° when ordering this backrest.</i>	No Charge	

# C350 Corpus 3G

## Backrest Options

### CHOOSE BACKREST TYPE:

Part Number	Description	Price	HCPCS Code
<input type="checkbox"/> I10669	Solid Back 18"W x 20"T - Corpus 3G Seat <i>Power Recline will be limited to 150° when ordering this backrest.</i>	No Charge	
<input type="checkbox"/> I10883	Solid Back 20"W x 20"T - Corpus 3G Seat <i>Power Recline will be limited to 150° when ordering this backrest.</i>	No Charge	

### PLEASE NOTE:

*All taller backrests listed below include: lateral wedges, lumbar pads, a UT rail mounted at the top for lateral adjustment of the headrest, and split backrest feature with height adjustment in 1" increments. Please indicate below where you want the height set at so that a properly sized backrest cushion can be provided.*

<input type="checkbox"/> I105110-99-0	Corpus 3G Ergo Back -14"W x 23 - 28"T (C3GERGBCA)	1,058.00	E2620
<input type="checkbox"/> I105111-99-0	Corpus 3G Ergo Back -16"W x 23 - 28"T (C3GERGBCA)	1,058.00	E2620
<input type="checkbox"/> I105112-99-0	Corpus 3G Ergo Back -18"W x 23 - 28"T (C3GERGBCA)	1,058.00	E2620
<input type="checkbox"/> I105853-99-0	Corpus 3G Ergo Back -20"W x 23 - 28"T (C3GERGBCA)	1,058.00	E2620
<input type="checkbox"/> I10668	Solid Back 16"W x 23 - 28"T - Corpus 3G Seat	No Charge	
<input type="checkbox"/> I10670	Solid Back 18"W x 23 - 28"T - Corpus 3G Seat	No Charge	
<input type="checkbox"/> I10884	Solid Back 20"W x 23 - 28"T - Corpus 3G Seat	No Charge	
	Set Backrest at 23" _____		
	Set Backrest at 24" _____		
	Set Backrest at 25" _____		
	Set Backrest at 26" _____		
	Set Backrest at 27" _____		
	Set Backrest at 28" _____		

### CHOOSE BACKREST CUSHION:

Part Number	Description	Price	HCPCS Code
<input type="checkbox"/> I10659	Leatherette Cushion - Corpus Ergo Back	No Charge	
<input type="checkbox"/> I10660	Black Mesh Cushion - Corpus Ergo Back	42.00	K0108
<input type="checkbox"/> I105117-99-0	Omit Backrest Cushion	No Charge	

### OPTIONAL ITEMS:

Part Number	Description	Price	HCPCS Code
<input type="checkbox"/> I105782-99-0	Tubular Backrest Frame - Corpus (1" Canes) <i>Provides 1" diameter, vertical canes for mounting aftermarket backrests. Frame can be adjusted to support backrest widths from 14" to 20". Canes are attached with 2" long UniTrack linkage disks that can rotate to provide space for backrest hardware. The tops of the canes can be adjusted up to a height of 27" above the seat pan. Depending on the aftermarket backrest, several inches of seat depth may be lost when mounting this system.</i>	493.00	
<input type="checkbox"/> I106072-99-0	Direct Backrest Frame - Corpus <i>Used to mount ROHO Agility, Jay 3, Matrx Elite, Varilite Icon and Ride Java backrests ranging in width from 14" to 20". Slotted brackets and 1" UniTrack linkage disks provide a close, adjustable mount. Backrest pitch can be tuned by rotating the linkage disks. Some seat depth may be lost due to the position and orientation of the backrest.</i>	493.00	K0108
<input type="checkbox"/> I10725	Split Backrest Cushion - Corpus <i>Custom Option.</i>	333.00	E2617
<input type="checkbox"/> I10786	Tailored Backrest Cushion - Corpus <i>Sides of cushion are trimmed to provide additional clearance for lateral trunk supports. Must also select a backrest cushion (I10659 or I10660) when ordering this custom option.</i>	125.00	E2617

## Seat Sizing

### CHOOSE SEAT WIDTH:

Part Number	Description	Price	HCPCS Code
<input type="checkbox"/> I107023-99-0	Set Seat Width at 17": _____	No Charge	
<input type="checkbox"/> I107024-99-0	Set Seat Width at 19": _____	No Charge	
<input type="checkbox"/> I107025-99-0	Set Seat Width at 21": _____	No Charge	
<input type="checkbox"/> I107026-99-0	Set Seat Width at 23": _____	No Charge	

### CHOOSE SEAT DEPTH:

# C350 Corpus 3G

## Seat Sizing

### CHOOSE SEAT DEPTH:

Part Number	Description	Price	HCPCS Code
	PLEASE NOTE: <i>There are three different seat depth kits, each providing three seat depth settings: (1) 14" - 16", (2) 17" - 19" and (3) 20" - 22". Please select the required seat depth and the appropriate kit will be provided.</i>		
<input type="checkbox"/> I106232-99-0	Seat Depth 14"	No Charge	
<input type="checkbox"/> I106233-99-0	Seat Depth 15"	No Charge	
<input type="checkbox"/> I106239-99-0	Seat Depth 16"	No Charge	
<input type="checkbox"/> I106235-99-0	Seat Depth 17"	No Charge	
<input type="checkbox"/> I106237-99-0	Seat Depth 18"	No Charge	
<input type="checkbox"/> I106241-99-0	Seat Depth 19"	No Charge	
<input type="checkbox"/> I106242-99-0	Seat Depth 20"	No Charge	
<input type="checkbox"/> I106243-99-0	Seat Depth 21"	No Charge	
<input type="checkbox"/> I106244-99-0	Seat Depth 22"	No Charge	

### CHOOSE SEAT CUSHION:

Part Number	Description	Price	HCPCS Code
	NOTE: <i>A Corpus II Ergonomic Seat Cushion will be provided for select seat sizes.</i>		
<input type="checkbox"/> I10651	Corpus Ergo Cushion - Leatherette (C3GERGSCA)	553.00	E2605
<input type="checkbox"/> I10652	Corpus Ergo Cushion - Black Mesh (C3GERGSCA)	553.00	E2605
<input type="checkbox"/> I10965	Corpus Ergo Air Cushion (C3GRHOSCA)	586.00	E2622
	<i>Features ROHO Dry Floatation technology. The removable ROHO cell pad provides adjustable skin protection while the contoured foam and thigh support wedges help maintain proper positioning and stability. Includes removable stretch fabric cover. Available in widths of 17", 19" and 21" and depths of 17", 19" and 21", as well as 17" x 23".</i>		
<input type="checkbox"/> I11113	Incontinence Cover for Ergo Cushion	161.00	E2619
	<i>Made from a stretch Dartex fabric and suitable for all Ergo cushions including the Corpus Ergo Air. The cover is cut to the contours of the Ergo cushion and is designed to be used as an outer layer of protection. Cover size will be matched to cushion size. Sent as a part.</i>		
<input type="checkbox"/> I11166	ROHO QUADTRO SELECT HIGH PROFILE < 20"W	586.00	E2624
	<i>Registered in PDAC under ROHO model numbers QS98C - QS1014C.</i>		
<input type="checkbox"/> I11167	ROHO QUADTRO SELECT HIGH PROFILE > 20"W	628.00	E2625,E2624
	<i>Registered in PDAC under ROHO model numbers QS118C - QS1414C.</i>		
<input type="checkbox"/> I11168	ROHO HIGH PROFILE Single Comp < 20"W	586.00	E2622
	<i>Registered in PDAC under ROHO model numbers 1R98C - 1R1014C.</i>		
<input type="checkbox"/> I11169	ROHO HIGH PROFILE Single Comp > 20"W	628.00	E2623,E2622
	<i>Registered in PDAC under ROHO model numbers 1R118C - 1R1414C.</i>		
<input type="checkbox"/> I11170	ROHO HIGH PROFILE Sensor Ready < 20"W	586.00	E2622
	<i>Single compartment cushion designed to work with Smart Check system. Must also order I107505-99-0. Registered in PDAC under ROHO model numbers 1R98H-CA-SR - 1R1012H-CA-SR. Not available in all seat widths and depths.</i>		
<input type="checkbox"/> I11171	ROHO HIGH PROFILE Sensor Ready > 20"W	628.00	E2623,E2622
	<i>Single compartment cushion designed to work with Smart Check system. Must also order I107505-99-0. Registered in PDAC under ROHO model numbers 1R119H-CA-SR - 1R1411-CA-SR. Not available in all seat widths and depths.</i>		
<input type="checkbox"/> I107505-99-0	ROHO Smart Check	254.00	
	<i>Sensing device that allows users to setup and maintain proper inflation pressure in ROHO Sensor Ready cushions. Only available when ordered with a ROHO Sensor Ready cushion.</i>		
<input type="checkbox"/> I11173	Incontinence Cover for ROHO Cushion	161.00	E2619
	<i>Made from a stretch Dartex fabric and suitable for all ROHO seat cushions. Cover size will be matched to the cushion size. Sent as a part.</i>		
<input type="checkbox"/> I105118-99-0	Omit Seat Cushion	No Charge	

## Armrest Options

### CHOOSE ARMBAR KIT AND MOUNTING POSITION:

Part Number	Description	Price	HCPCS Code
-------------	-------------	-------	------------



# C350 Corpus 3G

## Armrest Options

### CHOOSE ARMBAR KIT AND MOUNTING POSITION:

Part Number	Description	Price	HCPSC Code
PLEASE NOTE:			
<i>Narrow Armbar Kit can be adjusted from 15"W to 19"W between the armrest pads without any additional parts. Wide Armbar Kit can be adjusted from 19"W to 23"W between the armrest pads without any additional parts. Please select both Armbar Kit and Width Setting below.</i>			
<input type="checkbox"/> I106543-99-0	Narrow Armbar Kit, Inf Adjustable 15"W - 19"W	No Charge	
<input type="checkbox"/> I106544-99-0	Wide Armrest Kit, Inf Adjustable 19"W - 23"W	No Charge	
<i>Required when selecting a 20" wide backrest.</i>			
Set the Width Between the Armrests to: _____			
<i>Must be wider than backrest width to ensure proper fit. Please round to nearest 1".</i>			

### CHOOSE ARMREST HEIGHT MOUNTING POSITION:

Part Number	Description	Price	HCPSC Code
<input type="checkbox"/> I105069-99-0	Height Adjustable Armrests - NC	No Charge	
<i>No charge when ordering power tilt and/or recline.</i>			
<input type="checkbox"/> I105081-99-0	Height Adjustable Armrests	346.00	
SET ARMREST HEIGHT AT: _____			
<i>Corpus Armrests are infinitely adjustable from 7" - 15", measured from the seat pan to the top of the armrest pad. No additional parts are needed to get this full range of adjustment. In one inch increments (ex: 8", 9", etc.), please indicate where you would like the armrests set.</i>			

### CHOOSE ARMREST SWIVEL LOCKED OR ENABLED:

Part Number	Description	Price	HCPSC Code
ARMREST MEDIAL SWIVEL NOTE:			
<i>Every Corpus Armrest can flip backwards for transfers. The armrests also include a feature that allows the armrest pads to swivel medially towards the center of the seating system. This swivel feature can be locked out by simply moving a bolt located under the armrest bracket. Please indicate below whether each armrest should have medial swivel enabled or locked.</i>			
<i>Note: Medial swivel must be locked on Corpus VS or when ordering arm trough adapters.</i>			
<input type="checkbox"/> I10645	Left Armrest Bracket - Medial Swivel Locked	No Charge	
<input type="checkbox"/> I10646	Right Armrest Bracket - Medial Swivel Locked	No Charge	
<input type="checkbox"/> I105073-99-0	Left Armrest Bracket - Medial Swivel Enabled	No Charge	
<input type="checkbox"/> I105072-99-0	Right Armrest Bracket - Medial Swivel Enabled	No Charge	

### CHOOSE ARMREST PAD RIGHT:

Part Number	Description	Price	HCPSC Code
<input type="checkbox"/> I105100-99-0	4"x10" Corpus Armrest Pad Right - Leatherette	No Charge	
<input type="checkbox"/> I105104-99-0	4"x13" Corpus Armrest Pad Right - Leatherette	No Charge	
<input type="checkbox"/> I105105-99-0	4"x16" Corpus Armrest Pad Right - Leatherette	No Charge	
<input type="checkbox"/> I105999-99-0	4"x18" Corpus HD Armrest Pad Right - Leatherette	No Charge	
<input type="checkbox"/> I107212-99-0	4"x10" Corpus Armrest Pad Right - Black Fabric	No Charge	
<input type="checkbox"/> I106047-99-0	4"x13" Corpus Armrest Pad Right - Black Fabric	No Charge	
<input type="checkbox"/> I107214-99-0	4"x16" Corpus Armrest Pad Right - Black Fabric	No Charge	
<input type="checkbox"/> I105995-99-0	4"x18" Corpus HD Armrest Pad Right - Black Fabric	No Charge	
<input type="checkbox"/> I106486-99-0	4"x10" Corpus Arm w/ Stealth Gel Pad - Right	169.00	
<input type="checkbox"/> I106488-99-0	4"x13" Corpus Arm w/ Stealth Gel Pad - Right	174.00	
<input type="checkbox"/> I106490-99-0	4"x16" Corpus Arm w/ Stealth Gel Pad - Right	179.00	
<input type="checkbox"/> I106281-99-0	4"x18" Corpus HD Arm w/ Stealth Gel Pad - Right	184.00	

### CHOOSE ARMREST PAD LEFT:

Part Number	Description	Price	HCPSC Code
<input type="checkbox"/> I105098-99-0	4"x10" Corpus Armrest Pad Left - Leatherette	No Charge	
<input type="checkbox"/> I105102-99-0	4"x13" Corpus Armrest Pad Left - Leatherette	No Charge	
<input type="checkbox"/> I105107-99-0	4"x16" Corpus Armrest Pad Left - Leatherette	No Charge	
<input type="checkbox"/> I105998-99-0	4"x18" Corpus HD Armrest Pad Left - Leatherette	No Charge	
<input type="checkbox"/> I107213-99-0	4"x10" Corpus Armrest Pad Left - Black Fabric	No Charge	
<input type="checkbox"/> I106046-99-0	4"x13" Corpus Armrest Pad Left - Black Fabric	No Charge	

# C350 Corpus 3G

## Armrest Options

### CHOOSE ARMREST PAD LEFT:

Part Number	Description	Price	HCPCS Code
<input type="checkbox"/> I107215-99-0	4"x16" Corpus Armrest Pad Left - Black Fabric	No Charge	
<input type="checkbox"/> I105994-99-0	4"x18" Corpus HD Armrest Pad Left - Black Fabric	No Charge	
<input type="checkbox"/> I106487-99-0	4"x10" Corpus Arm w/ Stealth Gel Pad - Left	169.00	
<input type="checkbox"/> I106489-99-0	4"x13" Corpus Arm w/ Stealth Gel Pad - Left	174.00	
<input type="checkbox"/> I106491-99-0	4"x16" Corpus Arm w/ Stealth Gel Pad - Left	179.00	
<input type="checkbox"/> I106303-99-0	4"x18" Corpus HD Arm w/ Stealth Gel Pad - Left	184.00	

### OPTIONAL ITEMS:

Part Number	Description	Price	HCPCS Code
<input type="checkbox"/> I105463-99-0	Corpus Arm Trough Adapter 1" - Pair <i>For aftermarket armrest pads. Please indicate below if they should be installed or on chair sent as parts. Please Note: Medial Swivel must be must be locked (I10645 &amp; I10646) when selecting this option.</i>	333.00	
<input type="checkbox"/> I105462-99-0	Corpus Arm Trough Adapter 7/8" - Pair <i>For aftermarket armrest pads. Please indicate below if they should be installed or on chair sent as parts. Please Note: Medial Swivel must be must be locked (I10645 &amp; I10646) when selecting this option.</i> Install Arm Trough Adapters: _____ Send Arm Trough Adapters as Parts: _____	333.00	
<input type="checkbox"/> I107339-99-0	Armrest Pouch - Right	236.00	
<input type="checkbox"/> I107338-99-0	Armrest Pouch - Left	236.00	
<input type="checkbox"/> I107840-99-0	UT Armrest Side Support - Right Outer <i>Includes 2"T x 6.5"W tapered pad.</i>	340.00	
<input type="checkbox"/> I107863-99-0	UT Armrest Side Support - Left Outer <i>Includes 2"T x 6.5"W tapered pad.</i>	340.00	
<input type="checkbox"/> I107864-99-0	UT Armrest Side Support - Right Inner <i>Includes 2"T x 6.5"W tapered pad. May require an increase to the width between the armrests to ensure armrests will flip-up.</i>	340.00	
<input type="checkbox"/> I107865-99-0	UT Armrest Side Support - Left Inner <i>Includes 2"T x 6.5"W tapered pad. May require an increase to the width between the armrests to ensure armrests will flip-up.</i>	340.00	
<input type="checkbox"/> I107839-99-0	UT Elbow Support - Right <i>Must also select PN I10881.</i>	168.00	
<input type="checkbox"/> I107838-99-0	UT Elbow Support - Left <i>Must also select PN I10882.</i>	168.00	
<input type="checkbox"/> I10881	Adjustable Elbow Support Hardware - Right	247.00	E1028
<input type="checkbox"/> I10882	Adjustable Elbow Support Hardware - Left	247.00	E1028

## Legrest Options

### CHOOSE LEGREST TYPE:

Part Number	Description	Price	HCPCS Code
<input type="checkbox"/> I105080-99-0	Corpus Manual Elevating Legrest	No Charge	
<input type="checkbox"/> I105079-99-0	Corpus Power Elevating Legrest - 170° <i>Includes mechanical articulation. HCPCS Determined by Pay Source.</i>	3,355.00	E1012

### PLEASE NOTE:

*If choosing the Power Extending Legrest or the Power Transfer Legrest, footplate height can be adjusted from 12"-16". For the Power Transfer Legrest, there is no pre-set position within the range. Footplate height is determined solely by user adjustment.*

<input type="checkbox"/> I10889	Power Extending Legrest <i>Provides 4" of independent legrest extension. Footplate height can be adjusted from 12"-16" before power extension. Restricts ability to adjust each footplate independently. Compatible with Manual and Power Elevating Legrests.</i>	3,473.00	K0108
<input type="checkbox"/> I105078-99-0	Power Transfer Legrest	4,048.00	K0108

# C350 Corpus 3G

## Legrest Options

### ADDITIONAL POWER OPTIONS:

Part Number	Description	Price	HCPCS Code
<input type="checkbox"/> I105078-99-0	Power Transfer Legrest <i>A single input command extends the footplate to the ground and then raises the seat for easier transfer. Footplate can be extended up to 27" from seat pan to footplate. Compatible with Manual and Power Elevating Legrests. Seat elevator is required for full functionality. Please Note: Calf support pads are only available as a custom option. Customization may also be required when ordered with a 5° anterior tilt mounting plate.</i>	4,048.00	K0108

### CHOOSE FOOTPLATE STYLE:

Part Number	Description	Price	HCPCS Code
<input type="checkbox"/> I10711	2pc Wide Corpus Footplates - 9"D x 7.5"W <i>Not compatible with Power Transfer Legrest. Note: Depending on legrest angle, the footplates can come in contact with the caster wheels.</i>	261.00	K0108
<input type="checkbox"/> I10712	2pc Narrow Corpus Footplates - 9"D x 6"W <i>Not compatible with Power Transfer Legrest.</i>	No Charge	
<input type="checkbox"/> I106760-99-0	1pc Narrow Corpus Footpl - 9.75"D x 12.5"W <i>Footplate is 12.5"W in front, 8"W in back.</i>	261.00	K0108
<input type="checkbox"/> I105793-99-0	1pc Narrow Power Transfer Footplate - 10"D x 13"W <i>Footplate is 13"W in front, 8"W in back.</i>	No Charge	

### FOOTPLATE HEIGHT POSITION:

Part Number	Description	Price	HCPCS Code
	SET FOOTPLATE HEIGHT AT: _____ <i>This is measured from the SEAT PAN to the top of the FOOTPLATE, when the legrest assembly is at 90°. DO NOT INCLUDE CUSHION THICKNESS IN YOUR MEASUREMENT. In one inch increments, please indicate desired footplate height. NOTE: Corpus footplate height range is 6"-17" and Corpus VS height range is 9"-17". Height may need to be reduced for shipping. Calf Pads are not available at heights less than 9".</i>		

### OPTIONAL ITEMS:

Part Number	Description	Price	HCPCS Code
<input type="checkbox"/> I105252-99-0	UT Calf Support Kit (B) w Mesh Cover & Adj Hardware <i>3.5"W x 4.5"T Pad.</i>	No Charge	
<input type="checkbox"/> I105258-99-0	UT Calf Support Kit (J) w Mesh Cover & Adj Hardware <i>5"W x 6.5"T Pad.</i>	No Charge	
<input type="checkbox"/> I105250-99-0	UT Calf Support Kit (A) w Mesh Cover & Adj Hardware <i>5.5"W x 6"T Pad.</i>	No Charge	
<input type="checkbox"/> I105256-99-0	UT Calf Support Kit (H) w Mesh Cover & Adj Hardware <i>7"W x 7"T Curved Pad.</i>	373.00	
<input type="checkbox"/> I105254-99-0	UT Calf Support Kit (D) w Mesh Cover & Adj Hardware <i>8"W x 8"T Pad.</i>	373.00	
<input type="checkbox"/> I106530-99-0	Stealth Gel Pads - Narrow Footplates (Pair)	370.00	
<input type="checkbox"/> I106225-99-0	Stealth Gel Pads - Wide Footplates (Pair)	381.00	

## Accessories

<input type="checkbox"/> I104518-99-0	Bodypoint Non-Padded Hip Belt - Medium	No Charge	
<input type="checkbox"/> I104520-99-0	Bodypoint Padded Hip Belt - Medium	149.00	E0978
<input type="checkbox"/> I104514-99-0	Bodypoint Non-Padded Hip Belt - Large	126.00	E0978
<input type="checkbox"/> I104516-99-0	Bodypoint Padded Hip Belt - Large	198.00	E0978
<input type="checkbox"/> I100519-99-0	Positioning Belt Retractable	295.00	K0108
<input type="checkbox"/> I100600-99-0	Positioning Belt Retractable Long	306.00	K0108

### HEADREST OPTIONS:

Part Number	Description	Price	HCPCS Code
<input type="checkbox"/> I107470-99-0	Permobil 7"W Headrest w/ Link Hardware <i>Must also select I107480-99-0 or I107481-99-0.</i>	385.00	E0955

# C350 Corpus 3G

## Accessories

### HEADREST OPTIONS:

Part Number	Description	Price	HCPCS Code
<input type="checkbox"/> I107469-99-0	Permobil 10"W Headrest w/ Link Hardware <i>Must also select I107480-99-0 or I107481-99-0.</i>	385.00	E0955
<input type="checkbox"/> I107468-99-0	Permobil 14"W Headrest w/ Link Hardware <i>Must also select I107480-99-0 or I107481-99-0.</i>	385.00	E0955
<input type="checkbox"/> I107480-99-0	Adj/Removable Headrest Link Hardware - Standard	265.00	E1028
<input type="checkbox"/> I107481-99-0	Adj/Removable Headrest Link Hardware - Extended <i>Includes an additional 3" link to extend headrest pad range.</i>	285.00	E1028
<input type="checkbox"/> I106331-99-0	Permobil 10"W Headrest, Leatherette <i>Must also select PN I10552.</i>	352.00	E0955
<input type="checkbox"/> I106333-99-0	Permobil Ergo Shape Headrest, Black Fabric <i>Must also select PN I10552.</i>	370.00	E0955
<input type="checkbox"/> I10552	Adj/Removable Headrest Hardware	272.00	E1028
<input type="checkbox"/> I107499-99-0	Universal Headrest Adapter <i>Please Note: Recline is limited to 160° or less (depending on backrest height) when ordering a Universal Headrest Adapter on a 23" or 24" tall Corpus backrest.</i>	258.00	

### TRUNK SUPPORT OPTIONS:

Part Number	Description	Price	HCPCS Code
<input type="checkbox"/> I107488-99-0	Permobil 3" x 5" Lateral Supports - Pair <i>Must also select I107490-99-0. Only compatible with Corpus backrest shells.</i>	405.00	E0956
<input type="checkbox"/> I107489-99-0	Permobil 5" x 7" Lateral Supports - Pair <i>Must also select I107490-99-0. Only compatible with Corpus backrest shells.</i>	405.00	E0956
<input type="checkbox"/> I107490-99-0	Adj/Removable Lateral Link Hardware - Pair	555.00	E1028

### OPTIONAL ITEMS:

Part Number	Description	Price	HCPCS Code
<input type="checkbox"/> I103070-99-0	UT Amp. Support Pad (H) w/ Mesh Cover - (1 Each) <i>7"W x 7"T Curved Pad. Must also select PN I10553.</i>	160.00	E1020
<input type="checkbox"/> I103071-99-0	UT Amp. Support Pad (D) w/ Mesh Cover - (1 Each) <i>8"W x 8"T Flat Pad. Must also select PN I10553.</i>	160.00	E1020
<input type="checkbox"/> I10553	Adj Removable - Amp Support Hardware (1 Each) <i>Please indicate if you would like the amputation support mounted on the right or left side of the seat frame: Mount Amp Suppt on Left: _____ Mount Amp Suppt on Right: _____</i>	272.00	E1028
<input type="checkbox"/> I106104-99-0	UT Hip Suppt (B) w/ Mesh Cover - High (Pair) <i>4.5"W x 3.5"T Pad w/ 4" Tall Mounting Bracket. Must also select PN I10550.</i>	280.00	E0956
<input type="checkbox"/> I106103-99-0	UT Hip Suppt (A) w/ Mesh Cover - High (Pair) <i>6"W x 5.5"T Pad w/ 4" Tall Mounting Bracket. Must also select PN I10550.</i>	280.00	E0956
<input type="checkbox"/> I103094-99-0	UT Thigh/Hip Suppt (J) w/ Mesh Cover - High (Pair) <i>6.5"W x 5"T Pad w/ 4" Tall Mounting Bracket. Must also select PN I10550.</i>	280.00	E0956
<input type="checkbox"/> I103095-99-0	UT Thigh/Hip Suppt (C) w/ Mesh Cover - High (Pair) <i>8"W x 3.5"T Pad w/ 4" Tall Mounting Bracket. Must also select PN I10550.</i>	280.00	E0956
<input type="checkbox"/> I102977-99-0	UT Thigh/Hip Suppt (G) w/ Mesh Cover - High (Pair) <i>12"W x 3.5"T Pad w/ 4" Tall Mounting Bracket. Must also select PN I10550.</i>	280.00	E0956
<input type="checkbox"/> I10550	Adj Removable - Thigh Support Hardware (Pair)	535.00	E1028
<input type="checkbox"/> I106762-99-0	Abductor Kit Corpus - Leatherette <i>Must also select PN I10555.</i>	376.00	E0957
<input type="checkbox"/> I10555	Adj Removable - Abductor Hardware (1 Each)	272.00	E1028
<input type="checkbox"/> I105260-99-0	Transfer Handles, Corpus Seat	488.00	K0108
<input type="checkbox"/> I105292-99-0	Push Handles, Corpus Seat	426.00	
<input type="checkbox"/> I106210-99-0	RAM Self-Leveling Cup Holder	83.00	
<input type="checkbox"/> I106806-99-0	RAM X-Grip Phone Holder	83.00	
<input type="checkbox"/> I107100-99-0	RAM X-Grip 10" Tablet Holder w/Swing Away Mount <i>Compatible with most 10" screen tablets with or without a protective case. Includes adjustable swing away mount (I103476-99-0). Customization may be required if ordered with an Omni.</i>	415.00	E1028



# C350 Corpus 3G

## Accessories

### OPTIONAL ITEMS:

Part Number	Description	Price	HCPCS Code
<input type="checkbox"/> I107222-99-0	USB Charger <i>Provides up to 1.5A (at 5V) to charge most smartphones and tablets. The charger features an isolated circuit design for improved safety; and by pulling power from the ICS, prevents battery drainage when the chair is powered off. Includes mounting hardware and ICS connector. Please select armrest for mounting, keeping in mind that charger may need to be relocated. Mount on left armrest: _____ Mount on right armrest: _____</i>	202.00	
<input type="checkbox"/> I105805-99-0	Medical Necessities Bag Hooks <i>Please Note: Recline is limited to 165° or less (depending on backrest height) when ordering with 23" or 24" tall Corpus backrest.</i>	98.00	K0108
<input type="checkbox"/> I102059-99-0	Medical Necessities Bag <i>Includes hardware to attach bag to backrest. Please Note: Recline is limited to 165° or less (depending on backrest height) when ordering with a 23" or 24" tall Corpus backrest.</i>	284.00	
<input type="checkbox"/> I105840-99-0	Bodypoint Monoflex Medium Chest Support Belt Kit <i>18.25"L x 4"W Pad Size. Item will be shipped as a part.</i>	182.00	E0960
<input type="checkbox"/> I105839-99-0	Bodypoint Monoflex Large Chest Support Belt Kit <i>21.5"L x 4.5"W Pad Size. Item will be shipped as a part.</i>	182.00	E0960
<input type="checkbox"/> I105850-99-0	Bodypoint Monoflex X-Lrge Chest Support Belt Kit <i>24.5"L x 5"W Pad Size. Item will be shipped as a part.</i>	182.00	E0960
<input type="checkbox"/> I10843	Upper Extremity Support Small, Angle Adjustable <i>Must also select PN I10556. Mount on Left: _____ Mount on Right: _____</i>	300.00	E0950
<input type="checkbox"/> I105221-99-0	Upper Extremity Support Angle Adjustable. 14" W <i>Must also select PN I10556.</i>	337.00	E0950
<input type="checkbox"/> I105222-99-0	Upper Extremity Support Angle Adjustable. 16" W <i>Must also select PN I10556.</i>	353.00	E0950
<input type="checkbox"/> I105220-99-0	Upper Extremity Support Angle Adjustable. 18" W <i>Must also select PN I10556.</i>	370.00	E0950
<input type="checkbox"/> I105262-99-0	Upper Extremity Support Joystick Cutout 16" W <i>Must also select PN I10556.</i>	389.00	E0950
<input type="checkbox"/> I105263-99-0	Upper Extremity Support Joystick Cutout 18" W <i>Must also select PN I10556.</i>	389.00	E0950
<input type="checkbox"/> I106863-99-0	Upper Extremity Support Angle Adjust - UniTrack <i>Side mounted tray that attaches to the UniTrack Accessory Rail. Dimensions are: 19.5"W x 13"D. Can only be ordered with chairs that offer the UniTrack System. Must also select PN I10556.</i>	365.00	E0950
<input type="checkbox"/> I10556	Adj Removable - Tray Hardware (1 Each)	272.00	E1028
<input type="checkbox"/> I105224-99-0	Oxygen Holder, Corpus	267.00	E2208
<input type="checkbox"/> I105223-99-0	Crutch Holder, Corpus	267.00	E2207

## Alternative Drive Systems

### ALTERNATIVE PROPORTIONAL JOYSTICK SYSTEMS:

Part Number	Description	Price	HCPCS Code
<input type="checkbox"/> I10751	Mini-Joystick for Omni / IOM (9-pin connector) <i>Operating Force = 0.5 oz (13 g) Neutral to Full Deflection = 0.2" (5mm)</i>	3,187.00	E2312
<input type="checkbox"/> I103893-99-0	Compact Joystick for Omni / IOM (9-pin connector) <i>Operating Force = 7.8 oz (220 g) Neutral to Full Deflection = 1.1" (28mm) Joystick includes two 1/8" mono jacks for Mode and On/Off. Includes Standard Joystick Knob and 1.75" (42mm) Foam Ball Knob.</i>	1,559.00	E2373
<input type="checkbox"/> I103894-99-0	Compact Joystick LITE for Omni / IOM (9-pin conn) <i>Operating Force = 1.8 oz (50 g) Neutral to Full Deflection = 0.7" (17mm) Joystick includes two 1/8" mono jacks for Mode and On/Off. Includes Standard Joystick Knob and 1.75" (42mm) Foam Ball Knob.</i>	1,763.00	E2373

### ALT PROPORTIONAL JOYSTICK MOUNTING HARDWARE:

Part Number	Description	Price	HCPCS Code
-------------	-------------	-------	------------

# C350 Corpus 3G

## Alternative Drive Systems

### ALT PROPORTIONAL JOYSTICK MOUNTING HARDWARE:

Part Number	Description	Price	HCPCS Code
	<b>NOTE:</b> <i>Please see the PERMOfix Mounting Solutions, within the Alt Drive Accessories section, for additional mounting options.</i>		
<input type="checkbox"/> I104690-99-0	Alt JS Swing Away Mount Kit - Armrest <i>Mount in front of Left Arm: _____ Mount in front of Right Arm: _____</i>	273.00	E1028
<input type="checkbox"/> I107135-99-0	Alt JS / Omni Combo Swing Away Mount Kit - Armrest <i>Used to mount Omni and Alternative Proportional Joystick to the same armrest. Includes Swing Away Omni/Tablet Mount, with a longer mounting rod, two PERMOfix modules, one 3" rod and one 4" rod.</i> <i>Mount in front of Left Arm: _____ Mount in front of Right Arm: _____</i>	516.00	E1028
<input type="checkbox"/> I104695-99-0	Alt JS Adjustable Mount Kit - Midline <i>Mount on Left Arm: _____ Mount on Right Arm: _____</i>	801.00	E1028
<input type="checkbox"/> I104143-99-0	Alt JS Adjustable Mount Kit - Chin Boom <i>Swing-Away / Removable Chin Boom mounted on right of backrest shell. Custom charge may be applied if ordering lower backrest heights.</i>	851.00	E1028
<input type="checkbox"/> I10752	Alt JS Adjustable Mount Kit - Haims Harness, Lrg <i>For neck width up to 7½".</i>	605.00	E1028
<input type="checkbox"/> I10753	Alt JS Adjustable Mount Kit - Haims Harness, Sml <i>For neck width up to 5½".</i>	605.00	E1028

### ALT PROPORTIONAL JOYSTICK ACCESSORIES:

Part Number	Description	Price	HCPCS Code
<input type="checkbox"/> I104904-99-0	Dual Switches for Compact Joystick - Fixed Mount <i>Includes two pushbutton switches mounted on rigid arms. Switches can be used for Mode &amp; On/Off. Switch operation force = 4.5oz (130 g)</i> <i>Mount at front of Joystick: _____ Mount at rear of Joystick: _____</i> <i>Mount at left of Joystick: _____ Mount at right of Joystick: _____</i>	253.00	K0108
<input type="checkbox"/> I104641-99-0	Dual Switches for Compact Joystick -Gooseneck Mount <i>Includes two pushbutton switches mounted on adjustable gooseneck arms. Switches are positioned on Left and Right sides of Compact Joystick. Switches can be used for Mode &amp; On/Off. Switch operation force = 4.5oz (130 g)</i>	456.00	K0108

### SIP & PUFF SYSTEMS:

Part Number	Description	Price	HCPCS Code
<input type="checkbox"/> I10573	Sip & Puff System for Omni <i>Requires Omni Control Module. Cannot use this part with an Input Module. Includes breath tube kit.</i>	1,705.00	E2325,E2326
<input type="checkbox"/> I10574	Sip & Puff System with Fiber Optic Opti-Stop <i>Requires Omni Control Module. Cannot use this part with an Input Module. Includes breath tube kit.</i>	2,269.00	E2325,E2326
<input type="checkbox"/> I106807-99-0	Sip & Puff Kit for Total Control Head Array <i>For use with Total Control Head Array. Requires Omni Control Module. Cannot be used with an Input Module. Includes breath tube kit.</i>	879.00	E2325
<input type="checkbox"/> I10762	Adjustable Mounting Hrdw - Sip & Puff <i>Required when selecting PN I10573, I10574 or I106807-99-0.</i>	273.00	E1028

### PROPORTIONAL FOOT CONTROL:

Part Number	Description	Price	HCPCS Code
<input type="checkbox"/> I10767	Proportional Foot Control <i>Mount on Left: _____ Mount on Right: _____</i>	3,249.00	
<input type="checkbox"/> I10766	Adjustable Mounting Hrdw - Foot Control <i>Required when selecting PN I10767.</i>	273.00	E1028

## Total Control Head Array

### CONFIGURE HEAD ARRAY:

Part Number	Description	Price	HCPCS Code
-------------	-------------	-------	------------

# C350 Corpus 3G

## Total Control Head Array

### CONFIGURE HEAD ARRAY:

Part Number	Description	Price	HCPCS Code
<input type="checkbox"/> I106805-99-0	Permobil Head Array - 3-Switch System <i>Features an occipital pad sensor and two stainless steel arms configured with proximity sensors. The arms are mounted to the upper attachment points on the head array. Longer arms are also included to extend the sensor adjustment range. ESP recommended. Choose Omni Control Module (ID51154IR) separately.</i>	4,750.00	E2330
<input type="checkbox"/> I11072	Program Head Array as a 5-Switch System <i>Select additional arms/sensors below and indicate the intended use of each sensor in the Notes &amp; Comments section. Choose Omni Control Module (ID51154IR) or Input Module (I10581) separately.</i>	No Charge	
CHOOSE OCCIPITAL PAD SIZE (REQUIRED):			
<input type="checkbox"/> I105920-99-0	Small Occipital Pad for Head Array - 4.25"W x 3"H	No Charge	
<input type="checkbox"/> I105921-99-0	Large Occipital Pad for Head Array - 6"W x 4"H	No Charge	
CHOOSE OCCIPITAL MOUNT (REQUIRED):			
<input type="checkbox"/> I10847	Fixed Mount for Head Array Occipital Pad	No Charge	
<input type="checkbox"/> I105936-99-0	Adj. Swivel Mount for Head Array Occipital Pad <i>Enables Occipital Pad to be angled up/down and left/right independent of the lateral arms.</i>	426.00	K0108
CHOOSE BACKREST MOUNTING HARDWARE (REQUIRED):			
<input type="checkbox"/> I105927-99-0	Stealth TWB480 Backrest Mounting Hdwr with Links <i>Please select the length of the vertical stem below: Standard length (10"): _____ Short length (6"): _____</i>	284.00	E1028
<input type="checkbox"/> I105941-99-0	Stealth TWB480 Backrest Mnt Hdwr with Flipback	459.00	E1028
<input type="checkbox"/> I105925-99-0	Stealth HMO475 MultiAxis Backrest Mounting Hdwr	252.00	E1028
<input type="checkbox"/> I106319-99-0	Stealth TWB480 Backrest Mnt Hdwr w SubOcptl Mount	493.00	E1028
<input type="checkbox"/> I106282-99-0	Stealth 600 SubOccipital Pad - 8.5" W x 2.75" H <i>Must select when PN I106319-99-0 is ordered.</i>	101.00	K0108
SELECT ADDITIONAL LEFT ARM PACKAGE (OPTIONAL):			
<input type="checkbox"/> I10976	Add Left Lateral Arm w Proximity Sensor <i>Choose Lateral Arm Adjustment Range: Short (5.5" - 10"): _____ Long (5.5" - 14"): _____</i>	612.00	K0108
<input type="checkbox"/> I10977	Add Left Lateral Arm w Egg Switch <i>Choose Lateral Arm Adjustment Range: Short (5.5" - 10"): _____ Long (5.5" - 14"): _____</i>	500.00	K0108
<input type="checkbox"/> I10978	Add Left Lateral Arm w Pad (No Sensor or Switch) <i>Choose Lateral Arm Adjustment Range: Short (5.5" - 10"): _____ Long (5.5" - 14"): _____</i>	408.00	K0108
SELECT ADDITIONAL RIGHT ARM PACKAGE (OPTIONAL):			
<input type="checkbox"/> I10979	Add Right Lateral Arm w Proximity Sensor <i>Choose Lateral Arm Adjustment Range: Short (5.5" - 10"): _____ Long (5.5" - 14"): _____</i>	612.00	K0108
<input type="checkbox"/> I10980	Add Right Lateral Arm w Egg Switch <i>Choose Lateral Arm Adjustment Range: Short (5.5" - 10"): _____ Long (5.5" - 14"): _____</i>	500.00	K0108
<input type="checkbox"/> I10981	Add Right Lateral Arm w Pad (No Sensor or Switch) <i>Choose Lateral Arm Adjustment Range: Short (5.5" - 10"): _____ Long (5.5" - 14"): _____</i>	408.00	K0108

## Alt Drive Systems Setup

### CHOOSE OMNI OR INPUT MODULE (REQUIRED):

Part Number	Description	Price	HCPCS Code
-------------	-------------	-------	------------

# C350 Corpus 3G

## Alt Drive Systems Setup

### CHOOSE OMNI OR INPUT MODULE (REQUIRED):

Part Number	Description	Price	HCPCS Code
<input type="checkbox"/> ID51154IR	Omni Control Module for R-net w/ Infrared <i>Must also select mounting hardware, either I103476-99-0 or I107135-99-0. Please specify which type of alternative drive input you will be using if not ordering it from Permobil. Input Device: _____</i>	1,772.00	
<input type="checkbox"/> I103476-99-0	Adjustable Swing-Away Hardware - Omni/Tablet <i>Mount in front of Left Arm: _____ Mount in front of Right Arm: _____</i>	273.00	E1028
<input type="checkbox"/> I10581	Input Module (IOM) for R-net Control System <i>This device allows a proportional or switched drive input to be connected to an R-net Control System without the use of an Omni. It must be ordered in conjunction with an R-net Remote Joystick. It does NOT support any of the following input devices: Sip-n-Puff, 3-Switch / 4-switch Head Array or RIM control.</i>	1,254.00	K0108

### CHOOSE MOUNTING POSITION OF SEAT CONTROL SWITCHBOX:

Part Number	Description	Price	HCPCS Code
<input type="checkbox"/> I10080	Mount Seat Control Switchbox on Backrest	No Charge	

## Alt Drive Accessories

### PERMOFIX MOUNTING SOLUTIONS:

Part Number	Description	Price	HCPCS Code
<input type="checkbox"/> I106655-99-0	PERMOfix UniTrack Fixed Mount Kit <i>Includes two base joints, two 4" rods and one 8" rod.</i>	275.00	K0108
<input type="checkbox"/> I106656-99-0	PERMOfix UniTrack Removable Mount Kit <i>Includes quick release mount and thigh support hardware, two base joints and three different length rods.</i>	483.00	E1028
<input type="checkbox"/> I106657-99-0	PERMOfix UniTrack Swing Away Mount Kit <i>Includes PERMOfix Swing Away module, two base joints, one 2" rod and one 4" rod.</i>	691.00	E1028
<input type="checkbox"/> I105907-99-0	Proximity Sensor + Swivel Assembly <i>Only available with the Total Control Head Array. Swivel assembly allows sensor to be mounted with PERMOfix or directly to lateral arm of Total Control Head Array. Since Proximity Sensors require power, they CANNOT be used as On/Off switches.</i>	349.00	K0108
<input type="checkbox"/> I105919-99-0	Egg Switch + Swivel Mount Assembly <i>Swivel assembly allows switch to be mounted with PERMOfix or directly to lateral arm of Total Control Head Array.</i>	236.00	K0108
<input type="checkbox"/> I105912-99-0	Solid Oval Pad + Swivel Assembly <i>Swivel assembly allows pad to be mounted with PERMOfix or directly to lateral arm of Total Control Head Array.</i>	142.00	K0108
<input type="checkbox"/> I106675-99-0	Oval Plate + Swivel Assembly - No Switch <i>Blank plate provides mounting option for Micro Lite, Buddy Button or other switch options. Swivel assembly allows pad to be mounted with PERMOfix or directly to lateral arm of Total Control Head Array.</i>	125.00	K0108

### OPTIONAL ITEMS:

Part Number	Description	Price	HCPCS Code
<input type="checkbox"/> I101769-99-0	Micro Lite Switch Qty: _____	190.00	
<input type="checkbox"/> I101767-99-0	Buddy Button Qty: _____	161.00	
<input type="checkbox"/> I106293-99-0	Egg Switch Black Qty: _____	161.00	K0108
<input type="checkbox"/> I101771-99-0	Mini Cup Switch Qty: _____	161.00	
<input type="checkbox"/> I107500-99-0	Proximity Switch Qty: _____	637.00	
<input type="checkbox"/> I107495-99-0	Mini Joystick 5 Switch OC for Omni/IOM (9-pin)	690.00	E2312
<input type="checkbox"/> I107498-99-0	Touch Contact with CR2032 Battery Qty: _____	392.00	
<input type="checkbox"/> I10772	Dome Switch, Backrest Mnt w Gooseneck Qty: _____ <i>Mount on Left: _____ Mount on Right: _____ Mounts at top of backrest shell. Custom charge may be applied if ordering with short backrest heights Kit Includes: One PERMOfix joint with Flat Plate, two PERMOfix joints with 45° swivel adjustment, one 6" rod, one 8" rod and one 9" gooseneck switch with 50" cable and mounting clamp.</i>	504.00	K0108



# C350 Corpus 3G

## Alt Drive Accessories

### OPTIONAL ITEMS:

Part Number	Description	Price	HCPCS Code
<input type="checkbox"/> I10846	Dome Switch, Backrest Swg Awy Mnt w Gsneck Qty:_____ Mount on Left: _____ Mount on Right: _____ w 90° Swing-Away Function. Mounts at top of back. Custom charge applies if ordered with short back. Kit Includes: One PERMOfix joint with Flat Plate, one PERMOfix joint with 45° swivel adjustment, two 6" rods, one 90° Swing-Away joint, one 9" gooseneck switch w 50" cable and mounting clamp.	745.00	K0108
<input type="checkbox"/> I105935-99-0	3.5mm Ext. Cable for Sensor/Switch, 0.3m Long	32.00	K0108
<input type="checkbox"/> I105931-99-0	3.5mm Ext. Cable for Sensor/Switch, 1m Long	35.00	K0108
<input type="checkbox"/> I105932-99-0	3.5mm Ext. Cable for Sensor/Switch, 2m Long	40.00	K0108
<input type="checkbox"/> I105933-99-0	3.5mm Ext. Cable for Sensor/Switch, 3m Long	44.00	K0108
<input type="checkbox"/> I105934-99-0	3.5mm Y-splitter Cable for Sensor/Switch Allows two Proximity Sensors/Mechanical Switches to connect to the same direction jack.	35.00	K0108

Additional Order Instructions (for Permobil):

Notes & Comments (not for Permobil):

### Chair Order Policy:

You are highly encouraged to contact Permobil Sales with ANY questions regarding proper completion of our order forms at 1-800-736-0925. The order cannot be processed without the following:

- A client name or code, the client's measurements, the client's weight
- The order form must also be signed by a person authorized by your company to acknowledge the items selected. Permobil, Inc. is not responsible for configuration or size discrepancies resulting from customer errors on the order form.
- Order confirmations will be provided by Permobil summarizing the items selected. Please review this carefully.
- Once the order is shipped, any changes are subject to the returns policy and may be prohibited.

### Order Acknowledgement:

I, \_\_\_\_\_, am an agent of the medical equipment provider named on this order form and I have the authority to contract for the purchase of powered wheelchairs and related parts on behalf of said provider. I acknowledge that I have reviewed this order and that it is complete and accurate to the best of my knowledge. I further acknowledge that any changes to this order after submission of this order form are subject to the returns policy and may carry additional charges.

Signed: \_\_\_\_\_

Printed Name: \_\_\_\_\_

Title: \_\_\_\_\_

Date: \_\_\_\_\_

NOTE: All specifications and prices are subject to change without notice. Please note that prices displayed are only valid if ordered with the wheelchair. Any parts ordered for service or future alterations will carry different pricing and warranty. *The HCPCS codes provided are not intended to be billing or legal advice, rather our interpretation of the code definitions. Use of the codes does not ensure coverage or payment for the item. For coverage information, verify the policy of the appropriate payer.*