

# K450 MX



**\* Indicates a Required Field**

### Dealer Information

\* ATP: \_\_\_\_\_

Dealer Code: \_\_\_\_\_

\* Dealer Name: \_\_\_\_\_

Address: \_\_\_\_\_

\* City: \_\_\_\_\_

\* State/Zip: \_\_\_\_\_

\* Phone#: \_\_\_\_\_

\* Fax#: \_\_\_\_\_

PO#: \_\_\_\_\_

Email Address: \_\_\_\_\_

### Client Information

*Permobil recommends that the client is evaluated by a certified rehab specialist.*

\* First Name: \_\_\_\_\_

\* Last Name: \_\_\_\_\_

Diagnosis: \_\_\_\_\_

Funding Source: \_\_\_\_\_

Client Age: \_\_\_\_\_

### Client Measurements

\* Weight: \_\_\_\_\_

\* Height: \_\_\_\_\_

A) Top of Shoulders: \_\_\_\_\_

B) Chest Depth: \_\_\_\_\_

C) Chest Width: \_\_\_\_\_

D) Seat Depth: \_\_\_\_\_


E) Top of Head: \_\_\_\_\_

F) Elbow to Hand: \_\_\_\_\_

G) Seat Pan to Elbow: \_\_\_\_\_

H) Hip Width: \_\_\_\_\_

\*I) Knee to Foot: \_\_\_\_\_



Three column pricing on order/quote:

Display HCPCS Codes on order/quote:

Please send order/quote to fax#: **(800) 231-3256**

Email to: [sales@permobil.com](mailto:sales@permobil.com)

**Permobil Inc.**  
300 Duke Dr.  
Lebanon, TN 37090  
Tel: (800) 736-0925  
Fax: (800) 231-3256  
[www.permobil.com](http://www.permobil.com)

Prices effective **March 1, 2017**

# K450 MX

## Chair Model and Colors

### CHOOSE BASE:

Part Number	Description	Price	HCPCS Code
<input checked="" type="checkbox"/> I104524-99-0	K450 Base w/ R-net 5.2mph - K450 PS2 <i>Price Includes: K450 Power Base Set Up For Multiple Power Seat Functions, Tie-Down Hardware for Strap Systems, Charger, Flat Free Tires, and Standard Positioning Belt. Must Choose Flat Free or Pneumatic Tires Below. 5.0mph, Weight Capacity = 125lbs, STFH = 3"-25".</i>	14,735.00	K0891

### CHOOSE DRIVE TIRES:

Part Number	Description	Price	HCPCS Code
<input type="checkbox"/> I106204-99-0	Flat Free Casters - K450	No Charge	
<input type="checkbox"/> I106205-99-0	Flat Free Drive Tires - K450	No Charge	
<input type="checkbox"/> I106206-99-0	Pneumatic Casters - K450	No Charge	
<input type="checkbox"/> I106207-99-0	Pneumatic Drive Tires - K450	No Charge	

### CHOOSE SHROUD COLOR:

Part Number	Description	Price	HCPCS Code
<input type="checkbox"/> I104535-99-0	Shroud for K450 - Onyx Black	No Charge	
<input type="checkbox"/> I104536-99-0	Shroud for K450 - Platinum Silver	No Charge	
<input type="checkbox"/> I104537-99-0	Shroud for K450 - Pop Star Pink	No Charge	
<input type="checkbox"/> I104539-99-0	Shroud for K450 - Extreme Yellow	No Charge	
<input type="checkbox"/> I106367-99-0	Shroud for K450 - Volcano Red	No Charge	
<input type="checkbox"/> I106368-99-0	Shroud for K450 - Cobalt Blue	No Charge	
<input type="checkbox"/> I107234-99-0	Shroud for K450 - Purple Rain	No Charge	

## Base Options

### MUST CHOOSE SEAT ELEVATOR:

Part Number	Description	Price	HCPCS Code
<input checked="" type="checkbox"/> I10510	Power Adjustable Seat Height - 8" Travel <i>This item must be selected.</i>	3,145.00	E2300

### MUST CHOOSE BATTERIES:

Part Number	Description	Price	HCPCS Code
<input checked="" type="checkbox"/> I107032-99-0	Batteries, Grp 34 (60Ah) Sealed Gel, x2, Installed	932.00	E2359

### OPTIONAL ITEMS:

Part Number	Description	Price	HCPCS Code
<input type="checkbox"/> I104550-99-0	Anti-Tippers for K450	467.00	
<input type="checkbox"/> I104540-99-0	Light Kit - K450 R-net	1,195.00	
<input type="checkbox"/> I104553-99-0	WC 19 Transit Kit - K450 <i>This kit includes additional anchor points for tie down straps, seat stabilizing bars and hardware, and a WC-19 Approved integrated lap and shoulder belt. This kit has been approved according to ANSI/RESNA WC19. An ICS Switchbox must be installed on the chair in order to access WC19 mode.</i>	1,864.00	
<input type="checkbox"/> I104554-99-0	Fixed Ventilator Tray - K450 <i>Tray is 9.75"D x 13"W and increases the overall length of a K450 by approximately 11". The ventilator tray is designed to accommodate an LP10 style ventilator and one U1 Battery. A custom charge will be applied if ordering vent tray with any of the following items: WC-19 Transit Kit (I104553-99-0), Attendant Control (ID50882), Push Handles (I104576-99-0) or Bag Hooks (I105805-99-0). Please Note: Must also select Anti-Tippers (I104550-99-0).</i>	617.00	E1029
<input type="checkbox"/> ID51113	Programming Key for R-net - Permobil Version	511.00	

## Joystick Options

### CHOOSE JOYSTICK TYPE:

Part Number	Description	Price	HCPCS Code
<input type="checkbox"/> I105246-99-0	R-net LED Remote Joystick <b>PLEASE NOTE:</b> <i>For additional input devices, please refer to the Alternative Drive Systems Section. This Joystick does NOT have a Display Screen and does NOT contain mono jacks. Limited to five profiles with this Joystick Module.</i>	No Charge	

# K450 MX

## Joystick Options

### CHOOSE JOYSTICK TYPE:

Part Number	Description	Price	HCPCS Code
<input type="checkbox"/> I107492-99-0	R-net Remote Color JS w/ Mono Jacks, No Lights	860.00	
<input type="checkbox"/> I107491-99-0	R-net Remote Color JS w/ Toggle Switches, No Lghts	1,082.00	
<input type="checkbox"/> I104923-99-0	R-net Remote Color JS w/ Mono Jacks, Lights	860.00	
<input type="checkbox"/> I107493-99-0	R-net Remote Color JS w/ Toggle Switches, Lights	1,082.00	
<input type="checkbox"/> I103814-99-0	Compact Advanced Joystick for R-net <i>Compact sized joystick that does not contain an LCD Display. Must order Expandable Controller with this part. Contains two mono jacks for mode &amp; on/off. Does NOT contain a wrist support pad. If you want a wrist support pad, please select I10621 below.</i>	823.00	E2373
<input type="checkbox"/> I10621	Add Wrist Pad to Compact Joystick Advanced	258.00	

### CHOOSE JOYSTICK HANDLE TYPE:

Part Number	Description	Price	HCPCS Code
<input type="checkbox"/> I10045	Cone Shaped Joystick Knob - Standard	No Charge	
<input type="checkbox"/> I100960-99-0	Chin Cup For Joystick Knob	123.00	E2324
<input type="checkbox"/> I100961-99-0	"T" Handle For Joystick Knob	174.00	E2323
<input type="checkbox"/> I100962-99-0	"Large Ball" For Joystick Knob	123.00	E2323
<input type="checkbox"/> I100963-99-0	"Softball" For Joystick Knob	162.00	E2323
<input type="checkbox"/> I100964-99-0	"Mushroom" For Joystick Knob	162.00	E2323
<input type="checkbox"/> I100965-99-0	"Stick" For Joystick Knob	162.00	E2323
<input type="checkbox"/> I105442-99-0	1.75" Black Foam Ball For Joystick Knob	123.00	E2323
<input type="checkbox"/> I105846-99-0	Bodypoint J/S Handle 3" U-Shaped w/ Flex-Shaft	162.00	E2323
<input type="checkbox"/> I105847-99-0	Bodypoint J/S Handle 4" U-Shaped w/ Flex-Shaft	162.00	E2323
<input type="checkbox"/> I105848-99-0	Bodypoint J/S Handle Dome Shaped, Rubber	162.00	E2323

### CHOOSE JOYSTICK MOUNT:

Part Number	Description	Price	HCPCS Code
<input type="checkbox"/> I10412	Joystick Mount Right - Fixed, R-net Only	No Charge	
<input type="checkbox"/> I10413	Joystick Mount Left - Fixed, R-net Only	No Charge	
<input type="checkbox"/> I102537-99-0	Joystick Mount Right - Swing Away	363.00	E1028
<input type="checkbox"/> I102538-99-0	Joystick Mount Left - Swing Away	363.00	E1028

### OPTIONAL ITEMS:

Part Number	Description	Price	HCPCS Code
<input type="checkbox"/> I10411	Expandable Controller - R-net <i>This item is required when ordering power tilt, power recline, and power elevating legrests together OR when ordering an Alternative Drive Control.</i>	679.00	E2377
<input type="checkbox"/> I10513	Harness for Expandable Controller <i>Must be ordered with any Expandable Controller (I10411 or I11068).</i>	602.00	E2313
<input type="checkbox"/> ID50882-A	Attendant Control for R-net Control System Only <i>On the K450, the mounting location for the Attendant Control is off of the Push Handle Bracket. If ordering an attendant control, you must also order I104576-99-0 (Push Handles - MX Seat).</i>	1,332.00	E2331
<input type="checkbox"/> I104571-99-0	Co-Pilot Attendant Control - K450	3,203.00	
<input type="checkbox"/> I107350-99-0	R-net Bluetooth iDevice Module <i>Allows input device commands (such as joystick nudges or switch clicks) to be used to operate an iPhone, iPad or iPod. By assigning each input to an iOS function, users can make phone calls, send texts and emails and access applications. Kit includes a Bluetooth iDevice Module, an R-net Bus Cable, and two mono jack cables for direct switch inputs (switches not included). Must be ordered with a Color JSM or Omni for mode display.</i>	934.00	K0108
<input type="checkbox"/> I107367-99-0	R-net Bluetooth Mouse Module <i>Enables mouse control through any R-net input device. It can operate the mouse on a computer and on select Android devices (with R-net for Android app). Kit includes a Bluetooth Mouse Module, an R-net Bus Cable, two mono jack cables for L/R click (switches not included) and a USB BT Adapter for user's computer. Note: BT Adapter may not work with all computers. Must be ordered with a Color JSM or Omni for mode display.</i>	934.00	K0108

# K450 MX

## Joystick Options

### OPTIONAL ITEMS:

Part Number	Description	Price	HCPCS Code
<input type="checkbox"/> I107368-99-0	Add a second Bluetooth Mouse Module <i>Allows R-net to simultaneously pair to a second compatible computer or Android device. Each module is accessible through a separate R-net mode. Kit includes a unique Bluetooth Mouse Module, an R-net Bus Cable, two mono jack cables for L/R click and a USB BT Adapter.</i>	934.00	K0108

## Seating System

### CHOOSE ONE OF THE FOLLOWING:

Part Number	Description	Price	HCPCS Code
<input checked="" type="checkbox"/> I104324-99-0	MX Power Seating System <i>Price Includes: MX Seating System, Adjustable Seat Frame, Power Seat to Floor System, Adjustable Backrest Angle, and One Piece Footplate. Weight Limit 125lbs.</i>	8,122.00	K0108

### CHOOSE ONE OF THE FOLLOWING:

Part Number	Description	Price	HCPCS Code
<input type="checkbox"/> I104529-99-0	No Tilt <i>Weight Limit 125lbs.</i>	No Charge	
<input type="checkbox"/> I104551-99-0	45° Power Tilt <i>Weight Limit 125lbs.</i>	6,526.00	E1002

## Seating System Options

### CHOOSE SEAT FUNCTION CONTROL TYPE:

Part Number	Description	Price	HCPCS Code
<input checked="" type="checkbox"/> I10585	Multiple Seat Function Control Kit for R-net <i>MUST be selected if ordering two or more seat functions. This device is used to operate two or more seat functions through the wheelchair electronics via the drive control and/or ICS switch box. Order with R-net electronics or remote joystick options. Please indicate below if you want push buttons or toggles on the ICS Switchbox.</i>	3,533.00	E2311

### CHOOSE PUSH BUTTONS OR TOGGLES:

Part Number	Description	Price	HCPCS Code
<input type="checkbox"/> I104135-99-0	Push Buttons on ICS Switchbox (3m)	No Charge	
<input type="checkbox"/> I104136-99-0	Toggle Switches on ICS Switchbox (3m)	316.00	

### CHOOSE MOUNTING POSITION OF ICS SWITCHBOX:

Part Number	Description	Price	HCPCS Code
<input type="checkbox"/> I10626	Mount ICS Switchbox onto R-net Joystick Bracket <i>ICS Switchbox will be mounted to the same bracket as the R-net Joystick Module. If this item is selected, the ICS Switchbox will swing away with the Joystick Module.</i>	No Charge	
<input type="checkbox"/> I104651-99-0	Mount ICS Switchbox onto Armrest <i>Mount on Left: _____ Mount on Right: _____</i>	No Charge	
<input type="checkbox"/> I11050	Mount ICS Switchbox on Backrest <i>On the K450, the backrest mounting location for the ICS Switchbox is on the Push Handle bracket; therefore, you must also order I104576-99-0 (Push Handles - MX). Mount on Left: _____ Mount on Right: _____</i>	152.00	

### OPTIONAL ITEMS:

Part Number	Description	Price	HCPCS Code
<input type="checkbox"/> I105409-99-0	Do Not Install ICS Switchbox, Send as Part	No Charge	
<input type="checkbox"/> I104196-99-0	ICS Alternative Switchbox for R-net Electronics <i>Seat function interface box that accommodates up to 8 momentary switches. Provides an additional method for accessing the seat functions. 1/8" mini-mono switches not included.</i>	1,235.00	K0108
<input type="checkbox"/> I106743-99-0	Remote Stop System - All Electronics <i>Consists of a Radio Frequency (RF) handheld transmitter and an RF receiver mounted on the wheelchair. Custom charges may be required when ordering with a PS/PS Jr seat, a low backrest height or a pediatric wheelchair. If ordering on a K450, push handles must also be selected.</i>	1,221.00	E2321
<input type="checkbox"/> I10515	Memory Seat Program for K450 MX <i>If this is selected, the chair will be programmed to allow the user to store and save up to 3 different memory seat positions through the ICS Switchbox. On the K450 MX, ALL seat functions are capable of recalling a saved seat position. Please refer to the K450 Owner's Manual for more information on saving and recalling a memory seat position.</i>	No Charge	

# K450 MX

## Seating System Options

### Backrest Options

#### CHOOSE ONE OF THE FOLLOWING:

Part Number	Description	Price	HCPCS Code
	<b>MX STYLE BACKRESTS:</b> <i>Use this section to select MX Style Backrests. Back heights measured from the seat pan to the top of the backrest shell with the backrest set at 90°. Please be sure to pick both a backrest shell size and a backrest cushion style (Leatherette or Black Mesh Fabric). Note: A 16"W MX Backrest is NOT available.</i>		
<input type="checkbox"/> I104321-99-0	MX Backrest - 10"W x 14"H <i>Backrest can be adjusted to either 14" or 16" high. If ordering a headrest, backrest will be mounted in the 16" high position unless client measurements require a 14" high backrest, in which case, custom modifications will be made.</i>	No Charge	
<input type="checkbox"/> I104322-99-0	MX Backrest - 12"W x 16"H <i>Backrest can be adjusted to either 16" or 18" high.</i>	No Charge	
<input type="checkbox"/> I104323-99-0	MX Backrest - 14"W x 20"H <i>Backrest can be adjusted to either 20" or 22" high.</i>	No Charge	

#### CHOOSE ONE OF THE FOLLOWING:

Part Number	Description	Price	HCPCS Code
	<b>Please Select Cushion Style Below:</b>		
<input type="checkbox"/> I104309-99-0	MX/MF Backrest Cushion, 10"W x 14"T - Leatherette	No Charge	
<input type="checkbox"/> I104311-99-0	MX/MF Backrest Cushion, 12"W x 16"T - Leatherette	No Charge	
<input type="checkbox"/> I104313-99-0	MX/MF Backrest Cushion, 14"W x 20"T - Leatherette	No Charge	
<input type="checkbox"/> I104308-99-0	MX/MF Backrest Cushion, 10"W x 14"T - Black Mesh	No Charge	
<input type="checkbox"/> I104310-99-0	MX/MF Backrest Cushion, 12"W x 16"T - Black Mesh	No Charge	
<input type="checkbox"/> I104312-99-0	MX/MF Backrest Cushion, 14"W x 20"T - Black Mesh	No Charge	

#### CHOOSE ONE OF THE FOLLOWING:

Part Number	Description	Price	HCPCS Code
	<b>PS JR TUBULAR STYLE BACKRESTS</b> <i>Use this section to select PS Jr Tubular Style Backrest Canes. Note: A 10"W Tubular Backrest is NOT available.</i>		
<input type="checkbox"/> I103717-99-0	PS Jr Back Canes, Height & Angle Adjustable <i>Must also select PS Jr Style Armrests if ordering this backrest. Choose PS Jr Style Backrest Type (Metal Back or Omit Back) below. Backrest width should match seat width. Back height is adjustable from 16.5"T - 22.5"T, measured from the seat pan to the top of the back.</i>	No Charge	

#### CHOOSE ONE OF THE FOLLOWING:

Part Number	Description	Price	HCPCS Code
<input type="checkbox"/> I103719-99-0	PS Jr Solid Backrest - 12"W	No Charge	
<input type="checkbox"/> I103721-99-0	PS Jr Solid Backrest - 14"W	No Charge	
<input type="checkbox"/> I103723-99-0	PS Jr Solid Backrest - 16"W	No Charge	
<input type="checkbox"/> I10449	PS Jr - Omit Backrest	No Charge	

## Seat Sizing

#### CHOOSE ONE OF THE FOLLOWING:

Part Number	Description	Price	HCPCS Code
<input type="checkbox"/> I104317-99-0	MX Seat Size 10"W x 10"D <i>Seat depth can be adjusted to 10", 11", or 12". If ordering this seat size, the lowest allowed seat position when using the seat to floor function will be approximately 3.5" from the floor to the top of the seat pan.</i>	No Charge	
<input type="checkbox"/> I104318-99-0	MX Seat Size 12"W x 12"D <i>Seat depth can be adjusted to 12", 13", or 14". If ordering this seat size with MX Arms, the lowest allowed seat position when using the seat to floor function will be approximately 3.5" from the floor to the top of the seat pan. If ordering Flip Back Arms, the lowest allowed seat position will be 8".</i>	No Charge	

# K450 MX

## Seat Sizing

### CHOOSE ONE OF THE FOLLOWING:

Part Number	Description	Price	HCPCS Code
<input type="checkbox"/> I104319-99-0	MX Seat Size 14"W x 14"D <i>Seat depth can be adjusted to 14", 15", or 16". If ordering this seat size with MX Arms, the lowest allowed seat position when using the seat to floor function will be approximately 3.5" from the floor to the top of the seat pan. If ordering Flip Back Arms, the lowest allowed seat position will be 8".</i>	No Charge	
<input type="checkbox"/> I104320-99-0	MX Seat Size 16"W x 16"D <i>Seat depth can be adjusted to 16", 17", or 18". If ordering this seat size with MX Arms, the lowest allowed seat position when using the seat to floor function will be approximately 10" from the floor to the top of the seat pan. If ordering Flip Back Arms, the lowest allowed seat position will be 16".</i>	No Charge	

### OPTIONAL ITEMS:

Part Number	Description	Price	HCPCS Code
	Seat Cushions Below Are Approximately 2" Thick:		
<input type="checkbox"/> I306384-99-90-0	Leatherette Seat Cush MX - 10"W x 10"D (250x250)	572.00	
<input type="checkbox"/> I306385-99-90-0	Leatherette Seat Cush MX - 12"W x 12"D (300x300)	572.00	
<input type="checkbox"/> I306386-99-90-0	Leatherette Seat Cush MX - 14"W x 14"D (350x350)	572.00	
<input type="checkbox"/> I104703-99-0	Leatherette Seat Cush MX - 16"W x 16"D (400x400)	572.00	
<input type="checkbox"/> I104305-99-0	Black Mesh Seat Cushion MX/MF - 10"W x 10"D	576.00	
<input type="checkbox"/> I104306-99-0	Black Mesh Seat Cushion MX/MF - 12"W x 12"D	576.00	
<input type="checkbox"/> I104307-99-0	Black Mesh Seat Cushion MX/MF - 14"W x 14"D	576.00	
<input type="checkbox"/> I104702-99-0	Black Mesh Seat Cush MX - 16"W x 16"D (400x400)	572.00	

## Armrest Options

### CHOOSE ARMREST STYLE:

Part Number	Description	Price	HCPCS Code
<input type="checkbox"/> I10512	Height Adjustable Armrest Assembly MX Seat <i>If you are not selecting Tilt, you MUST order this item.</i>	334.00	E0973
<input type="checkbox"/> I10514	Height Adjustable Armrest Assembly MX Seat - NC <i>No Charge if ordering Power Tilt (I104551-99-0)</i>	No Charge	
<input type="checkbox"/> I10824	Flip Back, Height Adjustable Armrests <i>Mounted off the backrest canes. Includes: 45° of angle adjustability, rotational adjustment, 3" offset adjustment both medially and laterally, and height adjustment of 6" - 14" from the seat pan. In one inch increments, please indicate where you would like the armrests set.</i> <i>SET ARMRESTS AT (6" - 14"): _____</i>	376.00	E0973
<input type="checkbox"/> I10825	Flip Back, Height Adjustable Armrests - NC <i>No charge when ordering power tilt. In one inch increments, please indicate where you would like the armrests set.</i> <i>SET ARMRESTS AT (6" - 14"): _____</i>	No Charge	

### MX ARMREST PADS:

Part Number	Description	Price	HCPCS Code
<input type="checkbox"/> I104303-99-0	MX Armrest Pad Short Leatherette 8"L x 3.5"W	No Charge	
<input type="checkbox"/> I104304-99-0	MX Armrest Pad Long 12.5"Lx3.5"W	No Charge	

### FLIP BACK ARMREST PAD RIGHT:

Part Number	Description	Price	HCPCS Code
<input type="checkbox"/> I105100-99-0	4"x10" Corpus Armrest Pad Right - Leatherette	No Charge	
<input type="checkbox"/> I105104-99-0	4"x13" Corpus Armrest Pad Right - Leatherette	No Charge	
<input type="checkbox"/> I107212-99-0	4"x10" Corpus Armrest Pad Right - Black Fabric	No Charge	
<input type="checkbox"/> I106047-99-0	4"x13" Corpus Armrest Pad Right - Black Fabric	No Charge	
<input type="checkbox"/> I106486-99-0	4"x10" Corpus Arm w/ Stealth Gel Pad - Right	169.00	
<input type="checkbox"/> I106488-99-0	4"x13" Corpus Arm w/ Stealth Gel Pad - Right	174.00	
<input type="checkbox"/> I104951-99-0	3.5"x8" Unitrack Armrest Pad Right (C)	No Charge	
<input type="checkbox"/> I104950-99-0	3.5"x12.5" Unitrack Armrest Pad Right (G)	No Charge	

### FLIP BACK ARMREST PAD LEFT:

Part Number	Description	Price	HCPCS Code
-------------	-------------	-------	------------

# K450 MX

## Armrest Options

### FLIP BACK ARMREST PAD LEFT:

Part Number	Description	Price	HCPCS Code
<input type="checkbox"/> I105098-99-0	4"x10" Corpus Armrest Pad Left - Leatherette	No Charge	
<input type="checkbox"/> I105102-99-0	4"x13" Corpus Armrest Pad Left - Leatherette	No Charge	
<input type="checkbox"/> I107213-99-0	4"x10" Corpus Armrest Pad Left - Black Fabric	No Charge	
<input type="checkbox"/> I106046-99-0	4"x13" Corpus Armrest Pad Left - Black Fabric	No Charge	
<input type="checkbox"/> I106487-99-0	4"x10" Corpus Arm w/ Stealth Gel Pad - Left	169.00	
<input type="checkbox"/> I106489-99-0	4"x13" Corpus Arm w/ Stealth Gel Pad - Left	174.00	
<input type="checkbox"/> I104952-99-0	3.5"x8" Unitrack Armrest Pad Left (C)	No Charge	
<input type="checkbox"/> I104949-99-0	3.5"x12.5" Unitrack Armrest Pad Left (G)	No Charge	

## Legrest Options

### CHOOSE ONE OF THE FOLLOWING:

Part Number	Description	Price	HCPCS Code
<input type="checkbox"/> I104314-99-0	MX Footplate - Small (11.5" W x 8"D)	No Charge	
<input type="checkbox"/> I104315-99-0	MX Footplate - Large (14" W x 9.5"D)	No Charge	

### OPTIONAL ITEMS:

Part Number	Description	Price	HCPCS Code
<input type="checkbox"/> I103427-99-0	Footstrap Kit	194.00	
<input type="checkbox"/> I103104-99-0	Footstrap Kit - Padded	215.00	
<input type="checkbox"/> I104338-99-0	Small Anti-Skid Overlay for Footplate - K450	158.00	
<input type="checkbox"/> I104337-99-0	Large Anti-Skid Overlay for Footplate - K450	158.00	
<input type="checkbox"/> I107503-99-0	BP Ankle Huggers, Med Side Release, Adj (Pair)	280.00	

*This item will be sent as a part with the wheelchair. It will not be installed.*

## Accessories

<input type="checkbox"/> I104518-99-0	Bodypoint Non-Padded Hip Belt - Medium	No Charge	
<input type="checkbox"/> I105849-99-0	Bodypoint Padded Hip Belt - Extra Small	149.00	E0978
<input type="checkbox"/> I104520-99-0	Bodypoint Padded Hip Belt - Medium	149.00	E0978
<input type="checkbox"/> I100519-99-0	Positioning Belt Retractable	295.00	K0108

### HEADREST OPTIONS:

Part Number	Description	Price	HCPCS Code
<input type="checkbox"/> I100749-99-0	MF Headrest for Koala/PS Jr - Leatherette <i>Must also select PN I10552.</i>	332.00	E0955
<input type="checkbox"/> I104339-99-0	MF Headrest for MX Back - Leatherette <i>Must also select PN I10552. Also includes a UniTrack Bracket for Lateral Adjustment of the headrest. Please select appropriate UT Bracket below:</i> UT Bracket for Lat Head Adj (10"W MX Back): _____ UT Brack for Lat Head Adj (12-14"W MX Back): _____	332.00	E0955
<input type="checkbox"/> I10552	Adj/Removable Headrest Hardware	272.00	E1028
<input type="checkbox"/> I107499-99-0	Universal Headrest Adapter <i>Please Note: Recline is limited to 160° or less (depending on backrest height) when ordering a Universal Headrest Adapter on a 23" or 24" tall Corpus backrest.</i>	258.00	

### OPTIONAL ITEMS:

Part Number	Description	Price	HCPCS Code
<input type="checkbox"/> I106104-99-0	UT Hip Suppt (B) w/ Mesh Cover - High (Pair) <i>4.5"W x 3.5"T Pad w/ 4" Tall Mounting Bracket. Must also select PN I10550.</i>	280.00	E0956
<input type="checkbox"/> I106103-99-0	UT Hip Suppt (A) w/ Mesh Cover - High (Pair) <i>6"W x 5.5"T Pad w/ 4" Tall Mounting Bracket. Must also select PN I10550.</i>	280.00	E0956
<input type="checkbox"/> I10550	Adj Removable - Thigh Support Hardware (Pair)	535.00	E1028
<input type="checkbox"/> I104575-99-0	UT Trunk Suppt (B) for MX Back - 3.5"T x 4.5"D <i>Must also select the Swing Away Trunk Support Hardware (I10551).</i>	335.00	E0956
<input type="checkbox"/> I104574-99-0	UT Trunk Suppt (F) for MX Back - 3.5"T x 8.5"D <i>Must also select the Swing Away Trunk Support Hardware (I10551).</i>	335.00	E0956

# K450 MX

## Accessories

### OPTIONAL ITEMS:

Part Number	Description	Price	HCPCS Code
<input type="checkbox"/> I10551	Adj Removable - Trunk Support Hardware (Pair)	538.00	E1028
<input type="checkbox"/> I104576-99-0	Push Handles - K450 <i>Mounted off the chassis. Height and angle adjustable. Must be selected if ordering bag clips or an attendant control, or if mounting the ICS Switchbox on the backrest.</i>	435.00	
<input type="checkbox"/> I106210-99-0	RAM Self-Leveling Cup Holder	83.00	
<input type="checkbox"/> I106806-99-0	RAM X-Grip Phone Holder	83.00	
<input type="checkbox"/> I107100-99-0	RAM X-Grip 10" Tablet Holder w/Swing Away Mount <i>Compatible with most 10" screen tablets with or without a protective case. Includes adjustable swing away mount (I103476-99-0). Customization may be required if ordered with an Omni.</i>	415.00	E1028
<input type="checkbox"/> I107222-99-0	USB Charger <i>Provides up to 1.5A (at 5V) to charge most smartphones and tablets. The charger features an isolated circuit design for improved safety; and by pulling power from the ICS, prevents battery drainage when the chair is powered off. Includes mounting hardware and ICS connector. Please select armrest for mounting, keeping in mind that charger may need to be relocated. Mount on left armrest: _____ Mount on right armrest: _____</i>	202.00	
<input type="checkbox"/> I105805-99-0	Medical Necessities Bag Hooks <i>Please Note: Recline is limited to 165° or less (depending on backrest height) when ordering with 23" or 24" tall Corpus backrest.</i>	98.00	K0108
<input type="checkbox"/> I100598-99-0	Chest Harness, Position Belt w/Bracket <i>Chest harness kit includes a standard positioning belt.</i>	443.00	E0960
<input type="checkbox"/> ISH212	BP Shoulder Harness Medium, Front Pull <i>This item will be sent as a part with the wheelchair. It will not be installed.</i>	244.00	E0960
<input type="checkbox"/> I100578-99-0	Upper Extremity Support Small, Angle Adjustable <i>Must also select PN I10556.</i>	301.00	E0950
<input type="checkbox"/> I106863-99-0	Upper Extremity Support Angle Adjust - UniTrack <i>Side mounted tray that attaches to the UniTrack Accessory Rail. Dimensions are: 19.5"W x 13"D. Can only be ordered with chairs that offer the UniTrack System. Must also select PN I10556.</i>	365.00	E0950
<input type="checkbox"/> I10556	Adj Removable - Tray Hardware (1 Each)	272.00	E1028

## Stealth Accessories

### OPTIONAL ITEMS:

Part Number	Description	Price	HCPCS Code
<input type="checkbox"/> I106222-99-0	Stealth 3" x 5" Lateral Supports - PS Jr <i>Must also select I10723. Includes THIN Stealth hardware and specialized mounting bracket for Flip Back Armrest assembly. Only compatible with Flip Back Armrest hardware.</i>	415.00	E0956
<input type="checkbox"/> I10723	Adjustable Stealth Lateral Hardware (Pair)	562.00	E1028

## Alternative Drive Systems

### ALTERNATIVE PROPORTIONAL JOYSTICK SYSTEMS:

Part Number	Description	Price	HCPCS Code
<input type="checkbox"/> I10751	Mini-Joystick for Omni / IOM (9-pin connector) <i>Operating Force = 0.5 oz (13 g) Neutral to Full Deflection = 0.2" (5mm)</i>	3,187.00	E2312
<input type="checkbox"/> I103893-99-0	Compact Joystick for Omni / IOM (9-pin connector) <i>Operating Force = 7.8 oz (220 g) Neutral to Full Deflection = 1.1" (28mm) Joystick includes two 1/8" mono jacks for Mode and On/Off. Includes Standard Joystick Knob and 1.75" (42mm) Foam Ball Knob.</i>	1,559.00	E2373
<input type="checkbox"/> I103894-99-0	Compact Joystick LITE for Omni / IOM (9-pin conn) <i>Operating Force = 1.8 oz (50 g) Neutral to Full Deflection = 0.7" (17mm) Joystick includes two 1/8" mono jacks for Mode and On/Off. Includes Standard Joystick Knob and 1.75" (42mm) Foam Ball Knob.</i>	1,763.00	E2373

### ALT PROPORTIONAL JOYSTICK MOUNTING HARDWARE:

Part Number	Description	Price	HCPCS Code
	NOTE: <i>Please see the PERMOfix Mounting Solutions, within the Alt Drive Accessories section, for additional mounting options.</i>		



# K450 MX

## Alternative Drive Systems

### ALT PROPORTIONAL JOYSTICK MOUNTING HARDWARE:

Part Number	Description	Price	HCPCS Code
<input type="checkbox"/> I104690-99-0	Alt JS Swing Away Mount Kit - Armrest <i>Mount in front of Left Arm: _____ Mount in front of Right Arm: _____</i>	273.00	E1028
<input type="checkbox"/> I107135-99-0	Alt JS / Omni Combo Swing Away Mount Kit - Armrest <i>Used to mount Omni and Alternative Proportional Joystick to the same armrest. Includes Swing Away Omni/Tablet Mount, with a longer mounting rod, two PERMOfix modules, one 3" rod and one 4" rod. Mount in front of Left Arm: _____ Mount in front of Right Arm: _____</i>	516.00	E1028
<input type="checkbox"/> I104695-99-0	Alt JS Adjustable Mount Kit - Midline <i>Mount on Left Arm: _____ Mount on Right Arm: _____</i>	801.00	E1028
<input type="checkbox"/> I104143-99-0	Alt JS Adjustable Mount Kit - Chin Boom <i>Swing-Away / Removable Chin Boom mounted on right of backrest shell. Custom charge may be applied if ordering lower backrest heights.</i>	851.00	E1028
<input type="checkbox"/> I10752	Alt JS Adjustable Mount Kit - Haims Harness, Lrg <i>For neck width up to 7½".</i>	605.00	E1028
<input type="checkbox"/> I10753	Alt JS Adjustable Mount Kit - Haims Harness, Sml <i>For neck width up to 5½".</i>	605.00	E1028

### ALT PROPORTIONAL JOYSTICK ACCESSORIES:

Part Number	Description	Price	HCPCS Code
<input type="checkbox"/> I104904-99-0	Dual Switches for Compact Joystick - Fixed Mount <i>Includes two pushbutton switches mounted on rigid arms. Switches can be used for Mode &amp; On/Off. Switch operation force = 4.5oz (130 g) Mount at front of Joystick: _____ Mount at rear of Joystick: _____ Mount at left of Joystick: _____ Mount at right of Joystick: _____</i>	253.00	K0108
<input type="checkbox"/> I104641-99-0	Dual Switches for Compact Joystick -Gooseneck Mount <i>Includes two pushbutton switches mounted on adjustable gooseneck arms. Switches are positioned on Left and Right sides of Compact Joystick. Switches can be used for Mode &amp; On/Off. Switch operation force = 4.5oz (130 g)</i>	456.00	K0108

### SIP & PUFF SYSTEMS:

Part Number	Description	Price	HCPCS Code
<input type="checkbox"/> I10573	Sip & Puff System for Omni <i>Requires Omni Control Module. Cannot use this part with an Input Module. Includes breath tube kit.</i>	1,705.00	E2325,E2326
<input type="checkbox"/> I10574	Sip & Puff System with Fiber Optic Opti-Stop <i>Requires Omni Control Module. Cannot use this part with an Input Module. Includes breath tube kit.</i>	2,269.00	E2325,E2326
<input type="checkbox"/> I106807-99-0	Sip & Puff Kit for Total Control Head Array <i>For use with Total Control Head Array. Requires Omni Control Module. Cannot be used with an Input Module. Includes breath tube kit.</i>	879.00	E2325
<input type="checkbox"/> I10762	Adjustable Mounting Hrdw - Sip & Puff <i>Required when selecting PN I10573, I10574 or I106807-99-0.</i>	273.00	E1028

### PROPORTIONAL FOOT CONTROL:

Part Number	Description	Price	HCPCS Code
<input type="checkbox"/> I10767	Proportional Foot Control <i>Mount on Left: _____ Mount on Right: _____</i>	3,249.00	
<input type="checkbox"/> I10766	Adjustable Mounting Hrdw - Foot Control <i>Required when selecting PN I10767.</i>	273.00	E1028

## Total Control Head Array

### CONFIGURE HEAD ARRAY:

Part Number	Description	Price	HCPCS Code
<input type="checkbox"/> I106805-99-0	Permobil Head Array - 3-Switch System <i>Features an occipital pad sensor and two stainless steel arms configured with proximity sensors. The arms are mounted to the upper attachment points on the head array. Longer arms are also included to extend the sensor adjustment range. ESP recommended. Choose Omni Control Module (ID51154IR) separately.</i>	4,750.00	E2330

# K450 MX

## Total Control Head Array

### CONFIGURE HEAD ARRAY:

Part Number	Description	Price	HCPCS Code
<input type="checkbox"/> I11072	Program Head Array as a 5-Switch System <i>Select additional arms/sensors below and indicate the intended use of each sensor in the Notes &amp; Comments section. Choose Omni Control Module (ID51154IR) or Input Module (I10581) separately.</i>	No Charge	
CHOOSE OCCIPITAL PAD SIZE (REQUIRED):			
<input type="checkbox"/> I105920-99-0	Small Occipital Pad for Head Array - 4.25"W x 3"H	No Charge	
<input type="checkbox"/> I105921-99-0	Large Occipital Pad for Head Array - 6"W x 4"H	No Charge	
CHOOSE OCCIPITAL MOUNT (REQUIRED):			
<input type="checkbox"/> I10847	Fixed Mount for Head Array Occipital Pad	No Charge	
<input type="checkbox"/> I105936-99-0	Adj. Swivel Mount for Head Array Occipital Pad <i>Enables Occipital Pad to be angled up/down and left/right independent of the lateral arms.</i>	426.00	K0108
CHOOSE BACKREST MOUNTING HARDWARE (REQUIRED):			
<input type="checkbox"/> I105927-99-0	Stealth TWB480 Backrest Mounting Hdwr with Links <i>Please select the length of the vertical stem below: Standard length (10"): _____ Short length (6"): _____</i>	284.00	E1028
<input type="checkbox"/> I105941-99-0	Stealth TWB480 Backrest Mnt Hdwr with Flipback	459.00	E1028
<input type="checkbox"/> I105925-99-0	Stealth HMO475 MultiAxis Backrest Mounting Hdwr	252.00	E1028
<input type="checkbox"/> I106319-99-0	Stealth TWB480 Backrest Mnt Hdwr w SubOcptl Mount	493.00	E1028
<input type="checkbox"/> I106282-99-0	Stealth 600 SubOccipital Pad - 8.5" W x 2.75" H <i>Must select when PN I106319-99-0 is ordered.</i>	101.00	K0108
SELECT ADDITIONAL LEFT ARM PACKAGE (OPTIONAL):			
<input type="checkbox"/> I10976	Add Left Lateral Arm w Proximity Sensor <i>Choose Lateral Arm Adjustment Range: Short (5.5" - 10"): _____ Long (5.5" - 14"): _____</i>	612.00	K0108
<input type="checkbox"/> I10977	Add Left Lateral Arm w Egg Switch <i>Choose Lateral Arm Adjustment Range: Short (5.5" - 10"): _____ Long (5.5" - 14"): _____</i>	500.00	K0108
<input type="checkbox"/> I10978	Add Left Lateral Arm w Pad (No Sensor or Switch) <i>Choose Lateral Arm Adjustment Range: Short (5.5" - 10"): _____ Long (5.5" - 14"): _____</i>	408.00	K0108
SELECT ADDITIONAL RIGHT ARM PACKAGE (OPTIONAL):			
<input type="checkbox"/> I10979	Add Right Lateral Arm w Proximity Sensor <i>Choose Lateral Arm Adjustment Range: Short (5.5" - 10"): _____ Long (5.5" - 14"): _____</i>	612.00	K0108
<input type="checkbox"/> I10980	Add Right Lateral Arm w Egg Switch <i>Choose Lateral Arm Adjustment Range: Short (5.5" - 10"): _____ Long (5.5" - 14"): _____</i>	500.00	K0108
<input type="checkbox"/> I10981	Add Right Lateral Arm w Pad (No Sensor or Switch) <i>Choose Lateral Arm Adjustment Range: Short (5.5" - 10"): _____ Long (5.5" - 14"): _____</i>	408.00	K0108

## Alt Drive Systems Setup

### CHOOSE OMNI OR INPUT MODULE (REQUIRED):

Part Number	Description	Price	HCPCS Code
<input type="checkbox"/> ID51154IR	Omni Control Module for R-net w/ Infrared <i>Must also select mounting hardware, either I103476-99-0 or I107135-99-0. Please specify which type of alternative drive input you will be using if not ordering it from Permobil. Input Device: _____</i>	1,772.00	
<input type="checkbox"/> I103476-99-0	Adjustable Swing-Away Hardware - Omni/Tablet <i>Mount in front of Left Arm: _____ Mount in front of Right Arm: _____</i>	273.00	E1028

# K450 MX

## Alt Drive Systems Setup

### CHOOSE OMNI OR INPUT MODULE (REQUIRED):

Part Number	Description	Price	HCPCS Code
<input type="checkbox"/> I10581	Input Module (IOM) for R-net Control System <i>This device allows a proportional or switched drive input to be connected to an R-net Control System without the use of an Omni. It must be ordered in conjunction with an R-net Remote Joystick. It does NOT support any of the following input devices: Sip-n-Puff, 3-Switch / 4-switch Head Array or RIM control.</i>	1,254.00	K0108

### CHOOSE MOUNTING POSITION OF SEAT CONTROL SWITCHBOX:

Part Number	Description	Price	HCPCS Code
<input type="checkbox"/> I10080	Mount Seat Control Switchbox on Backrest	No Charge	

## Alt Drive Accessories

### PERMOFIX MOUNTING SOLUTIONS:

Part Number	Description	Price	HCPCS Code
<input type="checkbox"/> I106655-99-0	PERMOfix UniTrack Fixed Mount Kit <i>Includes two base joints, two 4" rods and one 8" rod.</i>	275.00	K0108
<input type="checkbox"/> I106656-99-0	PERMOfix UniTrack Removable Mount Kit <i>Includes quick release mount and thigh support hardware, two base joints and three different length rods.</i>	483.00	E1028
<input type="checkbox"/> I106657-99-0	PERMOfix UniTrack Swing Away Mount Kit <i>Includes PERMOfix Swing Away module, two base joints, one 2" rod and one 4" rod.</i>	691.00	E1028
<input type="checkbox"/> I105907-99-0	Proximity Sensor + Swivel Assembly <i>Only available with the Total Control Head Array. Swivel assembly allows sensor to be mounted with PERMOfix or directly to lateral arm of Total Control Head Array. Since Proximity Sensors require power, they CANNOT be used as On/Off switches.</i>	349.00	K0108
<input type="checkbox"/> I105919-99-0	Egg Switch + Swivel Mount Assembly <i>Swivel assembly allows switch to be mounted with PERMOfix or directly to lateral arm of Total Control Head Array.</i>	236.00	K0108
<input type="checkbox"/> I105912-99-0	Solid Oval Pad + Swivel Assembly <i>Swivel assembly allows pad to be mounted with PERMOfix or directly to lateral arm of Total Control Head Array.</i>	142.00	K0108
<input type="checkbox"/> I106675-99-0	Oval Plate + Swivel Assembly - No Switch <i>Blank plate provides mounting option for Micro Lite, Buddy Button or other switch options. Swivel assembly allows pad to be mounted with PERMOfix or directly to lateral arm of Total Control Head Array.</i>	125.00	K0108

### OPTIONAL ITEMS:

Part Number	Description	Price	HCPCS Code
<input type="checkbox"/> I101769-99-0	Micro Lite Switch Qty: _____	190.00	
<input type="checkbox"/> I101767-99-0	Buddy Button Qty: _____	161.00	
<input type="checkbox"/> I106293-99-0	Egg Switch Black Qty: _____	161.00	K0108
<input type="checkbox"/> I101771-99-0	Mini Cup Switch Qty: _____	161.00	
<input type="checkbox"/> I107500-99-0	Proximity Switch Qty: _____	637.00	
<input type="checkbox"/> I107495-99-0	Mini Joystick 5 Switch OC for Omni/IOM (9-pin)	690.00	E2312
<input type="checkbox"/> I107498-99-0	Touch Contact with CR2032 Battery Qty: _____	392.00	
<input type="checkbox"/> I10772	Dome Switch, Backrest Mnt w Gooseneck Qty: _____ <i>Mount on Left: _____ Mount on Right: _____</i> <i>Mounts at top of backrest shell. Custom charge may be applied if ordering with short backrest heights Kit Includes: One PERMOfix joint with Flat Plate, two PERMOfix joints with 45° swivel adjustment, one 6" rod, one 8" rod and one 9" gooseneck switch with 50" cable and mounting clamp.</i>	504.00	K0108
<input type="checkbox"/> I10846	Dome Switch, Backrest Swg Awy Mnt w Gsneck Qty: _____ <i>Mount on Left: _____ Mount on Right: _____</i> <i>w 90° Swing-Away Function. Mounts at top of back. Custom charge applies if ordered with short back. Kit Includes: One PERMOfix joint with Flat Plate, one PERMOfix joint with 45° swivel adjustment, two 6" rods, one 90° Swing-Away joint, one 9" gooseneck switch w 50" cable and mounting clamp.</i>	745.00	K0108
<input type="checkbox"/> I105935-99-0	3.5mm Ext. Cable for Sensor/Switch, 0.3m Long	32.00	K0108
<input type="checkbox"/> I105931-99-0	3.5mm Ext. Cable for Sensor/Switch, 1m Long	35.00	K0108

# K450 MX

## Alt Drive Accessories

### OPTIONAL ITEMS:

Part Number	Description	Price	HCPCS Code
<input type="checkbox"/> I105932-99-0	3.5mm Ext. Cable for Sensor/Switch, 2m Long	40.00	K0108
<input type="checkbox"/> I105933-99-0	3.5mm Ext. Cable for Sensor/Switch, 3m Long	44.00	K0108
<input type="checkbox"/> I105934-99-0	3.5mm Y-splitter Cable for Sensor/Switch	35.00	K0108

*Allows two Proximity Sensors/Mechanical Switches to connect to the same direction jack.*

Additional Order Instructions (for Permobil):

Notes & Comments (not for Permobil):

### Chair Order Policy:

You are highly encouraged to contact Permobil Sales with ANY questions regarding proper completion of our order forms at 1-800-736-0925. The order cannot be processed without the following:

- A client name or code, the client's measurements, the client's weight
- The order form must also be signed by a person authorized by your company to acknowledge the items selected. Permobil, Inc. is not responsible for configuration or size discrepancies resulting from customer errors on the order form.
- Order confirmations will be provided by Permobil summarizing the items selected. Please review this carefully.
- Once the order is shipped, any changes are subject to the returns policy and may be prohibited.

### Order Acknowledgement:

I, \_\_\_\_\_, am an agent of the medical equipment provider named on this order form and I have the authority to contract for the purchase of powered wheelchairs and related parts on behalf of said provider. I acknowledge that I have reviewed this order and that it is complete and accurate to the best of my knowledge. I further acknowledge that any changes to this order after submission of this order form are subject to the returns policy and may carry additional charges.

Signed: \_\_\_\_\_

Printed Name: \_\_\_\_\_

Title: \_\_\_\_\_

Date: \_\_\_\_\_

NOTE: All specifications and prices are subject to change without notice. Please note that prices displayed are only valid if ordered with the wheelchair. Any parts ordered for service or future alterations will carry different pricing and warranty. *The HCPCS codes provided are not intended to be billing or legal advice, rather our interpretation of the code definitions. Use of the codes does not ensure coverage or payment for the item. For coverage information, verify the policy of the appropriate payer.*