



ZR

HEAD TURNING
STYLE.



With the ZR's minimalist mono-tube frame, less truly is more. The clean aesthetic of the ZR turns heads while the titanium frame makes sure you have the strength and vibration-dampening ride to go with your beauty. Roll with style and strength.



TECHNICAL SPECIFICATIONS.

Frame Features	Titanium Mono-Tube Rigid TiFit Limited Adjustability
Frame Options	Black Edition SuperLite Edition Heavy Duty Frame Power Adaptable Frame Reinforced Frame
Transport Weight*	9.3 lb.
Weight Limit	265 lb. Standard 300 lb. Heavy Duty
Total Chair Width*	23"
Seat Width	10" - 20"
Front Seat Width (Seat Taper)	3 Options
Seat Depth	10" - 20"
Ergonomic Seat	7 Options
Front Seat Height	16" - 21"
Rear Seat Height	12.5" - 20"
Camber	0°, 2°, 4°, 6°, 8°, or 12°
Front Frame Angle	70°, 75°, 80°, 85°, or 90°
Footrest Width Options	Standard Taper or V-Front End
Center of Gravity Adjustment	2.25"
Front Wheels	23 Options
Rear Wheels	67 Options
Handrims	15 Options
Convenience	Transport Bracket Option Folding Angle Adjustable Seat Back
Aesthetics	6 Titanium Frame Finishes 16 Frame Colors 6 Color Anodize Packages 7 Ultrasuede Colors 8 Spinergy Wheel Colors
Warranty	Lifetime on Frame One Year on TiLite Components

*Weights and widths based on 16" x 16" chair with 0° camber, minimum rear wheel spacing and all standard features. Rear wheels are not included in transport weight. For reference, a pair of standard 24" rear wheels weighs 7.2 lb.



TIFIT. PERSONAL ERGONOMICS.

The difference between a TiFit® chair and all other chairs is that a TiFit chair is a wheeled prosthetic: Every frame is made to measure. The optimized fit of a TiFit chair brings unparalleled performance and function. You will feel the difference in everything you do.

CUSTOM CONFIGURATION.

Custom Configuration allows you to add options and components to your chair to meet your needs and to optimize your riding experience. This includes choices of wheels, handrims, tires, wheel locks, color, tattoos, and more. TiLite offers more options than any other manufacturer!



TiLite | The Ultimate Ride™ | © 2015 TiLite, LLC.
800.545.2266 | tilite.com | MS0022_RevB_ZR
TiLite wheelchairs are shown with available options.