

F3 Corpus

permobil



** Indicates a Required Field*

Dealer Information

* Contact: _____

Dealer Code: _____

* Dealer Name: _____

Address: _____

* City: _____

* State: _____

* Phone#: _____

* Fax#: _____

PO#: _____

Email Address: _____

Client Information

Permobil recommends that the client is evaluated by a certified rehab specialist.

* First Name: _____

* Last Name: _____

Diagnosis: _____

Funding Source: _____

Client Age: _____

Client Measurements

* Weight: _____

* Height: _____

A) Top of Shoulders: _____

B) Chest Depth: _____

C) Chest Width: _____

D) Seat Depth: _____


E) Top of Head: _____

F) Elbow to Hand: _____

G) Seat Pan to Elbow: _____

H) Hip Width: _____

* I) Knee to Foot: _____



Permobil Australia Pty Ltd.
Unit 3A / 6 Boundary Rd.
Northmead,
NSW 2152
Australia

www.permobil.com

Prices effective: April 1, 2017



Chair Model and Colors

CHOOSE BASE:

Part Number	Description	Price AUD
<input checked="" type="checkbox"/> I107157-99-0	F3 Base Corpus <i>Price Includes: FWD F3 Power Base, Four-Wheel Independent Suspension System, ESP (Enhanced Steering Performance), Tie-Down Hardware for Strap Systems, Charger, Flat Free or Pneumatic Tires and a Standard Positioning Belt. Max Speed = 10kmh, Weight Capacity = 136kg, STFH = 17.5".</i>	

CHOOSE SPEED PACKAGE:

Part Number	Description	Price AUD
<input type="checkbox"/> I107224-99-0	Set Max Speed to 10 kmh <i>Anti-Tippers are required when selecting this speed package.</i>	No Charge
<input type="checkbox"/> I107226-99-0	Set Max Speed to 8 kmh	No Charge
<input type="checkbox"/> I107330-99-0	Anti-Tippers for F3 / F5 <i>Required for user weight above 100kg or max speed above 8kmh.</i>	

CHOOSE DRIVE TIRES:

Part Number	Description	Price AUD
<input type="checkbox"/> I106999-99-0	Flat Free Tire 3"W x 14"T Black - Standard Tread <i>Designed for a smooth ride on indoor and outdoor surfaces. Includes flat free standard tread caster wheels.</i>	No Charge
<input type="checkbox"/> I107005-99-0	Pneumatic Tire 3"W x 14"T Black - Standard Tread <i>Designed for a smooth ride on indoor and outdoor surfaces. Includes flat free standard tread caster wheels.</i>	No Charge

OPTIONAL DRIVE TIRES:

Part Number	Description	Price AUD
<input type="checkbox"/> I107289-99-0	Pneumatic/Snow Tire 3"W x 14"T Black - Aggr Tread <i>Snow tires are designed for use in cold weather conditions. The nylon compound in these tires remains flexible at low temperatures, but wears more easily. Therefore, the snow tires (sold as wheel assemblies) are not recommended for year-long use and must be ordered along with standard drive wheels. Does not include caster wheels.</i>	

CHOOSE SHROUD/RIM COLOR:

Part Number	Description	Price AUD
<input type="checkbox"/> I107202-99-0	F3 Color Kit - Sunburst Orange	No Charge
<input type="checkbox"/> I107197-99-0	F3 Color Kit - Galactic Green	No Charge
<input type="checkbox"/> I107198-99-0	F3 Color Kit - Midnight Black	No Charge
<input type="checkbox"/> I107199-99-0	F3 Color Kit - Ocean Blue	No Charge
<input type="checkbox"/> I107200-99-0	F3 Color Kit - Power Pink	No Charge
<input type="checkbox"/> I107201-99-0	F3 Color Kit - Radiant Red	No Charge
<input type="checkbox"/> I107599-99-0	F3 Color Kit - Carbon Fiber	
<input type="checkbox"/> I107602-99-0	F3 Color Kit - Mossy Oak Break-Up Country	

OPTIONAL ITEMS:

Part Number	Description	Price AUD
<input checked="" type="checkbox"/> I10901	Australian Freight & Handling Fee	

Base Options

MUST CHOOSE EITHER SEAT ELEVATOR OR FIXED SEAT TUBE:

Part Number	Description	Price AUD
<input type="checkbox"/> I11110	Power Adjustable Seat Height - 12" Travel <i>Must also select Power Tilt when selecting this option.</i>	
<input type="checkbox"/> I11117	Seat Tube Fixed	No Charge



Base Options

SET SEAT HEIGHT:

Part Number	Description	Price AUD
<input type="checkbox"/> I107435-99-0	Set Seat Height to 17.5"	No Charge
<input type="checkbox"/> I107433-99-0	Set Seat Height at 18.5" <i>If ordering Power Adjustable Seat Height, 1" of elevator travel is used to set seat at 18.5".</i>	No Charge
<input type="checkbox"/> I107434-99-0	Set Seat Height at 19.5" <i>If ordering Power Adjustable Seat Height, 2" of elevator travel is used to set seat at 19.5".</i>	No Charge

MUST CHOOSE BATTERIES:

Part Number	Description	Price AUD
<input type="checkbox"/> I107031-99-0	Batteries, Grp 24 (72Ah) Sealed Gel, x2, Installed	
<input type="checkbox"/> I107032-99-0	Batteries, Grp 34 (60Ah) Sealed Gel, x2, Installed	

OPTIONAL ITEMS:

Part Number	Description	Price AUD
<input type="checkbox"/> I107459-99-0	Essentials Box with Reflectors <i>Please Note: The total backrest angle will be limited to 150° to avoid contact between backrest components and the Essentials Box.</i>	
<input type="checkbox"/> I107204-99-0	Essentials Box with Lights <i>Please Note: The total backrest angle will be limited to 150° to avoid contact between backrest components and the Essentials Box. Must also select Light Kit.</i>	
<input type="checkbox"/> I106981-99-0	LED Light Kit for F3 / F5	
<input type="checkbox"/> I10715	Compact Laptop Ventilator Tray - Corpus Seat <i>Accommodates Pulmonetics LTV Series and similar sized Laptop Ventilators. The tray is 13"T, adjustable in width from 8" - 13", and adjustable in depth from 2" - 6.5". Overall backrest angle limited to 150° when mounting this ventilator tray. Custom Option.</i>	
<input type="checkbox"/> I10787	Battery Tray for Laptop Ventilator Tray <i>Measures 11 3/8" W x 6 1/2" H x 1 5/8" D. Accommodates Pulmonetic Systems Lithium-Ion Battery Pack, CareFusion SprintPack Battery System, or Richardson PowerTech Vent Power Center. Custom Option.</i>	
<input type="checkbox"/> I10788	O2 Holder for Laptop Ventilator Tray <i>Mounts to Left or Right Side of Laptop Vent Tray. O2 Holder is 7" Tall. Width infinitely adjustable from 3.5" - 7". Depth infinitely adjustable from 3.75" - 5". Custom Option.</i>	
<input type="checkbox"/> ID51113	Programming Key for R-net - Permobil Version	

Joystick Options

PLEASE NOTE:

For additional input devices, refer to Alternative Drive Systems section.

<input type="checkbox"/> I105246-99-0	R-net LED Remote Joystick Module <i>This Joystick does NOT have a Display Screen and does NOT contain mono jacks. Limited to five profiles with this Joystick Module.</i>	No Charge
<input type="checkbox"/> I107778-99-0	Permobil Joystick Module <i>Designed exclusively for Permobil, this R-net color joystick features: control of up to 4 Bluetooth devices, Infrared control, programmable buttons and a programmable switch jack; all in a compact, durable Aluminum housing.</i>	
<input type="checkbox"/> I107829-99-0	3.5mm Y-Splitter Cable for PJSM <i>For use with the PJSM Mode/Profile jack. Allows 2 mechanical switches to control 4 separate functions. Requires PC Programming to assign the short hold and long hold commands to each switch. Switches sold separately.</i>	
<input type="checkbox"/> I105700-99-0	R-net Remote Color JS w/ Lights, Heavy Duty Gimbal <i>Heavy-Duty modification enables gimbal to hold up to extreme operating forces. Reinforced housing is attached to a standard Color LCD R-net Joystick. 3" long joystick knob has 3/4" diameter. Not compatible with alternative joystick handles.</i>	
<input type="checkbox"/> I103814-99-0	Compact Advanced Joystick for R-net <i>Compact sized joystick that does not contain an LCD Display. Contains two mono jacks for mode & on/off. Does NOT contain a wrist support pad. If you want a wrist support pad, please select I10621 below.</i>	



Joystick Options

		Price AUD
<input type="checkbox"/>	I10621 Add Wrist Pad to Compact Joystick Advanced	
CHOOSE JOYSTICK HANDLE TYPE:		
Part Number	Description	Price AUD
<input type="checkbox"/>	I10045 Cone Shaped Joystick Knob - Standard	No Charge
<input type="checkbox"/>	I100960-99-0 Chin Cup For Joystick Knob	
<input type="checkbox"/>	I100961-99-0 "T" Handle For Joystick Knob	
<input type="checkbox"/>	I100962-99-0 "Large Ball" For Joystick Knob	
<input type="checkbox"/>	I100963-99-0 "Softball" For Joystick Knob	
<input type="checkbox"/>	I100964-99-0 "Mushroom" For Joystick Knob	
<input type="checkbox"/>	I100965-99-0 "Stick" For Joystick Knob	
<input type="checkbox"/>	I105442-99-0 1.75" Black Foam Ball For Joystick Knob	
<input type="checkbox"/>	I105846-99-0 Bodypoint J/S Handle 3" U-Shaped w/ Flex-Shaft	
<input type="checkbox"/>	I105847-99-0 Bodypoint J/S Handle 4" U-Shaped w/ Flex-Shaft	
<input type="checkbox"/>	I105848-99-0 Bodypoint J/S Handle Dome Shaped, Rubber	
CHOOSE JOYSTICK MOUNT:		
Part Number	Description	Price AUD
<input type="checkbox"/>	I10664 Joystick Mount Right - Fixed	No Charge
<input type="checkbox"/>	I10665 Joystick Mount Left - Fixed	No Charge
<input type="checkbox"/>	I102537-99-0 Joystick Mount Right - Swing Away	
<input type="checkbox"/>	I102538-99-0 Joystick Mount Left - Swing Away	
<input type="checkbox"/>	I103820-99-0 Retractable Joystick Mount Right	
<input type="checkbox"/>	I103819-99-0 Retractable Joystick Mount Left	
<input type="checkbox"/>	I105210-99-0 Mount Standard Joystick on Corpus Back Hinge <i>Standard R-net Joystick is mounted to backrest hinge. Can be chosen instead of PN I105408-99-0. Mount on Left: _____ Mount on Right: _____</i>	
<input type="checkbox"/>	I105225-99-0 Co-Pilot Attendant Control R-net - Corpus	
<input type="checkbox"/>	I105408-99-0 Attendant Control R-net	
<input type="checkbox"/>	I107350-99-0 R-net Bluetooth iDevice Module <i>Allows input device commands (such as joystick nudges or switch clicks) to be used to operate an iPhone, iPad or iPod. By assigning each input to an iOS function, users can make phone calls, send texts and emails and access applications. Kit includes a Bluetooth iDevice Module, an R-net Bus Cable, and two mono jack cables for direct switch inputs (switches not included). Must be ordered with a Color JSM or Omni for mode display.</i>	
<input type="checkbox"/>	I107367-99-0 R-net Bluetooth Mouse Module <i>Enables mouse control through any R-net input device. It can operate the mouse on a computer and on select Android devices (with R-net for Android app). Kit includes a Bluetooth Mouse Module, an R-net Bus Cable, two mono jack cables for L/R click (switches not included) and a USB BT Adapter for user's computer. Note: BT Adapter may not work with all computers. Must be ordered with a Color JSM or Omni for mode display.</i>	
<input type="checkbox"/>	I107368-99-0 Add a second Bluetooth Mouse Module <i>Allows R-net to simultaneously pair to a second compatible computer or Android device. Each module is accessible through a separate R-net mode. Kit includes a unique Bluetooth Mouse Module, an R-net Bus Cable, two mono jack cables for L/R click and a USB BT Adapter.</i>	
<input type="checkbox"/>	I10582 Output Module (IOM) for R-net Control System <i>This allows a user to control an external device (ex: Communication Device, ECU, Door Opener, etc.) through the drive input (JSM, Head Array, etc.). Each module can provide up to 4 switch closures. May require specialty cable or wiring. Up to three output modules can be installed on one wheelchair.</i>	

Seating System



Seating System

CHOOSE POSTERIOR TILT/RECLINE PACKAGE:

Part Number	Description
	PLEASE NOTE: <i>A half speed limitation is imposed when the Total Backrest Angle (TBA) > 150°. TBA = Tilt Angle + Backrest Recline Angle. Example 1: 45° Tilt + 100° Recline = 145° TBA (Full Speed). Example 2: 35° Tilt + 120° Recline = 155° TBA (Half Speed).</i>
<input type="checkbox"/> I11130	Corpus Seat w/ No Tilt & Manual Recline <i>Manual Recline angle limited to 120°.</i>
<input type="checkbox"/> I11131	Corpus Seat w/ No Tilt & 175° Power Recline
<input type="checkbox"/> I11132	Corpus Seat w/ 50° Tilt & Manual Recline <i>Manual Recline angle limited to 120°.</i>
<input type="checkbox"/> I11133	Corpus Seat w/ 50° Tilt & 175° Power Recline

ACTIVE REACH PACKAGES:

Part Number	Description	Price AUD
	PLEASE NOTE: <i>Active Reach should be recommended by a licensed clinical professional and/or ATP. While operating Active Reach, positioning equipment such as lap belts and chest supports should be considered to meet the user's specific safety needs. Note: A low speed limitation is imposed from 5° to 10° of Active Reach, an extra low speed limitation is imposed from 10° to 20° of Active Reach, and a drive lockout is imposed above 20° of Active Reach.</i>	
<input type="checkbox"/> I11141	5° Active Reach <i>When ordered with a Power Adjustable Seat Height, tilt range will be -5° to +50°. When ordered without a Power Adjustable Seat Height, tilt range will be -5° to +45°. This option will raise the minimum seat to floor height by about 1" when the chair is at 0° of tilt and ordered with Power Tilt and a Fixed Seat Tube.</i>	No Charge
<input type="checkbox"/> I107090-99-0	10° Active Reach <i>Please Note: This option is only available on configurations ordered with Power Adjustable Seat Height and Power Tilt. A chest support belt (see Accessories section) is strongly recommended when selecting this option.</i>	No Charge
<input type="checkbox"/> I107091-99-0	20° Active Reach <i>Please Note: This option is only available on configurations ordered with Power Adjustable Seat Height, Power Tilt, Power Recline and Power Legrest Elevation. Must select Corpus VS Power Elevating Legrest, Corpus VS Footplates and Corpus VS Knee Support from VS Legrest Options section. A chest support belt (see Accessories section) is strongly recommended when selecting this option. Not compatible with retractable positioning belts.</i>	
<input type="checkbox"/> I11100	30° Active Reach w/Transfer <i>Please Note: This option is only available on configurations ordered with Power Adjustable Seat Height, Power Tilt, Power Recline and Power Legrest Elevation. Must select Corpus VS Power Elevating Legrest, Corpus VS Footplates and Corpus VS Knee Support from VS Legrest Options section. A chest support belt (see Accessories section) is strongly recommended when selecting this option. Not compatible with retractable positioning belts.</i>	

Seating System Options

<input checked="" type="checkbox"/> I10585	Seat Function Control Kit for R-net <i>Please indicate below if you want push buttons or toggles on the ICS Switchbox.</i>
--	---

CHOOSE PUSH BUTTONS OR TOGGLES:

Part Number	Description	Price AUD
<input type="checkbox"/> I102798-99-0	Push Buttons on ICS Switchbox	No Charge
<input type="checkbox"/> I102799-99-0	Toggle Switches on ICS Switchbox	
<input type="checkbox"/> I105409-99-0	Do Not Install ICS Switchbox, Send as Part	No Charge

CHOOSE MOUNTING POSITION OF ICS SWITCHBOX:

Part Number	Description	Price AUD
<input type="checkbox"/> I10626	Mount ICS Switchbox onto R-net Joystick Bracket <i>ICS Switchbox will be mounted to the same bracket as the R-net Joystick Module. If this item is selected, the ICS Switchbox will swing away with the Joystick Module.</i>	No Charge



Seating System Options

CHOOSE MOUNTING POSITION OF ICS SWITCHBOX:

Part Number	Description	Price AUD
<input type="checkbox"/> I104651-99-0	Mount ICS Switchbox onto Armrest <i>Mount on Left: _____ Mount on Right: _____</i>	No Charge
<input type="checkbox"/> I106256-99-0	Mount ICS Switchbox on Backrest <i>Mount on Left: _____ Mount on Right: _____</i>	

OPTIONAL ITEMS:

Part Number	Description	Price AUD
<input type="checkbox"/> I107705-99-0	Virtual Seating Coach Package <i>The Virtual Seating Coach (VSC) Package includes a Bluetooth dongle, a USB charger, a RAM X-grip Phone Holder and access to VSC coaching and regimen management. The VSC app can be downloaded for free from the App Store or Google Play. Requires an iPhone (iOS 9 or higher) or an Android phone (OS 4. 1. 2 or higher).</i> Mount VSC to Left Armrest: _____ <i>To provide access, USB Charger and Phone Holder may require relocation.</i> Mount VSC to Right Armrest: _____ <i>To provide access, USB Charger and Phone Holder may require relocation.</i>	
<input type="checkbox"/> I104196-99-0	ICS Alternative Switchbox for R-net Electronics <i>Seat function interface box that accommodates up to 8 momentary switches. Provides an additional method for accessing the seat functions. 1/8" mini-mono switches not included.</i>	
<input type="checkbox"/> I107088-99-0	Remote Stop System - F/M Series <i>Consists of a Radio Frequency (RF) handheld transmitter and an RF receiver mounted on the wheelchair.</i>	
<input type="checkbox"/> I10750	Memory Seat Program, Corpus w/ R-net	
<input type="checkbox"/> I106209-99-0	Legrest Movement Linked to Recline Electronically <i>Power Legrest will raise when Power Recline moves backward and will lower when back moves forward. If Power Legrest is operated, Power Recline will NOT automatically reposition. Feature works through ICS Switchbox and Driver Control.</i>	
<input type="checkbox"/> I10898	Enable Independent Repositioning Mode <i>IRM allows for sequential movement of power seat functions into a posterior position (35° tilt, 115° recline, 135° legrest) and returns to upright position (5° tilt, 100° recline, 110° legrest). This can be accessed using a standard drive control or a single switch. On the Permobil Joystick Module, IRM is accessible through the mode/profile jack (long switch hold). Please see the Alt Drive Accessories section for switch options.</i>	No Charge

Backrest Options

CHOOSE FIXED OR SLIDING BACKREST KIT:

Part Number	Description	Price AUD
<input type="checkbox"/> I105067-99-0	Sliding Backrest Kit - 5" Mechanical Shear	No Charge
<input type="checkbox"/> I105066-99-0	Fixed Backrest Kit - No Mechanical Shear	

CHOOSE BACKREST TYPE:

Part Number	Description
	BACKREST HEIGHT NOTE: <i>Backrest heights are measured from the seat pan to the top of the backrest shell with the backrest set at 90°.</i>
<input type="checkbox"/> I106545-99-0	Corpus Ergo Back - 14"W x 20"T (C3GERGBCA) <i>Power Recline will be limited to 150° when ordering this backrest.</i>
<input type="checkbox"/> I106546-99-0	Corpus Ergo Back - 16"W x 20"T (C3GERGBCA) <i>Power Recline will be limited to 150° when ordering this backrest.</i>
<input type="checkbox"/> I106547-99-0	Corpus Ergo Back - 18"W x 20"T (C3GERGBCA) <i>Power Recline will be limited to 150° when ordering this backrest.</i>
<input type="checkbox"/> I106548-99-0	Corpus Ergo Back - 20"W x 20"T (C3GERGBCA) <i>Power Recline will be limited to 150° when ordering this backrest.</i>
	PLEASE NOTE:



Backrest Options

CHOOSE BACKREST TYPE:

Part Number	Description
PLEASE NOTE:	
<i>All taller backrests listed below include: lateral wedges, lumbar pads, a UT rail mounted at the top for lateral adjustment of the headrest, and split backrest feature with height adjustment in 1" increments. Please indicate below where you want the height set at so that a properly sized backrest cushion can be provided.</i>	
<input type="checkbox"/> I106532-99-0	Corpus Ergo Back - 14"W x 23" - 28"T (C3GERGBCA)
<input type="checkbox"/> I106533-99-0	Corpus Ergo Back - 16"W x 23" - 28"T (C3GERGBCA)
<input type="checkbox"/> I106534-99-0	Corpus Ergo Back - 18"W x 23" - 28"T (C3GERGBCA)
<input type="checkbox"/> I106535-99-0	Corpus Ergo Back - 20"W x 23" - 28"T (C3GERGBCA)
	Set Backrest at 23" _____
	Set Backrest at 24" _____
	Set Backrest at 25" _____
	Set Backrest at 26" _____
	Set Backrest at 27" _____
	Set Backrest at 28" _____

CHOOSE BACKREST CUSHION:

Part Number	Description	Price AUD
<input type="checkbox"/> I10659	Leatherette Cushion - Corpus Ergo Back	No Charge
<input type="checkbox"/> I10660	Black Mesh Cushion - Corpus Ergo Back	
<input type="checkbox"/> I105117-99-0	Omit Backrest Cushion	No Charge

OPTIONAL ITEMS:

Part Number	Description	Price AUD
<input type="checkbox"/> I105782-99-0	Tubular Backrest Frame - Corpus (1" Canes) <i>Provides 1" diameter, vertical canes for mounting aftermarket backrests. Frame can be adjusted to support backrest widths from 14" to 20". Canes are attached with 2" long UniTrack linkage disks that can rotate to provide space for backrest hardware. The tops of the canes can be adjusted up to a height of 27" above the seat pan. Depending on the aftermarket backrest, several inches of seat depth may be lost when mounting this system.</i>	
<input type="checkbox"/> I106072-99-0	Direct Backrest Frame - Corpus <i>Used to mount ROHO Agility, Jay 3, Matrx Elite, Varilite Icon and Ride Java backrests ranging in width from 14" to 20". Slotted brackets and 1" UniTrack linkage disks provide a close, adjustable mount. Backrest pitch can be tuned by rotating the linkage disks. Some seat depth may be lost due to the position and orientation of the backrest.</i>	
<input type="checkbox"/> I10725	Split Backrest Cushion - Corpus <i>Custom Option.</i>	
<input type="checkbox"/> I10786	Tailored Backrest Cushion - Corpus <i>Sides of cushion are trimmed to provide additional clearance for lateral trunk supports. Must also select a backrest cushion (I10659 or I10660) when ordering this custom option.</i>	

Seat Sizing

CHOOSE SEAT WIDTH:

Part Number	Description	Price AUD
<input type="checkbox"/> I107023-99-0	Set Seat Width at 17": _____	No Charge
<input type="checkbox"/> I107024-99-0	Set Seat Width at 19": _____	No Charge
<input type="checkbox"/> I107025-99-0	Set Seat Width at 21": _____	No Charge
<input type="checkbox"/> I107026-99-0	Set Seat Width at 23": _____	No Charge

CHOOSE SEAT DEPTH:

Part Number	Description
-------------	-------------



Seat Sizing

CHOOSE SEAT DEPTH:

Part Number	Description	
	PLEASE NOTE:	
	<i>There are three seat depth kits, each providing three seat depth settings: (1) 14" - 16", (2) 17" - 19" and (3) 20" - 22". Please select the required seat depth and the appropriate kit will be provided.</i>	
<input type="checkbox"/> I107013-99-0	Seat Depth 14"	No Charge
<input type="checkbox"/> I107014-99-0	Seat Depth 15"	No Charge
<input type="checkbox"/> I107015-99-0	Seat Depth 16"	No Charge
<input type="checkbox"/> I107016-99-0	Seat Depth 17"	No Charge
<input type="checkbox"/> I107017-99-0	Seat Depth 18"	No Charge
<input type="checkbox"/> I107018-99-0	Seat Depth 19"	No Charge
<input type="checkbox"/> I107019-99-0	Seat Depth 20"	No Charge
<input type="checkbox"/> I107020-99-0	Seat Depth 21"	No Charge
<input type="checkbox"/> I107021-99-0	Seat Depth 22"	No Charge

CHOOSE SEAT CUSHION:

Part Number	Description	
	NOTE:	
	<i>See supplemental order form for ROHO seat cushion options. A Corpus II Ergonomic Seat Cushion will be provided for select seat sizes.</i>	
<input type="checkbox"/> I10651	Corpus Ergo Cushion - Leatherette (C3GERGSCA)	
<input type="checkbox"/> I10652	Corpus Ergo Cushion - Black Mesh (C3GERGSCA)	
<input type="checkbox"/> I10965	Corpus Ergo Air Cushion (C3GRHOSCA)	
	<i>Features ROHO Dry Floatation technology. The removable ROHO cell pad provides adjustable skin protection while the contoured foam and thigh support wedges help maintain proper positioning and stability. Includes removable stretch fabric cover. Available in widths of 17", 19" and 21" and depths of 17", 19" and 21", as well as 17" x 23".</i>	
<input type="checkbox"/> I11113	Incontinence Cover for Ergo Cushion	
	<i>Made from a stretch Dartex fabric and suitable for all Ergo cushions including the Corpus Ergo Air. The cover is cut to the contours of the Ergo cushion and is designed to be used as an outer layer of protection. Cover size will be matched to cushion size. Sent as a part.</i>	
<input type="checkbox"/> I105118-99-0	Omit Seat Cushion	No Charge

Armrest Options

CHOOSE ARMBAR KIT AND MOUNTING POSITION:

Part Number	Description	
	PLEASE NOTE:	
	<i>Narrow Armbar Kit can be adjusted from 15"W to 19"W between the armrest pads without any additional parts. Wide Armbar Kit can be adjusted from 19"W to 23"W between the armrest pads without any additional parts. Please select both Armbar Kit and Width Setting below.</i>	
<input type="checkbox"/> I106543-99-0	Narrow Armbar Kit, Inf Adjustable 15"W - 19"W	No Charge
<input type="checkbox"/> I106544-99-0	Wide Armrest Kit, Inf Adjustable 19"W - 23"W	No Charge
	<i>Required when selecting a 20" wide backrest.</i>	
	Set the Width Between the Armrests to: _____	
	<i>Must be wider than backrest width to ensure proper fit. Please round to nearest 1".</i>	

CHOOSE ARMREST HEIGHT MOUNTING POSITION:

Part Number	Description	Price AUD
<input type="checkbox"/> I105069-99-0	Height Adjustable Armrests - NC	No Charge
	<i>No charge when ordering power tilt and/or recline.</i>	
<input type="checkbox"/> I105081-99-0	Height Adjustable Armrests	



Armrest Options

CHOOSE ARMREST HEIGHT MOUNTING POSITION:

Part Number **Description**

SET ARMREST HEIGHT AT: _____

Corpus Armrests are infinitely adjustable from 7" - 15", measured from the seat pan to the top of the armrest pad. No additional parts are needed to get this full range of adjustment. In one inch increments (ex: 8", 9", etc.), please indicate where you would like the armrests set.

CHOOSE ARMREST SWIVEL LOCKED OR ENABLED:

Part Number **Description**

ARMREST MEDIAL SWIVEL NOTE:

Every Corpus Armrest can flip backwards for transfers. The armrests also include a feature that allows the armrest pads to swivel medially towards the center of the seating system. This swivel feature can be locked out by simply moving a bolt located under the armrest bracket. Please indicate below whether each armrest should have medial swivel enabled or locked. Note: Medial swivel must be locked on Corpus VS or when ordering arm trough adapters.

<input type="checkbox"/>	I10645	Left Armrest Bracket - Medial Swivel Locked	No Charge
<input type="checkbox"/>	I10646	Right Armrest Bracket - Medial Swivel Locked	No Charge
<input type="checkbox"/>	I105073-99-0	Left Armrest Bracket - Medial Swivel Enabled	No Charge
<input type="checkbox"/>	I105072-99-0	Right Armrest Bracket - Medial Swivel Enabled	No Charge

CHOOSE ARMREST PAD RIGHT:

Part Number **Description**

Price AUD

<input type="checkbox"/>	I105100-99-0	4"x10" Corpus Armrest Pad Right - Leatherette	No Charge
<input type="checkbox"/>	I105104-99-0	4"x13" Corpus Armrest Pad Right - Leatherette	No Charge
<input type="checkbox"/>	I105105-99-0	4"x16" Corpus Armrest Pad Right - Leatherette	No Charge
<input type="checkbox"/>	I105999-99-0	4"x18" Corpus Armrest Pad Right - Leatherette	No Charge
<input type="checkbox"/>	I107212-99-0	4"x10" Corpus Armrest Pad Right - Black Fabric	No Charge
<input type="checkbox"/>	I106047-99-0	4"x13" Corpus Armrest Pad Right - Black Fabric	No Charge
<input type="checkbox"/>	I107214-99-0	4"x16" Corpus Armrest Pad Right - Black Fabric	No Charge
<input type="checkbox"/>	I105995-99-0	4"x18" Corpus Armrest Pad Right - Black Fabric	No Charge
<input type="checkbox"/>	I106486-99-0	4"x10" Corpus Arm w/ Stealth Gel Pad - Right	
<input type="checkbox"/>	I106488-99-0	4"x13" Corpus Arm w/ Stealth Gel Pad - Right	
<input type="checkbox"/>	I106490-99-0	4"x16" Corpus Arm w/ Stealth Gel Pad - Right	
<input type="checkbox"/>	I106281-99-0	4"x18" Corpus Arm w/ Stealth Gel Pad - Right	

CHOOSE ARMREST PAD LEFT:

Part Number **Description**

Price AUD

<input type="checkbox"/>	I105098-99-0	4"x10" Corpus Armrest Pad Left - Leatherette	No Charge
<input type="checkbox"/>	I105102-99-0	4"x13" Corpus Armrest Pad Left - Leatherette	No Charge
<input type="checkbox"/>	I105107-99-0	4"x16" Corpus Armrest Pad Left - Leatherette	No Charge
<input type="checkbox"/>	I105998-99-0	4"x18" Corpus Armrest Pad Left - Leatherette	No Charge
<input type="checkbox"/>	I107213-99-0	4"x10" Corpus Armrest Pad Left - Black Fabric	No Charge
<input type="checkbox"/>	I106046-99-0	4"x13" Corpus Armrest Pad Left - Black Fabric	No Charge
<input type="checkbox"/>	I107215-99-0	4"x16" Corpus Armrest Pad Left - Black Fabric	No Charge
<input type="checkbox"/>	I105994-99-0	4"x18" Corpus Armrest Pad Left - Black Fabric	No Charge
<input type="checkbox"/>	I106487-99-0	4"x10" Corpus Arm w/ Stealth Gel Pad - Left	
<input type="checkbox"/>	I106489-99-0	4"x13" Corpus Arm w/ Stealth Gel Pad - Left	
<input type="checkbox"/>	I106491-99-0	4"x16" Corpus Arm w/ Stealth Gel Pad - Left	
<input type="checkbox"/>	I106303-99-0	4"x18" Corpus Arm w/ Stealth Gel Pad - Left	

OPTIONAL ITEMS:

Part Number **Description**

Price AUD



Armrest Options

OPTIONAL ITEMS:

Part Number	Description	Price AUD
<input type="checkbox"/> I105463-99-0	Corpus Arm Trough Adapter 1" - Pair <i>For aftermarket armrest pads. Please indicate below if they should be installed or on chair sent as parts. Please Note: Medial Swivel must be must be locked (I10645 & I10646) when selecting this option.</i>	
<input type="checkbox"/> I105462-99-0	Corpus Arm Trough Adpater 7/8" - Pair <i>For aftermarket armrest pads. Please indicate below if they should be installed or on chair sent as parts. Please Note: Medial Swivel must be must be locked (I10645 & I10646) when selecting this option.</i> Install Arm Trough Adapters: _____ Send Arm Trough Adapters as Parts: _____	
<input type="checkbox"/> I107339-99-0	Armrest Pouch - Right	
<input type="checkbox"/> I107338-99-0	Armrest Pouch - Left	
<input type="checkbox"/> I107840-99-0	UT Armrest Side Support - Right Outer <i>Includes 2"T x 6.5"W tapered pad.</i>	
<input type="checkbox"/> I107863-99-0	UT Armrest Side Support - Left Outer <i>Includes 2"T x 6.5"W tapered pad.</i>	
<input type="checkbox"/> I107864-99-0	UT Armrest Side Support - Right Inner <i>Includes 2"T x 6.5"W tapered pad. May require an increase to the width between the armrests to ensure armrests will flip-up.</i>	
<input type="checkbox"/> I107865-99-0	UT Armrest Side Support - Left Inner <i>Includes 2"T x 6.5"W tapered pad. May require an increase to the width between the armrests to ensure armrests will flip-up.</i>	
<input type="checkbox"/> I107839-99-0	UT Elbow Support - Right <i>Must also select PN I10881.</i>	
<input type="checkbox"/> I107838-99-0	UT Elbow Support - Left <i>Must also select PN I10882.</i>	
<input type="checkbox"/> I10881	Adjustable Elbow Support Hardware - Right	
<input type="checkbox"/> I10882	Adjustable Elbow Support Hardware - Left	

Legrest Options

CORPUS LEGREST OPTIONS:

Part Number	Description	Price AUD
<input type="checkbox"/> I105080-99-0	CHOOSE LEGREST TYPE: Corpus Manual Elevating Legrest	No Charge
<input type="checkbox"/> I105079-99-0	Corpus Power Elevating Legrest - 170° <i>Includes mechanical articulation.</i>	
<input type="checkbox"/> I105437-99-0	CHOOSE FOOTPLATE STYLE: 2pc Wide Corpus Footplates - 9"D x 7.5"W	No Charge
<input type="checkbox"/> I105436-99-0	2pc Narrow Corpus Footplates - 9"D x 6"W	
<input type="checkbox"/> I106760-99-0	1pc Narrow Corpus Footpl - 9.75"D x 12.5"W <i>Footplate is 12.5"W in front, 8"W in back.</i>	
<input type="checkbox"/> I105252-99-0	CALF SUPPORT OPTIONS: UT Calf Support Kit (B) w Mesh Cover & Adj Hardware <i>3.5"W x 4.5"T Pad.</i>	No Charge
<input type="checkbox"/> I105258-99-0	UT Calf Support Kit (J) w Mesh Cover & Adj Hardware <i>5"W x 6.5"T Pad.</i>	No Charge
<input type="checkbox"/> I105250-99-0	UT Calf Support Kit (A) w Mesh Cover & Adj Hardware <i>5.5"W x 6"T Pad.</i>	No Charge



Legrest Options

CORPUS LEGREST OPTIONS:

Part Number	Description	Price AUD
<input type="checkbox"/> I105256-99-0	UT Calf Support Kit (H) w Mesh Cover & Adj Hardware 7"W x 7"T Curved Pad.	
<input type="checkbox"/> I105254-99-0	UT Calf Support Kit (D) w Mesh Cover & Adj Hardware 8"W x 8"T Pad.	

CORPUS VS LEGREST OPTIONS (FOR SELECT ACTIVE REACH PACKAGES):

Part Number	Description	Price AUD
CORPUS VS LEGREST AND KNEE SUPPORT:		
<i>NOTE: The Corpus VS Legrest is only available with 20°, 30° and 45° of Active Reach.</i>		
<input type="checkbox"/> I107186-99-0	Corpus VS Power Elevating Legrest - 180° <i>Includes power articulation.</i>	
<input type="checkbox"/> I107038-99-0	Knee Support - Corpus VS <i>Select when ordering Corpus VS Legrest. Must also select I11129.</i>	
<input type="checkbox"/> I11129	Adj Removable Knee Support Hardware (1 Each)	
CHOOSE FOOTPLATE STYLE:		
<input type="checkbox"/> I106969-99-0	1pc Corpus VS Footplate - 9"D x 16.5"W	No Charge
<input type="checkbox"/> I106971-99-0	2pc Wide Corpus VS Footplates - 9"D x 7.5"W	
<input type="checkbox"/> I106967-99-0	2pc Narrow Corpus VS Footplates - 9"D x 6"W	
CALF SUPPORT OPTIONS:		
<input type="checkbox"/> I107211-99-0	UT Calf Support Kit (J) w Mesh Cover & Adj Hardware 5"W x 6.5"T Pad. <i>Compatible with footplate heights of 10" - 18".</i>	No Charge
<input type="checkbox"/> I106190-99-0	UT Calf Support Kit (A) w Mesh Cover & Adj Hardware 5.5"W x 6"T Pad. <i>Compatible with footplate heights of 10" - 18".</i>	No Charge
<input type="checkbox"/> I106188-99-0	UT Calf Support Kit (H) w Mesh Cover & Adj Hardware 7"W x 7"T Curved Pad. <i>Compatible with footplate heights of 11" - 18".</i>	
<input type="checkbox"/> I106186-99-0	UT Calf Support Kit (D) w Mesh Cover & Adj Hardware 8"W x 8"T Pad. <i>Compatible with footplate heights of 12" - 18".</i>	

FOOTPLATE HEIGHT POSITION:

Part Number	Description
	SET FOOTPLATE HEIGHT AT: _____ <i>This is measured from the SEAT PAN to the top of the FOOTPLATE, when the legrest assembly is at 90°. DO NOT INCLUDE CUSHION THICKNESS IN YOUR MEASUREMENT. In one inch increments, please indicate desired footplate height. NOTE: Corpus footplate height range is 6"-17" and Corpus VS height range is 9"-17". Height may need to be reduced for shipping. Calf Pads are not available at heights less than 9".</i>

OPTIONAL ITEMS:

Part Number	Description	Price AUD
<input type="checkbox"/> I106530-99-0	Stealth Gel Pads - Narrow Footplates (Pair)	
<input type="checkbox"/> I106225-99-0	Stealth Gel Pads - Wide Footplates (Pair)	

Accessories

MUST CHOOSE POSITIONING BELT:

Part Number	Description	Price AUD
<input checked="" type="checkbox"/> I104518-99-0	NOTE: <i>See supplemental order form for more Bodypoint options.</i> Bodypoint Non-Padded Hip Belt - Medium	No Charge



Accessories

CHOOSE ONE OF THE FOLLOWING:

Part Number	Description	Price AUD
<input type="checkbox"/> I107470-99-0	Permobil 7"W Headrest w/ Link Hardware <i>Must also select I107480-99-0 or I107481-99-0.</i>	
<input type="checkbox"/> I107469-99-0	Permobil 10"W Headrest w/ Link Hardware <i>Must also select I107480-99-0 or I107481-99-0.</i>	
<input type="checkbox"/> I107468-99-0	Permobil 14"W Headrest w/ Link Hardware <i>Must also select I107480-99-0 or I107481-99-0.</i>	
<input type="checkbox"/> I107480-99-0	Adj/Removable Headrest Link Hardware - Standard	
<input type="checkbox"/> I107481-99-0	Adj/Removable Headrest Link Hardware - Extended <i>Includes an additional 3" link to extend headrest pad range.</i>	
<input type="checkbox"/> I106331-99-0	Permobil 10"W Headrest, Leatherette <i>Must also select PN I10552.</i>	
<input type="checkbox"/> I106333-99-0	Permobil Ergo Shape Headrest, Black Fabric <i>Must also select PN I10552.</i>	
<input type="checkbox"/> I10552	Adj/Removable Headrest Hardware	
<input type="checkbox"/> I107499-99-0	Universal Headrest Adapter <i>Please Note: Recline is limited to 160° or less (depending on backrest height) when ordering a Universal Headrest Adapter on a 23" or 24" tall Corpus backrest.</i>	

CHOOSE ONE OF THE FOLLOWING:

Part Number	Description	Price AUD
<input type="checkbox"/> I107488-99-0	Permobil 3" x 5" Lateral Supports - Pair <i>Must also select I107490-99-0. Only compatible with Corpus backrest shells.</i>	
<input type="checkbox"/> I107489-99-0	Permobil 5" x 7" Lateral Supports - Pair <i>Must also select I107490-99-0. Only compatible with Corpus backrest shells.</i>	
<input type="checkbox"/> I107490-99-0	Adj/Removable Lateral Link Hardware - Pair	

OPTIONAL ITEMS:

Part Number	Description	Price AUD
<input type="checkbox"/> I103070-99-0	UT Amp. Support Pad (H) w/ Mesh Cover - (1 Each) <i>7"W x 7"T Curved Pad. Must also select PN I10553.</i>	
<input type="checkbox"/> I103071-99-0	UT Amp. Support Pad (D) w/ Mesh Cover - (1 Each) <i>8"W x 8"T Flat Pad. Must also select PN I10553.</i>	
<input type="checkbox"/> I10553	Adj Removable - Amp Support Hardware (1 Each) <i>Please indicate if you would like the amputation support mounted on the right or left side of the seat frame: Mount Amp Suppt on Left: _____ Mount Amp Suppt on Right: _____</i>	
<input type="checkbox"/> I106104-99-0	UT Hip Suppt (B) w/ Mesh Cover - High (Pair) <i>4.5"W x 3.5"T Pad w/ 4" Tall Mounting Bracket. Must also select PN I10550.</i>	
<input type="checkbox"/> I106103-99-0	UT Hip Suppt (A) w/ Mesh Cover - High (Pair) <i>6"W x 5.5"T Pad w/ 4" Tall Mounting Bracket. Must also select PN I10550.</i>	
<input type="checkbox"/> I103094-99-0	UT Thigh/Hip Suppt (J) w/ Mesh Cover - High (Pair) <i>6.5"W x 5"T Pad w/ 4" Tall Mounting Bracket. Must also select PN I10550.</i>	
<input type="checkbox"/> I103095-99-0	UT Thigh/Hip Suppt (C) w/ Mesh Cover - High (Pair) <i>8"W x 3.5"T Pad w/ 4" Tall Mounting Bracket. Must also select PN I10550.</i>	
<input type="checkbox"/> I102977-99-0	UT Thigh/Hip Suppt (G) w/ Mesh Cover - High (Pair) <i>12"W x 3.5"T Pad w/ 4" Tall Mounting Bracket. Must also select PN I10550.</i>	
<input type="checkbox"/> I10550	Adj Removable - Thigh Support Hardware (Pair)	
<input type="checkbox"/> I106762-99-0	Abductor Kit Corpus - Leatherette <i>Must also select PN I10555.</i>	
<input type="checkbox"/> I10555	Adj Removable - Abductor Hardware (1 Each)	
<input type="checkbox"/> I105260-99-0	Transfer Handles, Corpus Seat	
<input type="checkbox"/> I105292-99-0	Push Handles, Corpus Seat	



Accessories

OPTIONAL ITEMS:

Part Number	Description	Price AUD
<input type="checkbox"/> I106210-99-0	RAM Self-Leveling Cup Holder	
<input type="checkbox"/> I106806-99-0	RAM X-Grip Phone Holder	
<input type="checkbox"/> I107100-99-0	RAM X-Grip 10" Tablet Holder w/Swing Away Mount <i>Compatible with most 10" screen tablets with or without a protective case. Includes adjustable swing away mount (I103476-99-0). Note: Customization may be required if ordered along with an Omni.</i>	
<input type="checkbox"/> I107222-99-0	USB Charger <i>Provides up to 1.5A (at 5V) to charge most smartphones and tablets. The charger features an isolated circuit design for improved safety; and by pulling power from the ICS, prevents battery drainage when the chair is powered off. Includes mounting hardware and ICS connector. Please select armrest for mounting, keeping in mind that charger may need to be relocated. Mount on left armrest: _____ Mount on right armrest: _____</i>	
<input type="checkbox"/> I105805-99-0	Medical Necessities Bag Hooks <i>Please Note: Recline is limited to 165° or less (depending on backrest height) when ordering with 23" or 24" tall Corpus backrest.</i>	
<input type="checkbox"/> I102059-99-0	Medical Necessities Bag <i>Includes hardware to attach bag to backrest. Please Note: Recline is limited to 165° or less (depending on backrest height) when ordering with a 23" or 24" tall Corpus backrest.</i>	
<input type="checkbox"/> I105840-99-0	Bodypoint Monoflex Medium Chest Support Belt Kit <i>18.25"L x 4"W Pad Size. Item will be shipped as a part.</i>	
<input type="checkbox"/> I105839-99-0	Bodypoint Monoflex Large Chest Support Belt Kit <i>21.5"L x 4.5"W Pad Size. Item will be shipped as a part.</i>	
<input type="checkbox"/> I105850-99-0	Bodypoint Monoflex X-Lrge Chest Support Belt Kit <i>24.5"L x 5"W Pad Size. Item will be shipped as a part.</i>	
<input type="checkbox"/> I10843	Tray, Small, Angle Adjustable <i>Must also select PN I10556. Mount on Left: _____ Mount on Right: _____</i>	
<input type="checkbox"/> I105221-99-0	Tray, Angle Adjustable, 14" W <i>Must also select PN I10556.</i>	
<input type="checkbox"/> I105222-99-0	Tray, Angle Adjustable, 16" W <i>Must also select PN I10556.</i>	
<input type="checkbox"/> I105220-99-0	Tray, Angle Adjustable, 18" W <i>Must also select PN I10556.</i>	
<input type="checkbox"/> I105262-99-0	Tray with Joystick Cutout, 16" W <i>Must also select PN I10556.</i>	
<input type="checkbox"/> I105263-99-0	Tray with Joystick Cutout, 18" W <i>Must also select PN I10556.</i>	
<input type="checkbox"/> I106863-99-0	Tray, UniTrack, Angle Adjustable <i>Side mounted tray that attaches to the UniTrack Accessory Rail. Dimensions are: 19.5"W x 13"D. Can only be ordered with chairs that offer the UniTrack System. Must also select PN I10556.</i>	
<input type="checkbox"/> I10556	Adj Removable - Tray Hardware (1 Each)	
<input type="checkbox"/> I105224-99-0	Oxygen Holder, Corpus	
<input type="checkbox"/> I105223-99-0	Crutch Holder, Corpus	

Alternative Drive Systems

ALTERNATIVE PROPORTIONAL JOYSTICK SYSTEMS:

Part Number	Description	Price AUD
<input type="checkbox"/> I10751	Mini-Joystick for Omni / IOM (9-pin connector) <i>Operating Force = 13 g Neutral to Full Deflection = 5mm</i>	



Alternative Drive Systems

ALTERNATIVE PROPORTIONAL JOYSTICK SYSTEMS:

Part Number	Description	Price AUD
<input type="checkbox"/> I103893-99-0	Compact Joystick for Omni / IOM (9-pin connector) <i>Operating Force = 220 g Neutral to Full Deflection = 28mm Joystick includes two 3.5mm mono jacks for Mode and On/Off. Includes Standard Joystick Knob and 1.75" (42mm) Foam Ball Knob.</i>	
<input type="checkbox"/> I103894-99-0	Compact Joystick LITE for Omni / IOM (9-pin conn) <i>Operating Force = 50 g Neutral to Full Deflection = 17mm Joystick includes two 3.5mm mono jacks for Mode and On/Off. Includes Standard Joystick Knob and 1.75" (42mm) Foam Ball Knob.</i>	

ALT PROPORTIONAL JOYSTICK MOUNTING HARDWARE:

Part Number	Description
	NOTE: <i>Please see the PERMOfix Mounting Solutions, within the Alt Drive Accessories section, for additional mounting options.</i>
<input type="checkbox"/> I104690-99-0	Alt. JS Swing Away Mount Kit - Armrest <i>Mount in front of Left Arm: _____ Mount in front of Right Arm: _____</i>
<input type="checkbox"/> I107135-99-0	Alt JS/Omni Combo Swing Away Mount Kit - Armrest <i>Used to mount Omni and Alternative Proportional Joystick to the same armrest. Includes Swing Away Omni/Tablet Mount, with a longer mounting rod, two PERMOfix modules, one 3" rod and one 4" rod.</i> <i>Mount in front of Left Arm: _____ Mount in front of Right Arm: _____</i>
<input type="checkbox"/> I104695-99-0	Alt. JS Adjustable Mount Kit - Midline <i>Mount on Left Arm: _____ Mount on Right Arm: _____</i>
<input type="checkbox"/> I104143-99-0	Alt. JS Adjustable Mount Kit - Chin Boom <i>Swing-Away / Removable Chin Boom mounted on right of backrest shell. Custom charge may be applied if ordering lower backrest heights.</i>
<input type="checkbox"/> I10752	Alt. JS Adjustable Mount Kit - Haims Harness, Lrg <i>For neck width up to 19cm.</i>
<input type="checkbox"/> I10753	Alt. JS Adjustable Mount Kit - Haims Harness, Sml <i>For neck width up to 14cm.</i>

ALT PROPORTIONAL JOYSTICK ACCESSORIES:

Part Number	Description	Price AUD
<input type="checkbox"/> I104904-99-0	Dual Switches for Compact Joystick - Fixed Mount <i>Includes two pushbutton switches mounted on rigid Switches can be used for Mode & On/Off. Switch operation force = 130 g.</i> <i>Mount at front of Joystick: _____ Mount at rear of Joystick: _____</i> <i>Mount at left of Joystick: _____ Mount at right of Joystick: _____</i>	
<input type="checkbox"/> I104641-99-0	Dual Switches for Compact Joystick - Gooseneck Mnt <i>Includes two pushbutton switches mounted on adjustable gooseneck arms. Switches are positioned on Left and Right sides of Compact Joystick. Switches can be used for Mode & On/Off. Switch operation force = 130 g</i>	

SIP & PUFF SYSTEMS:

Part Number	Description	Price AUD
<input type="checkbox"/> I10573	Sip & Puff System for Omni <i>Requires Omni Control Module. Cannot use this part with an Input Module. Includes breath tube kit.</i>	
<input type="checkbox"/> I10574	Sip & Puff System with Fiber Optic Opti-Stop <i>Requires Omni Control Module. Cannot use this part with an Input Module. Includes breath tube kit.</i>	
<input type="checkbox"/> I106807-99-0	Sip & Puff Kit for Total Control Head Array <i>For use with Total Control Head Array. Requires Omni Control Module. Cannot be used with an Input Module. Includes breath tube kit.</i>	
<input type="checkbox"/> I10762	Adjustable Mounting Hrdw - Sip & Puff <i>Required when selecting PN I10573, I10574 or I106807-99-0.</i>	



Alternative Drive Systems

PROPORTIONAL FOOT CONTROL:

Part Number	Description	Price AUD
<input type="checkbox"/> I10767	Proportional Foot Control Mount on Left: _____ Mount on Right: _____	
<input type="checkbox"/> I10766	Adjustable Mounting Hrdw - Foot Control Required when selecting PN I10767.	

Total Control Head Array

CONFIGURE HEAD ARRAY:

Part Number	Description	Price AUD
<input type="checkbox"/> I106805-99-0	Permobil Head Array - 3-Switch System <i>Features an occipital pad sensor and two stainless steel arms configured with proximity sensors. The arms are mounted to the upper attachment points on the head array. Longer arms are also included to extend the sensor adjustment range. ESP recommended. Choose Omni Control Module (ID51154IR) separately.</i>	
<input type="checkbox"/> I11072	Program Head Array as a 5-Switch System <i>Select additional arms/sensors below and indicate the intended use of each sensor in the Notes & Comments section. Choose Omni Control Module (ID51154IR) or Input Module (I10581) separately.</i>	No Charge
CHOOSE OCCIPITAL PAD SIZE (REQUIRED):		
<input type="checkbox"/> I105920-99-0	Small Occipital Pad for Head Array - 4.25"W x 3"H	No Charge
<input type="checkbox"/> I105921-99-0	Large Occipital Pad for Head Array - 6"W x 4"H	No Charge
CHOOSE OCCIPITAL PAD MOUNT (REQUIRED):		
<input type="checkbox"/> I10847	Fixed Mount for Occipital Pad - Head Array	No Charge
<input type="checkbox"/> I105936-99-0	Adj. Swivel Mount for Head Array Occipital Pad <i>Enables Occipital Pad to be angled up/down and left/right independent of the lateral arms.</i>	No Charge
CHOOSE BACKREST MOUNTING HARDWARE (REQUIRED):		
<input type="checkbox"/> I105927-99-0	Stealth TWB480 Backrest Mounting Hdwr with Links <i>Please select the length of the vertical stem below: Standard length (10"): _____ Short length (6"): _____</i>	
<input type="checkbox"/> I105941-99-0	Stealth TWB480 Backrest Mtng Hdwr with Flipback	
<input type="checkbox"/> I105925-99-0	Stealth HMO475 Backr Mtng Hdwr w 17mm Ball	
<input type="checkbox"/> I106319-99-0	Stealth TWB480 Backrest Mnt Hdwr w SubOcptl Mount	
<input type="checkbox"/> I106282-99-0	Stealth 600 SubOccipital Pad - 8.5"W x 2.75"H <i>Must select when PN I106319-99-0 is ordered.</i>	
SELECT ADDITIONAL LEFT ARM PACKAGE (OPTIONAL):		
<input type="checkbox"/> I10976	Add Left Lateral Arm w Proximity Sensor <i>Choose Lateral Arm Adjustment Range: Short (5.5" - 10"): _____ Long (5.5" - 14"): _____</i>	
<input type="checkbox"/> I10977	Add Left Lateral Arm w Egg Switch <i>Choose Lateral Arm Adjustment Range: Short (5.5" - 10"): _____ Long (5.5" - 14"): _____</i>	
<input type="checkbox"/> I10978	Add Left Lateral Arm w Pad (No Sensor or Switch) <i>Choose Lateral Arm Adjustment Range: Short (5.5" - 10"): _____ Long (5.5" - 14"): _____</i>	
SELECT ADDITIONAL RIGHT ARM PACKAGE (OPTIONAL):		
<input type="checkbox"/> I10979	Add Right Lateral Arm w Proximity Sensor <i>Choose Lateral Arm Adjustment Range: Short (5.5" - 10"): _____ Long (5.5" - 14"): _____</i>	



Total Control Head Array

CONFIGURE HEAD ARRAY:

Part Number	Description	Price AUD
<input type="checkbox"/> I10980	Add Right Lateral Arm w Egg Switch <i>Choose Lateral Arm Adjustment Range: Short (5.5" - 10"): _____ Long (5.5" - 14"): _____</i>	
<input type="checkbox"/> I10981	Add Right Lateral Arm w Pad (No Sensor or Switch) <i>Choose Lateral Arm Adjustment Range: Short (5.5" - 10"): _____ Long (5.5" - 14"): _____</i>	

Alt Drive Systems Setup

CHOOSE OMNI OR INPUT MODULE (REQUIRED):

Part Number	Description	Price AUD
<input type="checkbox"/> ID51154IR	Omni Control Module for R-net w/ Infrared <i>Must also select mounting hardware, either I103476-99-0 or I107135-99-0. Please specify which type of alternative drive input you will be using if not ordering it from Permobil. Input Device: _____</i>	
<input type="checkbox"/> I103476-99-0	Adjustable Swing-Away Hardware - Omni/Tablet <i>Mount in front of Left Arm: _____ Mount in front of Right Arm: _____</i>	
<input type="checkbox"/> I10581	Input Module (IOM) for R-net Control System <i>This device allows a proportional or switched drive input to be connected to an R-net Control System without the use of an Omni. It must be ordered in conjunction with an R-net Remote Joystick. It does NOT support any of the following input devices: Sip-n-Puff, 3-Switch Head Array or RIM control.</i>	

CHOOSE MOUNTING POSITION OF SEAT CONTROL SWITCHBOX:

Part Number	Description	Price AUD
<input type="checkbox"/> I10080	Mount Seat Control Switchbox on Backrest	No Charge

Alt Drive Accessories

PERMOFIX MOUNTING SOLUTIONS:

Part Number	Description	Price AUD
<input type="checkbox"/> I106655-99-0	PERMOfix UniTrack Fixed Mount Kit <i>Includes two base joints, two 4" rods and one 8" rod.</i>	
<input type="checkbox"/> I106656-99-0	PERMOfix UniTrack Removable Mount Kit <i>Includes quick release mount and thigh support hardware, two base joints and three different length rods.</i>	
<input type="checkbox"/> I106657-99-0	PERMOfix UniTrack Swing Away Mount Kit <i>Includes PERMOfix Swing Away module, two base joints, one 2" rod and one 4" rod.</i>	
<input type="checkbox"/> I105907-99-0	Proximity Sensor + Swivel Assm - Not Mounted <i>Only available with the Total Control Head Array. Swivel assembly allows sensor to be mounted with PERMOfix or directly to lateral arm of Total Control Head Array. Since Proximity Sensors require power, they CANNOT be used as On/Off switches.</i>	
<input type="checkbox"/> I105919-99-0	Egg Switch + Swivel Mount Assm - Not Mounted <i>Swivel assembly allows switch to be mounted with PERMOfix or directly to lateral arm of Total Control Head Array.</i>	
<input type="checkbox"/> I105912-99-0	Solid Lateral Pad+Swivel Joint - Not Mounted <i>Swivel assembly allows pad to be mounted with PERMOfix or directly to lateral arm of Total Control Head Array.</i>	
<input type="checkbox"/> I106675-99-0	Oval Plate + Swivel Assm - No Switch - Not Mounted <i>Blank plate provides mounting option for Micro Lite, Buddy Button or other switch options. Swivel assembly allows pad to be mounted with PERMOfix or directly to lateral arm of Total Control Head Array.</i>	

OPTIONAL ITEMS:

Part Number	Description	Price AUD
-------------	-------------	-----------



Alt Drive Accessories

OPTIONAL ITEMS:

Part Number	Description	Price AUD
<input type="checkbox"/> I101769-99-0	Micro Lite Switch Qty: _____	
<input type="checkbox"/> I101767-99-0	Buddy Button Qty: _____	
<input type="checkbox"/> I106293-99-0	Egg Switch Black Qty: _____	
<input type="checkbox"/> I101771-99-0	Mini Cup Switch Qty: _____	
<input type="checkbox"/> I107500-99-0	Proximity Switch Qty: _____	
<input type="checkbox"/> I107495-99-0	Mini Joystick 5 Switch OC for Omni/IOM (9-pin)	
<input type="checkbox"/> I107498-99-0	Touch Contact with CR2032 Battery Qty: _____	
<input type="checkbox"/> I10772	Dome Switch, Backrest Mnt w Gooseneck Qty: _____ <i>Mount on Left: _____ Mount on Right: _____</i> <i>Mounts at top of backrest shell. Custom charge may be applied if ordering with short backrest heights. Kit Includes: One PERMOfix joint with Flat Plate, two PERMOfix joints with 45° swivel adjustment, one 150mm rod, one 200mm rod and one 225mm gooseneck switch with 1.2m cable and mounting clamp.</i>	
<input type="checkbox"/> I10846	Dome Switch, Backrest Swg Awy Mnt w Gsneck Qty: _____ <i>Mount on Left: _____ Mount on Right: _____</i> <i>w 90° Swing-Away Function. Mounts at top of back. Custom charge applies if ordered with short back. Kit Includes: One PERMOfix joint with Flat Plate, one PERMOfix joint with 45° swivel adjustment, two 6" rods, one 90° Swing-Away joint, one 9"</i>	
<input type="checkbox"/> I105935-99-0	3.5mm Ext. Cable for Sensor/Switch, 0.3m Long	
<input type="checkbox"/> I105931-99-0	3.5mm Ext. Cable for Sensor/Switch, 1m Long	
<input type="checkbox"/> I105932-99-0	3.5mm Ext. Cable for Sensor/Switch, 2m Long	
<input type="checkbox"/> I105933-99-0	3.5mm Ext. Cable for Sensor/Switch, 3m Long	
<input type="checkbox"/> I105934-99-0	3.5mm Y-splitter Cable for Sensor/Switch <i>Allows two Proximity Sensors/Mechanical Switches to connect to the same direction jack.</i>	

F3 / F5 Upgrade Packages

Additional Order Instructions:

Chair Order Policy:

You are highly encouraged to contact Permobil Sales with ANY questions regarding proper completion of our order forms. The order cannot be processed without the following:

- A client name or code, the client’s measurements, the client’s weight
- The order form must also be signed by a person authorized by your company to acknowledge the items selected. Permobil is not responsible for configuration or size discrepancies resulting from customer errors on the order form.
- Order confirmations will be provided by Permobil summarizing the items selected. Please review this carefully.
- Once the order is shipped, any changes are subject to the returns policy and may be prohibited.

Order Acknowledgement:

I, _____, am an agent of the medical equipment provider named on this order form and I have the authority to contract for the purchase of powered wheelchairs and related parts on behalf of said provider. I acknowledge that I have reviewed this order and that it is complete and accurate to the best of my knowledge. I further acknowledge that any changes to this order after submission of this order form are subject to the returns policy and may carry additional charges.

Signed: _____

Printed Name: _____

Title: _____

Date: _____

NOTE: All specifications and prices are subject to change without notice. Please note that prices displayed are only valid if ordered with the wheelchair. Any parts ordered for service or future alterations will carry different pricing and warranty.