

# F5 Corpus

# permobil



## Dealer Information

\* **Contact:** \_\_\_\_\_

Dealer Code: \_\_\_\_\_

\* **Dealer Name:** \_\_\_\_\_

Address: \_\_\_\_\_

\* **City:** \_\_\_\_\_

\* **Province:** \_\_\_\_\_

\* **Phone#:** \_\_\_\_\_

\* **Fax#:** \_\_\_\_\_

PO#: \_\_\_\_\_

Email Address: \_\_\_\_\_

## Client Information

*Permobil recommends that the client is evaluated by a certified rehab specialist.*

\* **First Name:** \_\_\_\_\_

\* **Last Name:** \_\_\_\_\_

Diagnosis: \_\_\_\_\_

Funding Source: \_\_\_\_\_

Client Age: \_\_\_\_\_

## Client Measurements

\* **Weight:** \_\_\_\_\_

\* **Height:** \_\_\_\_\_

A) Top of Shoulders: \_\_\_\_\_

B) Chest Depth: \_\_\_\_\_

C) Chest Width: \_\_\_\_\_

D) Seat Depth: \_\_\_\_\_


E) Top of Head: \_\_\_\_\_

F) Elbow to Hand: \_\_\_\_\_

G) Seat Pan to Elbow: \_\_\_\_\_

H) Hip Width: \_\_\_\_\_

\* **I) Knee to Foot:** \_\_\_\_\_



Display ADP Codes on order/quote:

Three column pricing on order/quote:

Please send order/quote to fax#: **(800) 231-3256**

Email to: [sales@permobil.com](mailto:sales@permobil.com)

**Permobil Ltd.**  
12 - 75 Mary St.  
Aurora, Ontario  
L4G 1G3  
Tel: (800) 736-0925  
Fax: (800) 231-3256  
[www.permobil.com](http://www.permobil.com)

Prices effective March 1, 2017 and payable in Canadian dollars.

**Chair Model and Colors**

**CHOOSE BASE:**

<input checked="" type="checkbox"/>	323539	F5 Base - Corpus	14,079.00	WP0001257
<i>Price Includes: FWD F5 Power Base Set Up For Multiple Power Seat Functions, Four-Wheel Independent Suspension System, Tie-Down Hardware for Strap Systems, Charger, Flat Free or Pneumatic Tires and a Standard Positioning Belt. Weight Capacity = 300lbs, STFH = 17.5".</i>				

**CHOOSE SPEED PACKAGE:**

Part Number	Description	MSRP - CA	ADP Code
<input type="checkbox"/>	920065 Set Max Speed to 7.5 MPH <i>Anti-Tippers are required when selecting this speed package.</i>	No Charge	
<input type="checkbox"/>	920064 Set Max Speed to 6 MPH	No Charge	
<input type="checkbox"/>	323070 Anti-Tippers for F3/F5 Corpus <i>Anti-Tippers are required for user weight above 220lbs or max speed above 6.0mph for the F5 or above 5mph for the F3.</i>	471.00	

**CHOOSE DRIVE TIRES:**

Part Number	Description	MSRP - CA	ADP Code
<input type="checkbox"/>	324117 Flat Free Tire 3"W x 14"T Black - Standard Tread <i>Designed for a smooth ride on indoor and outdoor surfaces. Includes flat free standard tread caster wheels.</i>	No Charge	
<input type="checkbox"/>	324118 Pneumatic Tire 3"W x 14"T Black - Standard Tread <i>Designed for a smooth ride on indoor and outdoor surfaces. Includes flat free standard tread caster wheels.</i>	No Charge	
<input type="checkbox"/>	325401 Flat Free Tire 3"W x 14"T Black - Aggressive Tread <i>Designed for improved performance on outdoor surfaces. Includes flat free aggressive tread caster wheels.</i>	No Charge	

**OPTIONAL DRIVE TIRES:**

Part Number	Description	MSRP - CA	ADP Code
<input type="checkbox"/>	325402 Pneumatic/Snow Tire 3"W x 14"T Black - Aggressive Tread <i>Snow tires are designed for use in cold weather conditions. The nylon compound in these tires remains flexible at low temperatures, but wears more easily. Therefore, the snow tires (sold as wheel assemblies) are not recommended for year-long use and must be ordered along with standard drive wheels. Does not include caster wheels. Price shown is for 1(ea). A quantity of 2 will be provided.</i>	263.00	

**CHOOSE SHROUD/RIM COLOR:**

Part Number	Description	MSRP - CA	ADP Code
<input type="checkbox"/>	323622063 F5 Color Kit - Sunburst Orange	No Charge	
<input type="checkbox"/>	323622067 F5 Color Kit - Galactic Green	No Charge	
<input type="checkbox"/>	323622060 F5 Color Kit - Midnight Black	No Charge	
<input type="checkbox"/>	323622066 F5 Color Kit - Ocean Blue	No Charge	
<input type="checkbox"/>	323622061 F5 Color Kit - Power Pink	No Charge	
<input type="checkbox"/>	323622062 F5 Color Kit - Radiant Red	No Charge	
<input type="checkbox"/>	323622079 F5 Color Kit - Phantom Grey	No Charge	
<input type="checkbox"/>	323622074 F5 Color Kit - Titan Blue	No Charge	
<input type="checkbox"/>	323622081 F5 Color Kit - Carbon Fiber	No Charge	
<input type="checkbox"/>	323622082 F5 Color Kit - Mossy Oak Break-Up Country	No Charge	

**Base Options**

**MUST CHOOSE EITHER SEAT ELEVATOR OR FIXED SEAT TUBE:**

Part Number	Description	MSRP - CA	ADP Code
<input type="checkbox"/>	324425 Power Adjustable Seat Height - 14" Travel <i>Must also select Power Tilt when selecting this option.</i>	3,441.00	

## Base Options

### MUST CHOOSE EITHER SEAT ELEVATOR OR FIXED SEAT TUBE:

Part Number	Description	MSRP - CA	ADP Code
<input type="checkbox"/> 324107	Seat Tube Fixed	No Charge	

### SET SEAT HEIGHT:

Part Number	Description	MSRP - CA	ADP Code
<input type="checkbox"/>	Set Seat Height to 17.5"	No Charge	
<input type="checkbox"/>	Set Seat Height at 18.5" <i>If ordering Power Adjustable Seat Height, 1" of elevator travel is used to set seat at 18.5".</i>	No Charge	
<input type="checkbox"/>	Set Seat Height at 19.5" <i>If ordering Power Adjustable Seat Height, 2" of elevator travel is used to set seat at 19.5".</i>	No Charge	

### MUST CHOOSE BATTERIES:

Part Number	Description	MSRP - CA	ADP Code
<input checked="" type="checkbox"/> 602130	Batteries Grp 24 (72Ah) Sealed Gel Installed <i>Price is for 1(ea). A quantity of 2 will be provided.</i>	483.00	

### OPTIONAL ITEMS:

Part Number	Description	MSRP - CA	ADP Code
<input type="checkbox"/> 325969	Essentials Box for Chassis w/ Reflectors <i>Please Note: The total backrest angle will be limited to 150° to avoid contact between backrest components and the Essentials Box.</i>	536.00	
<input type="checkbox"/> 326641	Essentials Box for Chassis w/ Lights <i>Please Note: The total backrest angle will be limited to 150° to avoid contact between backrest components and the Essentials Box. Must also select the Light Kit.</i>	892.00	
<input type="checkbox"/> 323619	LED Light Kit for F3 / F5	1,233.00	
<input checked="" type="checkbox"/> 330415	Enhanced Steering Perform Unit (ESP) R-net	2,158.00	
<input type="checkbox"/> 320273	Compact Laptop Ventilator Tray - Corpus Seat <i>Mounts to rear of Corpus backrest hinge. Accommodates Pulmonetic Systems / CareFusion LTV Series and similar sized Laptop Ventilators. The tray is 330mm Tall, adjustable Width between 200mm - 330mm, and adjustable Depth between 50mm - 165mm. Overall backrest angle limited to 150° when mounting this ventilator tray. Custom Option.</i>	1,682.00	
<input type="checkbox"/> 320399	Battery Tray for Laptop Ventilator Tray <i>Mounts to rear of Laptop Ventilator Tray. Measures 290mm W x 170mm H x 42mm D. Accommodates Pulmonetic Systems Lithium-Ion Battery Pack, CareFusion SprintPack Battery System, or Richardson PowerTech Vent Power Center. Custom Option.</i>	269.00	
<input type="checkbox"/> 320390	O2 Holder for Laptop Ventilator Tray <i>Mounts to Left or Right Side of Laptop Vent Tray. O2 Holder is 175mm Tall. Width infinitely adjustable 90mm - 170mm. Depth infinitely adjustable 95mm - 125mm. Custom Option.</i>	380.00	WEPX
<input type="checkbox"/> 613399	Programming Key for R-net - Permobil Version	524.00	

## Joystick Options

### CHOOSE JOYSTICK TYPE:

Part Number	Description	MSRP - CA	ADP Code
PLEASE NOTE: <i>For additional input devices, refer to Alternative Drive Systems Section.</i>			
<input type="checkbox"/> 615264	Permobil Joystick Module <i>Designed exclusively for Permobil, this R-net color joystick features: control of up to 4 Bluetooth devices, Infrared control, programmable buttons and a programmable switch jack; all in a compact, durable Aluminum housing.</i>	1,185.00	
<input type="checkbox"/> 920147	Permobil Joystick Module NC <i>Designed exclusively for Permobil, this R-net color joystick features: control of up to 4 Bluetooth devices, Infrared control, programmable buttons and a programmable switch jack; all in a compact, durable Aluminum housing.</i>	No Charge	
<input type="checkbox"/> 615351	3.5mm Y-Splitter Cable for PJSM <i>For use with the PJSM Mode/Profile jack. Allows 2 mechanical switches to control 4 separate functions. Requires PC Programming to assign the short hold and long hold commands to each switch. Switches sold separately.</i>	42.00	

## Joystick Options

## CHOOSE JOYSTICK TYPE:

Part Number	Description	MSRP - CA	ADP Code
<input type="checkbox"/> 613977	R-net Remote Color JS w/ Lights Heavy Duty Gimble <i>Heavy-Duty modification enables gimbal to hold up to extreme operating forces. Reinforced housing is attached to a standard Colour LCD R-net Joystick. 3" long joystick knob has 3/4" diameter. Not compatible with alternative joystick handles.</i>	2,042.00	WEX1
<input type="checkbox"/> 613411	Compact Advanced Joystick for R-net <i>Compact sized joystick that does not contain an LCD Display. Contains two mono jacks for mode &amp; on/off. Does NOT contain a wrist support pad. If you want a wrist support pad, please select 614934 below.</i>	1,990.00	WEX1
<input type="checkbox"/> 614934	Wrist Pad for Compact Joystick Advanced	253.00	

## CHOOSE JOYSTICK HANDLE TYPE:

Part Number	Description	MSRP - CA	ADP Code
<input type="checkbox"/> 603086	Cone Shaped Joystick Knob - Standard	No Charge	
<input type="checkbox"/> 602154	Chin Cup For Joystick Knob	119.00	
<input type="checkbox"/> 602151	T Handle For Joystick Knob	176.00	
<input type="checkbox"/> 602152	Large Ball For Joystick Knob	119.00	
<input type="checkbox"/> 602155	Softball For Joystick Knob	160.00	
<input type="checkbox"/> 602153	Mushroom For Joystick Knob	160.00	
<input type="checkbox"/> 602156	Stick For Joystick Knob	160.00	
<input type="checkbox"/> 613703	1.75" Black Foam Ball For Joystick Knob	119.00	
<input type="checkbox"/> 613970	Bodypoint J/S Handle 3" U-Shaped w/ Flex-Shaft	153.00	
<input type="checkbox"/> 613971	Bodypoint J/S Handle 4" U-Shaped w/ Flex-Shaft	153.00	
<input type="checkbox"/> 613972	Bodypoint J/S Handle Dome Shaped Rubber	153.00	

Part Number	Description	MSRP - CA	ADP Code
<input type="checkbox"/> 317769	Joystick Mount - Fixed <i>Mount on Left: _____ Mount on Right: _____</i>	No Charge	
<input type="checkbox"/> 317404	Retractable Joystick Mount <i>Mount on Left: _____ Mount on Right: _____</i>	452.00	
<input type="checkbox"/> 319940	Mount Standard Joystick on Corpus Back Hinge <i>Standard R-net Joystick is mounted to backrest hinge. Can be chosen instead of PN 323725. Mount on Left: _____ Mount on Right: _____</i>	314.00	
<input type="checkbox"/> 318238	Midline Control Mount <i>Mount Midline on Left Armrest: _____ Mount Midline on Right Armrest: _____</i>	392.00	WE3D
<input type="checkbox"/> 614290	Stealth Gatlin Midline Mount - 8 Position <i>Features eight locking positions each 45° apart. Adjustments are made by pulling out on mounting arm and rotating to desired angle. Joystick bracket sold separately (select from the brackets listed below). Please indicate if you want it mounted on left or right armrest: Mount Midline on Left Armrest: _____ Mount Midline on Right Armrest: _____</i>	690.00	WE3D
<input type="checkbox"/> 614291	Stealth Gatlin Midline Mount - 2 Position <i>Designed for easy swing away. Mount can be set to swing from horizontal to 95° or from horizontal to 180°. Joystick bracket sold separately (select from the brackets listed below). Please indicate if you want it mounted on left or right armrest: Mount Midline on Left Armrest: _____ Mount Midline on Right Armrest: _____</i>	690.00	WE3D
<input type="checkbox"/> 614292	Stealth Gatlin Bracket for R-net/Compact JS Adv	283.00	
<input type="checkbox"/> 614295	Stealth Gatlin Bracket for Compact Joystick	283.00	
<input type="checkbox"/> 614296	Stealth Gatlin Bracket for Mini Joystick	283.00	
<input type="checkbox"/> 614294	Stealth Gatlin Swivel Bracket for Omni	322.00	
<input type="checkbox"/> 613649	PERMOfix Base Joint with 9/16"(14mm) Through-Mount <i>Includes a shim that allows the PERMOFix Joint to be mounted on the 1/2" Stealth Gatlin arm. Provides an adjustable mount for a Mini Joystick or Compact Joystick.</i>	114.00	
<input type="checkbox"/> 614280	Stealth Height Adjustable Mount <i>Replaces standard UT panel bracket. Provides 1.75" of height adjustability to all Permobil joystick mounts and Stealth Gatlin Midline mounts.</i>	252.00	

**Joystick Options**

OPTIONAL ITEMS:

<input type="checkbox"/>	323702	Co-Pilot Attendant Control R-net - Corpus	5,034.00
<input type="checkbox"/>	323725	Attendant Control R-net	1,307.00
<input type="checkbox"/>	325446	R-net Bluetooth iDevice Module	1,055.00
		<i>Allows input device commands (such as joystick nudges or switch clicks) to be used to operate an iPhone, iPad or iPod. By assigning each input to an iOS function, users can make phone calls, send texts and emails and access applications. Kit includes a Bluetooth iDevice Module, an R-net Bus Cable, and two mono jack cables for direct switch inputs (switches not included). Must be ordered with a Color JSM or Omni for mode display.</i>	
<input type="checkbox"/>	325447	R-net Bluetooth Mouse Module	1,055.00
		<i>Enables mouse control through any R-net input device. It can operate the mouse on a computer and on select Android devices (with R-net for Android app). Kit includes a Bluetooth Mouse Module, an R-net Bus Cable, two mono jack cables for L/R click (switches not included) and a USB BT Adapter for user's computer. Note: BT Adapter may not work with all computers. Must be ordered with a Color JSM or Omni for mode display.</i>	
<input type="checkbox"/>	325468	Add a second Bluetooth Mouse Module	1,055.00
		<i>Allows R-net to simultaneously pair to a second compatible computer or Android device. Each module is accessible through a separate R-net mode. Kit includes a unique Bluetooth Mouse Module, an R-net Bus Cable, two mono jack cables for L/R click and a USB BT Adapter.</i>	
<input type="checkbox"/>	613774	Output Module (IOM) for R-net Control System	1,749.00
		<i>This allows a user to control an external device (ex: Communication Device, ECU, Door Opener, etc.) through the drive input (JSM, Head Array, etc.). Each module can provide up to 4 switch closures. May require specialty cable or wiring. Up to three output modules can be installed on one wheelchair.</i>	

**Seating System**

CHOOSE POSTERIOR TILT/RECLINE PACKAGE:

Part Number	Description	MSRP - CA	ADP Code
	PLEASE NOTE:		
	<i>A half speed limitation is imposed when the Total Backrest Angle (TBA) &gt; 150°. TBA = Tilt Angle + Backrest Recline Angle. Example 1: 45° Tilt + 100° Recline = 145° TBA (Full Speed). Example 2: 35° Tilt + 120° Recline = 155° TBA (Half Speed).</i>		
<input checked="" type="checkbox"/>	324095	Corpus Seat Frame	1,422.00
<input type="checkbox"/>	920121	Fixed Tilt - AP Unit	No Charge
<input type="checkbox"/>	920122	50° Power Tilt - AP Unit	7,069.00
<input type="checkbox"/>	324241	120° Manual Recline - Corpus	No Charge
<input type="checkbox"/>	325178	175° Power Recline - F5	8,690.00
<input type="checkbox"/>	920131	50° Power Tilt & 175° Power Recline - F5	12,732.00
	<i>Select this option when both Power Tilt and Power Recline are needed.</i>		

ANTERIOR TILT PACKAGES:

Part Number	Description	MSRP - CA	ADP Code
	PLEASE NOTE:		
	<i>Active Reach should be recommended by a licensed clinical professional and/or ATP. While operating Active Reach, positioning equipment such as lap belts and chest supports should be considered to meet the user's specific safety needs. Note: A low speed limitation is imposed from 5° to 10° of Active Reach, an extra low speed limitation is imposed from 10° to 20° of Active Reach and a drive lockout is imposed above 20° of Active Reach.</i>		
<input type="checkbox"/>	325128	5° Active Reach	No Charge
	<i>When ordered with a Power Adjustable Seat Height, the tilt range will be -5° to +50°. When ordered without a Power Adjustable Seat Height, the tilt range will be -5° to +45°.</i>		
<input type="checkbox"/>	325129	10° Active Reach	No Charge
	<i>Please Note: This option is only available on configurations ordered with a Power Adjustable Seat Height and Power Tilt. A chest support belt (see accessories section) is strongly recommended when selecting this option.</i>		
<input type="checkbox"/>	325130	20° Active Reach	2,958.00



## Seating System

## ANTERIOR TILT PACKAGES:

Part Number	Description	MSRP - CA	ADP Code
<input type="checkbox"/> 325130	20° Active Reach <i>Please Note: This option is only available on configurations ordered with Power Adjustable Seat Height, Power Tilt, Power Recline and Power Legrest Elevation. Must select a Corpus VS Footplate and Corpus VS Knee Support and hardware from VS Legrest Options section. A chest support belt (see Accessories section) is strongly recommended when selecting this option. Not compatible with retractable positioning belts.</i>	2,958.00	
<input type="checkbox"/> 325131	45° Active Reach <i>Please Note: This option is only available on configurations ordered with Power Adjustable Seat Height, Power Tilt, Power Recline and Power Legrest Elevation. Must select a Corpus VS Footplate and Corpus VS Knee Support and hardware from VS Legrest Options section. A chest support belt (see Accessories section) is strongly recommended when selecting this option. Not compatible with retractable positioning belts.</i>	3,326.00	

## Seating System Options

## CHOOSE SEAT FUNCTION CONTROL TYPE:

Part Number	Description	MSRP - CA	ADP Code
<input type="checkbox"/> 920019	Single Seat Function Control Switch for R-net <i>Please choose either Push Buttons or Toggle Switches below.</i>	No Charge	
<input type="checkbox"/> 920017	Single Seat Function Control Kit for R-net <i>For operation of a single seat function through the wheelchair electronics via the drive control and/or ICS Switchbox. Order with R-net Controller or remote joystick options. Please indicate below if you want push buttons or toggles on the ICS Switchbox.</i>	1,661.00	
<input type="checkbox"/> 920018	Multiple Seat Function Control Kit for R-net <i>MUST be selected if ordering two or more seat functions. This device is used to operate two or more seat functions through the wheelchair electronics via the drive control and/or ICS switch box. Order with R-net electronics or remote joystick options. Please indicate below if you want push buttons or toggles on the ICS Switchbox.</i>	2,768.00	SEC1P1005

## CHOOSE PUSH BUTTONS OR TOGGLES:

Part Number	Description	MSRP - CA	ADP Code
<input type="checkbox"/> 920135	Push Buttons on ICS Switchbox (1.4m)	No Charge	
<input type="checkbox"/> 920137	Toggle Switches on ICS Switchbox (1.4m)	296.00	

## CHOOSE MOUNTING POSITION OF ICS SWITCHBOX:

Part Number	Description	MSRP - CA	ADP Code
<input type="checkbox"/> 316437	Mount ICS Switchbox onto R-net Joystick Bracket <i>ICS Switchbox will be mounted to the same bracket as the R-net Joystick Module. If this item is selected, the ICS Switchbox will swing away with the Joystick Module.</i>	No Charge	
<input type="checkbox"/> 315973	Mount ICS Switchbox onto UT Armrest <i>Mount on Left: _____ Mount on Right: _____</i>	No Charge	
<input type="checkbox"/> 321798	Mount ICS Switchbox on Backrest <i>Mount on Left: _____ Mount on Right: _____</i>	166.00	

## OPTIONAL ITEMS:

Part Number	Description	MSRP - CA	ADP Code
<input type="checkbox"/> 920124	Virtual Seating Coach Package w/ Phone Holder <i>The Virtual Seating Coach (VSC) Package includes a Bluetooth dongle, a USB charger, a RAM X-grip Phone Holder and access to VSC coaching and regimen management. The VSC app can be downloaded for free from the App Store or Google Play. Requires an iPhone (iOS 9 or higher) or an Android phone (OS 4. 1. 2 or higher).</i> Mount VSC to Left Armrest: _____ <i>To provide access, USB Charger and Phone Holder may require relocation.</i> Mount VSC to Right Armrest: _____ <i>To provide access, USB Charger and Phone Holder may require relocation.</i>	2,261.00	
<input type="checkbox"/> 317909	ICS Alternative Switchbox for R-net Electronics <i>Seat function interface box that accommodates up to 8 momentary switches. Provides an additional method for accessing the seat functions. 1/8" mini-mono switches not included.</i>	1,307.00	

**Seating System Options**

**OPTIONAL ITEMS:**

Part Number	Description	MSRP - CA	ADP Code
<input type="checkbox"/> 324664	Remote Stop System - F3 / F5 <i>Consists of a Radio Frequency (RF) handheld transmitter and an RF receiver mounted on the wheelchair.</i>	1,439.00	
<input type="checkbox"/> 330365	Memory Seat Program Corpus w/ R-Net (Recall & Save)	121.00	
<input type="checkbox"/> 322020	Legrest Movement Linked to Recline Electronically <i>Power Legrest will raise when Power Recline moves backward and will lower when back moves forward. If Power Legrest is operated, Power Recline will NOT automatically reposition. Feature works through ICS Switchbox and Driver Control.</i>	49.00	
<input type="checkbox"/> 321953	Enable Independent Repositioning Mode <i>IRM allows for sequential movement of power seat functions into a posterior position (35° tilt, 115° recline, 135° legrest) and returns to upright position (5° tilt, 100° recline, 110° legrest). This can be accessed using a standard drive control or a single switch. On the Permobil Joystick Module, IRM is accessible through the mode/profile jack (long switch hold). Please see the Alt Drive Accessories section for switch options.</i>	No Charge	

**Backrest Options**

**CHOOSE FIXED OR SLIDING BACKREST KIT:**

Part Number	Description	MSRP - CA	ADP Code
<input type="checkbox"/> 318608	Sliding Backrest Kit - 5" Mechanical Shear	No Charge	
<input type="checkbox"/> 319640	Fixed Backrest Kit - No Mechanical Shear	53.00	

**CHOOSE BACKREST TYPE:**

Part Number	Description	MSRP - CA	ADP Code
<b>BACKREST HEIGHT NOTE:</b> <i>Backrest heights are measured from the seat pan to the top of the backrest shell with the backrest set at 90°. Ergo Backs include lateral wedges and lumbar pads.</i>			
<input type="checkbox"/> 325309	Corpus Ergo Back - 14"W x 20"T <i>Power Recline will be limited to 150° when ordering this backrest.</i>	1,072.00	
<input type="checkbox"/> 325310	Corpus Ergo Back - 16"W x 20"T <i>Power Recline will be limited to 150° when ordering this backrest.</i>	1,072.00	
<input type="checkbox"/> 325311	Corpus Ergo Back - 18"W x 20"T <i>Power Recline will be limited to 150° when ordering this backrest.</i>	1,072.00	
<input type="checkbox"/> 325312	Corpus Ergo Back - 20"W x 20"T <i>Power Recline will be limited to 150° when ordering this backrest.</i>	1,072.00	
<b>PLEASE NOTE:</b> <i>All taller backrests listed below include: lateral wedges, lumbar pads, a UT rail mounted at the top for lateral adjustment of the headrest, and split backrest feature with height adjustment in 1" increments. Please indicate below where you want the height set at so that a properly sized backrest cushion can be provided.</i>			
<input type="checkbox"/> 325313	Corpus Ergo Back - 14"W x 23" - 28"T	1,072.00	
<input type="checkbox"/> 325314	Corpus Ergo Back - 16"W x 23" - 28"T	1,072.00	
<input type="checkbox"/> 325315	Corpus Ergo Back - 18"W x 23" - 28"T	1,072.00	
<input type="checkbox"/> 325316	Corpus Ergo Back - 20"W x 23" - 28"T	1,072.00	
	Set Backrest at 23" _____		
	Set Backrest at 24" _____		
	Set Backrest at 25" _____		
	Set Backrest at 26" _____		
	Set Backrest at 27" _____		
	Set Backrest at 28" _____		

**CHOOSE BACKREST CUSHION:**

Part Number	Description	MSRP - CA	ADP Code
<input type="checkbox"/>	Leatherette Cushion - Corpus Ergo Back	No Charge	
<input type="checkbox"/>	Black Mesh Cushion - Corpus Ergo Back	53.00	
<input type="checkbox"/>	Omit Backrest Cushion	No Charge	

**Backrest Options**

**OPTIONAL ITEMS:**

Part Number	Description	MSRP - CA	ADP Code
<input type="checkbox"/> 320523	Tubular Backrest Frame - Corpus (1" Canes) <i>Provides 1" diameter, vertical canes for mounting aftermarket backrests. Frame can be adjusted to support backrest widths from 14" to 20". Canes are attached with 2" long UniTrack linkage disks that can rotate to provide space for backrest hardware. The tops of the canes can be adjusted up to a height of 27" above the seat pan. Depending on the aftermarket backrest, several inches of seat depth may be lost when mounting this system.</i>	477.00	
<input type="checkbox"/> 321190	Direct Backrest Frame - Corpus <i>Used to mount ROHO Agility, Jay 3, Matrx Elite, Varilite Icon and Ride Java backrests ranging in width from 14" to 20". Slotted brackets and 1" UniTrack linkage disks provide a close, adjustable mount. Backrest pitch can be tuned by rotating the linkage disks. Some seat depth may be lost due to the position and orientation of the backrest.</i>	477.00	
<input type="checkbox"/> 920068	Split Backrest Cushion - Corpus <i>Custom Option.</i>	355.00	
<input type="checkbox"/> 920069	Tailored Backrest Cushion - Corpus <i>Sides of cushion are trimmed to provide additional clearance for lateral trunk supports. Must also select a backrest cushion when ordering this custom option.</i>	111.00	

**Seat Sizing**

**CHOOSE SEAT WIDTH:**

Part Number	Description	MSRP - CA	ADP Code
<input type="checkbox"/> 319197	Corpus Seat Width Kit - 17"W-19"W <i>Set Seat Width at 17"W ____ Set Seat Width at 19"W ____</i>	No Charge	
<input type="checkbox"/> 320231	Corpus Seat Width Kit - 21"W-23"W <i>Set Seat Width at 21"W ____ Set Seat Width at 23"W ____</i>	No Charge	

**CHOOSE SEAT DEPTH:**

Part Number	Description	MSRP - CA	ADP Code
	PLEASE NOTE: <i>There are three seat depth kits, each providing three seat depth settings: (1) 14" - 16", (2) 17" - 19" and (3) 20" - 22". Please select the required seat depth and the appropriate kit will be provided.</i>		
<input type="checkbox"/> 321435	Seat Depth Kit 14" - 16" <i>Set Seat Depth at 14"D ____ 15"D ____ 16"D ____</i>	No Charge	
<input type="checkbox"/> 321441	Seat Depth Kit 17" - 19" <i>Set Seat Depth at: 17"D ____ 18"D ____ 19"D ____</i>	No Charge	
<input type="checkbox"/> 321445	Seat Depth Kit 20" - 22" <i>Set Seat Depth at: 20"D ____ 21"D ____ 22"D ____</i>	No Charge	

**CHOOSE SEAT CUSHION:**

Part Number	Description	MSRP - CA	ADP Code
	NOTE: <i>A Corpus II Ergonomic Seat Cushion will be provided for select seat sizes.</i>		
<input type="checkbox"/>	Corpus Ergo Cushion - Leatherette (C3GERGSCA) <i>Part number will be provided on the quote / order confirmation.</i>	587.00	SES1P0005
<input type="checkbox"/>	Corpus Ergo Cushion - Black Mesh (C3GERGSCA) <i>Part number will be provided on the quote / order confirmation.</i>	587.00	SES1P0005
<input type="checkbox"/>	Corpus Ergo Air Cushion (C3GRHOSCA) <i>Features ROHO Dry Floatation technology. The removable ROHO cell pad provides adjustable skin protection while the contoured foam and thigh support wedges help maintain proper positioning and stability. Includes removable stretch fabric cover. Available in widths of 17", 19" and 21" and depths of 17", 19" and 21", as well as 17" x 23". Part number will be provided on the quote / order confirmation.</i>	623.00	SE0001267



## Seat Sizing

## CHOOSE SEAT CUSHION:

Part Number	Description	MSRP - CA	ADP Code
<input type="checkbox"/>	Incontinence Cover for Ergo Cushion <i>Made from a stretch Dartex fabric and suitable for all Ergo cushions including the Corpus Ergo Air. The cover is cut to the contours of the Ergo cushion and is designed to be used as an outer layer of protection. Cover size will be matched to cushion size. Sent as a part. Part Number will be provided on the quote / order confirmation based on the Seat Width and Seat Depth selected.</i>	186.00	SEMCF106L
<input type="checkbox"/>	Omit Seat Cushion	No Charge	

## Armrest Options

## CHOOSE ARMBAR KIT AND MOUNTING POSITION:

Part Number	Description	MSRP - CA	ADP Code
	PLEASE NOTE: <i>Narrow Armbar Kit can be adjusted from 15"W to 19"W between the armrest pads without any additional parts. Wide Armbar Kit can be adjusted from 19"W to 23"W between the armrest pads without any additional parts. Please select both Armbar Kit and Width Setting below.</i>		
<input type="checkbox"/> 325306	Narrow Armbar Kit Inf Adjustable 15"W - 19"W	No Charge	
<input type="checkbox"/> 325308	Wide Armrest Kit Inf Adjustable 19"W - 23"W <i>Required when selecting a 20" wide backrest.</i> Set the Width Between the Armrests to: _____ <i>Must be wider than backrest width to ensure proper fit. Please round to nearest 1".</i>	No Charge	

## CHOOSE ARMREST HEIGHT MOUNTING POSITION:

Part Number	Description	MSRP - CA	ADP Code
<input checked="" type="checkbox"/> 920037	Height Adjustable Armrests - NC SET ARMREST HEIGHT AT: _____ <i>Corpus Armrests are infinitely adjustable from 7" - 15", measured from the seat pan to the top of the armrest pad. No additional parts are needed to get this full range of adjustment. In one inch increments (ex: 8", 9", etc.), please indicate where you would like the armrests set.</i>	No Charge	

## CHOOSE ARMREST SWIVEL LOCKED OR ENABLED:

Part Number	Description	MSRP - CA	ADP Code
	ARMREST MEDIAL SWIVEL NOTE: <i>Every Corpus Armrest can flip backwards for transfers. The armrests also include a feature that allows the armrest pads to swivel medially towards the center of the seating system. This swivel feature can be locked out by simply moving a bolt located under the armrest bracket. Please indicate below whether each armrest should have medial swivel enabled or locked.</i> <i>Note: Medial swivel must be locked on Corpus VS or when ordering arm trough adapters.</i>		
<input type="checkbox"/> 319345	Armrest Bracket - Medial Swivel Assembly Locked <i>Left Armrest Locked: ____ Right Armrest Locked: ____</i>	No Charge	
<input type="checkbox"/> 319347	Armrest Bracket - Medial Swivel Assembly Enabled <i>Left Armrest Enabled: ____ Right Armrest Enabled: ____</i>	No Charge	

## CHOOSE ARMREST PAD RIGHT:

Part Number	Description	MSRP - CA	ADP Code
<input type="checkbox"/> 324280	4"x10" Corpus Armrest Pad - Leatherette	No Charge	
<input type="checkbox"/> 324282	4"x13" Corpus Armrest Pad - Leatherette	No Charge	
<input type="checkbox"/> 324283	4"x16" Corpus Armrest Pad - Leatherette	No Charge	
<input type="checkbox"/> 324299	4"x18" Corpus HD Armrest Pad - Leatherette	No Charge	
<input type="checkbox"/> 326504	4"x10" Corpus Armrest Pad - Black Fabric	No Charge	
<input type="checkbox"/> 324306	4"x13" Corpus Armrest Pad - Black Fabric	No Charge	
<input type="checkbox"/> 325254	4"x16" Corpus Armrest Pad - Black Fabric	No Charge	
<input type="checkbox"/> 324297	4"x18" Corpus HD Armrest Pad - Black Fabric	No Charge	
<input type="checkbox"/> 326523	4"x10" Corpus Arm w/ Stealth Gel Pad	191.00	
<input type="checkbox"/> 326524	4"x13" Corpus Arm w/ Stealth Gel Pad	197.00	
<input type="checkbox"/> 326525	4"x16" Corpus Arm w/ Stealth Gel Pad	202.00	

**Armrest Options**

**CHOOSE ARMREST PAD RIGHT:**

Part Number	Description	MSRP - CA	ADP Code
<input type="checkbox"/> 324309	4"x18" Corpus HD Arm w/ Stealth Gel Pad	210.00	

**CHOOSE ARMREST PAD LEFT:**

Part Number	Description	MSRP - CA	ADP Code
<input type="checkbox"/> 324280	4"x10" Corpus Armrest Pad - Leatherette	No Charge	
<input type="checkbox"/> 324282	4"x13" Corpus Armrest Pad - Leatherette	No Charge	
<input type="checkbox"/> 324283	4"x16" Corpus Armrest Pad - Leatherette	No Charge	
<input type="checkbox"/> 324299	4"x18" Corpus HD Armrest Pad - Leatherette	No Charge	
<input type="checkbox"/> 326504	4"x10" Corpus Armrest Pad - Black Fabric	No Charge	
<input type="checkbox"/> 324306	4"x13" Corpus Armrest Pad - Black Fabric	No Charge	
<input type="checkbox"/> 325254	4"x16" Corpus Armrest Pad - Black Fabric	No Charge	
<input type="checkbox"/> 324297	4"x18" Corpus HD Armrest Pad - Black Fabric	No Charge	
<input type="checkbox"/> 326523	4"x10" Corpus Arm w/ Stealth Gel Pad	191.00	
<input type="checkbox"/> 326524	4"x13" Corpus Arm w/ Stealth Gel Pad	197.00	
<input type="checkbox"/> 326525	4"x16" Corpus Arm w/ Stealth Gel Pad	202.00	
<input type="checkbox"/> 324309	4"x18" Corpus HD Arm w/ Stealth Gel Pad	210.00	

**OPTIONAL ITEMS:**

Part Number	Description	MSRP - CA	ADP Code
<input type="checkbox"/> 320094	Corpus Arm Trough Adapter 1" <i>For aftermarket armrest pads. Please indicate below if they should be installed or on chair sent as parts. Please Note: Medial Swivel must be must be locked when selecting this option. Price shown is for 1(ea). A quantity of 2 will be provided.</i>	174.00	
<input type="checkbox"/> 320095	Corpus Arm Trough Adapter 7/8" <i>For aftermarket armrest pads. Please indicate below if they should be installed or on chair sent as parts. Please Note: Medial Swivel must be must be locked when selecting this option. Price shown is for 1(ea). A quantity of 2 will be provided.</i> Install Arm Trough Adapters: _____ Send Arm Trough Adapters as Parts: _____	174.00	
<input type="checkbox"/> 325819	Armrest Pouch <i>Armrest Pouch Left: _____ Armrest Pouch Right: _____</i>	252.00	
<input type="checkbox"/> 303666	Corpus Armrest Side Support <i>Custom Option. Leatherette Pad (2"T x 8.5"L). Not compatible with Stealth Elbow Support, except when ordered with an 18" armrest. Left Inner Side Support: _____ Right Inner Side Support: _____ Left Outer Side Support: _____ Right Outer Side Support: _____</i>	259.00	SEMCF601L, SEMCF602L
<input type="checkbox"/> 614346	Stealth 3" x 10" Elbow Support - Right <i>Must also select PN 920030.</i>	168.00	
<input type="checkbox"/> 614345	Stealth 3" x 10" Elbow Support - Left <i>Must also select PN 920030.</i>	168.00	
<input type="checkbox"/> 920030	Adjustable Elbow Support Hardware	283.00	
<input type="checkbox"/> 920030	Adjustable Elbow Support Hardware	283.00	

**Legrest Options**

**CORPUS LEGREST OPTIONS:**

Part Number	Description	MSRP - CA	ADP Code
<input type="checkbox"/> 319556	CHOOSE LEGREST TYPE: Corpus Manual Elevating Legrest <i>Charge for Manual Legrest Elevation: _____ Do not Charge for Manual Legrest Elevation: _____</i>	797.00	
<input type="checkbox"/> 318467	Corpus Power Elevating Legrest - 170° <i>Includes mechanical articulation.</i>	3,549.00	

CHOOSE FOOTPLATE STYLE:

## Legrest Options

## CORPUS LEGREST OPTIONS:

Part Number	Description	MSRP - CA	ADP Code
<input type="checkbox"/> 319179	2pc Wide Corpus Footplates - 9"D x 7.5"W	No Charge	
<input type="checkbox"/> 319636	2pc Narrow Corpus Footplates - 9"D x 6"W	379.00	
<input type="checkbox"/> 322004	1pc Narrow Corpus Footpl - 9.75"D x 12.5"W <i>Footplate is 12.5"W in front, 8"W in back.</i>	379.00	
<b>CALF SUPPORT OPTIONS:</b>			
<input type="checkbox"/> 323706	UT Calf Support Kit (B) w Mesh Cover & Adj Hardware <i>3.5"W x 4.5"T Pad.</i>	No Charge	
<input type="checkbox"/> 323711	UT Calf Support Kit (J) w Mesh Cover & Adj Hardware <i>5"W x 6.5"T Pad.</i>	No Charge	
<input type="checkbox"/> 323704	UT Calf Support Kit (A) w Mesh Cover & Adj Hardware <i>5.5"W x 6"T Pad.</i>	No Charge	
<input type="checkbox"/> 323709	UT Calf Support Kit (H) w Mesh Cover & Adj Hardware <i>7"W x 7"T Curved Pad.</i>	392.00	
<input type="checkbox"/> 323707	UT Calf Support Kit (D) w Mesh Cover & Adj Hardware <i>8"W x 8"T Pad.</i>	392.00	

## CORPUS VS LEGREST OPTIONS (FOR SELECT ANTERIOR TILT PACKAGES):

Part Number	Description	MSRP - CA	ADP Code
<b>CORPUS VS LEGREST AND KNEE SUPPORT:</b>			
<i>NOTE: The Corpus VS Legrest is only available with 20°, 30° and 45° of Active Reach.</i>			
<input type="checkbox"/> 325010	Corpus VS Power Elevating Legrest - 180° <i>Includes power articulation.</i>	3,645.00	
<input type="checkbox"/> 323856	Knee Support/Abductor - Corpus VS <i>Select when ordering Corpus VS Legrest. Must also select 920115.</i>	209.00	
<input type="checkbox"/> 920115	Adj Removable Knee Support Hardware	305.00	
<b>CHOOSE FOOTPLATE STYLE:</b>			
<input type="checkbox"/> 323855	1pc Corpus VS Footplate - 9"D x 16.5"W	No Charge	
<input type="checkbox"/> 323851	2pc Wide Corpus VS Footplates - 9"D x 7.5"W	379.00	
<input type="checkbox"/> 323852	2pc Narrow Corpus VS Footplates - 9"D x 6"W	379.00	
<b>CALF SUPPORT OPTIONS:</b>			
<input type="checkbox"/> 325126	UT Calf Support Kit (J) w Mesh Cover & Adj Hardware <i>5"W x 6.5"T Pad. Compatible with footplate heights of 10" - 18".</i>	No Charge	
<input type="checkbox"/> 323829	UT Calf Support Kit (A) w Mesh Cover & Adj Hardware <i>5.5"W x 6"T Pad. Compatible with footplate heights of 10" - 18".</i>	No Charge	
<input type="checkbox"/> 323828	UT Calf Support Kit (H) w Mesh Cover & Adj Hardware <i>7"W x 7"T Curved Pad. Compatible with footplate heights of 11" - 18".</i>	420.00	
<input type="checkbox"/> 323827	UT Calf Support Kit (D) w Mesh Cover & Adj Hardware <i>8"W x 8"T Pad. Compatible with footplate heights of 12" - 18".</i>	420.00	

## FOOTPLATE HEIGHT POSITION:

Part Number	Description	MSRP - CA	ADP Code
	SET FOOTPLATE HEIGHT AT: _____ <i>This is measured from the SEAT PAN to the top of the FOOTPLATE, when the legrest assembly is at 90°. DO NOT INCLUDE CUSHION THICKNESS IN YOUR MEASUREMENT. In one inch increments, please indicate desired footplate height. NOTE: Corpus footplate height range is 6"-17" and Corpus VS height range is 9"-17". Height may need to be reduced for shipping. Calf Pads are not available at heights less than 9" (10" on VS Legrest).</i>		

## OPTIONAL ITEMS:

Part Number	Description	MSRP - CA	ADP Code
<input type="checkbox"/> 614446	Stealth Gel Pads - Narrow Footplates (Pair)	399.00	

## Legrest Options

## OPTIONAL ITEMS:

Part Number	Description	MSRP - CA	ADP Code
<input type="checkbox"/> 614388	Stealth Gel Pads - Wide Footplates (Pair)	411.00	

## Accessories

<input type="checkbox"/> 330501	Bodypoint Non-Padded Hip Belt - Medium	No Charge	
<input type="checkbox"/> 330503	Bodypoint Padded Hip Belt - Medium	144.00	SERND0030
<input type="checkbox"/> 330502	Bodypoint Non-Padded Hip Belt - Large	116.00	SERND0030
<input type="checkbox"/> 330504	Bodypoint Padded Hip Belt - Large	197.00	SERND0030
<input type="checkbox"/> 305003	Positioning Belt Retractable	298.00	SERND0030
<input type="checkbox"/> 310189	Positioning Belt Retractable Long	323.00	SERND0030

## CHOOSE ONE OF THE FOLLOWING:

Part Number	Description	MSRP - CA	ADP Code
<input type="checkbox"/> 325711	Permobil 7"W Headrest for Link Hardware <i>Must also select 329886 or 329887.</i>	439.00	SEH1P0005
<input type="checkbox"/> 325712	Permobil 10"W Headrest for Link Hardware <i>Must also select 329886 or 329887.</i>	439.00	SEH1P0005
<input type="checkbox"/> 325713	Permobil 14"W Headrest for Link Hardware <i>Must also select 329886 or 329887.</i>	439.00	SEH1P0005
<input type="checkbox"/> 329886	Adj/Removable Headrest Link Hardware - Standard	305.00	
<input type="checkbox"/> 329887	Adj/Removable Headrest Link Hardware - Extended <i>Includes an additional 3" link to extend headrest pad range.</i>	329.00	
<input type="checkbox"/> 322147	Permobil 10"W Headrest Leatherette <i>Must also select PN 920047.</i>	621.00	SEH1P0005
<input type="checkbox"/> 324310	Permobil Ergo Shape Headrest Black Fabric <i>Must also select PN 920047.</i>	513.00	SEH1P0005
<input type="checkbox"/> 920047	Adjustable Removable Headrest Hardware	133.00	SEHND0020
<input type="checkbox"/> 614838	Universal Headrest Adapter <i>Please Note: Recline is limited to 160° or less (depending on backrest height) when ordering a Universal Headrest Adapter on a 23" or 24" tall Corpus backrest.</i>	250.00	

## CHOOSE ONE OF THE FOLLOWING:

Part Number	Description	MSRP - CA	ADP Code
<input type="checkbox"/> 329547	Permobil 3" x 5" Lateral Supports (Pair) <i>Must also select Lateral Link Hardware. Only compatible with Corpus backrest shells.</i>	459.00	SELND0010
<input type="checkbox"/> 329548	Permobil 5" x 7" Lateral Supports (Pair) <i>Must also select Lateral Link Hardware. Only compatible with Corpus backrest shells.</i>	459.00	SELND0010
<input type="checkbox"/> 329549	Adjustable Removable Lateral Link Hardware - (Pair)	629.00	SEICF203L

## OPTIONAL ITEMS:

Part Number	Description	MSRP - CA	ADP Code
<input type="checkbox"/> 323480	UT Amp. Support Pad (H) w/ Mesh Cover - (1 Each) <i>7"W x 7"T Curved Pad. Must also select PN 920023.</i>	373.00	SEFND0005
<input type="checkbox"/> 323481	UT Amp. Support Pad (D) w/ Mesh Cover - (1 Each) <i>8"W x 8"T Flat Pad. Must also select PN 920023.</i>	373.00	SEFND0005
<input type="checkbox"/> 920023	Adj Removable - Amp Support Hardware <i>Please indicate if you would like the amputation support mounted on the right or left side of the seat frame: Mount Amp Suppt on Left: _____ Mount Amp Suppt on Right: _____</i>	133.00	
<input type="checkbox"/> 324587	UT Thigh/Hip Support Pad (B) w/ Mesh Cover - Tall <i>4.5"W x 3.5"T Pad w/ 4" Tall Mounting Bracket. Must also select PN 920024. Price shown is for 1(ea). A quantity of 2 will be provided.</i>	247.00	
<input type="checkbox"/> 324694	UT Hip Suppt (A) w/ Mesh Cover - Tall <i>6"W x 5.5"T Pad w/ 4" Tall Mounting Bracket. Must also select PN 920024. Price shown is for 1(ea). A quantity of 2 will be provided.</i>	247.00	

## Accessories

## OPTIONAL ITEMS:

Part Number	Description	MSRP - CA	ADP Code
<input type="checkbox"/> 324589	UT Thigh/Hip Support Pad (J) w/ Mesh Cover - Tall <i>6.5"W x 5"T Pad w/ 4" Tall Mounting Bracket. Must also select PN 920024. Price shown is for 1(ea). A quantity of 2 will be provided.</i>	247.00	SELND0015
<input type="checkbox"/> 324588	UT Thigh/Hip Support Pad (C) w/ Mesh Cover - Tall <i>8"W x 3.5"T Pad w/ 4" Tall Mounting Bracket. Must also select PN 920024. Price shown is for 1(ea). A quantity of 2 will be provided.</i>	247.00	SELND0015
<input type="checkbox"/> 324570	UT Thigh/Hip Support Pad (G) w/ Mesh Cover - Tall <i>12"W x 3.5"T Pad w/ 4" Tall Mounting Bracket. Must also select PN 920024. Price shown is for 1(ea). A quantity of 2 will be provided.</i>	247.00	SELND0015
<input type="checkbox"/> 920024	Adj Removable - Thigh Support Hardware <i>Price shown is for 1(ea). A quantity of 2 will be provided.</i>	67.00	
<input type="checkbox"/> 920025	Adj Removable - Hip Support Hardware <i>Price shown is for 1(ea). A quantity of 2 will be provided.</i>	67.00	
<input type="checkbox"/> 320381	Abductor Kit Corpus - Leatherette <i>Must also select PN 920022.</i>	369.00	
<input type="checkbox"/> 920022	Adj Removable - Abductor Hardware	133.00	
<input type="checkbox"/> 323712	Transfer Handles Corpus Seat	513.00	
<input type="checkbox"/> 319397	Push Handles Corpus Seat	452.00	
<input type="checkbox"/> 324167	RAM Self-Leveling Cup Holder	88.00	
<input type="checkbox"/> 324166	RAM X-Grip Phone Holder	88.00	
<input type="checkbox"/> 324139	RAM X-Grip 10" Tablet Holder w/Swing Away Mount <i>Compatible with most 10" screen tablets with or without a protective case. Includes adjustable swing away mount. Note: Customization may be required if ordered along with an Omni.</i>	489.00	
<input type="checkbox"/> 325191	USB Charger <i>Provides up to 1.5A (at 5V) to charge most smartphones and tablets. The charger features an isolated circuit design for improved safety; and by pulling power from the ICS, prevents battery drainage when the chair is powered off. Includes mounting hardware and ICS connector. Please select armrest for mounting, keeping in mind that charger may need to be relocated. Mount on left armrest: _____ Mount on right armrest: _____</i>	209.00	
<input type="checkbox"/> 315797	Medical Necessities Bag <i>Includes hardware to attach bag to backrest. Please Note: Recline is limited to 165° or less (depending on backrest height) when ordering with a 23" or 24" tall Corpus backrest.</i>	272.00	
<input type="checkbox"/> 320316	Medical Necessities Bag Hooks <i>Please Note: Recline is limited to 165° or less (depending on backrest height) when ordering with 23" or 24" tall Corpus backrest.</i>	96.00	
<input type="checkbox"/> 330507	Bodypoint Monoflex Medium Chest Support Belt Kit <i>18.25"L x 4"W Pad Size. Item will be shipped as a part.</i>	187.00	SERND0001
<input type="checkbox"/> 330508	Bodypoint Monoflex Large Chest Support Belt Kit <i>21.5"L x 4.5"W Pad Size. Item will be shipped as a part.</i>	187.00	SERND0001
<input type="checkbox"/> 330509	Bodypoint Monoflex X-Lrge Chest Support Belt Kit <i>24.5"L x 5"W Pad Size. Item will be shipped as a part.</i>	187.00	SERND0001
<input type="checkbox"/> 320549	Upper Extremity Support Small C3G <i>Must also select PN 920029. Mount on Left: _____ Mount on Right: _____</i>	455.00	SETND0015
<input type="checkbox"/> 319559	Upper Extremity Support Angle Adjustable. 14" W <i>Must also select PN 920029.</i>	455.00	SETND0010
<input type="checkbox"/> 319558	Upper Extremity Support Angle Adjustable. 16" W <i>Must also select PN 920029.</i>	455.00	SETND0010
<input type="checkbox"/> 319560	Upper Extremity Support Angle Adjustable. 18" W <i>Must also select PN 920029.</i>	455.00	SETND0010
<input type="checkbox"/> 323714	Upper Extremity Support Joystick Cutout 16" W <i>Must also select PN 920029.</i>	455.00	SETND0010
<input type="checkbox"/> 323715	Upper Extremity Support Joystick Cutout 18" W <i>Must also select PN 920029.</i>	455.00	SETND0010



**Accessories**

**OPTIONAL ITEMS:**

Part Number	Description	MSRP - CA	ADP Code
<input type="checkbox"/> 321272	Upper Extremity Support Angle Adjust - UniTrack <i>Must also select PN 920029.</i>	455.00	
<input type="checkbox"/> 920029	Adj Removable - Tray Hardware	133.00	
<input type="checkbox"/> 323701	Oxygen Holder Corpus	331.00	WEPX
<input type="checkbox"/> 323700	Crutch Holder Corpus	344.00	

**Alternative Drive Systems**

**ALTERNATIVE PROPORTIONAL JOYSTICK SYSTEMS:**

Part Number	Description	MSRP - CA	ADP Code
<input type="checkbox"/> 613244	Mini-Joystick for Omni / IOM (9-pin connector) <i>Operating Force = 13 g Neutral to Full Deflection = 5mm</i>	3,031.00	
<input type="checkbox"/> 613440	Compact Joystick for Omni / IOM (9-pin connector) <i>Operating Force = 220 g Neutral to Full Deflection = 28mm Joystick includes two 3.5mm mono jacks for Mode and On/Off. Includes Standard Joystick Knob and 1.75" (42mm) Foam Ball Knob.</i>	1,479.00	
<input type="checkbox"/> 613441	Compact Joystick LITE for Omni / IOM (9-pin conn) <i>Operating Force = 50 g Neutral to Full Deflection = 17mm Joystick includes two 3.5mm mono jacks for Mode and On/Off. Includes Standard Joystick Knob and 1.75" (42mm) Foam Ball Knob.</i>	1,692.00	

**ALT PROPORTIONAL JOYSTICK MOUNTING HARDWARE:**

Part Number	Description	MSRP - CA	ADP Code
<b>NOTE:</b> <i>Please see the PERMOfix Mounting Solutions, within the Alt Drive Accessories section, for additional mounting options.</i>			
<input type="checkbox"/> 330322	Alt JS Swing Away Mount Kit - Armrest <i>Mount in front of Left Arm: _____ Mount in front of Right Arm: _____</i>	287.00	
<input type="checkbox"/> 324915	Alt JS / Omni Combo Swing Away Mount Kit - Armrest <i>Used to mount Omni and Alternative Proportional Joystick to the same armrest. Includes Swing Away Omni/Tablet Mount, with a longer mounting rod, two PERMOfix modules, one 3" rod and one 4" rod.</i> <i>Mount in front of Left Arm: _____ Mount in front of Right Arm: _____</i>	582.00	
<input type="checkbox"/> 323688	Alternative JS Adjustable Mount Kit - Midline <i>Mount on Left Arm: _____ Mount on Right Arm: _____</i>	873.00	
<input type="checkbox"/> 323674	Alt JS Adjustable Mount Kit - Chin Boom <i>Swing-Away / Removable Chin Boom mounted on right of backrest shell. Custom charge may be applied if ordering lower backrest heights.</i>	931.00	
<input type="checkbox"/> 614931	Alt JS Adjustable Mount Kit - Haims Harness Lrg <i>For neck width up to 19cm.</i>	665.00	
<input type="checkbox"/> 614932	Alt JS Adjustable Mount Kit - Haims Harness Sml <i>For neck width up to 14cm.</i>	665.00	

**ALT PROPORTIONAL JOYSTICK ACCESSORIES:**

Part Number	Description	MSRP - CA	ADP Code
<input type="checkbox"/> 613741	Dual Switches for Compact Joystick - Fixed Mount <i>Includes two pushbutton switches mounted on rigid Switches can be used for Mode &amp; On/Off. Switch operation force = 130 g</i> <i>Mount at front of Joystick: _____ Mount at rear of Joystick: _____</i> <i>Mount at left of Joystick: _____ Mount at right of Joystick: _____</i>	276.00	
<input type="checkbox"/> 613597	Dual Switches for Compact Joystick -Gooseneck Mount <i>Includes two pushbutton switches mounted on adjustable gooseneck arms. Switches are positioned on Left and Right sides of Compact Joystick. Switches can be used for Mode &amp; On/Off. Switch operation force = 130 g</i>	502.00	

**SIP & PUFF SYSTEMS:**

Part Number	Description	MSRP - CA	ADP Code
-------------	-------------	-----------	----------

## Alternative Drive Systems

### SIP & PUFF SYSTEMS:

Part Number	Description	MSRP - CA	ADP Code
<input type="checkbox"/> 614928	Sip & Puff System for Omni <i>Requires Omni Control Module. Cannot use this part with an Input Module. Includes breath tube kit.</i>	1,778.00	WEX5
<input type="checkbox"/> 329842	Sip & Puff System with Fiber Optic Opti-Stop <i>Requires Omni Control Module. Cannot use this part with an Input Module. Includes breath tube kit.</i>	2,366.00	WEX5
<input type="checkbox"/> 614559	Sip & Puff Kit for Total Control Head Array <i>For use with Total Control Head Array. Requires Omni Control Module. Cannot be used with an Input Module. Includes breath tube kit.</i>	982.00	
<input type="checkbox"/> 920073	Adjustable Mounting Hrdw - Sip & Puff <i>Required when selecting PN 614928, 329842 or 614559.</i>	287.00	

### ALTERNATIVE PROPORTIONAL CONTROL SYSTEMS:

Part Number	Description	MSRP - CA	ADP Code
<input type="checkbox"/> 614487	Proportional Foot Control <i>Mount on Left: _____ Mount on Right: _____</i>	3,121.00	
<input type="checkbox"/> 920074	Adjustable Mounting Hrdw - Foot Control <i>Required when selecting PN 614487.</i>	287.00	

## Total Control Head Array

### CONFIGURE HEAD ARRAY:

Part Number	Description	MSRP - CA	ADP Code
<input type="checkbox"/> 330548	Permobil Head Array - 3-Switch System <i>Features an occipital pad sensor and two stainless steel arms configured with proximity sensors. The arms are mounted to the upper attachment points on the head array. Longer arms are also included to extend the sensor adjustment range. ESP recommended. Choose Omni Control Module (613581) separately.</i>	5,485.00	WEX4
<input type="checkbox"/> 319505	Program R-net Omni for 5-Switch Input Device <i>Select additional arms/sensors below and indicate the intended use of each sensor in the Notes &amp; Comments section. Choose Omni Control Module (613581) or Input Module (920021) separately.</i>	No Charge	
CHOOSE OCCIPITAL PAD SIZE (REQUIRED):			
<input type="checkbox"/> 614108	Small Occipital Pad for Head Array - 4.25"W x 3"H	No Charge	
<input type="checkbox"/> 614109	Large Occipital Pad for Head Array - 6"W x 4"H	No Charge	
CHOOSE OCCIPITAL MOUNT (REQUIRED):			
<input type="checkbox"/>	Fixed Mount for Head Array Occipital Pad	No Charge	
<input type="checkbox"/> 614097	Adj. Swivel Mount for Head Array Occipital Pad <i>Enables Occipital Pad to be angled up/down and left/right independent of the lateral arms.</i>	491.00	
CHOOSE BACKREST MOUNTING HARDWARE (REQUIRED):			
<input type="checkbox"/> 614236	Stealth TWB480 Backrest Mounting Hdwr - 10" stem	305.00	
<input type="checkbox"/> 614259	Stealth TWB480 Backrest Mounting Hdwr - 6" stem	305.00	
<input type="checkbox"/> 614234	Stealth TWB480 Backrest Mnt Hdwr with Flipback	504.00	
<input type="checkbox"/> 614235	Stealth HMO475 MultiAxis Backrest Mounting Hdwr	258.00	
<input type="checkbox"/> 614260	Stealth TWB480 Backrest Mnt Hdwr w SubOccpl Mount	548.00	
<input type="checkbox"/> 614258	Stealth 600 SubOccipital Pad - 8.5" W x 2.75" H <i>Must select when PN 614260 is ordered.</i>	111.00	
SELECT ADDITIONAL LEFT ARM PACKAGE (OPTIONAL):			
<input type="checkbox"/>	Add Left Lateral Arm w Proximity Sensor <i>Choose Lateral Arm Adjustment Range: Short (5.5" - 10"): _____ Long (5.5" - 14"): _____</i>	404.00	

## Total Control Head Array

### CONFIGURE HEAD ARRAY:

Part Number	Description	MSRP - CA	ADP Code
<input type="checkbox"/>	Add Left Lateral Arm w Egg Switch <i>Choose Lateral Arm Adjustment Range: Short (5.5" - 10"): _____ Long (5.5" - 14"): _____</i>	269.00	
<input type="checkbox"/>	Add Left Lateral Arm w Pad (No Sensor or Switch) <i>Choose Lateral Arm Adjustment Range: Short (5.5" - 10"): _____ Long (5.5" - 14"): _____</i>	164.00	
<b>SELECT ADDITIONAL RIGHT ARM PACKAGE (OPTIONAL):</b>			
<input type="checkbox"/>	Add Right Lateral Arm w Proximity Sensor <i>Choose Lateral Arm Adjustment Range: Short (5.5" - 10"): _____ Long (5.5" - 14"): _____</i>	404.00	
<input type="checkbox"/>	Add Right Lateral Arm w Egg Switch <i>Choose Lateral Arm Adjustment Range: Short (5.5" - 10"): _____ Long (5.5" - 14"): _____</i>	269.00	
<input type="checkbox"/>	Add Right Lateral Arm w Pad (No Sensor or Switch) <i>Choose Lateral Arm Adjustment Range: Short (5.5" - 10"): _____ Long (5.5" - 14"): _____</i>	164.00	

## Alt Drive Systems Setup

### CHOOSE OMNI OR INPUT MODULE (REQUIRED):

Part Number	Description	MSRP - CA	ADP Code
<input type="checkbox"/> 613581	Omni Control Module for R-net w/ Infrared <i>Must also select mounting hardware, either 315315 or 324915. Please specify which type of alternative drive input you will be using if not ordering it from Permobil. Device: _____</i>	2,478.00	
<input type="checkbox"/> 315315	Adjustable Swing-Away Hardware - Omni/Tablet <i>Mount in front of Left Arm: _____ Mount in front of Right Arm: _____</i>	282.00	
<input type="checkbox"/> 920021	Input Module (IOM) for R-net Control System <i>This device allows a proportional or switched drive input to be connected to an R-net Control System without the use of an Omni. It must be ordered in conjunction with an R-net Remote Joystick. It does NOT support any of the following input devices: Sip-n-Puff, 3-Switch Head Array or RIM control.</i>	1,749.00	

### CHOOSE MOUNTING POSITION OF SEAT CONTROL SWITCHBOX:

Part Number	Description	MSRP - CA	ADP Code
<input type="checkbox"/>	Mount Seat Control Switchbox on Backrest	No Charge	

## Alt Drive Accessories

### PERMOFIX MOUNTING SOLUTIONS:

Part Number	Description	MSRP - CA	ADP Code
<input type="checkbox"/> 323835	PERMOfix UniTrack Fixed Mount Kit <i>Includes two base joints, two 4" rods and one 8" rod.</i>	309.00	
<input type="checkbox"/> 323836	PERMOfix UniTrack Removable Mount Kit <i>Includes quick release mount and thigh support hardware, two base joints and three different length rods.</i>	543.00	
<input type="checkbox"/> 330092	PERMOfix UniTrack Swing Away Mount Kit <i>Includes PERMOfix Swing Away module, two base joints, one 2" rod and one 4" rod.</i>	783.00	
<input type="checkbox"/> 614103	Proximity Sensor + Swivel Assembly <i>Only available with the Total Control Head Array. Swivel assembly allows sensor to be mounted with PERMOfix or directly to lateral arm of Total Control Head Array. Since Proximity Sensors require power, they CANNOT be used as On/Off switches.</i>	404.00	
<input type="checkbox"/> 614107	Egg Switch + Swivel Mount Assembly <i>Swivel assembly allows switch to be mounted with PERMOfix or directly to lateral arm of Total Control Head Array.</i>	269.00	
<input type="checkbox"/> 614105	Solid Oval Pad + Swivel Assembly <i>Swivel assembly allows pad to be mounted with PERMOfix or directly to lateral arm of Total Control Head Array.</i>	164.00	

## Alt Drive Accessories

## PERMOFIX MOUNTING SOLUTIONS:

Part Number	Description	MSRP - CA	ADP Code
<input type="checkbox"/> 614113	Oval Plate + Swivel Assembly - No Switch <i>Blank plate provides mounting option for Micro Lite, Buddy Button or other switch options. Swivel assembly allows pad to be mounted with PERMOfix or directly to lateral arm of Total Control Head Array.</i>	142.00	

## OPTIONAL ITEMS:

Part Number	Description	MSRP - CA	ADP Code
<input type="checkbox"/> 613306	Micro Lite Switch	178.00	
<input type="checkbox"/> 613303	Buddy Button - Red	134.00	
<input type="checkbox"/> 614356	Egg Switch Black	148.00	
<input type="checkbox"/> 613308	Mini Cup Switch	131.00	
<input type="checkbox"/> 614842	Proximity Switch	621.00	
<input type="checkbox"/> 614837	Mini Joystick 5 Switch OC for Omni/IOM (9-pin)	1,055.00	
<input type="checkbox"/> 613921	Touch Contact with CR2032 Battery	299.00	
<input type="checkbox"/> 330019	Dome Switch Backrest Mount w Gooseneck <i>Mount on Left: _____ Mount on Right: _____ Mounts at top of backrest shell. Custom charge may be applied if ordering with short backrest heights. Kit Includes: One PERMOfix joint with Flat Plate, two PERMOfix joints with 45° swivel adjustment, one 150mm rod, one 200mm rod and one 225mm gooseneck switch with 1.2m cable and mounting clamp.</i>	553.00	
<input type="checkbox"/> 330021	Dome Switch Backrest Swing Away Mount w Gooseneck <i>Mount on Left: _____ Mount on Right: _____ w 90° Swing-Away Function. Mounts at top of back. Custom charge applies if ordered with short back. Kit Includes: One PERMOfix joint with Flat Plate, one PERMOfix joint with 45° swivel adjustment, two 6" rods, one 90° Swing-Away joint, one 9" gooseneck switch w 50" cable and mounting clamp.</i>	830.00	
<input type="checkbox"/> 311814	3.5mm Ext. Cable for Sensor/Switch 0.3m Long	38.00	
<input type="checkbox"/> 311810	3.5mm Ext. Cable for Sensor/Switch 1m Long	42.00	
<input type="checkbox"/> 311811	3.5mm Ext. Cable for Sensor/Switch 2m Long	46.00	
<input type="checkbox"/> 311812	3.5mm Ext. Cable for Sensor/Switch 3m Long	50.00	
<input type="checkbox"/> 311813	3.5mm Y-splitter Cable for Sensor/Switch <i>Allows two Proximity Sensors/Mechanical Switches to connect to the same direction jack.</i>	42.00	

Notes & Additional Instructions:

**Chair Order Policy:**

You are highly encouraged to contact Permobil Sales with ANY questions regarding proper completion of our order forms at 1-800-736-0925. The order cannot be processed without the following:

- A client name or code, the client's measurements, the client's weight
- The order form must also be signed by a person authorized by your company to acknowledge the items selected. Permobil, LTD. is not responsible for configuration or size discrepancies resulting from customer errors on the order form.
- Order confirmations will be provided by Permobil summarizing the items selected. Please review this carefully.
- Once the order is shipped, any changes are subject to the returns policy and may be prohibited.

**Order Acknowledgement:**

I, \_\_\_\_\_, am an agent of the medical equipment provider named on this order form and I have the authority to contract for the purchase of powered wheelchairs and related parts on behalf of said provider. I acknowledge that I have reviewed this order and that it is complete and accurate to the best of my knowledge. I further acknowledge that any changes to this order after submission of this order form are subject to the returns policy and may carry additional charges.

Signed: \_\_\_\_\_

Printed Name: \_\_\_\_\_

Title: \_\_\_\_\_

Date: \_\_\_\_\_

NOTE: All specifications and prices are subject to change without notice. Please note that prices displayed are only valid if ordered with the wheelchair. Any parts ordered for service or future alterations will carry different pricing and warranty. *The ADP codes provided are not intended to be billing or legal advice, rather our interpretation of the code definitions. Use of the codes does not ensure coverage or payment for the item. For coverage information, verify the policy of the appropriate payer.*