


M300 Corpus 3G

permobil



Display ADP Codes on order/quote:

Three column pricing on order/quote:

Please send order/quote to fax#: (800) 231-3256

Email to: sales@permobil.com

Dealer Information

* Contact: _____

Dealer Code: _____

* Dealer Name: _____

Address: _____

* City: _____

* Province: _____

* Phone#: _____

* Fax#: _____

PO#: _____

Email Address: _____

Client Information

Permobil recommends that the client is evaluated by a certified rehab specialist.

* First Name: _____

* Last Name: _____

Diagnosis: _____

Funding Source: _____

Client Age: _____

Client Measurements

* Weight: _____

* Height: _____

A) Top of Shoulders: _____

B) Chest Depth: _____

C) Chest Width: _____

D) Seat Depth: _____

E) Top of Head: _____

F) Elbow to Hand: _____

G) Seat Pan to Elbow: _____

H) Hip Width: _____

*I) Knee to Foot: _____

Permobil Ltd.

12 - 75 Mary St.

Aurora, Ontario

L4G 1G3

Tel: (800) 736-0925

Fax: (800) 231-3256

www.permobil.com

Prices effective March 1, 2017 and payable in Canadian dollars.

Chair Model and Colors

CHOOSE BASE:

<input checked="" type="checkbox"/>	325253	M300 Power Base (6 mph) <i>Price Includes: M300 Power Base, Shock Absorbing Suspension System, R-net Controller, Tie-Down Hardware For Strap Systems, Charger, Flat Free Tires, and Standard Positioning Belt. Must Choose Flat Free or Pneumatic Tires Below. Weight Capacity = 300lbs. STFH = 17.5"</i>	8,069.00	WP0000987
<input type="checkbox"/>	920016	7.5mph High Speed Upgrade Package - M300 <i>Price Includes 7.5mph High Speed Motor Package.</i>	2,514.00	

CHOOSE WHEEL PACKAGE:

Part Number	Description	MSRP - CA	ADP Code
<input type="checkbox"/>	Grey Wheel Package - MWD <i>Includes grey tires on silver casters and drive wheel rims. Hub cap color is selected below. Please select caster size and drive wheel type: 6" Casters: _____ 8" Casters: _____ Flat free drive tires: _____ Pneumatic drive tires: _____</i>	No Charge	
<input type="checkbox"/>	Black Wheel Package - MWD <i>Includes black tires on silver casters and drive wheel rims. Hub cap color is selected below. Please select caster size and drive wheel type: 6" Casters: _____ 8" Casters: _____ Flat free drive tires: _____ Pneumatic drive tires: _____</i>	No Charge	
<input type="checkbox"/>	Blackout Wheel Package - MWD <i>Custom Option. Includes black tires and custom matte black casters, drive wheels and hub caps. Please select caster size and drive wheel type: 6" Casters: _____ 8" Casters: _____ Flat free drive tires: _____ Pneumatic drive tires: _____</i>	392.00	

OPTIONAL DRIVE TIRES:

Part Number	Description	MSRP - CA	ADP Code
<input type="checkbox"/>	317788 Tri-Spoke Split Rim w/ Snow Tire Pneumatic 3x8" <i>Snow tires are designed for use in cold weather conditions. The nylon compound in these tires remains flexible at low temperatures, but wears more easily. Therefore, the snow tires (sold as wheel assemblies) are not recommended for year-long use and must be ordered along with standard drive wheels. Not available with black tires or black rims. Price shown is for 1(ea). A quantity of 2 will be provided.</i>	263.00	

CHOOSE SHROUD COLOR:

Part Number	Description	MSRP - CA	ADP Code
<input type="checkbox"/>	319900071 M300 Shroud - Onyx Black	No Charge	
<input type="checkbox"/>	319900072 M300 Shroud - Platinum Silver	No Charge	
<input type="checkbox"/>	319900055 M300 Shroud - Sublime Green	No Charge	
<input type="checkbox"/>	319900077 M300 Shroud - Cobalt Blue	No Charge	
<input type="checkbox"/>	319900079 M300 Shroud - Phantom Grey	No Charge	
<input type="checkbox"/>	319900073 M300 Shroud - Pop Star Pink	No Charge	
<input type="checkbox"/>	319900078 M300 Shroud - Purple Rain	No Charge	
<input type="checkbox"/>	319900080 M300 Shroud - Volcano Red	No Charge	

CHOOSE HUBCAP COLOR:

Part Number	Description	MSRP - CA	ADP Code
<input type="checkbox"/>	Match Hub Cap Color to Shroud: _____ <i>Note: Hub caps will be matched to shroud color if no selection is made.</i>		
<input type="checkbox"/>	Solid Color Hub Caps - Pair <i>Platinum Silver: _____ Onyx Black: _____ Cobalt Blue: _____ Phantom Gray: _____ Pop Star Pink: _____ Volcano Red: _____ Purple Rain: _____ Titan Blue: _____ Sublime Green: _____ Crystal Blue: _____ Vintage Wine: _____ Extreme Yellow: _____ Chili Pepper Red: _____</i>	No Charge	

Base Options

MUST CHOOSE EITHER SEAT ELEVATOR OR FIXED SEAT TUBE:

Part Number	Description	MSRP - CA	ADP Code
-------------	-------------	-----------	----------

Base Options

MUST CHOOSE EITHER SEAT ELEVATOR OR FIXED SEAT TUBE:

Part Number	Description	MSRP - CA	ADP Code
<input type="checkbox"/> 323012	Power Adjustable Seat Height - 8" Travel	3,008.00	
<input type="checkbox"/> 320240	Seat Tube Fixed	No Charge	

MUST CHOOSE BATTERIES:

Part Number	Description	MSRP - CA	ADP Code
<input type="checkbox"/> 602130	Batteries Grp 24 (72Ah) Sealed Gel Installed <i>Price is for 1(ea). A quantity of 2 will be provided.</i>	483.00	
<input type="checkbox"/> 603073	Batteries Grp 34 (60Ah) Sealed Gel Installed <i>Price is for 1(ea). A quantity of 2 will be provided.</i>	476.00	

OPTIONAL ITEMS:

Part Number	Description	MSRP - CA	ADP Code
<input type="checkbox"/>	Seat Stabilizer Option	71.00	
<input type="checkbox"/> 320252	Enhanced Steering Perform Unit (ESP) R-net <i>Recommended item with switched drive inputs for improved tracking.</i>	1,508.00	WEPQ
<input type="checkbox"/> 320178	LED Light Kit for M300 / M400 w/ R-net Electronics	1,153.00	
<input type="checkbox"/> 320273	Compact Laptop Ventilator Tray - Corpus Seat <i>Mounts to rear of Corpus backrest hinge. Accommodates Pulmonetic Systems / CareFusion LTV Series and similar sized Laptop Ventilators. The tray is 330mm Tall, adjustable Width between 200mm - 330mm, and adjustable Depth between 50mm - 165mm. Overall backrest angle limited to 150° when mounting this ventilator tray. Custom Option.</i>	1,682.00	
<input type="checkbox"/> 320399	Battery Tray for Laptop Ventilator Tray <i>Mounts to rear of Laptop Ventilator Tray. Measures 290mm W x 170mm H x 42mm D. Accommodates Pulmonetic Systems Lithium-Ion Battery Pack, CareFusion SprintPack Battery System, or Richardson PowerTech Vent Power Center. Custom Option.</i>	269.00	
<input type="checkbox"/> 320390	O2 Holder for Laptop Ventilator Tray <i>Mounts to Left or Right Side of Laptop Vent Tray. O2 Holder is 175mm Tall. Width infinitely adjustable 90mm - 170mm. Depth infinitely adjustable 95mm - 125mm. Custom Option.</i>	380.00	WEPX
<input type="checkbox"/> 613399	Programming Key for R-net - Permobil Version	524.00	
<input type="checkbox"/> 614390	Stealth BPAP Vent Tray - MWD Base <i>Designed for BPAP systems, the tray has a 12"W x 7"D surface that can support up to 40 lbs. A 1" lip and adjustable strap help secure contents. The tray mounts to the rear tie-down loops, is tapered to fit between the rear swing arms and extends 1" beyond the trailing casters.</i>	766.00	

Joystick Options

CHOOSE JOYSTICK TYPE:

Part Number	Description	MSRP - CA	ADP Code
	PLEASE NOTE: <i>For additional input devices, refer to Alternative Drive Systems Section.</i>		
<input type="checkbox"/> 615264	Permobil Joystick Module <i>Designed exclusively for Permobil, this R-net color joystick features: control of up to 4 Bluetooth devices, Infrared control, programmable buttons and a programmable switch jack; all in a compact, durable Aluminum housing.</i>	1,185.00	
<input type="checkbox"/> 920147	Permobil Joystick Module NC	No Charge	
<input type="checkbox"/> 615351	3.5mm Y-Splitter Cable for PJSM <i>For use with the PJSM Mode/Profile jack. Allows 2 mechanical switches to control 4 separate functions. Requires PC Programming to assign the short hold and long hold commands to each switch. Switches sold separately.</i>	42.00	
<input type="checkbox"/> 613977	R-net Remote Color JS w/ Lights Heavy Duty Gimble <i>Heavy-Duty modification enables gimbal to hold up to extreme operating forces. Reinforced housing is attached to a standard Colour LCD R-net Joystick. 3" long joystick knob has 3/4" diameter. Not compatible with alternative joystick handles.</i>	2,042.00	WEX1
<input type="checkbox"/> 613411	Compact Advanced Joystick for R-net <i>Compact sized joystick that does not contain an LCD Display. Contains two mono jacks for mode & on/off. Does NOT contain a wrist support pad. If you want a wrist support pad, please select 614934 below.</i>	1,990.00	WEX1
<input type="checkbox"/> 614934	Wrist Pad for Compact Joystick Advanced	253.00	

Joystick Options

CHOOSE JOYSTICK HANDLE TYPE:

Part Number	Description	MSRP - CA	ADP Code
<input type="checkbox"/> 603086	Cone Shaped Joystick Knob - Standard	No Charge	
<input type="checkbox"/> 602154	Chin Cup For Joystick Knob	119.00	
<input type="checkbox"/> 602151	T Handle For Joystick Knob	176.00	
<input type="checkbox"/> 602152	Large Ball For Joystick Knob	119.00	
<input type="checkbox"/> 602155	Softball For Joystick Knob	160.00	
<input type="checkbox"/> 602153	Mushroom For Joystick Knob	160.00	
<input type="checkbox"/> 602156	Stick For Joystick Knob	160.00	
<input type="checkbox"/> 613703	1.75" Black Foam Ball For Joystick Knob	119.00	
<input type="checkbox"/> 613970	Bodypoint J/S Handle 3" U-Shaped w/ Flex-Shaft	153.00	
<input type="checkbox"/> 613971	Bodypoint J/S Handle 4" U-Shaped w/ Flex-Shaft	153.00	
<input type="checkbox"/> 613972	Bodypoint J/S Handle Dome Shaped Rubber	153.00	

Part Number	Description	MSRP - CA	ADP Code
<input type="checkbox"/> 317769	Joystick Mount - Fixed <i>Mount on Left: _____ Mount on Right: _____</i>	No Charge	
<input type="checkbox"/> 317404	Retractable Joystick Mount <i>Mount on Left: _____ Mount on Right: _____</i>	452.00	
<input type="checkbox"/> 319940	Mount Standard Joystick on Corpus Back Hinge <i>Standard R-net Joystick is mounted to backrest hinge. Can be chosen instead of PN 323725. Mount on Left: _____ Mount on Right: _____</i>	314.00	
<input type="checkbox"/> 318238	Midline Control Mount <i>Mount Midline on Left Armrest: _____ Mount Midline on Right Armrest: _____</i>	392.00	WE3D
<input type="checkbox"/> 614290	Stealth Gatlin Midline Mount - 8 Position <i>Features eight locking positions each 45° apart. Adjustments are made by pulling out on mounting arm and rotating to desired angle. Joystick bracket sold separately (select from the brackets listed below). Please indicate if you want it mounted on left or right armrest: Mount Midline on Left Armrest: _____ Mount Midline on Right Armrest: _____</i>	690.00	WE3D
<input type="checkbox"/> 614291	Stealth Gatlin Midline Mount - 2 Position <i>Designed for easy swing away. Mount can be set to swing from horizontal to 95° or from horizontal to 180°. Joystick bracket sold separately (select from the brackets listed below). Please indicate if you want it mounted on left or right armrest: Mount Midline on Left Armrest: _____ Mount Midline on Right Armrest: _____</i>	690.00	WE3D
<input type="checkbox"/> 614292	Stealth Gatlin Bracket for R-net/Compact JS Adv	283.00	
<input type="checkbox"/> 614295	Stealth Gatlin Bracket for Compact Joystick	283.00	
<input type="checkbox"/> 614296	Stealth Gatlin Bracket for Mini Joystick	283.00	
<input type="checkbox"/> 614294	Stealth Gatlin Swivel Bracket for Omni	322.00	
<input type="checkbox"/> 613649	PERMOFix Base Joint with 9/16"(14mm) Through-Mount <i>Includes a shim that allows the PERMOFix Joint to be mounted on the 1/2" Stealth Gatlin arm. Provides an adjustable mount for a Mini Joystick or Compact Joystick.</i>	114.00	
<input type="checkbox"/> 614280	Stealth Height Adjustable Mount <i>Replaces standard UT panel bracket. Provides 1.75" of height adjustability to all Permobil joystick mounts and Stealth Gatlin Midline mounts.</i>	252.00	

OPTIONAL ITEMS:

<input type="checkbox"/> 323702	Co-Pilot Attendant Control R-net - Corpus	5,034.00	
<input type="checkbox"/> 323725	Attendant Control R-net	1,307.00	
<input type="checkbox"/> 325446	R-net Bluetooth iDevice Module <i>Allows input device commands (such as joystick nudges or switch clicks) to be used to operate an iPhone, iPad or iPod. By assigning each input to an iOS function, users can make phone calls, send texts and emails and access applications. Kit includes a Bluetooth iDevice Module, an R-net Bus Cable, and two mono jack cables for direct switch inputs (switches not included). Must be ordered with a Color JSM or Omni for mode display.</i>	1,055.00	

Joystick Options

Part Number	Description	MSRP - CA	ADP Code
<input type="checkbox"/> 325447	R-net Bluetooth Mouse Module <i>Enables mouse control through any R-net input device. It can operate the mouse on a computer and on select Android devices (with R-net for Android app). Kit includes a Bluetooth Mouse Module, an R-net Bus Cable, two mono jack cables for L/R click (switches not included) and a USB BT Adapter for user's computer. Note: BT Adapter may not work with all computers. Must be ordered with a Color JSM or Omni for mode display.</i>	1,055.00	
<input type="checkbox"/> 325468	Add a second Bluetooth Mouse Module <i>Allows R-net to simultaneously pair to a second compatible computer or Android device. Each module is accessible through a separate R-net mode. Kit includes a unique Bluetooth Mouse Module, an R-net Bus Cable, two mono jack cables for L/R click and a USB BT Adapter.</i>	1,055.00	
<input type="checkbox"/> 613774	Output Module (IOM) for R-net Control System <i>This allows a user to control an external device (ex: Communication Device, ECU, Door Opener, etc.) through the drive input (JSM, Head Array, etc.). Each module can provide up to 4 switch closures. May require specialty cable or wiring. Up to three output modules can be installed on one wheelchair.</i>	1,749.00	

Seating System

CHOOSE SEAT MOUNTING PLATE TYPE:

Part Number	Description	MSRP - CA	ADP Code
<input type="checkbox"/> 319130	0° Mount Plate - Corpus 3G Seat <i>Range of Tilt will be 0° to +50°.</i>	No Charge	
<input type="checkbox"/> 319471	5° Anterior Tilt Mounting Plate <i>Range of Tilt will be -5° to +45°. This option will raise the minimum seat to floor height by about 1" when the seat is at 0° of tilt.</i>	344.00	

CHOOSE SEATING SYSTEM:

Part Number	Description	MSRP - CA	ADP Code
	PLEASE NOTE: <i>Weight Limit 300lbs. For the Corpus 3G Seat on the M300/M400, there is a half speed limitation when the Total Backrest Angle (TBA) > 150°. TBA = Tilt Angle + Backrest Recline Angle. Example 1: 45° Tilt + 100° Recline = 145° TBA (Full Speed) Example 2: 35° Tilt + 120° Recline = 155° TBA (Half Speed)</i>		
<input checked="" type="checkbox"/> 325116	Seat Frame - Corpus 3G	1,293.00	
<input type="checkbox"/> 318469	Fixed Tilt - Corpus 3G	No Charge	
<input type="checkbox"/> 318470	50° CG Power Tilt	6,531.00	
<input type="checkbox"/> 324241	120° Manual Recline - Corpus	No Charge	
<input type="checkbox"/> 920132	175° Power Recline - M300 C3G	7,431.00	
<input type="checkbox"/> 920035	50° Power Tilt & 175° Power Recline - Corpus 3G <i>Select this option when both Power Tilt and Power Recline are needed.</i>	11,598.00	

Seating System Options

CHOOSE SEAT FUNCTION CONTROL TYPE:

Part Number	Description	MSRP - CA	ADP Code
<input type="checkbox"/> 920019	Single Seat Function Control Switch for R-net <i>Please choose either Push Buttons or Toggle Switches below.</i>	No Charge	
<input type="checkbox"/> 920017	Single Seat Function Control Kit for R-net <i>For operation of a single seat function through the wheelchair electronics via the drive control and/or ICS Switchbox. Order with R-net Controller or remote joystick options. Please indicate below if you want push buttons or toggles on the ICS Switchbox.</i>	1,661.00	
<input type="checkbox"/> 920018	Multiple Seat Function Control Kit for R-net <i>MUST be selected if ordering two or more seat functions. This device is used to operate two or more seat functions through the wheelchair electronics via the drive control and/or ICS switch box. Order with R-net electronics or remote joystick options. Please indicate below if you want push buttons or toggles on the ICS Switchbox.</i>	2,768.00	SEC1P1005

CHOOSE PUSH BUTTONS OR TOGGLES:

Part Number	Description	MSRP - CA	ADP Code
-------------	-------------	-----------	----------

Seating System Options

CHOOSE PUSH BUTTONS OR TOGGLES:

Part Number	Description	MSRP - CA	ADP Code
<input type="checkbox"/> 920135	Push Buttons on ICS Switchbox (1.4m)	No Charge	
<input type="checkbox"/> 920137	Toggle Switches on ICS Switchbox (1.4m)	296.00	

CHOOSE MOUNTING POSITION OF ICS SWITCHBOX:

Part Number	Description	MSRP - CA	ADP Code
<input type="checkbox"/> 316437	Mount ICS Switchbox onto R-net Joystick Bracket <i>ICS Switchbox will be mounted to the same bracket as the R-net Joystick Module. If this item is selected, the ICS Switchbox will swing away with the Joystick Module.</i>	No Charge	
<input type="checkbox"/> 315973	Mount ICS Switchbox onto UT Armrest <i>Mount on Left: _____ Mount on Right: _____</i>	No Charge	
<input type="checkbox"/> 321798	Mount ICS Switchbox on Backrest <i>Mount on Left: _____ Mount on Right: _____</i>	166.00	

OPTIONAL ITEMS:

Part Number	Description	MSRP - CA	ADP Code
<input type="checkbox"/> 920124	Virtual Seating Coach Package w/ Phone Holder <i>The Virtual Seating Coach (VSC) Package includes a Bluetooth dongle, a USB charger, a RAM X-grip Phone Holder and access to VSC coaching and regimen management. The VSC app can be downloaded for free from the App Store or Google Play. Requires an iPhone (iOS 9 or higher) or an Android phone (OS 4. 1. 2 or higher).</i> Mount VSC to Left Armrest: _____ <i>To provide access, USB Charger and Phone Holder may require relocation.</i> Mount VSC to Right Armrest: _____ <i>To provide access, USB Charger and Phone Holder may require relocation.</i>	2,261.00	
<input type="checkbox"/> 317909	ICS Alternative Switchbox for R-net Electronics <i>Seat function interface box that accommodates up to 8 momentary switches. Provides an additional method for accessing the seat functions. 1/8" mini-mono switches not included.</i>	1,307.00	
<input type="checkbox"/> 614512	Remote Stop System - All Electronics <i>Consists of a Radio Frequency (RF) handheld transmitter and an RF receiver mounted on the wheelchair. Custom charges may be required when ordering with a PS/PS Jr seat, a low backrest height or a pediatric wheelchair. If ordering on a K450, push handles must also be selected.</i>	1,247.00	
<input type="checkbox"/> 330365	Memory Seat Program Corpus w/ R-Net (Recall & Save)	121.00	
<input type="checkbox"/> 322020	Legrest Movement Linked to Recline Electronically <i>Power Legrest will raise when Power Recline moves backward and will lower when back moves forward. If Power Legrest is operated, Power Recline will NOT automatically reposition. Feature works through ICS Switchbox and Driver Control.</i>	49.00	
<input type="checkbox"/> 321953	Enable Independent Repositioning Mode <i>IRM allows for sequential movement of power seat functions into a posterior position (35° tilt, 115° recline, 135° legrest) and returns to upright position (5° tilt, 100° recline, 110° legrest). This can be accessed using a standard drive control or a single switch. On the Permobil Joystick Module, IRM is accessible through the mode/profile jack (long switch hold). Please see the Alt Drive Accessories section for switch options.</i>	No Charge	

Backrest Options

CHOOSE FIXED OR SLIDING BACKREST KIT:

Part Number	Description	MSRP - CA	ADP Code
<input type="checkbox"/> 318608	Sliding Backrest Kit - 5" Mechanical Shear	No Charge	
<input type="checkbox"/> 319640	Fixed Backrest Kit - No Mechanical Shear	53.00	

CHOOSE BACKREST TYPE:

Part Number	Description	MSRP - CA	ADP Code
<input type="checkbox"/> 325309	Corpus Ergo Back - 14"W x 20"T <i>BACKREST HEIGHT NOTE: Backrest heights are measured from the seat pan to the top of the backrest shell with the backrest set at 90°. Ergo Backs include lateral wedges and lumbar pads.</i> <i>Power Recline will be limited to 150° when ordering this backrest.</i>	1,072.00	

Backrest Options

CHOOSE BACKREST TYPE:

Part Number	Description	MSRP - CA	ADP Code
<input type="checkbox"/> 325310	Corpus Ergo Back - 16"W x 20" <i>Power Recline will be limited to 150° when ordering this backrest.</i>	1,072.00	
<input type="checkbox"/> 325311	Corpus Ergo Back - 18"W x 20" <i>Power Recline will be limited to 150° when ordering this backrest.</i>	1,072.00	
<input type="checkbox"/> 325312	Corpus Ergo Back - 20"W x 20" <i>Power Recline will be limited to 150° when ordering this backrest.</i>	1,072.00	

PLEASE NOTE:

All taller backrests listed below include: lateral wedges, lumbar pads, a UT rail mounted at the top for lateral adjustment of the headrest, and split backrest feature with height adjustment in 1" increments. Please indicate below where you want the height set at so that a properly sized backrest cushion can be provided.

<input type="checkbox"/> 325313	Corpus Ergo Back - 14"W x 23" - 28"T	1,072.00	
<input type="checkbox"/> 325314	Corpus Ergo Back - 16"W x 23" - 28"T	1,072.00	
<input type="checkbox"/> 325315	Corpus Ergo Back - 18"W x 23" - 28"T	1,072.00	
<input type="checkbox"/> 325316	Corpus Ergo Back - 20"W x 23" - 28"T	1,072.00	
	Set Backrest at 23" _____		
	Set Backrest at 24" _____		
	Set Backrest at 25" _____		
	Set Backrest at 26" _____		
	Set Backrest at 27" _____		
	Set Backrest at 28" _____		

CHOOSE BACKREST CUSHION:

Part Number	Description	MSRP - CA	ADP Code
<input type="checkbox"/>	Leatherette Cushion - Corpus Ergo Back	No Charge	
<input type="checkbox"/>	Black Mesh Cushion - Corpus Ergo Back	53.00	
<input type="checkbox"/>	Omit Backrest Cushion	No Charge	

OPTIONAL ITEMS:

Part Number	Description	MSRP - CA	ADP Code
<input type="checkbox"/> 320523	Tubular Backrest Frame - Corpus (1" Canes) <i>Provides 1" diameter, vertical canes for mounting aftermarket backrests. Frame can be adjusted to support backrest widths from 14" to 20". Canes are attached with 2" long UniTrack linkage disks that can rotate to provide space for backrest hardware. The tops of the canes can be adjusted up to a height of 27" above the seat pan. Depending on the aftermarket backrest, several inches of seat depth may be lost when mounting this system.</i>	477.00	
<input type="checkbox"/> 321190	Direct Backrest Frame - Corpus <i>Used to mount ROHO Agility, Jay 3, Matrix Elite, Varilite Icon and Ride Java backrests ranging in width from 14" to 20". Slotted brackets and 1" UniTrack linkage disks provide a close, adjustable mount. Backrest pitch can be tuned by rotating the linkage disks. Some seat depth may be lost due to the position and orientation of the backrest.</i>	477.00	
<input type="checkbox"/> 920068	Split Backrest Cushion - Corpus <i>Custom Option.</i>	355.00	
<input type="checkbox"/> 920069	Tailored Backrest Cushion - Corpus <i>Sides of cushion are trimmed to provide additional clearance for lateral trunk supports. Must also select a backrest cushion when ordering this custom option.</i>	111.00	

Seat Sizing

CHOOSE SEAT WIDTH:

Part Number	Description	MSRP - CA	ADP Code
<input type="checkbox"/> 319197	Corpus Seat Width Kit - 17"W-19"W <i>Set Seat Width at 17"W _____ Set Seat Width at 19"W _____</i>	No Charge	
<input type="checkbox"/> 320231	Corpus Seat Width Kit - 21"W-23"W <i>Set Seat Width at 21"W _____ Set Seat Width at 23"W _____</i>	No Charge	

Seat Sizing

CHOOSE SEAT DEPTH:

Part Number	Description	MSRP - CA	ADP Code
	PLEASE NOTE: <i>There are three seat depth kits, each providing three seat depth settings: (1) 14" - 16", (2) 17" - 19" and (3) 20" - 22". Please select the required seat depth and the appropriate kit will be provided.</i>		
<input type="checkbox"/> 321435	Seat Depth Kit 14" - 16" <i>Set Seat Depth at: 14"D ____ 15"D ____ 16"D ____</i>	No Charge	
<input type="checkbox"/> 321441	Seat Depth Kit 17" - 19" <i>Set Seat Depth at: 17"D ____ 18"D ____ 19"D ____</i>	No Charge	
<input type="checkbox"/> 321445	Seat Depth Kit 20" - 22" <i>Set Seat Depth at: 20"D ____ 21"D ____ 22"D ____</i>	No Charge	

CHOOSE SEAT CUSHION:

Part Number	Description	MSRP - CA	ADP Code
	NOTE: <i>A Corpus II Ergonomic Seat Cushion will be provided for select seat sizes.</i>		
<input type="checkbox"/>	Corpus Ergo Cushion - Leatherette (C3GERGSCA) <i>Part number will be provided on the quote / order confirmation.</i>	587.00	SES1P0005
<input type="checkbox"/>	Corpus Ergo Cushion - Black Mesh (C3GERGSCA) <i>Part number will be provided on the quote / order confirmation.</i>	587.00	SES1P0005
<input type="checkbox"/>	Corpus Ergo Air Cushion (C3GRHOSCA) <i>Features ROHO Dry Floatation technology. The removable ROHO cell pad provides adjustable skin protection while the contoured foam and thigh support wedges help maintain proper positioning and stability. Includes removable stretch fabric cover. Available in widths of 17", 19" and 21" and depths of 17", 19" and 21", as well as 17" x 23". Part number will be provided on the quote / order confirmation.</i>	623.00	SE0001267
<input type="checkbox"/>	Incontinence Cover for Ergo Cushion <i>Made from a stretch Dartex fabric and suitable for all Ergo cushions including the Corpus Ergo Air. The cover is cut to the contours of the Ergo cushion and is designed to be used as an outer layer of protection. Cover size will be matched to cushion size. Sent as a part. Part Number will be provided on the quote / order confirmation based on the Seat Width and Seat Depth selected.</i>	186.00	SEMCF106L
<input type="checkbox"/>	Omit Seat Cushion	No Charge	

Armrest Options

CHOOSE ARMBAR KIT AND MOUNTING POSITION:

Part Number	Description	MSRP - CA	ADP Code
	PLEASE NOTE: <i>Narrow Armbar Kit can be adjusted from 15"W to 19"W between the armrest pads without any additional parts. Wide Armbar Kit can be adjusted from 19"W to 23"W between the armrest pads without any additional parts. Please select both Armbar Kit and Width Setting below.</i>		
<input type="checkbox"/> 325306	Narrow Armbar Kit Inf Adjustable 15"W - 19"W	No Charge	
<input type="checkbox"/> 325308	Wide Armrest Kit Inf Adjustable 19"W - 23"W <i>Required when selecting a 20" wide backrest.</i>	No Charge	
	Set the Width Between the Armrests to: _____ <i>Must be wider than backrest width to ensure proper fit. Please round to nearest 1".</i>		

CHOOSE ARMREST HEIGHT MOUNTING POSITION:

Part Number	Description	MSRP - CA	ADP Code
<input type="checkbox"/> 920037	Height Adjustable Armrests - NC SET ARMREST HEIGHT AT: _____ <i>Corpus Armrests are infinitely adjustable from 7" - 15", measured from the seat pan to the top of the armrest pad. No additional parts are needed to get this full range of adjustment. In one inch increments (ex: 8", 9", etc.), please indicate where you would like the armrests set.</i>	No Charge	

CHOOSE ARMREST SWIVEL LOCKED OR ENABLED:

Part Number	Description	MSRP - CA	ADP Code
-------------	-------------	-----------	----------

Armrest Options

CHOOSE ARMREST SWIVEL LOCKED OR ENABLED:

Part Number	Description	MSRP - CA	ADP Code
ARMREST MEDIAL SWIVEL NOTE:			
<i>Every Corpus Armrest can flip backwards for transfers. The armrests also include a feature that allows the armrest pads to swivel medially towards the center of the seating system. This swivel feature can be locked out by simply moving a bolt located under the armrest bracket. Please indicate below whether each armrest should have medial swivel enabled or locked. Note: Medial swivel must be locked on Corpus VS or when ordering arm trough adapters.</i>			
<input type="checkbox"/> 319345	Armrest Bracket - Medial Swivel Assembly Locked <i>Left Armrest Locked: _____ Right Armrest Locked: _____</i>	No Charge	
<input type="checkbox"/> 319347	Armrest Bracket - Medial Swivel Assembly Enabled <i>Left Armrest Enabled: _____ Right Armrest Enabled: _____</i>	No Charge	

CHOOSE ARMREST PAD RIGHT:

Part Number	Description	MSRP - CA	ADP Code
<input type="checkbox"/> 324280	4"x10" Corpus Armrest Pad - Leatherette	No Charge	
<input type="checkbox"/> 324282	4"x13" Corpus Armrest Pad - Leatherette	No Charge	
<input type="checkbox"/> 324283	4"x16" Corpus Armrest Pad - Leatherette	No Charge	
<input type="checkbox"/> 324299	4"x18" Corpus HD Armrest Pad - Leatherette	No Charge	
<input type="checkbox"/> 326504	4"x10" Corpus Armrest Pad - Black Fabric	No Charge	
<input type="checkbox"/> 324306	4"x13" Corpus Armrest Pad - Black Fabric	No Charge	
<input type="checkbox"/> 325254	4"x16" Corpus Armrest Pad - Black Fabric	No Charge	
<input type="checkbox"/> 324297	4"x18" Corpus HD Armrest Pad - Black Fabric	No Charge	
<input type="checkbox"/> 326523	4"x10" Corpus Arm w/ Stealth Gel Pad	191.00	
<input type="checkbox"/> 326524	4"x13" Corpus Arm w/ Stealth Gel Pad	197.00	
<input type="checkbox"/> 326525	4"x16" Corpus Arm w/ Stealth Gel Pad	202.00	
<input type="checkbox"/> 324309	4"x18" Corpus HD Arm w/ Stealth Gel Pad	210.00	

CHOOSE ARMREST PAD LEFT:

Part Number	Description	MSRP - CA	ADP Code
<input type="checkbox"/> 324280	4"x10" Corpus Armrest Pad - Leatherette	No Charge	
<input type="checkbox"/> 324282	4"x13" Corpus Armrest Pad - Leatherette	No Charge	
<input type="checkbox"/> 324283	4"x16" Corpus Armrest Pad - Leatherette	No Charge	
<input type="checkbox"/> 324299	4"x18" Corpus HD Armrest Pad - Leatherette	No Charge	
<input type="checkbox"/> 326504	4"x10" Corpus Armrest Pad - Black Fabric	No Charge	
<input type="checkbox"/> 324306	4"x13" Corpus Armrest Pad - Black Fabric	No Charge	
<input type="checkbox"/> 325254	4"x16" Corpus Armrest Pad - Black Fabric	No Charge	
<input type="checkbox"/> 324297	4"x18" Corpus HD Armrest Pad - Black Fabric	No Charge	
<input type="checkbox"/> 326523	4"x10" Corpus Arm w/ Stealth Gel Pad	191.00	
<input type="checkbox"/> 326524	4"x13" Corpus Arm w/ Stealth Gel Pad	197.00	
<input type="checkbox"/> 326525	4"x16" Corpus Arm w/ Stealth Gel Pad	202.00	
<input type="checkbox"/> 324309	4"x18" Corpus HD Arm w/ Stealth Gel Pad	210.00	

OPTIONAL ITEMS:

Part Number	Description	MSRP - CA	ADP Code
<input type="checkbox"/> 320094	Corpus Arm Trough Adapter 1" <i>For aftermarket armrest pads. Please indicate below if they should be installed or on chair sent as parts. Please Note: Medial Swivel must be must be locked when selecting this option. Price shown is for 1(ea). A quantity of 2 will be provided.</i>	174.00	
<input type="checkbox"/> 320095	Corpus Arm Trough Adapter 7/8" <i>For aftermarket armrest pads. Please indicate below if they should be installed or on chair sent as parts. Please Note: Medial Swivel must be must be locked when selecting this option. Price shown is for 1(ea). A quantity of 2 will be provided.</i>	174.00	
Install Arm Trough Adapters: _____			
Send Arm Trough Adapters as Parts: _____			

Armrest Options

OPTIONAL ITEMS:

Part Number	Description	MSRP - CA	ADP Code
<input type="checkbox"/> 325819	Armrest Pouch <i>Armrest Pouch Left: _____ Armrest Pouch Right: _____</i>	252.00	
<input type="checkbox"/> 303666	Corpus Armrest Side Support <i>Custom Option. Leatherette Pad (2"T x 8.5"L). Not compatible with Stealth Elbow Support, except when ordered with an 18" armrest. Left Inner Side Support: _____ Right Inner Side Support: _____ Left Outer Side Support: _____ Right Outer Side Support: _____</i>	259.00	SEMCF601L, SEMCF602L
<input type="checkbox"/> 614346	Stealth 3" x 10" Elbow Support - Right <i>Must also select PN 920030.</i>	168.00	
<input type="checkbox"/> 614345	Stealth 3" x 10" Elbow Support - Left <i>Must also select PN 920030.</i>	168.00	
<input type="checkbox"/> 920030	Adjustable Elbow Support Hardware	283.00	
<input type="checkbox"/> 920030	Adjustable Elbow Support Hardware	283.00	

Legrest Options

CHOOSE LEGREST TYPE:

Part Number	Description	MSRP - CA	ADP Code
<input type="checkbox"/> 319556	Corpus Manual Elevating Legrest <i>Charge for Manual Legrest Elevation: _____ Do not Charge for Manual Legrest Elevation: _____</i>	797.00	
<input type="checkbox"/> 320166	Power Elevating Center Mount Legrests - 170° <i>Includes mechanical articulation.</i>	3,549.00	

ADDITIONAL POWER OPTIONS:

Part Number	Description	MSRP - CA	ADP Code
PLEASE NOTE: <i>If choosing the Power Extending Legrest or the Power Transfer Legrest, footplate height can be adjusted from 12"-16". For the Power Transfer Legrest, there is no pre-set position within the range. Footplate height is determined solely by user adjustment.</i>			
<input type="checkbox"/> 320905	Power Extending Legrest w Wide Footplates - C3G <i>Provides 4" of independent legrest extension. Footplate height can be adjusted from 12"-16" before power extension. Restricts ability to adjust each footplate independently. Compatible with Manual and Power Elevating Legrests. Includes 2 Piece Wide Corpus Footplates 9"D x 7.5"W.</i>	3,856.00	
<input type="checkbox"/> 320912	Power Extending Legrest w Narrow Footplates - C3G <i>Provides 4" of independent legrest extension. Footplate height can be adjusted from 12"-16" before power extension. Restricts ability to adjust each footplate independently. Compatible with Manual and Power Elevating Legrests. Includes 2 Piece Narrow Corpus Footplates 9"D x 6"W.</i>	3,856.00	
<input type="checkbox"/> 320898	Power Transfer Legrest <i>A single input command extends the footplate to the ground and then raises the seat for easier transfer. Footplate can be extended up to 27" from seat pan to footplate. Compatible with Manual and Power Elevating Legrests. Seat elevator is required for full functionality. Please Note: Calf support pads are only available as a custom option. Customization may also be required when ordered with a 5° anterior tilt mounting plate.</i>	4,242.00	

CHOOSE FOOTPLATE STYLE:

Part Number	Description	MSRP - CA	ADP Code
<input type="checkbox"/> 319179	2pc Wide Corpus Footplates - 9"D x 7.5"W	379.00	
<input type="checkbox"/> 319636	2pc Narrow Corpus Footplates - 9"D x 6"W	No Charge	
<input type="checkbox"/> 322004	1pc Narrow Corpus Footpl - 9.75"D x 12.5"W <i>Footplate is 12.5"W in front, 8"W in back.</i>	379.00	
<input type="checkbox"/> 321010	1pc Narrow Power Transfer Footplate - 10"D x 13"W <i>Footplate is 13"W in front, 8"W in back.</i>	No Charge	

Legrest Options

FOOTPLATE HEIGHT POSITION:

Part Number	Description	MSRP - CA	ADP Code
	SET FOOTPLATE HEIGHT AT:_____		

This is measured from the SEAT PAN to the top of the FOOTPLATE, when the legrest assembly is at 90°. DO NOT INCLUDE CUSHION THICKNESS IN YOUR MEASUREMENT. In one inch increments, please indicate desired footplate height. NOTE: Corpus footplate height range is 6"-17" and Corpus VS height range is 9"-17". Height may need to be reduced for shipping. Calf Pads are not available at heights less than 9" (10" on VS Legrest).

OPTIONAL ITEMS:

Part Number	Description	MSRP - CA	ADP Code
<input type="checkbox"/> 323706	UT Calf Support Kit (B) w Mesh Cover & Adj Hardware 3.5"W x 4.5"T Pad.	No Charge	
<input type="checkbox"/> 323711	UT Calf Support Kit (J) w Mesh Cover & Adj Hardware 5"W x 6.5"T Pad.	No Charge	
<input type="checkbox"/> 323704	UT Calf Support Kit (A) w Mesh Cover & Adj Hardware 5.5"W x 6"T Pad.	No Charge	
<input type="checkbox"/> 323709	UT Calf Support Kit (H) w Mesh Cover & Adj Hardware 7"W x 7"T Curved Pad.	392.00	
<input type="checkbox"/> 323707	UT Calf Support Kit (D) w Mesh Cover & Adj Hardware 8"W x 8"T Pad.	392.00	
<input type="checkbox"/> 614446	Stealth Gel Pads - Narrow Footplates (Pair)	399.00	
<input type="checkbox"/> 614388	Stealth Gel Pads - Wide Footplates (Pair)	411.00	

Accessories

<input type="checkbox"/> 330501	Bodypoint Non-Padded Hip Belt - Medium	No Charge	
<input type="checkbox"/> 330503	Bodypoint Padded Hip Belt - Medium	144.00	SERND0030
<input type="checkbox"/> 330502	Bodypoint Non-Padded Hip Belt - Large	116.00	SERND0030
<input type="checkbox"/> 330504	Bodypoint Padded Hip Belt - Large	197.00	SERND0030
<input type="checkbox"/> 305003	Positioning Belt Retractable	298.00	SERND0030
<input type="checkbox"/> 310189	Positioning Belt Retractable Long	323.00	SERND0030

CHOOSE ONE OF THE FOLLOWING:

Part Number	Description	MSRP - CA	ADP Code
<input type="checkbox"/> 325711	Permobil 7"W Headrest for Link Hardware <i>Must also select 329886 or 329887.</i>	439.00	SEH1P0005
<input type="checkbox"/> 325712	Permobil 10"W Headrest for Link Hardware <i>Must also select 329886 or 329887.</i>	439.00	SEH1P0005
<input type="checkbox"/> 325713	Permobil 14"W Headrest for Link Hardware <i>Must also select 329886 or 329887.</i>	439.00	SEH1P0005
<input type="checkbox"/> 329886	Adj/Removable Headrest Link Hardware - Standard	305.00	
<input type="checkbox"/> 329887	Adj/Removable Headrest Link Hardware - Extended <i>Includes an additional 3" link to extend headrest pad range.</i>	329.00	
<input type="checkbox"/> 322147	Permobil 10"W Headrest Leatherette <i>Must also select PN 920047.</i>	621.00	SEH1P0005
<input type="checkbox"/> 324310	Permobil Ergo Shape Headrest Black Fabric <i>Must also select PN 920047.</i>	513.00	SEH1P0005
<input type="checkbox"/> 920047	Adjustable Removable Headrest Hardware	133.00	SEHND0020
<input type="checkbox"/> 614838	Universal Headrest Adapter <i>Please Note: Recline is limited to 160° or less (depending on backrest height) when ordering a Universal Headrest Adapter on a 23" or 24" tall Corpus backrest.</i>	250.00	

CHOOSE ONE OF THE FOLLOWING:

Part Number	Description	MSRP - CA	ADP Code
<input type="checkbox"/> 329547	Permobil 3" x 5" Lateral Supports (Pair) <i>Must also select Lateral Link Hardware. Only compatible with Corpus backrest shells.</i>	459.00	SELND0010

Accessories

CHOOSE ONE OF THE FOLLOWING:

Part Number	Description	MSRP - CA	ADP Code
<input type="checkbox"/> 329548	Permobil 5" x 7" Lateral Supports (Pair) <i>Must also select Lateral Link Hardware. Only compatible with Corpus backrest shells.</i>	459.00	SELND0010
<input type="checkbox"/> 329549	Adjustable Removable Lateral Link Hardware - (Pair)	629.00	SEICF203L

OPTIONAL ITEMS:

Part Number	Description	MSRP - CA	ADP Code
<input type="checkbox"/> 323480	UT Amp. Support Pad (H) w/ Mesh Cover - (1 Each) <i>7"W x 7"T Curved Pad. Must also select PN 920023.</i>	373.00	SEFND0005
<input type="checkbox"/> 323481	UT Amp. Support Pad (D) w/ Mesh Cover - (1 Each) <i>8"W x 8"T Flat Pad. Must also select PN 920023.</i>	373.00	SEFND0005
<input type="checkbox"/> 920023	Adj Removable - Amp Support Hardware <i>Please indicate if you would like the amputation support mounted on the right or left side of the seat frame: Mount Amp Suppt on Left: _____ Mount Amp Suppt on Right: _____</i>	133.00	
<input type="checkbox"/> 324587	UT Thigh/Hip Support Pad (B) w/ Mesh Cover - Tall <i>4.5"W x 3.5"T Pad w/ 4" Tall Mounting Bracket. Must also select PN 920024. Price shown is for 1(ea). A quantity of 2 will be provided.</i>	247.00	
<input type="checkbox"/> 324694	UT Hip Suppt (A) w/ Mesh Cover - Tall <i>6"W x 5.5"T Pad w/ 4" Tall Mounting Bracket. Must also select PN 920024. Price shown is for 1(ea). A quantity of 2 will be provided.</i>	247.00	
<input type="checkbox"/> 324589	UT Thigh/Hip Support Pad (J) w/ Mesh Cover - Tall <i>6.5"W x 5"T Pad w/ 4" Tall Mounting Bracket. Must also select PN 920024. Price shown is for 1(ea). A quantity of 2 will be provided.</i>	247.00	SELND0015
<input type="checkbox"/> 324588	UT Thigh/Hip Support Pad (C) w/ Mesh Cover - Tall <i>8"W x 3.5"T Pad w/ 4" Tall Mounting Bracket. Must also select PN 920024. Price shown is for 1(ea). A quantity of 2 will be provided.</i>	247.00	SELND0015
<input type="checkbox"/> 324570	UT Thigh/Hip Support Pad (G) w/ Mesh Cover - Tall <i>12"W x 3.5"T Pad w/ 4" Tall Mounting Bracket. Must also select PN 920024. Price shown is for 1(ea). A quantity of 2 will be provided.</i>	247.00	SELND0015
<input type="checkbox"/> 920024	Adj Removable - Thigh Support Hardware <i>Price shown is for 1(ea). A quantity of 2 will be provided.</i>	67.00	
<input type="checkbox"/> 920025	Adj Removable - Hip Support Hardware <i>Price shown is for 1(ea). A quantity of 2 will be provided.</i>	67.00	
<input type="checkbox"/> 320381	Abductor Kit Corpus - Leatherette <i>Must also select PN 920022.</i>	369.00	
<input type="checkbox"/> 920022	Adj Removable - Abductor Hardware	133.00	
<input type="checkbox"/> 323712	Transfer Handles Corpus Seat	513.00	
<input type="checkbox"/> 319397	Push Handles Corpus Seat	452.00	
<input type="checkbox"/> 324167	RAM Self-Leveling Cup Holder	88.00	
<input type="checkbox"/> 324166	RAM X-Grip Phone Holder	88.00	
<input type="checkbox"/> 324139	RAM X-Grip 10" Tablet Holder w/Swing Away Mount <i>Compatible with most 10" screen tablets with or without a protective case. Includes adjustable swing away mount. Note: Customization may be required if ordered along with an Omni.</i>	489.00	
<input type="checkbox"/> 325191	USB Charger <i>Provides up to 1.5A (at 5V) to charge most smartphones and tablets. The charger features an isolated circuit design for improved safety; and by pulling power from the ICS, prevents battery drainage when the chair is powered off. Includes mounting hardware and ICS connector. Please select armrest for mounting, keeping in mind that charger may need to be relocated. Mount on left armrest: _____ Mount on right armrest: _____</i>	209.00	
<input type="checkbox"/> 315797	Medical Necessities Bag <i>Includes hardware to attach bag to backrest. Please Note: Recline is limited to 165° or less (depending on backrest height) when ordering with a 23" or 24" tall Corpus backrest.</i>	272.00	
<input type="checkbox"/> 320316	Medical Necessities Bag Hooks <i>Please Note: Recline is limited to 165° or less (depending on backrest height) when ordering with 23" or 24" tall Corpus backrest.</i>	96.00	

Accessories

OPTIONAL ITEMS:

Part Number	Description	MSRP - CA	ADP Code
<input type="checkbox"/> 330507	Bodypoint Monoflex Medium Chest Support Belt Kit <i>18.25"L x 4"W Pad Size. Item will be shipped as a part.</i>	187.00	SERND0001
<input type="checkbox"/> 330508	Bodypoint Monoflex Large Chest Support Belt Kit <i>21.5"L x 4.5"W Pad Size. Item will be shipped as a part.</i>	187.00	SERND0001
<input type="checkbox"/> 330509	Bodypoint Monoflex X-Lrge Chest Support Belt Kit <i>24.5"L x 5"W Pad Size. Item will be shipped as a part.</i>	187.00	SERND0001
<input type="checkbox"/> 320549	Upper Extremity Support Small C3G <i>Must also select PN 920029.</i>	455.00	SETND0015
<input type="checkbox"/> 319559	Upper Extremity Support Angle Adjustable. 14" W <i>Mount on Left: _____ Mount on Right: _____ Must also select PN 920029.</i>	455.00	SETND0010
<input type="checkbox"/> 319558	Upper Extremity Support Angle Adjustable. 16" W <i>Must also select PN 920029.</i>	455.00	SETND0010
<input type="checkbox"/> 319560	Upper Extremity Support Angle Adjustable. 18" W <i>Must also select PN 920029.</i>	455.00	SETND0010
<input type="checkbox"/> 323714	Upper Extremity Support Joystick Cutout 16" W <i>Must also select PN 920029.</i>	455.00	SETND0010
<input type="checkbox"/> 323715	Upper Extremity Support Joystick Cutout 18" W <i>Must also select PN 920029.</i>	455.00	SETND0010
<input type="checkbox"/> 321272	Upper Extremity Support Angle Adjust - UniTrack <i>Must also select PN 920029.</i>	455.00	
<input type="checkbox"/> 920029	Adj Removable - Tray Hardware	133.00	
<input type="checkbox"/> 323701	Oxygen Holder Corpus	331.00	WEPX
<input type="checkbox"/> 323700	Crutch Holder Corpus	344.00	

Alternative Drive Systems

ALTERNATIVE PROPORTIONAL JOYSTICK SYSTEMS:

Part Number	Description	MSRP - CA	ADP Code
<input type="checkbox"/> 613244	Mini-Joystick for Omni / IOM (9-pin connector) <i>Operating Force = 13 g Neutral to Full Deflection = 5mm</i>	3,031.00	
<input type="checkbox"/> 613440	Compact Joystick for Omni / IOM (9-pin connector) <i>Operating Force = 220 g Neutral to Full Deflection = 28mm Joystick includes two 3.5mm mono jacks for Mode and On/Off. Includes Standard Joystick Knob and 1.75" (42mm) Foam Ball Knob.</i>	1,479.00	
<input type="checkbox"/> 613441	Compact Joystick LITE for Omni / IOM (9-pin conn) <i>Operating Force = 50 g Neutral to Full Deflection = 17mm Joystick includes two 3.5mm mono jacks for Mode and On/Off. Includes Standard Joystick Knob and 1.75" (42mm) Foam Ball Knob.</i>	1,692.00	

ALT PROPORTIONAL JOYSTICK MOUNTING HARDWARE:

Part Number	Description	MSRP - CA	ADP Code
	NOTE: <i>Please see the PERMOfix Mounting Solutions, within the Alt Drive Accessories section, for additional mounting options.</i>		
<input type="checkbox"/> 330322	Alt JS Swing Away Mount Kit - Armrest <i>Mount in front of Left Arm: _____ Mount in front of Right Arm: _____</i>	287.00	
<input type="checkbox"/> 324915	Alt JS / Omni Combo Swing Away Mount Kit - Armrest <i>Used to mount Omni and Alternative Proportional Joystick to the same armrest. Includes Swing Away Omni/Tablet Mount, with a longer mounting rod, two PERMOfix modules, one 3" rod and one 4" rod.</i> <i>Mount in front of Left Arm: _____ Mount in front of Right Arm: _____</i>	582.00	
<input type="checkbox"/> 323688	Alternative JS Adjustable Mount Kit - Midline <i>Mount on Left Arm: _____ Mount on Right Arm: _____</i>	873.00	

Alternative Drive Systems

ALT PROPORTIONAL JOYSTICK MOUNTING HARDWARE:

Part Number	Description	MSRP - CA	ADP Code
<input type="checkbox"/> 323674	Alt JS Adjustable Mount Kit - Chin Boom <i>Swing-Away / Removable Chin Boom mounted on right of backrest shell. Custom charge may be applied if ordering lower backrest heights.</i>	931.00	
<input type="checkbox"/> 614931	Alt JS Adjustable Mount Kit - Haims Harness Lrg <i>For neck width up to 19cm.</i>	665.00	
<input type="checkbox"/> 614932	Alt JS Adjustable Mount Kit - Haims Harness Sml <i>For neck width up to 14cm.</i>	665.00	

ALT PROPORTIONAL JOYSTICK ACCESSORIES:

Part Number	Description	MSRP - CA	ADP Code
<input type="checkbox"/> 613741	Dual Switches for Compact Joystick - Fixed Mount <i>Includes two pushbutton switches mounted on rigid Switches can be used for Mode & On/Off. Switch operation force = 130 g Mount at front of Joystick: _____ Mount at rear of Joystick: _____ Mount at left of Joystick: _____ Mount at right of Joystick: _____</i>	276.00	
<input type="checkbox"/> 613597	Dual Switches for Compact Joystick -Gooseneck Mount <i>Includes two pushbutton switches mounted on adjustable gooseneck arms. Switches are positioned on Left and Right sides of Compact Joystick. Switches can be used for Mode & On/Off. Switch operation force = 130 g</i>	502.00	

SIP & PUFF SYSTEMS:

Part Number	Description	MSRP - CA	ADP Code
<input type="checkbox"/> 614928	Sip & Puff System for Omni <i>Requires Omni Control Module. Cannot use this part with an Input Module. Includes breath tube kit.</i>	1,778.00	WEX5
<input type="checkbox"/> 329842	Sip & Puff System with Fiber Optic Opti-Stop <i>Requires Omni Control Module. Cannot use this part with an Input Module. Includes breath tube kit.</i>	2,366.00	WEX5
<input type="checkbox"/> 614559	Sip & Puff Kit for Total Control Head Array <i>For use with Total Control Head Array. Requires Omni Control Module. Cannot be used with an Input Module. Includes breath tube kit.</i>	982.00	
<input type="checkbox"/> 920073	Adjustable Mounting Hrdw - Sip & Puff <i>Required when selecting PN 614928, 329842 or 614559.</i>	287.00	

ALTERNATIVE PROPORTIONAL CONTROL SYSTEMS:

Part Number	Description	MSRP - CA	ADP Code
<input type="checkbox"/> 614487	Proportional Foot Control <i>Mount on Left: _____ Mount on Right: _____</i>	3,121.00	
<input type="checkbox"/> 920074	Adjustable Mounting Hrdw - Foot Control <i>Required when selecting PN 614487.</i>	287.00	

Total Control Head Array

CONFIGURE HEAD ARRAY:

Part Number	Description	MSRP - CA	ADP Code
<input type="checkbox"/> 330548	Permobil Head Array - 3-Switch System <i>Features an occipital pad sensor and two stainless steel arms configured with proximity sensors. The arms are mounted to the upper attachment points on the head array. Longer arms are also included to extend the sensor adjustment range. ESP recommended. Choose Omni Control Module (613581) separately.</i>	5,485.00	WEX4
<input type="checkbox"/> 319505	Program R-net Omni for 5-Switch Input Device <i>Select additional arms/sensors below and indicate the intended use of each sensor in the Notes & Comments section. Choose Omni Control Module (613581) or Input Module (920021) separately.</i>	No Charge	
CHOOSE OCCIPITAL PAD SIZE (REQUIRED):			
<input type="checkbox"/> 614108	Small Occipital Pad for Head Array - 4.25"W x 3"H	No Charge	

Total Control Head Array

CONFIGURE HEAD ARRAY:

Part Number	Description	MSRP - CA	ADP Code
<input type="checkbox"/> 614109	Large Occipital Pad for Head Array - 6"W x 4"H	No Charge	
CHOOSE OCCIPITAL MOUNT (REQUIRED):			
<input type="checkbox"/>	Fixed Mount for Head Array Occipital Pad	No Charge	
<input type="checkbox"/> 614097	Adj. Swivel Mount for Head Array Occipital Pad <i>Enables Occipital Pad to be angled up/down and left/right independent of the lateral arms.</i>	491.00	
CHOOSE BACKREST MOUNTING HARDWARE (REQUIRED):			
<input type="checkbox"/> 614236	Stealth TWB480 Backrest Mounting Hdwr - 10" stem	305.00	
<input type="checkbox"/> 614259	Stealth TWB480 Backrest Mounting Hdwr - 6" stem	305.00	
<input type="checkbox"/> 614234	Stealth TWB480 Backrest Mnt Hdwr with Flipback	504.00	
<input type="checkbox"/> 614235	Stealth HMO475 MultiAxis Backrest Mounting Hdwr	258.00	
<input type="checkbox"/> 614260	Stealth TWB480 Backrest Mnt Hdwr w SubOcptl Mount	548.00	
<input type="checkbox"/> 614258	Stealth 600 SubOccipital Pad - 8.5" W x 2.75" H <i>Must select when PN 614260 is ordered.</i>	111.00	
SELECT ADDITIONAL LEFT ARM PACKAGE (OPTIONAL):			
<input type="checkbox"/>	Add Left Lateral Arm w Proximity Sensor <i>Choose Lateral Arm Adjustment Range: Short (5.5" - 10"): _____ Long (5.5" - 14"): _____</i>	404.00	
<input type="checkbox"/>	Add Left Lateral Arm w Egg Switch <i>Choose Lateral Arm Adjustment Range: Short (5.5" - 10"): _____ Long (5.5" - 14"): _____</i>	269.00	
<input type="checkbox"/>	Add Left Lateral Arm w Pad (No Sensor or Switch) <i>Choose Lateral Arm Adjustment Range: Short (5.5" - 10"): _____ Long (5.5" - 14"): _____</i>	164.00	
SELECT ADDITIONAL RIGHT ARM PACKAGE (OPTIONAL):			
<input type="checkbox"/>	Add Right Lateral Arm w Proximity Sensor <i>Choose Lateral Arm Adjustment Range: Short (5.5" - 10"): _____ Long (5.5" - 14"): _____</i>	404.00	
<input type="checkbox"/>	Add Right Lateral Arm w Egg Switch <i>Choose Lateral Arm Adjustment Range: Short (5.5" - 10"): _____ Long (5.5" - 14"): _____</i>	269.00	
<input type="checkbox"/>	Add Right Lateral Arm w Pad (No Sensor or Switch) <i>Choose Lateral Arm Adjustment Range: Short (5.5" - 10"): _____ Long (5.5" - 14"): _____</i>	164.00	

Alt Drive Systems Setup

CHOOSE OMNI OR INPUT MODULE (REQUIRED):

Part Number	Description	MSRP - CA	ADP Code
<input type="checkbox"/> 613581	Omni Control Module for R-net w/ Infrared <i>Must also select mounting hardware, either 315315 or 324915. Please specify which type of alternative drive input you will be using if not ordering it from Permobil. Device: _____</i>	2,478.00	
<input type="checkbox"/> 315315	Adjustable Swing-Away Hardware - Omni/Tablet <i>Mount in front of Left Arm: _____ Mount in front of Right Arm: _____</i>	282.00	
<input type="checkbox"/> 920021	Input Module (IOM) for R-net Control System <i>This device allows a proportional or switched drive input to be connected to an R-net Control System without the use of an Omni. It must be ordered in conjunction with an R-net Remote Joystick. It does NOT support any of the following input devices: Sip-n-Puff, 3-Switch Head Array or RIM control.</i>	1,749.00	

CHOOSE MOUNTING POSITION OF SEAT CONTROL SWITCHBOX:

Part Number	Description	MSRP - CA	ADP Code
<input type="checkbox"/>	Mount Seat Control Switchbox on Backrest	No Charge	

Alt Drive Accessories

PERMOFIX MOUNTING SOLUTIONS:

Part Number	Description	MSRP - CA	ADP Code
<input type="checkbox"/> 323835	PERMOfix UniTrack Fixed Mount Kit <i>Includes two base joints, two 4" rods and one 8" rod.</i>	309.00	
<input type="checkbox"/> 323836	PERMOfix UniTrack Removable Mount Kit <i>Includes quick release mount and thigh support hardware, two base joints and three different length rods.</i>	543.00	
<input type="checkbox"/> 330092	PERMOfix UniTrack Swing Away Mount Kit <i>Includes PERMOfix Swing Away module, two base joints, one 2" rod and one 4" rod.</i>	783.00	
<input type="checkbox"/> 614103	Proximity Sensor + Swivel Assembly <i>Only available with the Total Control Head Array. Swivel assembly allows sensor to be mounted with PERMOfix or directly to lateral arm of Total Control Head Array. Since Proximity Sensors require power, they CANNOT be used as On/Off switches.</i>	404.00	
<input type="checkbox"/> 614107	Egg Switch + Swivel Mount Assembly <i>Swivel assembly allows switch to be mounted with PERMOfix or directly to lateral arm of Total Control Head Array.</i>	269.00	
<input type="checkbox"/> 614105	Solid Oval Pad + Swivel Assembly <i>Swivel assembly allows pad to be mounted with PERMOfix or directly to lateral arm of Total Control Head Array.</i>	164.00	
<input type="checkbox"/> 614113	Oval Plate + Swivel Assembly - No Switch <i>Blank plate provides mounting option for Micro Lite, Buddy Button or other switch options. Swivel assembly allows pad to be mounted with PERMOfix or directly to lateral arm of Total Control Head Array.</i>	142.00	

OPTIONAL ITEMS:

Part Number	Description	MSRP - CA	ADP Code
<input type="checkbox"/> 613306	Micro Lite Switch	178.00	
<input type="checkbox"/> 613303	Buddy Button - Red	134.00	
<input type="checkbox"/> 614356	Egg Switch Black	148.00	
<input type="checkbox"/> 613308	Mini Cup Switch	131.00	
<input type="checkbox"/> 614842	Proximity Switch	621.00	
<input type="checkbox"/> 614837	Mini Joystick 5 Switch OC for Omni/IOM (9-pin)	1,055.00	
<input type="checkbox"/> 613921	Touch Contact with CR2032 Battery	299.00	
<input type="checkbox"/> 330019	Dome Switch Backrest Mount w Gooseneck <i>Mount on Left: _____ Mount on Right: _____ Mounts at top of backrest shell. Custom charge may be applied if ordering with short backrest heights. Kit Includes: One PERMOfix joint with Flat Plate, two PERMOfix joints with 45° swivel adjustment, one 150mm rod, one 200mm rod and one 225mm gooseneck switch with 1.2m cable and mounting clamp.</i>	553.00	
<input type="checkbox"/> 330021	Dome Switch Backrest Swing Away Mount w Gooseneck <i>Mount on Left: _____ Mount on Right: _____ w 90° Swing-Away Function. Mounts at top of back. Custom charge applies if ordered with short back. Kit Includes: One PERMOfix joint with Flat Plate, one PERMOfix joint with 45° swivel adjustment, two 6" rods, one 90° Swing-Away joint, one 9" gooseneck switch w 50" cable and mounting clamp.</i>	830.00	
<input type="checkbox"/> 311814	3.5mm Ext. Cable for Sensor/Switch 0.3m Long	38.00	
<input type="checkbox"/> 311810	3.5mm Ext. Cable for Sensor/Switch 1m Long	42.00	
<input type="checkbox"/> 311811	3.5mm Ext. Cable for Sensor/Switch 2m Long	46.00	
<input type="checkbox"/> 311812	3.5mm Ext. Cable for Sensor/Switch 3m Long	50.00	
<input type="checkbox"/> 311813	3.5mm Y-splitter Cable for Sensor/Switch <i>Allows two Proximity Sensors/Mechanical Switches to connect to the same direction jack.</i>	42.00	

Notes & Additional Instructions:

Chair Order Policy:

You are highly encouraged to contact Permobil Sales with ANY questions regarding proper completion of our order forms at 1-800-736-0925. The order cannot be processed without the following:

- A client name or code, the client's measurements, the client's weight
- The order form must also be signed by a person authorized by your company to acknowledge the items selected. Permobil, LTD. is not responsible for configuration or size discrepancies resulting from customer errors on the order form.
- Order confirmations will be provided by Permobil summarizing the items selected. Please review this carefully.
- Once the order is shipped, any changes are subject to the returns policy and may be prohibited.

Order Acknowledgement:

I, _____, am an agent of the medical equipment provider named on this order form and I have the authority to contract for the purchase of powered wheelchairs and related parts on behalf of said provider. I acknowledge that I have reviewed this order and that it is complete and accurate to the best of my knowledge. I further acknowledge that any changes to this order after submission of this order form are subject to the returns policy and may carry additional charges.

Signed: _____

Printed Name: _____

Title: _____

Date: _____

NOTE: All specifications and prices are subject to change without notice. Please note that prices displayed are only valid if ordered with the wheelchair. Any parts ordered for service or future alterations will carry different pricing and warranty. *The ADP codes provided are not intended to be billing or legal advice, rather our interpretation of the code definitions. Use of the codes does not ensure coverage or payment for the item. For coverage information, verify the policy of the appropriate payer.*