

**permobil**

TILITE MANUAL WHEELCHAIRS

# AERO X

## MULTIFACETED SIMPLICITY.

The Aero X features an innovative modular design with clean modern lines and configurability. Featuring multiple configurations that include heavy-duty, hemi, and amputee options. Add superior rolling dynamics and full adjustability, and you get an Aero X that excels in both performance and price.



# AERO X



## TECHNICAL SPECIFICATIONS

Frame Features	Aluminum Box Frame Folding Modular Design Full Adjustability	Camber	0°, 2°, 4°, or 6°
Total Width Folded*	10.5"	Center of Gravity Adjustment	6"
Total Chair Width*	22.5"	Transit Options	Transit Tie-Down Transport Bracket
Seat Width	12" – 22"	Transport Weight*	Fixed Front - 15.6 lb. Swing Away - 18 lb.
Seat Depth	12" – 20"	Weight Limit	265 lb. Standard 350 lb. Heavy Duty
Front Seat Height	13.5" – 21.5"	HCPCS Code	K0005
Rear Seat Height	13" – 20"	Warranty	Five Years on Frame One Year on TiLite Components
Front Frame Angle	Fixed Front - 70°, 80°, or 85° Swing Away - 70°, 80°, 85°, or 90°	* Weights and widths based on 16" x 16" chair with 0° camber, minimum rear wheel spacing and all standard features. Rear wheels are not included in transport weight. For reference, a pair of standard 24" rear wheels weighs 7.2 lb.	
Seat Back Type	<i>Folding Seat Back</i> Adj. Height 11" – 20.5" Adj. Angle 80° – 101°  <i>Fixed Seat Back</i> Adj. Height 11" – 20"		

## MODULAR DESIGN

Optimize your ride with custom configuration. This chair's modular design allows you to add options including after market seating systems, switching fixed and swing away front ends, and growing chair width and / or depth to meet changing customer needs.

## CUSTOM CONFIGURATION.


Custom Configuration allows you to add options and components to your chair to meet your needs and to optimize your riding experience. This includes choices of wheels, handrims, tires, wheel locks, color, tattoos, and more. TiLite offers more options than any other manufacturer!

## COLOR OPTIONS


### HIGH GLOSS PAINT




### MATTE PAINT



### METALLIC PAINT



### CANDY PAINT



### PEARLESCENT PAINT

