

M300 PS Jr



** Indicates a Required Field*

Dealer Information

* Contact: _____

Dealer Code: _____

* Dealer Name: _____

Address: _____

* City: _____

* State: _____

* Phone#: _____

* Fax#: _____

PO#: _____

Email Address: _____

Client Information

Permobil recommends that the client is evaluated by a certified rehab specialist.

* First Name: _____

* Last Name: _____

Diagnosis: _____

Funding Source: _____

Client Age: _____

Client Measurements

* Weight: _____

* Height: _____

A) Top of Shoulders: _____

B) Chest Depth: _____

C) Chest Width: _____

D) Seat Depth: _____


E) Top of Head: _____

F) Elbow to Hand: _____

G) Seat Pan to Elbow: _____

H) Hip Width: _____

*I) Knee to Foot: _____





Chair Model and Colors

CHOOSE BASE:

Part Number	Description	Price AUD
<input checked="" type="checkbox"/> I10857	M300 Base PS Jr <i>Price Includes: M300 Power Base, Shock Absorbing Suspension System, R-net Controller, Tie-Down Hardware For Strap Systems, Charger, Flat Free Tires, and Standard Positioning Belt. Must Choose Flat Free or Pneumatic Tires Below. 10kmh, Weight Capacity = 75kg, STFH = 17".</i>	

CHOOSE WHEEL PACKAGE:

Part Number	Description	Price AUD
<input type="checkbox"/> I11054	Grey Tire Package - MWD <i>Includes grey tires on silver casters and drive wheel rims. Hub cap color is selected below. Please select caster size and drive wheel type: 6" Casters: _____ 8" Casters: _____ Flat free drive tires: _____ Pneumatic drive tires: _____</i>	No Charge
<input type="checkbox"/> I11055	Black Tire Package - MWD <i>Includes black tires on silver casters and drive wheel rims. Hub cap color is selected below. Please select caster size and drive wheel type: 6" Casters: _____ 8" Casters: _____ Flat free drive tires: _____ Pneumatic drive tires: _____</i>	No Charge
<input type="checkbox"/> I11056	Blackout Wheel Package - MWD <i>Custom Option. Includes black tires and custom matte black casters, drive wheels and hub caps. Please select caster size and drive wheel type: 6" Casters: _____ 8" Casters: _____ Flat free drive tires: _____ Pneumatic drive tires: _____</i>	

OPTIONAL DRIVE TIRES:

Part Number	Description	Price AUD
<input type="checkbox"/> I10782	Tri-Spoke Split Rim w/ Snow Tire Pneumatic 3x8" <i>Snow tires are designed for use in cold weather conditions. The nylon compound in these tires remains flexible at low temperatures, but wears more easily. Therefore, the snow tires (sold as wheel assemblies) are not recommended for year-long use and must be ordered along with standard drive wheels. Not available with black tires or black rims.</i>	

CHOOSE SHROUD COLOR:

Part Number	Description	Price AUD
<input type="checkbox"/> I105239-99-0	M300 Shroud - Onyx Black	No Charge
<input type="checkbox"/> I105243-99-0	M300 Shroud - Platinum Silver	No Charge
<input type="checkbox"/> I105242-99-0	M300 Shroud - Sublime Green	No Charge
<input type="checkbox"/> I106172-99-0	M300 Shroud - Cobalt Blue	No Charge
<input type="checkbox"/> I106174-99-0	M300 Shroud - Phantom Grey	No Charge
<input type="checkbox"/> I106161-99-0	M300 Shroud - Pop Star Pink	No Charge
<input type="checkbox"/> I106173-99-0	M300 Shroud - Purple Rain	No Charge
<input type="checkbox"/> I106162-99-0	M300 Shroud - Volcano Red	No Charge

CHOOSE HUBCAP COLOR:

Part Number	Description	Price AUD
<input type="checkbox"/> I10871	Match Hubcap Color to Shroud: _____ <i>Note: Hub caps will be matched to shroud color if no selection is made.</i> Solid Color Hub Cap - Pair <i>Platinum Silver: _____ Onyx Black: _____ Cobalt Blue: _____ Phantom Gray: _____ Pop Star Pink: _____ Volcano Red: _____ Purple Rain: _____ Titan Blue: _____ Sublime Green: _____ Crystal Blue: _____ Vintage Wine: _____ Extreme Yellow: _____ Chili Pepper Red: _____</i>	No Charge

CHOOSE ONE OF THE FOLLOWING:

Part Number	Description	Price AUD
<input checked="" type="checkbox"/> I10901	Australian Freight & Handling Fee	



Base Options

MUST CHOOSE EITHER SEAT ELEVATOR OR FIXED SEAT TUBE:

Part Number	Description	Price AUD
<input type="checkbox"/> I101994-99-0	Seat Tube Fixed	No Charge
<input type="checkbox"/> I106958-99-0	Power Adjustable Seat Height - 8" Travel	

MUST CHOOSE BATTERIES:

Part Number	Description	Price AUD
<input type="checkbox"/> I107031-99-0	Batteries, Grp 24 (72Ah) Sealed Gel, x2, Installed	
<input type="checkbox"/> I107032-99-0	Batteries, Grp 34 (60Ah) Sealed Gel, x2, Installed	

OPTIONAL ITEMS:

Part Number	Description	Price AUD
<input type="checkbox"/> I105428-99-0	Enhanced Steering Perform Unit (ESP), R-net <i>Recommended item with switched drive inputs for improved tracking.</i>	
<input type="checkbox"/> I105297-99-0	LED Light Kit for M300 / M400 w/ R-net Electronics	
<input type="checkbox"/> ID51113	Programming Key for R-net - Permobil Version	
<input type="checkbox"/> I106223-99-0	Stealth BPAP Vent Tray - MWD Base <i>Designed for BPAP systems, the tray has a 12"W x 7"D surface that can support up to 40 lbs. A 1" lip and adjustable strap help secure contents. The tray mounts to the rear tie-down loops, is tapered to fit between the rear swing arms and extends 1" beyond the trailing casters.</i>	

Joystick Options

PLEASE NOTE:

For additional input devices, refer to Alternative Drive Systems section.

<input type="checkbox"/> I105246-99-0	R-net LED Remote Joystick Module <i>This Joystick does NOT have a Display Screen and does NOT contain mono jacks. Limited to five profiles with this Joystick Module.</i>	No Charge
<input type="checkbox"/> I107491-99-0	R-net Remote Color JS w/ Toggle Switches, No Lights	
<input type="checkbox"/> I107492-99-0	R-net Remote Color JS w/ Mono Jacks, No Lights	
<input type="checkbox"/> I104923-99-0	R-net Remote Color JS w/ Mono Jacks, Lights	
<input type="checkbox"/> I107493-99-0	R-net Remote Color JS w/ Toggle Switches, Lights	
<input type="checkbox"/> I105700-99-0	R-net Remote Color JS w/ Lights, Heavy Duty Gimbal <i>Heavy-Duty modification enables gimbal to hold up to extreme operating forces. Reinforced housing is attached to a standard Color LCD R-net Joystick. 3" long joystick knob has 3/4" diameter. Not compatible with alternative joystick handles.</i>	
<input type="checkbox"/> I103814-99-0	Compact Advanced Joystick for R-net <i>Compact sized joystick that does not contain an LCD Display. Contains two mono jacks for mode & on/off. Does NOT contain a wrist support pad. If you want a wrist support pad, please select I10621 below.</i>	
<input type="checkbox"/> I10621	Add Wrist Pad to Compact Joystick Advanced	

CHOOSE JOYSTICK HANDLE TYPE:

Part Number	Description	Price AUD
<input type="checkbox"/> I10045	Cone Shaped Joystick Knob - Standard	No Charge
<input type="checkbox"/> I100960-99-0	Chin Cup For Joystick Knob	
<input type="checkbox"/> I100961-99-0	"T" Handle For Joystick Knob	
<input type="checkbox"/> I100962-99-0	"Large Ball" For Joystick Knob	
<input type="checkbox"/> I100963-99-0	"Softball" For Joystick Knob	
<input type="checkbox"/> I100964-99-0	"Mushroom" For Joystick Knob	
<input type="checkbox"/> I100965-99-0	"Stick" For Joystick Knob	
<input type="checkbox"/> I105442-99-0	1.75" Black Foam Ball For Joystick Knob	
<input type="checkbox"/> I105846-99-0	Bodypoint J/S Handle 3" U-Shaped w/ Flex-Shaft	
<input type="checkbox"/> I105847-99-0	Bodypoint J/S Handle 4" U-Shaped w/ Flex-Shaft	
<input type="checkbox"/> I105848-99-0	Bodypoint J/S Handle Dome Shaped, Rubber	



Joystick Options

CHOOSE JOYSTICK MOUNT:

Part Number	Description	Price AUD
<input type="checkbox"/> I10412	Joystick Mount Right - Fixed, R-net Only	No Charge
<input type="checkbox"/> I10413	Joystick Mount Left - Fixed, R-net Only	No Charge
<input type="checkbox"/> I102537-99-0	Joystick Mount Right - Swing Away	
<input type="checkbox"/> I102538-99-0	Joystick Mount Left - Swing Away	
<input type="checkbox"/> I108037-99-0	Slimline Retractable Joystick Mount - Right	
<input type="checkbox"/> I108038-99-0	Slimline Retractable Joystick Mount - Left	
<input type="checkbox"/> I11232	Adjustable Height Panel Bracket <i>Provides 1.77" of vertical adjustment below the armrest and 25° of angle adjustment. Compatible with all Permobil joystick mounts except for the M1 Fixed Joystick Mount.</i>	
<input type="checkbox"/> I103472-99-0	Mount Standard Joystick on Backrest Shell <i>Standard R-net Joystick is mounted to backrest. Can be chosen instead of PN ID50882. Mount on Left: _____ Mount on Right: _____</i>	
<input type="checkbox"/> ID50882	Attendant Control for R-net Control System Only	
<input type="checkbox"/> I107350-99-0	R-net Bluetooth iDevice Module <i>Allows input device commands (such as joystick nudges or switch clicks) to be used to operate an iPhone, iPad or iPod. By assigning each input to an iOS function, users can make phone calls, send texts and emails and access applications. Kit includes a Bluetooth iDevice Module, an R-net Bus Cable, and two mono jack cables for direct switch inputs (switches not included). Must be ordered with a Color JSM or Omni for mode display.</i>	
<input type="checkbox"/> I107367-99-0	R-net Bluetooth Mouse Module <i>Enables mouse control through any R-net input device. It can operate the mouse on a computer and on select Android devices (with R-net for Android app). Kit includes a Bluetooth Mouse Module, an R-net Bus Cable, two mono jack cables for L/R click (switches not included) and a USB BT Adapter for user's computer. Note: BT Adapter may not work with all computers. Must be ordered with a Color JSM or Omni for mode display.</i>	
<input type="checkbox"/> I107368-99-0	Add a second Bluetooth Mouse Module <i>Allows R-net to simultaneously pair to a second compatible computer or Android device. Each module is accessible through a separate R-net mode. Kit includes a unique Bluetooth Mouse Module, an R-net Bus Cable, two mono jack cables for L/R click and a USB BT Adapter.</i>	
<input type="checkbox"/> I10582	Output Module (IOM) for R-net Control System <i>This allows a user to control an external device (ex: Communication Device, ECU, Door Opener, etc.) through the drive input (JSM, Head Array, etc.). Each module can provide up to 4 switch closures. May require specialty cable or wiring. Up to three output modules can be installed on one wheelchair.</i>	

Seating System

CHOOSE ONE OF THE FOLLOWING:

Part Number	Description	Price AUD
<input type="checkbox"/> I103712-99-0	PS Junior Seat Frame - Only <i>Weight limit 75kg.</i>	No Charge
<input type="checkbox"/> I103715-99-0	PS Jr Seat w/ 45° Manual Tilt <i>Must be ordered with Stroller Handles (I103811-99-0). Weight limit 75kg.</i>	
<input type="checkbox"/> I103716-99-0	PS Junior Seat Frame with 45° Power Tilt <i>Weight limit 75kg.</i>	

OPTIONAL ITEMS:

Part Number	Description	Price AUD
<input type="checkbox"/> I103811-99-0	Stroller Handle for PS Jr Seating System <i>Angle adjustable. This item must be selected when ordering manual tilt (I103715-99-0).</i>	

Seating System Options

- I10585 Seat Function Control Kit for R-net



Seating System Options

CHOOSE SEAT FUNCTION CONTROL TYPE:

Part Number	Description	Price AUD
<input checked="" type="checkbox"/> I10585	Seat Function Control Kit for R-net <i>Please indicate below if you want push buttons or toggles on the ICS Switchbox.</i>	

CHOOSE PUSH BUTTONS OR TOGGLES:

Part Number	Description	Price AUD
<input type="checkbox"/> I102798-99-0	Push Buttons on ICS Switchbox	No Charge
<input type="checkbox"/> I102799-99-0	Toggle Switches on ICS Switchbox	
<input type="checkbox"/> I105409-99-0	Do Not Install ICS Switchbox, Send as Part	No Charge

CHOOSE MOUNTING POSITION OF ICS SWITCHBOX:

Part Number	Description	Price AUD
<input type="checkbox"/> I10626	Mount ICS Switchbox onto R-net Joystick Bracket <i>ICS Switchbox will be mounted to the same bracket as the R-net Joystick Module. If this item is selected, the ICS Switchbox will swing away with the Joystick Module.</i>	No Charge
<input type="checkbox"/> I104651-99-0	Mount ICS Switchbox onto Armrest <i>Mount on Left: _____ Mount on Right: _____</i>	No Charge
<input type="checkbox"/> I106256-99-0	Mount ICS Switchbox on Backrest <i>Mount on Left: _____ Mount on Right: _____</i>	

OPTIONAL ITEMS:

Part Number	Description	Price AUD
<input type="checkbox"/> I104196-99-0	ICS Alternative Switchbox for R-net Electronics <i>Seat function interface box that accommodates up to 8 momentary switches. Provides an additional method for accessing the seat functions. 1/8" mini-mono switches not included.</i>	
<input type="checkbox"/> I106743-99-0	Remote Stop System - All Electronics <i>Consists of a Radio Frequency (RF) handheld transmitter and an RF receiver mounted on the wheelchair. Custom charges may be required when ordering with a PS/PS Jr seat, a low backrest height or a pediatric wheelchair. If ordering on a K450, push handles must also be selected.</i>	

Backrest Options

CHOOSE ONE OF THE FOLLOWING:

Part Number	Description	Price AUD
<input checked="" type="checkbox"/> I103717-99-0	PS Jr Back Canes, Height & Angle Adjustable <i>Must also select PS Jr Style Armrests if ordering this backrest. Choose PS Jr Style Backrest Type (Metal Back or Omit Back) below. Backrest width should match seat width. Back height is adjustable from 16.5"T - 22.5"T, measured from the seat pan to the top of the back.</i>	No Charge

CHOOSE ONE OF THE FOLLOWING:

Part Number	Description	Price AUD
<input type="checkbox"/> I10448	PS Jr Curved Metal Backrest <i>Backrest width will match seat width. Back height is adjustable from 16.5"T - 22.5"T, measured from the seat pan to the top of the backrest. Please note: As back height is increased, a proportional gap is created between the seat pan and the bottom of the backrest.</i>	No Charge
<input type="checkbox"/> I10449	PS Jr - Omit Backrest	No Charge
<input type="checkbox"/> I104616-99-0	Backrest Adjustment Quick Lock, PS / PS Jr	

Seat Sizing

CHOOSE ONE OF THE FOLLOWING:

Part Number	Description	Price AUD
<input type="checkbox"/> I10432	PS Jr Seat Width 11"	No Charge
<input type="checkbox"/> I10433	PS Jr Seat Width 12"	No Charge
<input type="checkbox"/> I10434	PS Jr Seat Width 13"	No Charge



Seat Sizing

CHOOSE ONE OF THE FOLLOWING:

Part Number	Description	Price AUD
<input type="checkbox"/> I10435	PS Jr Seat Width 14"	No Charge
<input type="checkbox"/> I10436	PS Jr Seat Width 15"	No Charge
<input type="checkbox"/> I10437	PS Jr Seat Width 16"	No Charge

CHOOSE ONE OF THE FOLLOWING:

Part Number	Description	Price AUD
<input type="checkbox"/> I10985	SET SEAT DEPTH AT: _____ <i>Seat depth is adjustable from 5" to 18". Please round desired seat depth to nearest 0.5".</i>	No Charge

CHOOSE ONE OF THE FOLLOWING:

Part Number	Description	Price AUD
<input type="checkbox"/> I10276	Miniflex Seat Cushion Leatherette <i>Approximately 2" thick. Please select cushion size: 10"W x 10"D: _____ 12"W x 12"D: _____ 14"W x 14"D: _____ Custom sizes are available upon request.</i>	
<input type="checkbox"/> I10999	Ergonomic Seat Cushion Leatherette <i>Please select cushion size: 16.5"W x 15"D: _____ 17.5"W x 16.5"D: _____ 17.5"W x 19.5"D: _____ Custom sizes available upon request.</i>	

Armrest Options

CHOOSE ONE OF THE FOLLOWING:

Part Number	Description	Price AUD
<input type="checkbox"/> I10824	Flip Back, Height Adjustable Armrests <i>Mounted off the backrest canes. Includes: 45° of angle adjustability, rotational adjustment, 3" offset adjustment both medially and laterally, and and height adjustment of 6" - 14" from the seat pan. In one inch increments, please indicate where you would like the armrests set. SET ARMRESTS AT (6" - 14"): _____</i>	
<input type="checkbox"/> I10825	Flip Back, Height Adjustable Armrests - NC <i>No charge when ordering power tilt. In one inch increments, please indicate where you would like the armrests set. ARMRESTS AT (6" - 14"): _____</i>	No Charge

CHOOSE ONE OF THE FOLLOWING:

Part Number	Description	Price AUD
<input type="checkbox"/> I105100-99-0	4"x10" Corpus Armrest Pad Right - Leatherette	No Charge
<input type="checkbox"/> I105104-99-0	4"x13" Corpus Armrest Pad Right - Leatherette	No Charge
<input type="checkbox"/> I107212-99-0	4"x10" Corpus Armrest Pad Right - Black Fabric	No Charge
<input type="checkbox"/> I106047-99-0	4"x13" Corpus Armrest Pad Right - Black Fabric	No Charge
<input type="checkbox"/> I106486-99-0	4"x10" Corpus Arm w/ Stealth Gel Pad - Right	
<input type="checkbox"/> I106488-99-0	4"x13" Corpus Arm w/ Stealth Gel Pad - Right	
<input type="checkbox"/> I104951-99-0	3.5"x8" Unitrack Armrest Pad Right (C)	No Charge
<input type="checkbox"/> I104950-99-0	3.5"x12.5" Unitrack Armrest Pad Right (G)	No Charge

CHOOSE ONE OF THE FOLLOWING:

Part Number	Description	Price AUD
<input type="checkbox"/> I105098-99-0	4"x10" Corpus Armrest Pad Left - Leatherette	No Charge
<input type="checkbox"/> I105102-99-0	4"x13" Corpus Armrest Pad Left - Leatherette	No Charge
<input type="checkbox"/> I107213-99-0	4"x10" Corpus Armrest Pad Left - Black Fabric	No Charge
<input type="checkbox"/> I106046-99-0	4"x13" Corpus Armrest Pad Left - Black Fabric	No Charge
<input type="checkbox"/> I106487-99-0	4"x10" Corpus Arm w/ Stealth Gel Pad - Left	
<input type="checkbox"/> I106489-99-0	4"x13" Corpus Arm w/ Stealth Gel Pad - Left	
<input type="checkbox"/> I104952-99-0	3.5"x8" Unitrack Armrest Pad Left (C)	No Charge



Armrest Options

CHOOSE ONE OF THE FOLLOWING:

Part Number	Description	Price AUD
<input type="checkbox"/> I104949-99-0	3.5"x12.5" Unitrack Armrest Pad Left (G)	No Charge

OPTIONAL ITEMS:

Part Number	Description	Price AUD
<input type="checkbox"/> I107339-99-0	Armrest Pouch - Right	
<input type="checkbox"/> I107338-99-0	Armrest Pouch - Left	
<input type="checkbox"/> I107840-99-0	UT Armrest Side Support - Right Outer <i>Includes 2"T x 6.5"W tapered pad.</i>	
<input type="checkbox"/> I107863-99-0	UT Armrest Side Support - Left Outer <i>Includes 2"T x 6.5"W tapered pad.</i>	
<input type="checkbox"/> I107864-99-0	UT Armrest Side Support - Right Inner <i>Includes 2"T x 6.5"W tapered pad. May require an increase to the width between the armrests to ensure armrests will flip-up.</i>	
<input type="checkbox"/> I107865-99-0	UT Armrest Side Support - Left Inner <i>Includes 2"T x 6.5"W tapered pad. May require an increase to the width between the armrests to ensure armrests will flip-up.</i>	
<input type="checkbox"/> I107839-99-0	UT Elbow Support - Right <i>Must also select PN I10881.</i>	
<input type="checkbox"/> I107838-99-0	UT Elbow Support - Left <i>Must also select PN I10882.</i>	
<input type="checkbox"/> I10881	Adjustable Elbow Support Hardware - Right	
<input type="checkbox"/> I10882	Adjustable Elbow Support Hardware - Left	

Legrest Options

CHOOSE ONE OF THE FOLLOWING:

Part Number	Description	Price AUD
<input type="checkbox"/> I10586	70° Swing Away Legrests w Multi Axis Adj Footplate <i>Includes 5"W x 8"D Multi Axis Footplate. Range of adjustment from the seat pan to the top of each footplate is: 6 - 17".</i>	No Charge
<input type="checkbox"/> I10823	Aeron Legrest Assembly for PS Jr Seat - Narrow <i>90° Legrest Hanger with Manual Adjustment Legrest Elevation. Price includes 6"W x 9"D Corpus Style, Angle Adjustable, Two Piece Footplates. Range of adjustment from the seat pan to the top of the footplate is: 12 - 17" (with calf pads) and 6.5 - 17" (without calf pads).</i>	
<input type="checkbox"/> I11052	Aeron Legrest Assembly for PS Jr Seat - Narrow <i>90° Legrest Hanger with Manual Adjustment Legrest Elevation. Price includes 12.5"W x 10"D, Angle Adjustable, One Piece Footplate. Range of adjustment from the seat pan to the top of the footplate is: 12 - 17" (with calf pads) and 6.5 - 17" (without calf pads).</i>	

OPTIONAL ITEMS:

Part Number	Description	Price AUD
<input type="checkbox"/> I104684-99-0	UT Calf Support Kit (B) w Mesh Cover & Adj Hardware <i>3.5"W x 4.5"T Pad. Recommended pad to be used with lower leg lengths of 12" or less. Bracket designed for use with Aeron Center Mount Legrest Assembly. Cannot be ordered with Swing Away Legs (I10586 & I10587).</i>	No Charge
<input type="checkbox"/> I104686-99-0	UT Calf Support Kit (J) w Mesh Cover & Adj Hardware <i>5"W x 6.5"T Pad. Bracket designed for use with Aeron Center Mount Legrest Assembly. Cannot be ordered with Swing Away Legs (I10586 & I10587).</i>	No Charge
<input type="checkbox"/> I107503-99-0	BP Ankle Huggers, Med Side Release, Adj (Pair) <i>This item will be sent as a part with the wheelchair. It will not be installed.</i>	

Accessories



Accessories

MUST CHOOSE POSITIONING BELT:

Part Number	Description	
	NOTE:	
	<i>See supplemental order form for more Bodypoint options.</i>	
<input checked="" type="checkbox"/> I104518-99-0	Bodypoint Non-Padded Hip Belt - Medium	No Charge
<input type="checkbox"/> I105849-99-0	Bodypoint Padded Hip Belt - Extra Small	

CHOOSE ONE OF THE FOLLOWING:

Part Number	Description	Price AUD
<input type="checkbox"/> I100749-99-0	MF Headrest for Koala/PS Jr - Leatherette <i>Must also select PN I10552.</i>	
<input type="checkbox"/> I10552	Adj/Removable Headrest Hardware	

OPTIONAL ITEMS:

Part Number	Description	Price AUD
<input type="checkbox"/> I103077-99-0	UT Hip Suppt (B) w/ Mesh Cover - Low (Pair) <i>4.5"W x 3.5"T Pad w/ 2" Tall Mounting Bracket. Must also select PN I10550.</i>	
<input type="checkbox"/> I103079-99-0	UT Hip Suppt (A) w/ Mesh Cover - Low (Pair) <i>6"W x 5"T Pad w/ 2" Tall Mounting Bracket. Must also select PN I10550.</i>	
<input type="checkbox"/> I103095-99-0	UT Thigh/Hip Suppt (C) w/ Mesh Cover - High (Pair) <i>8"W x 3.5"T Pad w/ 4" Tall Mounting Bracket. Must also select PN I10550.</i>	
<input type="checkbox"/> I10550	Adj Removable - Thigh Support Hardware (Pair)	
<input type="checkbox"/> I107499-99-0	Universal Headrest Adapter <i>Please Note: Recline is limited to 160° or less (depending on backrest height) when ordering a Universal Headrest Adapter on a 23" or 24" tall Corpus backrest.</i>	
<input type="checkbox"/> I106210-99-0	RAM Self-Leveling Cup Holder	
<input type="checkbox"/> I106806-99-0	RAM X-Grip Phone Holder	
<input type="checkbox"/> I107100-99-0	RAM X-Grip 10" Tablet Holder w/Swing Away Mount <i>Compatible with most 10" screen tablets with or without a protective case. Includes adjustable swing away mount (I103476-99-0). Note: Customization may be required if ordered along with an Omni.</i>	
<input type="checkbox"/> I107222-99-0	USB Charger <i>Provides up to 1.5A (at 5V) to charge most smartphones and tablets. The charger features an isolated circuit design for improved safety; and by pulling power from the ICS, prevents battery drainage when the chair is powered off. Includes mounting hardware and ICS connector. Please select armrest for mounting, keeping in mind that charger may need to be relocated. Mount on left armrest: _____ Mount on right armrest: _____</i>	
<input type="checkbox"/> I105805-99-0	Medical Necessities Bag Hooks <i>Please Note: Recline is limited to 165° or less (depending on backrest height) when ordering with 23" or 24" tall Ergo backrest. Custom required if ordered with a ROHO AGILITY backrest.</i>	
<input type="checkbox"/> I102059-99-0	Medical Necessities Bag <i>Includes hardware to attach bag to backrest. Please Note: Recline is limited to 165° or less (depending on backrest height) when ordering with a 23" or 24" tall Ergo backrest. Custom required if ordered with a ROHO AGILITY backrest.</i>	
<input type="checkbox"/> ISH212	BP Shoulder Harness Medium, Front Pull <i>This item will be sent as a part with the wheelchair. It will not be installed.</i>	
<input type="checkbox"/> I100578-99-0	Small Angle Adjustable Tray <i>Must also select PN I10556.</i>	
<input type="checkbox"/> I105221-99-0	Angle Adjustable Tray - 14" W <i>Must also select PN I10556.</i>	
<input type="checkbox"/> I105222-99-0	Angle Adjustable Tray - 16" W <i>Must also select PN I10556.</i>	
<input type="checkbox"/> I106863-99-0	Angle Adjustable UniTrack Tray <i>Side mounted tray that attaches to the UniTrack Accessory Rail. Dimensions are: 19.5"W x 13"D. Can only be ordered with chairs that offer the UniTrack System. Must also select PN I10556.</i>	
<input type="checkbox"/> I10556	Adj Removable - Tray Hardware (1 Each)	



Accessories

OPTIONAL ITEMS:

Part Number	Description	Price AUD
<input type="checkbox"/> I103798-99-0	After Market Seating System Adapter Used to connect after market seating systems to the UniTrack Accessory Rail on the PS Jr Seat Frame. Includes Mounting Hardware.	

Stealth Accessories

OPTIONAL ITEMS:

Part Number	Description	Price AUD
<input type="checkbox"/> I106222-99-0	Stealth 3" x 5" Lateral Supports - PS Jr Must also select I10723. Includes THIN Stealth hardware and specialized mounting bracket for Flip Back Armrest assembly. Only compatible with Flip Back Armrest hardware.	
<input type="checkbox"/> I10723	Adj Removable Stealth Lateral Hardware (Pair)	

Alternative Drive Systems

ALTERNATIVE PROPORTIONAL JOYSTICK SYSTEMS:

Part Number	Description	Price AUD
<input type="checkbox"/> I10751	Mini-Joystick for Omni / IOM (9-pin connector) Operating Force = 13 g Neutral to Full Deflection = 5mm	
<input type="checkbox"/> I103893-99-0	Compact Joystick for Omni / IOM (9-pin connector) Operating Force = 220 g Neutral to Full Deflection = 28mm Joystick includes two 3.5mm mono jacks for Mode and On/Off. Includes Standard Joystick Knob and 1.75" (42mm) Foam Ball Knob.	
<input type="checkbox"/> I103894-99-0	Compact Joystick LITE for Omni / IOM (9-pin conn) Operating Force = 50 g Neutral to Full Deflection = 17mm Joystick includes two 3.5mm mono jacks for Mode and On/Off. Includes Standard Joystick Knob and 1.75" (42mm) Foam Ball Knob.	

ALT PROPORTIONAL JOYSTICK MOUNTING HARDWARE:

Part Number	Description
<input type="checkbox"/> I104690-99-0	Alt. JS Swing Away Mount Kit - Armrest NOTE: Please see the PERMOfix Mounting Solutions, within the Alt Drive Accessories section, for additional mounting options. Mount in front of Left Arm: _____ Mount in front of Right Arm: _____
<input type="checkbox"/> I107135-99-0	Alt JS/Omni Combo Swing Away Mount Kit - Armrest Used to mount Omni and Alternative Proportional Joystick to the same armrest. Includes Swing Away Omni/Tablet Mount, with a longer mounting rod, two PERMOfix modules, one 3" rod and one 4" rod. Mount in front of Left Arm: _____ Mount in front of Right Arm: _____
<input type="checkbox"/> I104695-99-0	Alt. JS Adjustable Mount Kit - Midline Mount on Left Arm: _____ Mount on Right Arm: _____
<input type="checkbox"/> I104143-99-0	Alt. JS Adjustable Mount Kit - Chin Boom Swing-Away / Removable Chin Boom mounted on right of backrest shell. Custom charge may be applied if ordering lower backrest heights.
<input type="checkbox"/> I10752	Alt. JS Adjustable Mount Kit - Haims Harness, Lrg For neck width up to 19cm.
<input type="checkbox"/> I10753	Alt. JS Adjustable Mount Kit - Haims Harness, Sml For neck width up to 14cm.

ALT PROPORTIONAL JOYSTICK ACCESSORIES:

Part Number	Description	Price AUD
-------------	-------------	-----------



Alternative Drive Systems

ALT PROPORTIONAL JOYSTICK ACCESSORIES:

Part Number	Description	Price AUD
<input type="checkbox"/> I104904-99-0	Dual Switches for Compact Joystick - Fixed Mount <i>Includes two pushbutton switches mounted on rigid Switches can be used for Mode & On/Off. Switch operation force = 130 g. Mount at front of Joystick: _____ Mount at rear of Joystick: _____ Mount at left of Joystick: _____ Mount at right of Joystick: _____</i>	
<input type="checkbox"/> I104641-99-0	Dual Switches for Compact Joystick - Gooseneck Mnt <i>Includes two pushbutton switches mounted on adjustable gooseneck arms. Switches are positioned on Left and Right sides of Compact Joystick. Switches can be used for Mode & On/Off. Switch operation force = 130 g</i>	

SIP & PUFF SYSTEMS:

Part Number	Description	Price AUD
<input type="checkbox"/> I10573	Sip & Puff System for Omni <i>Requires Omni Control Module. Cannot use this part with an Input Module. Includes breath tube kit.</i>	
<input type="checkbox"/> I10574	Sip & Puff System with Fiber Optic Opti-Stop <i>Requires Omni Control Module. Cannot use this part with an Input Module. Includes breath tube kit.</i>	
<input type="checkbox"/> I106807-99-0	Sip & Puff Kit for Total Control Head Array <i>For use with Total Control Head Array. Requires Omni Control Module. Cannot be used with an Input Module. Includes breath tube kit.</i>	
<input type="checkbox"/> I10762	Adjustable Mounting Hrdw - Sip & Puff <i>Required when selecting PN I10573, I10574 or I106807-99-0.</i>	

PROPORTIONAL FOOT CONTROL:

Part Number	Description	Price AUD
<input type="checkbox"/> I10767	Proportional Foot Control <i>Mount on Left: _____ Mount on Right: _____</i>	
<input type="checkbox"/> I10766	Adjustable Mounting Hrdw - Foot Control <i>Required when selecting PN I10767.</i>	

Total Control Head Array

CONFIGURE HEAD ARRAY:

Part Number	Description	Price AUD
<input type="checkbox"/> I106805-99-0	Permobil Head Array - 3-Switch System <i>Features an occipital pad sensor and two stainless steel arms configured with proximity sensors. The arms are mounted to the upper attachment points on the head array. Longer arms are also included to extend the sensor adjustment range. ESP recommended. Choose Omni Control Module (ID51154IR) separately.</i>	
<input type="checkbox"/> I11072	Program Head Array as a 5-Switch System <i>Select additional arms/sensors below and indicate the intended use of each sensor in the Notes & Comments section. Choose Omni Control Module (ID51154IR) or Input Module (I10581) separately.</i>	No Charge
<input type="checkbox"/> I105920-99-0	CHOOSE OCCIPITAL PAD SIZE (REQUIRED): Small Occipital Pad for Head Array - 4.25"W x 3"H	No Charge
<input type="checkbox"/> I105921-99-0	Large Occipital Pad for Head Array - 6"W x 4"H	No Charge
<input type="checkbox"/> I10847	CHOOSE OCCIPITAL PAD MOUNT (REQUIRED): Fixed Mount for Occipital Pad - Head Array	No Charge
<input type="checkbox"/> I105936-99-0	Adj. Swivel Mount for Head Array Occipital Pad <i>Enables Occipital Pad to be angled up/down and left/right independent of the lateral arms.</i>	No Charge
	CHOOSE BACKREST MOUNTING HARDWARE (REQUIRED):	



Total Control Head Array

CONFIGURE HEAD ARRAY:

Part Number	Description	Price AUD
<input type="checkbox"/> I105927-99-0	Stealth TWB480 Backrest Mounting Hdwr with Links <i>Please select the length of the vertical stem below: Standard length (10"): _____ Short length (6"): _____</i>	
<input type="checkbox"/> I105941-99-0	Stealth TWB480 Backrest Mtng Hdwr with Flipback	
<input type="checkbox"/> I105925-99-0	Stealth HMO475 Backr Mtng Hdwr w 17mm Ball	
<input type="checkbox"/> I106319-99-0	Stealth TWB480 Backrest Mnt Hdwr w SubOcptl Mount	
<input type="checkbox"/> I106282-99-0	Stealth 600 SubOccipital Pad - 8.5"W x 2.75"H <i>Must select when PN I106319-99-0 is ordered.</i>	
SELECT ADDITIONAL LEFT ARM PACKAGE (OPTIONAL):		
<input type="checkbox"/> I10976	Add Left Lateral Arm w Proximity Sensor <i>Choose Lateral Arm Adjustment Range: Short (5.5" - 10"): _____ Long (5.5" - 14"): _____</i>	
<input type="checkbox"/> I10977	Add Left Lateral Arm w Egg Switch <i>Choose Lateral Arm Adjustment Range: Short (5.5" - 10"): _____ Long (5.5" - 14"): _____</i>	
<input type="checkbox"/> I10978	Add Left Lateral Arm w Pad (No Sensor or Switch) <i>Choose Lateral Arm Adjustment Range: Short (5.5" - 10"): _____ Long (5.5" - 14"): _____</i>	
SELECT ADDITIONAL RIGHT ARM PACKAGE (OPTIONAL):		
<input type="checkbox"/> I10979	Add Right Lateral Arm w Proximity Sensor <i>Choose Lateral Arm Adjustment Range: Short (5.5" - 10"): _____ Long (5.5" - 14"): _____</i>	
<input type="checkbox"/> I10980	Add Right Lateral Arm w Egg Switch <i>Choose Lateral Arm Adjustment Range: Short (5.5" - 10"): _____ Long (5.5" - 14"): _____</i>	
<input type="checkbox"/> I10981	Add Right Lateral Arm w Pad (No Sensor or Switch) <i>Choose Lateral Arm Adjustment Range: Short (5.5" - 10"): _____ Long (5.5" - 14"): _____</i>	

Alt Drive Systems Setup

CHOOSE OMNI OR INPUT MODULE (REQUIRED):

Part Number	Description	Price AUD
<input type="checkbox"/> ID51154IR	Omni Control Module for R-net w/ Infrared <i>Must also select mounting hardware, either I103476-99-0 or I107135-99-0. Please specify which type of alternative drive input you will be using if not ordering it from Permobil. Input Device: _____</i>	
<input type="checkbox"/> I103476-99-0	Adjustable Swing-Away Hardware - Omni/Tablet <i>Mount in front of Left Arm: _____ Mount in front of Right Arm: _____</i>	
<input type="checkbox"/> I10581	Input Module (IOM) for R-net Control System <i>This device allows a proportional or switched drive input to be connected to an R-net Control System without the use of an Omni. It must be ordered in conjunction with an R-net Remote Joystick. It does NOT support any of the following input devices: Sip-n-Puff, 3-Switch Head Array or RIM control.</i>	

CHOOSE MOUNTING POSITION OF SEAT CONTROL SWITCHBOX:

Part Number	Description	Price AUD
<input type="checkbox"/> I10080	Mount Seat Control Switchbox on Backrest	No Charge

Alt Drive Accessories

PERMOFIX MOUNTING SOLUTIONS:

Part Number	Description	Price AUD
-------------	-------------	-----------



Alt Drive Accessories

PERMOFIX MOUNTING SOLUTIONS:

Part Number	Description	Price AUD
<input type="checkbox"/> I106655-99-0	PERMOfix UniTrack Fixed Mount Kit <i>Includes two base joints, two 4" rods and one 8" rod.</i>	
<input type="checkbox"/> I106656-99-0	PERMOfix UniTrack Removable Mount Kit <i>Includes quick release mount and thigh support hardware, two base joints and three different length rods.</i>	
<input type="checkbox"/> I106657-99-0	PERMOfix UniTrack Swing Away Mount Kit <i>Includes PERMOfix Swing Away module, two base joints, one 2" rod and one 4" rod.</i>	
<input type="checkbox"/> I105907-99-0	Proximity Sensor + Swivel Assm - Not Mounted <i>Only available with the Total Control Head Array. Swivel assembly allows sensor to be mounted with PERMOfix or directly to lateral arm of Total Control Head Array. Since Proximity Sensors require power, they CANNOT be used as On/Off switches.</i>	
<input type="checkbox"/> I105919-99-0	Egg Switch + Swivel Mount Assm - Not Mounted <i>Swivel assembly allows switch to be mounted with PERMOfix or directly to lateral arm of Total Control Head Array.</i>	
<input type="checkbox"/> I105912-99-0	Solid Lateral Pad+Swivel Joint - Not Mounted <i>Swivel assembly allows pad to be mounted with PERMOfix or directly to lateral arm of Total Control Head Array.</i>	
<input type="checkbox"/> I106675-99-0	Oval Plate + Swivel Assm - No Switch - Not Mounted <i>Blank plate provides mounting option for Micro Lite, Buddy Button or other switch options. Swivel assembly allows pad to be mounted with PERMOfix or directly to lateral arm of Total Control Head Array.</i>	

OPTIONAL ITEMS:

Part Number	Description	Price AUD
<input type="checkbox"/> I101769-99-0	Micro Lite Switch Qty: _____	
<input type="checkbox"/> I101767-99-0	Buddy Button Qty: _____	
<input type="checkbox"/> I106293-99-0	Egg Switch Black Qty: _____	
<input type="checkbox"/> I101771-99-0	Mini Cup Switch Qty: _____	
<input type="checkbox"/> I107500-99-0	Proximity Switch Qty: _____	
<input type="checkbox"/> I107495-99-0	Mini Joystick 5 Switch OC for Omni/IOM (9-pin)	
<input type="checkbox"/> I107498-99-0	Touch Contact with CR2032 Battery Qty: _____	
<input type="checkbox"/> I105935-99-0	3.5mm Ext. Cable for Sensor/Switch, 0.3m Long	
<input type="checkbox"/> I105931-99-0	3.5mm Ext. Cable for Sensor/Switch, 1m Long	
<input type="checkbox"/> I105932-99-0	3.5mm Ext. Cable for Sensor/Switch, 2m Long	
<input type="checkbox"/> I105933-99-0	3.5mm Ext. Cable for Sensor/Switch, 3m Long	
<input type="checkbox"/> I105934-99-0	3.5mm Y-splitter Cable for Sensor/Switch <i>Allows two Proximity Sensors/Mechanical Switches to connect to the same direction jack.</i>	

Additional Order Instructions:

Chair Order Policy:

You are highly encouraged to contact Permobil Sales with ANY questions regarding proper completion of our order forms. The order cannot be processed without the following:

- A client name or code, the client's measurements, the client's weight
- The order form must also be signed by a person authorized by your company to acknowledge the items selected. Permobil is not responsible for configuration or size discrepancies resulting from customer errors on the order form.
- Order confirmations will be provided by Permobil summarizing the items selected. Please review this carefully.
- Once the order is shipped, any changes are subject to the returns policy and may be prohibited.

Order Acknowledgement:

I, _____, am an agent of the medical equipment provider named on this order form and I have the authority to contract for the purchase of powered wheelchairs and related parts on behalf of said provider. I acknowledge that I have reviewed this order and that it is complete and accurate to the best of my knowledge. I further acknowledge that any changes to this order after submission of this order form are subject to the returns policy and may carry additional charges.

Signed: _____

Printed Name: _____

Title: _____

Date: _____

NOTE: All specifications and prices are subject to change without notice. Please note that prices displayed are only valid if ordered with the wheelchair. Any parts ordered for service or future alterations will carry different pricing and warranty.