

K450 MX



* Indicates a Required Field

Dealer Information

* Contact: _____

Dealer Code: _____

* Dealer Name: _____

Address: _____

* City: _____

* State: _____

* Phone#: _____

* Fax#: _____

PO#: _____

Email Address: _____

Client Information

Permobil recommends that the client is evaluated by a certified rehab specialist.

* First Name: _____

* Last Name: _____

Diagnosis: _____

Funding Source: _____

Client Age: _____

Client Measurements

* Weight: _____

* Height: _____

A) Top of Shoulders: _____

B) Chest Depth: _____

C) Chest Width: _____

D) Seat Depth: _____

E) Top of Head: _____

F) Elbow to Hand: _____

G) Seat Pan to Elbow: _____

H) Hip Width: _____

*I) Knee to Foot: _____





Chair Model and Colors

		Price AUD
<input checked="" type="checkbox"/>	I104524-99-0 K450 Base w/ R-net 8kmh <i>Price Includes: K450 Power Base, Tie-Down Hardware for Strap Systems, Charger, Flat Free Tires, and Standard Positioning Belt. Must Choose Flat Free or Pneumatic Tires Below. Max. Speed = 8 km/h, Weight Capacity = 60 kg, STFH = 75 mm - 635 mm.</i>	

CHOOSE DRIVE TIRES:

	Part Number	Description	Price AUD
<input type="checkbox"/>	I106204-99-0	Flat Free Casters - K450	No Charge
<input type="checkbox"/>	I106205-99-0	Flat Free Drive Tires - K450	No Charge
<input type="checkbox"/>	I106206-99-0	Pneumatic Casters - K450	No Charge
<input type="checkbox"/>	I106207-99-0	Pneumatic Drive Tires - K450	No Charge

CHOOSE SHROUD COLOR:

	Part Number	Description	Price AUD
<input type="checkbox"/>	I104535-99-0	Shroud for K450 - Onyx Black	No Charge
<input type="checkbox"/>	I104536-99-0	Shroud for K450 - Platinum Silver	No Charge
<input type="checkbox"/>	I104537-99-0	Shroud for K450 - Pop Star Pink	No Charge
<input type="checkbox"/>	I104539-99-0	Shroud for K450 - Extreme Yellow	No Charge
<input type="checkbox"/>	I106367-99-0	Shroud for K450 - Volcano Red	No Charge
<input type="checkbox"/>	I106368-99-0	Shroud for K450 - Cobalt Blue	No Charge
<input type="checkbox"/>	I107234-99-0	Shroud for K450 - Purple Rain	No Charge

CHOOSE ONE OF THE FOLLOWING:

	Part Number	Description	Price AUD
<input checked="" type="checkbox"/>	I10901	Australian Freight & Handling Fee	

Base Options

MUST CHOOSE SEAT ELEVATOR:

	Part Number	Description	Price AUD
<input checked="" type="checkbox"/>	I10510	Power Adjustable Seat Height - 8" Travel <i>This item must be selected.</i>	

MUST CHOOSE BATTERIES:

	Part Number	Description	Price AUD
<input checked="" type="checkbox"/>	I107032-99-0	Batteries, Grp 34 (60Ah) Sealed Gel, x2, Installed	

OPTIONAL ITEMS:

	Part Number	Description	Price AUD
<input type="checkbox"/>	I104550-99-0	Support Wheels for K450	
<input type="checkbox"/>	I104540-99-0	Light Kit - K450 R-net	
<input type="checkbox"/>	I104553-99-0	WC 19 Transit Kit - K450 <i>This kit includes additional anchor points for tie down straps, seat stabilizing bars and hardware, and a WC-19 Approved integrated lap and shoulder belt. This kit has been approved according to ANSI/RESNA WC19. Cannot be ordered with a 10"W Seat Size (I104317-99-0). An ICS Switchbox must be installed on the chair in order to access WC19 mode.</i>	
<input type="checkbox"/>	I104554-99-0	Fixed Ventilator Tray - K450 <i>Tray is 9.75"D x 13"W and increases the overall length of a K450 by approximately 11". The ventilator tray is designed to accommodate an LP10 style ventilator and one U1 Battery. A custom charge will need to be applied if ordering vent tray with any of the following items: WC-19 Transit Kit (I104553-99-0), Attendant Control (ID50882), Push Handles (I104576-99-0) or Bag Hooks (I105805-99-0). Please Note: Must also select Support Wheels (I104550-99-0).</i>	
<input type="checkbox"/>	ID51113	Programming Key for R-net - Permobil Version	



Joystick Options

PLEASE NOTE:

For additional input devices, refer to Alternative Drive Systems section.

I108319-99-0

Permobil Joystick Module

Designed exclusively for Permobil, this R-net color joystick features: control of up to 4 Bluetooth devices, Infrared control, programmable buttons and a programmable switch jack; all in a compact, durable Aluminum housing.

I105246-99-0

R-net LED Remote Joystick Module

No Charge

This Joystick does NOT have a Display Screen and does NOT contain mono jacks. Limited to five profiles with this Joystick Module.

CHOOSE JOYSTICK HANDLE TYPE:

Part Number	Description	Price AUD
<input type="checkbox"/> I10045	Cone Shaped Joystick Knob - Standard	No Charge
<input type="checkbox"/> I100960-99-0	Chin Cup For Joystick Knob	
<input type="checkbox"/> I100961-99-0	"T" Handle For Joystick Knob	
<input type="checkbox"/> I100962-99-0	"Large Ball" For Joystick Knob	
<input type="checkbox"/> I100963-99-0	"Softball" For Joystick Knob	
<input type="checkbox"/> I100964-99-0	"Mushroom" For Joystick Knob	
<input type="checkbox"/> I100965-99-0	"Stick" For Joystick Knob	
<input type="checkbox"/> I105442-99-0	1.75" Black Foam Ball For Joystick Knob	
<input type="checkbox"/> I105846-99-0	Bodypoint J/S Handle 3" U-Shaped w/ Flex-Shaft	
<input type="checkbox"/> I105847-99-0	Bodypoint J/S Handle 4" U-Shaped w/ Flex-Shaft	
<input type="checkbox"/> I105848-99-0	Bodypoint J/S Handle Dome Shaped, Rubber	

CHOOSE JOYSTICK MOUNT:

Part Number	Description	Price AUD
<input type="checkbox"/> I10412	Joystick Mount Right - Fixed, R-net Only	No Charge
<input type="checkbox"/> I10413	Joystick Mount Left - Fixed, R-net Only	No Charge
<input type="checkbox"/> I102537-99-0	Joystick Mount Right - Swing Away	
<input type="checkbox"/> I102538-99-0	Joystick Mount Left - Swing Away	
<input type="checkbox"/> ID50882-A	Attendant Control for R-net Control System Only <i>On the K450, the mounting location for the Attendant Control is off of the Push Handle Bracket. If ordering an attendant control, you must also order I104576-99-0 (Push Handles - MX Seat).</i>	
<input type="checkbox"/> I104571-99-0	Co-Pilot Attendant Control - K450	
<input type="checkbox"/> I107350-99-0	R-net Bluetooth iDevice Module <i>Allows input device commands (such as joystick nudges or switch clicks) to be used to operate an iPhone, iPad or iPod. By assigning each input to an iOS function, users can make phone calls, send texts and emails and access applications. Kit includes a Bluetooth iDevice Module, an R-net Bus Cable, and two mono jack cables for direct switch inputs (switches not included). Must be ordered with a Color JSM or Omni for mode display.</i>	
<input type="checkbox"/> I107367-99-0	R-net Bluetooth Mouse Module <i>Enables mouse control through any R-net input device. It can operate the mouse on a computer and on select Android devices (with R-net for Android app). Kit includes a Bluetooth Mouse Module, an R-net Bus Cable, two mono jack cables for L/R click (switches not included) and a USB BT Adapter for user's computer. Note: BT Adapter may not work with all computers. Must be ordered with a Color JSM or Omni for mode display.</i>	
<input type="checkbox"/> I107368-99-0	Add a second Bluetooth Mouse Module <i>Allows R-net to simultaneously pair to a second compatible computer or Android device. Each module is accessible through a separate R-net mode. Kit includes a unique Bluetooth Mouse Module, an R-net Bus Cable, two mono jack cables for L/R click and a USB BT Adapter.</i>	

Seating System

CHOOSE ONE OF THE FOLLOWING:

Part Number	Description	Price AUD
-------------	-------------	-----------



Seating System

CHOOSE ONE OF THE FOLLOWING:

Part Number	Description	Price AUD
<input checked="" type="checkbox"/> I104324-99-0	MX Power Seating System <i>Price Includes: MX Seating System, Adjustable Seat Frame, Power Seat to Floor System, Adjustable Backrest Angle, and One Piece Footplate. Weight Limit 57kg.</i>	

CHOOSE ONE OF THE FOLLOWING:

Part Number	Description	Price AUD
<input type="checkbox"/> I104529-99-0	No Tilt <i>Weight Limit 57kg.</i>	No Charge
<input type="checkbox"/> I104551-99-0	45° Power Tilt <i>Weight Limit 57kg.</i>	

Seating System Options

CHOOSE SEAT FUNCTION CONTROL TYPE:

Part Number	Description	Price AUD
<input checked="" type="checkbox"/> I10585	Seat Function Control Kit for R-net <i>Please indicate below if you want push buttons or toggles on the ICS Switchbox.</i>	

CHOOSE PUSH BUTTONS OR TOGGLES:

Part Number	Description	Price AUD
<input type="checkbox"/> I104135-99-0	Push Buttons on Seat Function Control Kit	No Charge
<input type="checkbox"/> I104136-99-0	Toggles on Seat Function Control Kit	

CHOOSE MOUNTING POSITION OF ICS SWITCHBOX:

Part Number	Description	Price AUD
<input type="checkbox"/> I10626	Mount ICS Switchbox onto R-net Joystick Bracket <i>ICS Switchbox will be mounted to the same bracket as the R-net Joystick Module. If this item is selected, the ICS Switchbox will swing away with the Joystick Module.</i>	No Charge
<input type="checkbox"/> I104651-99-0	Mount ICS Switchbox onto Armrest <i>Mount on Left: _____ Mount on Right: _____</i>	No Charge
<input type="checkbox"/> I11050	Mount ICS Switchbox on Backrest <i>On the K450, the backrest mounting location for the ICS Switchbox is on the Push Handle bracket; therefore, you must also order I104576-99-0 (Push Handles - MX). Mount on Left: _____ Mount on Right: _____</i>	

OPTIONAL ITEMS:

Part Number	Description	Price AUD
<input type="checkbox"/> I105409-99-0	Do Not Install ICS Switchbox, Send as Part	No Charge
<input type="checkbox"/> I104196-99-0	ICS Alternative Switchbox for R-net Electronics <i>Seat function interface box that accommodates up to 8 momentary switches. Provides an additional method for accessing the seat functions. 1/8" mini-mono switches not included.</i>	
<input type="checkbox"/> I106743-99-0	Remote Stop System - All Electronics <i>Consists of a Radio Frequency (RF) handheld transmitter and an RF receiver mounted on the wheelchair. Custom charges may be required when ordering with a PS/PS Jr seat, a low backrest height or a pediatric wheelchair. If ordering on a K450, push handles must also be selected.</i>	
<input type="checkbox"/> I10515	Memory Seat Program for K450 MX <i>If this is selected, the chair will be programmed to allow the user to store and save up to 3 different memory seat positions through the ICS Switchbox. On the K450 MX, ALL seat functions are capable of recalling a saved seat position. Please refer to the K450 Owner's Manual for more information on saving and recalling a memory seat position.</i>	No Charge

Backrest Options

CHOOSE ONE OF THE FOLLOWING:

Part Number	Description	Price AUD
-------------	-------------	-----------



Backrest Options

CHOOSE ONE OF THE FOLLOWING:

Part Number	Description	
MX STYLE BACKRESTS:		
<i>Use this section to select MX Style Backrests. Back heights measured from the seat pan to the top of the backrest shell with the backrest set at 90°. Please be sure to pick both a backrest shell size and a backrest cushion style (Leatherette or Black Mesh Fabric). Note: A 16"W MX Backrest is NOT available.</i>		
<input type="checkbox"/> I104321-99-0	MX Backrest - 10"W x 14"H <i>Backrest can be adjusted to either 14" or 16" high. If ordering a headrest, backrest will be mounted in the 16" high position unless client measurements require a 14" backrest height, in which case, custom modifications will be done.</i>	No Charge
<input type="checkbox"/> I104322-99-0	MX Backrest - 12"W x 16"H <i>Backrest can be adjusted to either 16" or 18" high.</i>	No Charge
<input type="checkbox"/> I104323-99-0	MX Backrest - 14"W x 20"H <i>Backrest can be adjusted to either 20" or 22" high.</i>	No Charge

CHOOSE ONE OF THE FOLLOWING:

Part Number	Description	
Please Select Cushion Style Below:		
<input type="checkbox"/> I104309-99-0	MX/MF Backrest Cushion, 10"W x 14"T - Leatherette	No Charge
<input type="checkbox"/> I104311-99-0	MX/MF Backrest Cushion, 12"W x 16"T - Leatherette	No Charge
<input type="checkbox"/> I104313-99-0	MX/MF Backrest Cushion, 14"W x 20"T - Leatherette	No Charge
<input type="checkbox"/> I104308-99-0	MX/MF Backrest Cushion, 10"W x 14"T - Black Mesh	No Charge
<input type="checkbox"/> I104310-99-0	MX/MF Backrest Cushion, 12"W x 16"T - Black Mesh	No Charge
<input type="checkbox"/> I104312-99-0	MX/MF Backrest Cushion, 14"W x 20"T - Black Mesh	No Charge

CHOOSE ONE OF THE FOLLOWING:

Part Number	Description	
PS JR TUBULAR STYLE BACKRESTS		
<i>Use this section to select PS Jr Tubular Style Backrest Canes. Note: A 10"W Tubular Backrest is NOT available.</i>		
<input type="checkbox"/> I103717-99-0	PS Jr Back Canes, Height & Angle Adjustable <i>Must also select PS Jr Style Armrests if ordering this backrest. Choose PS Jr Style Backrest Type (Metal Back or Omit Back) below. Backrest width should match seat width. Back height is adjustable from 16.5"T - 22.5"T, measured from the seat pan to the top of the back.</i>	No Charge

CHOOSE ONE OF THE FOLLOWING:

Part Number	Description	Price AUD
<input type="checkbox"/> I103719-99-0	PS Jr Solid Backrest - 12"W	No Charge
<input type="checkbox"/> I103721-99-0	PS Jr Solid Backrest - 14"W	No Charge
<input type="checkbox"/> I103723-99-0	PS Jr Solid Backrest - 16"W	No Charge
<input type="checkbox"/> I10449	PS Jr - Omit Backrest	No Charge

Seat Sizing

CHOOSE ONE OF THE FOLLOWING:

Part Number	Description	Price AUD
<input type="checkbox"/> I104317-99-0	MX Seat Size 10"W x 10"D <i>Seat depth can be adjusted to 10", 11", or 12". If ordering this seat size, the lowest allowed seat position when using the seat to floor function will be approximately 3.5" from the floor to the top of the seat pan.</i>	No Charge
<input type="checkbox"/> I104318-99-0	MX Seat Size 12"W x 12"D <i>Seat depth can be adjusted to 12", 13", or 14". If ordering this seat size with MX Arms, the lowest allowed seat position when using the seat to floor function will be approximately 3.5" from the floor to the top of the seat pan. If ordering Flip Back Arms, the lowest allowed seat position will be 8".</i>	No Charge



Seat Sizing

CHOOSE ONE OF THE FOLLOWING:

Part Number	Description	Price AUD
<input type="checkbox"/> I104319-99-0	MX Seat Size 14"W x 14"D <i>Seat depth can be adjusted to 14", 15", or 16". If ordering this seat size with MX Arms, the lowest allowed seat position when using the seat to floor function will be approximately 3.5" from the floor to the top of the seat pan. If ordering Flip Back Arms, the lowest allowed seat position will be 8".</i>	No Charge
<input type="checkbox"/> I104320-99-0	MX Seat Size 16"W x 16"D <i>Seat depth can be adjusted to 16", 17", or 18". If ordering this seat size with MX Arms, the lowest allowed seat position when using the seat to floor function will be approximately 10" from the floor to the top of the seat pan. If ordering Flip Back Arms, the lowest allowed seat position will be 16".</i>	No Charge

OPTIONAL ITEMS:

Part Number	Description
<input type="checkbox"/> I306384-99-90-0	Seat Cushions Below Are Approximately 2" Thick: Leatherette Seat Cush MX - 10"W x 10"D (250x250)
<input type="checkbox"/> I306385-99-90-0	Leatherette Seat Cush MX - 12"W x 12"D (300x300)
<input type="checkbox"/> I306386-99-90-0	Leatherette Seat Cush MX - 14"W x 14"D (350x350)
<input type="checkbox"/> I104703-99-0	Leatherette Seat Cush MX - 16"W x 16"D (400x400)
<input type="checkbox"/> I104305-99-0	Black Mesh Seat Cushion MX/MF - 10"W x 10"D
<input type="checkbox"/> I104306-99-0	Black Mesh Seat Cushion MX/MF - 12"W x 12"D
<input type="checkbox"/> I104307-99-0	Black Mesh Seat Cushion MX/MF - 14"W x 14"D
<input type="checkbox"/> I104702-99-0	Black Mesh Seat Cush MX - 16"W x 16"D (400x400)

Armrest Options

CHOOSE ARMREST STYLE:

Part Number	Description	Price AUD
<input type="checkbox"/> I10512	Height Adjustable Armrest Assembly MX Seat <i>If you are not selecting Tilt, you MUST order this item.</i>	
<input type="checkbox"/> I10514	Height Adjustable Armrest Assembly MX Seat - NC <i>No Charge if ordering Power Tilt (I104551-99-0)</i>	No Charge
<input type="checkbox"/> I10824	Flip Back, Height Adjustable Armrests <i>Mounted off the backrest canes. Includes: 45° of angle adjustability, rotational adjustment, 3" offset adjustment both medially and laterally, and height adjustment of 6" - 14" from the seat pan. In one inch increments, please indicate where you would like the armrests set.</i> SET ARMRESTS AT (6" - 14"): _____	
<input type="checkbox"/> I10825	Flip Back, Height Adjustable Armrests - NC <i>No charge when ordering power tilt. In one inch increments, please indicate where you would like the armrests set.</i> ARMRESTS AT (6" - 14"): _____ SET	No Charge

MX ARMREST PADS:

Part Number	Description	Price AUD
<input type="checkbox"/> I104303-99-0	MX Armrest Pad Short Leatherette 8"L x 3.5"W	No Charge
<input type="checkbox"/> I104304-99-0	MX Armrest Pad Long Leatherette 12.5"Lx3.5"W	No Charge

FLIP BACK ARMREST PAD RIGHT:

Part Number	Description	Price AUD
<input type="checkbox"/> I105100-99-0	4"x10" Corpus Armrest Pad Right - Leatherette	No Charge
<input type="checkbox"/> I105104-99-0	4"x13" Corpus Armrest Pad Right - Leatherette	No Charge
<input type="checkbox"/> I107212-99-0	4"x10" Corpus Armrest Pad Right - Black Fabric	No Charge
<input type="checkbox"/> I106047-99-0	4"x13" Corpus Armrest Pad Right - Black Fabric	No Charge
<input type="checkbox"/> I106486-99-0	4"x10" Corpus Arm w/ Stealth Gel Pad - Right	
<input type="checkbox"/> I106488-99-0	4"x13" Corpus Arm w/ Stealth Gel Pad - Right	



Armrest Options

FLIP BACK ARMREST PAD RIGHT:

Part Number	Description	Price AUD
<input type="checkbox"/> I104951-99-0	3.5"x8" Unitrack Armrest Pad Right (C)	No Charge
<input type="checkbox"/> I104950-99-0	3.5"x12.5" Unitrack Armrest Pad Right (G)	No Charge

FLIP BACK ARMREST PAD LEFT:

Part Number	Description	Price AUD
<input type="checkbox"/> I105098-99-0	4"x10" Corpus Armrest Pad Left - Leatherette	No Charge
<input type="checkbox"/> I105102-99-0	4"x13" Corpus Armrest Pad Left - Leatherette	No Charge
<input type="checkbox"/> I107213-99-0	4"x10" Corpus Armrest Pad Left - Black Fabric	No Charge
<input type="checkbox"/> I106046-99-0	4"x13" Corpus Armrest Pad Left - Black Fabric	No Charge
<input type="checkbox"/> I106487-99-0	4"x10" Corpus Arm w/ Stealth Gel Pad - Left	
<input type="checkbox"/> I106489-99-0	4"x13" Corpus Arm w/ Stealth Gel Pad - Left	
<input type="checkbox"/> I104952-99-0	3.5"x8" Unitrack Armrest Pad Left (C)	No Charge
<input type="checkbox"/> I104949-99-0	3.5"x12.5" Unitrack Armrest Pad Left (G)	No Charge

Legrest Options

CHOOSE ONE OF THE FOLLOWING:

Part Number	Description	Price AUD
<input type="checkbox"/> I104314-99-0	MX Footplate - Small (11.5" W x 8"D)	No Charge
<input type="checkbox"/> I104315-99-0	MX Footplate - Large (14" W x 9.5"D)	No Charge

OPTIONAL ITEMS:

Part Number	Description	Price AUD
<input type="checkbox"/> I103427-99-0	Footstrap Kit	
<input type="checkbox"/> I103104-99-0	Footstrap Kit - Padded	
<input type="checkbox"/> I104338-99-0	Small Anti-Skid Overlay for Footplate - K450	
<input type="checkbox"/> I104337-99-0	Large Anti-Skid Overlay for Footplate - K450	
<input type="checkbox"/> I107503-99-0	BP Ankle Huggers, Med Side Release, Adj (Pair)	

This item will be sent as a part with the wheelchair. It will not be installed.

Accessories

MUST CHOOSE POSITIONING BELT:

Part Number	Description	Price AUD
<input checked="" type="checkbox"/> I104518-99-0	Bodypoint Non-Padded Hip Belt - Medium	No Charge
<input type="checkbox"/> I105849-99-0	Bodypoint Padded Hip Belt - Extra Small	

CHOOSE ONE OF THE FOLLOWING:

Part Number	Description	Price AUD
<input type="checkbox"/> I100749-99-0	MF Headrest for Koala/PS Jr - Leatherette <i>Must also select PN I10552.</i>	
<input type="checkbox"/> I104339-99-0	MF Headrest for MX Back - Leatherette <i>Must also select PN I10552. Also includes a UT Bracket for lateral adjustment of the headrest. Please indicate below which size you need.</i> UT Brack for Lat Head Adj (10"W MX Back): ____ <i>Allows for Lateral Adjustment of the headrest.</i> UT Brack for Lat Head Adj(12-14"W MX Back): ____ <i>Allows for Lateral Adjustment of the headrest.</i>	
<input type="checkbox"/> I10552	Adj/Removable Headrest Hardware	



Accessories

CHOOSE ONE OF THE FOLLOWING:

Part Number	Description	Price AUD
<input type="checkbox"/> I107499-99-0	Universal Headrest Adapter <i>Please Note: Recline is limited to 160° or less (depending on backrest height) when ordering a Universal Headrest Adapter on a 23" or 24" tall Corpus backrest.</i>	

OPTIONAL ITEMS:

Part Number	Description	Price AUD
<input type="checkbox"/> I106104-99-0	UT Hip Suppt (B) w/ Mesh Cover - High (Pair) <i>4.5"W x 3.5"T Pad w/ 4" Tall Mounting Bracket. Must also select PN I10550.</i>	
<input type="checkbox"/> I106103-99-0	UT Hip Suppt (A) w/ Mesh Cover - High (Pair) <i>6"W x 5.5"T Pad w/ 4" Tall Mounting Bracket. Must also select PN I10550.</i>	
<input type="checkbox"/> I10550	Adj Removable - Thigh Support Hardware (Pair)	
<input type="checkbox"/> I104575-99-0	UT Trunk Suppt (B) for MX Back - 3.5"T x 4.5"D <i>Must also select the Swing Away Trunk Support Hardware (I10551).</i>	
<input type="checkbox"/> I104574-99-0	UT Trunk Suppt (F) for MX Back - 3.5"T x 8.5"D <i>Must also select the Swing Away Trunk Support Hardware (I10551).</i>	
<input type="checkbox"/> I10551	Adj Removable - Trunk Support Hardware (Pair)	
<input type="checkbox"/> I104576-99-0	Push Handles - K450 <i>Mounted off the chassis. Height and angle adjustable. Must be selected if ordering bag clips or an attendant control or if mounting the ICS Switchbox on the backrest.</i>	
<input type="checkbox"/> I106210-99-0	RAM Self-Leveling Cup Holder	
<input type="checkbox"/> I106806-99-0	RAM X-Grip Phone Holder	
<input type="checkbox"/> I107100-99-0	RAM X-Grip 10" Tablet Holder w/Swing Away Mount <i>Compatible with most 10" screen tablets with or without a protective case. Includes adjustable swing away mount (I103476-99-0). Note: Customization may be required if ordered along with an Omni.</i>	
<input type="checkbox"/> I107222-99-0	USB Charger <i>Provides up to 1.5A (at 5V) to charge most smartphones and tablets. The charger features an isolated circuit design for improved safety; and by pulling power from the ICS, prevents battery drainage when the chair is powered off. Includes mounting hardware and ICS connector. Please select armrest for mounting, keeping in mind that charger may need to be relocated. Mount on left armrest: _____ Mount on right armrest: _____</i>	
<input type="checkbox"/> I105805-99-0	Medical Necessities Bag Hooks <i>Please Note: Recline is limited to 165° or less (depending on backrest height) when ordering with 23" or 24" tall Ergo backrest. Custom required if ordered with a ROHO AGILITY backrest.</i>	
<input type="checkbox"/> I100598-99-0	Chest Harness <i>Chest harness kit includes a standard positioning belt.</i>	
<input type="checkbox"/> ISH212	BP Shoulder Harness Medium, Front Pull <i>This item will be sent as a part with the wheelchair. It will not be installed.</i>	
<input type="checkbox"/> I100578-99-0	Small Angle Adjustable Tray <i>Must also select PN I10556.</i>	
<input type="checkbox"/> I106863-99-0	Angle Adjustable UniTrack Tray <i>Side mounted tray that attaches to the UniTrack Accessory Rail. Dimensions are: 19.5"W x 13"D. Can only be ordered with chairs that offer the UniTrack System. Must also select PN I10556.</i>	
<input type="checkbox"/> I10556	Adj Removable - Tray Hardware (1 Each)	

Alternative Drive Systems

ALTERNATIVE PROPORTIONAL JOYSTICK SYSTEMS:

Part Number	Description	Price AUD
<input type="checkbox"/> I10751	Mini-Joystick for Omni - (9-pin connector) <i>Operating Force = 13 g Neutral to Full Deflection = 5mm</i>	



Alternative Drive Systems

ALTERNATIVE PROPORTIONAL JOYSTICK SYSTEMS:

Part Number	Description	Price AUD
<input type="checkbox"/> I103893-99-0	Compact Joystick for Omni - (9-pin connector) <i>Operating Force = 220 g Neutral to Full Deflection = 28mm Joystick includes two 3.5mm mono jacks for Mode and On/Off. Includes Standard Joystick Knob and 1.75" (42mm) Foam Ball Knob.</i>	
<input type="checkbox"/> I103894-99-0	Compact Joystick LITE for Omni - (9-pin connector) <i>Operating Force = 50 g Neutral to Full Deflection = 17mm Joystick includes two 3.5mm mono jacks for Mode and On/Off. Includes Standard Joystick Knob and 1.75" (42mm) Foam Ball Knob.</i>	

ALT PROPORTIONAL JOYSTICK MOUNTING HARDWARE:

Part Number	Description
	NOTE: <i>Please see the PERMOfix Mounting Solutions, within the Alt Drive Accessories section, for additional mounting options.</i>
<input type="checkbox"/> I104690-99-0	Alt. JS Swing Away Mount Kit - Armrest <i>Mount in front of Left Arm: _____ Mount in front of Right Arm: _____</i>
<input type="checkbox"/> I107135-99-0	Alt JS/Omni Combo Swing Away Mount Kit - Armrest <i>Used to mount Omni and Alternative Proportional Joystick to the same armrest. Includes Swing Away Omni/Tablet Mount, with a longer mounting rod, two PERMOfix modules, one 3" rod and one 4" rod.</i> <i>Mount in front of Left Arm: _____ Mount in front of Right Arm: _____</i>
<input type="checkbox"/> I104695-99-0	Alt. JS Adjustable Mount Kit - Midline <i>Mount on Left Arm: _____ Mount on Right Arm: _____</i>
<input type="checkbox"/> I104143-99-0	Alt. JS Adjustable Mount Kit - Chin Boom <i>Swing-Away / Removable Chin Boom mounted on right of backrest shell. Custom charge may be applied if ordering lower backrest heights.</i>
<input type="checkbox"/> I10752	Haims Harness for Chin Control Mount, Large <i>For neck width up to 19cm.</i>
<input type="checkbox"/> I10753	Haims Harness for Chin Control Mount, Small <i>For neck width up to 14cm.</i>

ALT PROPORTIONAL JOYSTICK ACCESSORIES:

Part Number	Description	Price AUD
<input type="checkbox"/> I104904-99-0	Dual Switches for Compact Joystick - Fixed Mount <i>Includes two pushbutton switches mounted on rigid Switches can be used for Mode & On/Off. Switch operation force = 130 g.</i> <i>Mount at front of Joystick: _____ Mount at rear of Joystick: _____</i> <i>Mount at left of Joystick: _____ Mount at right of Joystick: _____</i>	
<input type="checkbox"/> I104641-99-0	Dual Switches for Compact Joystick - Gooseneck Mnt <i>Includes two pushbutton switches mounted on adjustable gooseneck arms. Switches are positioned on Left and Right sides of Compact Joystick. Switches can be used for Mode & On/Off. Switch operation force = 130 g</i>	

SIP & PUFF SYSTEMS:

Part Number	Description	Price AUD
<input type="checkbox"/> I10573	Sip & Puff System for Omni <i>Requires Omni Control Module. Includes breath tube kit.</i>	
<input type="checkbox"/> I10574	Sip & Puff System with Fiber Optic Opti-Stop <i>Requires Omni Control Module. Includes breath tube kit.</i>	
<input type="checkbox"/> I106807-99-0	Sip & Puff Kit for Total Control Head Array <i>For use with Total Control Head Array. Requires Omni Control Module. Includes breath tube kit.</i>	
<input type="checkbox"/> I10762	Adjustable Mounting Hrdw - Sip & Puff <i>Required when selecting PN I10573, I10574 or I106807-99-0.</i>	

PROPORTIONAL FOOT CONTROL:

Part Number	Description	Price AUD
-------------	-------------	-----------



Alternative Drive Systems

PROPORTIONAL FOOT CONTROL:

Part Number	Description	Price AUD
<input type="checkbox"/> I10767	Proportional Foot Control Mount on Left: _____ Mount on Right: _____	
<input type="checkbox"/> I10766	Adjustable Mounting Hrdw - Foot Control Required when selecting PN I10767.	

Total Control Head Array

CONFIGURE HEAD ARRAY:

Part Number	Description	Price AUD
<input type="checkbox"/> I106805-99-0	Permobil Head Array - 3-Switch System Features an occipital pad sensor and two stainless steel arms configured with proximity sensors. The arms are mounted to the upper attachment points on the head array. Longer arms are also included to extend the sensor adjustment range. ESP recommended. Choose Omni Control Module (ID51154IR) separately.	
<input type="checkbox"/> I11072	Program Head Array as a 5-Switch System Select additional arms/sensors below and indicate the intended use of each sensor in the Notes & Comments section. Choose Omni Control Module (ID51154IR) separately.	No Charge

CHOOSE OCCIPITAL PAD SIZE (REQUIRED):

<input type="checkbox"/> I105920-99-0	Small Occipital Pad for Head Array - 4.25"W x 3"H	No Charge
<input type="checkbox"/> I105921-99-0	Large Occipital Pad for Head Array - 6"W x 4"H	No Charge

CHOOSE OCCIPITAL PAD MOUNT (REQUIRED):

<input type="checkbox"/> I10847	Fixed Mount for Occipital Pad - Head Array	No Charge
<input type="checkbox"/> I105936-99-0	Adj. Swivel Mount for Head Array Occipital Pad Enables Occipital Pad to be angled up/down and left/right independent of the lateral arms.	No Charge

CHOOSE HEAD ARRAY MOUNT HARDWARE (REQUIRED):

<input type="checkbox"/> I105927-99-0	Stealth TWB480 Backrest Mounting Hdwr with Links Please select the length of the vertical stem below: Standard length (10"): _____ Short length (6"): _____	
<input type="checkbox"/> I105941-99-0	Stealth TWB480 Backrest Mtng Hdwr with Flipback	
<input type="checkbox"/> I105925-99-0	Stealth HMO475 Backr Mtng Hdwr w 17mm Ball	
<input type="checkbox"/> I106319-99-0	Stealth TWB480 Backrest Mnt Hdwr w SubOcpl Mount	
<input type="checkbox"/> I106282-99-0	Stealth 600 SubOccipital Pad - 8.5"W x 2.75"H Must select when PN I106319-99-0 is ordered.	

SELECT ADDITIONAL LEFT ARM PACKAGE (OPTIONAL):

<input type="checkbox"/> I10976	Add Left Lateral Arm w Proximity Sensor Choose Lateral Arm Adjustment Range: Short (5.5" - 10"): _____ Long (5.5" - 14"): _____	
<input type="checkbox"/> I10977	Add Left Lateral Arm w Egg Switch Choose Lateral Arm Adjustment Range: Short (5.5" - 10"): _____ Long (5.5" - 14"): _____	
<input type="checkbox"/> I10978	Add Left Lateral Arm w Pad (No Sensor or Switch) Choose Lateral Arm Adjustment Range: Short (5.5" - 10"): _____ Long (5.5" - 14"): _____	

SELECT ADDITIONAL RIGHT ARM PACKAGE (OPTIONAL):

<input type="checkbox"/> I10979	Add Right Lateral Arm w Proximity Sensor Choose Lateral Arm Adjustment Range: Short (5.5" - 10"): _____ Long (5.5" - 14"): _____	
---------------------------------	----------------------------------------------------------------------------------------------------------------------------------------	--



Total Control Head Array

CONFIGURE HEAD ARRAY:

Part Number	Description	Price AUD
<input type="checkbox"/> I10980	Add Right Lateral Arm w Egg Switch <i>Choose Lateral Arm Adjustment Range: Short (5.5" - 10"): _____ Long (5.5" - 14"): _____</i>	
<input type="checkbox"/> I10981	Add Right Lateral Arm w Pad (No Sensor or Switch) <i>Choose Lateral Arm Adjustment Range: Short (5.5" - 10"): _____ Long (5.5" - 14"): _____</i>	

Alt Drive Systems Setup

CHOOSE OMNI (REQUIRED):

Part Number	Description	Price AUD
<input type="checkbox"/> ID51154IR	Omni Control Module for R-net w/ Infrared <i>Must also select mounting hardware, either I103476-99-0 or I107135-99-0. Please specify which type of alternative drive input you will be using if not ordering it from Permobil. Input Device: _____</i>	
<input type="checkbox"/> I103476-99-0	Adjustable Swing-Away Hardware - Omni/Tablet <i>Mount in front of Left Arm: _____ Mount in front of Right Arm: _____</i>	

CHOOSE MOUNTING POSITION OF SEAT CONTROL SWITCHBOX:

Part Number	Description	Price AUD
<input type="checkbox"/> I10080	Mount Seat Control Switchbox on Backrest	No Charge

Alt Drive Accessories

PERMOFIX MOUNTING SOLUTIONS:

Part Number	Description	Price AUD
<input type="checkbox"/> I106655-99-0	PERMOfix UniTrack Fixed Mount Kit <i>Includes two base joints, two 4" rods and one 8" rod.</i>	
<input type="checkbox"/> I106656-99-0	PERMOfix UniTrack Removable Mount Kit <i>Includes quick release mount and thigh support hardware, two base joints and three different length rods.</i>	
<input type="checkbox"/> I106657-99-0	PERMOfix UniTrack Swing Away Mount Kit <i>Includes PERMOfix Swing Away module, two base joints, one 2" rod and one 4" rod.</i>	
<input type="checkbox"/> I105907-99-0	Proximity Sensor + Swivel Assm - Not Mounted <i>Only available with the Total Control Head Array. Swivel assembly allows sensor to be mounted with PERMOfix or directly to lateral arm of Total Control Head Array. Since Proximity Sensors require power, they CANNOT be used as On/Off switches.</i>	
<input type="checkbox"/> I105919-99-0	Egg Switch + Swivel Mount Assm - Not Mounted <i>Swivel assembly allows switch to be mounted with PERMOfix or directly to lateral arm of Total Control Head Array.</i>	
<input type="checkbox"/> I105912-99-0	Solid Lateral Pad+Swivel Joint - Not Mounted <i>Swivel assembly allows pad to be mounted with PERMOfix or directly to lateral arm of Total Control Head Array.</i>	
<input type="checkbox"/> I106675-99-0	Oval Plate + Swivel Assm - No Switch - Not Mounted <i>Blank plate provides mounting option for Micro Lite, Buddy Button or other switch options. Swivel assembly allows pad to be mounted with PERMOfix or directly to lateral arm of Total Control Head Array.</i>	

OPTIONAL ITEMS:

Part Number	Description	Price AUD
<input type="checkbox"/> I101769-99-0	Micro Lite Switch Qty: _____	
<input type="checkbox"/> I101767-99-0	Buddy Button Qty: _____	
<input type="checkbox"/> I106293-99-0	Egg Switch Black Qty: _____	
<input type="checkbox"/> I101771-99-0	Mini Cup Switch Qty: _____	



Alt Drive Accessories

OPTIONAL ITEMS:

Part Number	Description	Price AUD
<input type="checkbox"/> I107500-99-0	Proximity Switch Qty: _____	
<input type="checkbox"/> I107495-99-0	Mini Joystick 5 Switch OC for Omni - (9-pin)	
<input type="checkbox"/> I107498-99-0	Touch Contact with CR2032 Battery Qty: _____	
<input type="checkbox"/> I10772	Dome Switch, Backrest Mnt w Gooseneck Qty: _____ <i>Mount on Left: _____ Mount on Right: _____</i> <i>Mounts at top of backrest shell. Custom charge may be applied if ordering with short backrest heights. Kit Includes: One PERMOfix joint with Flat Plate, two PERMOfix joints with 45° swivel adjustment, one 150mm rod, one 200mm rod and one 225mm gooseneck switch with 1.2m cable and mounting clamp.</i>	
<input type="checkbox"/> I10846	Dome Switch, Backrest Swg Awy Mnt w Gsneck Qty: _____ <i>Mount on Left: _____ Mount on Right: _____</i> <i>w 90° Swing-Away Function. Mounts at top of back. Custom charge applies if ordered with short back. Kit Includes: One PERMOfix joint with Flat Plate, one PERMOfix joint with 45° swivel adjustment, two 6" rods, one 90° Swing-Away joint, one 9"</i>	
<input type="checkbox"/> I105935-99-0	3.5mm Ext. Cable for Sensor/Switch, 0.3m Long	
<input type="checkbox"/> I105931-99-0	3.5mm Ext. Cable for Sensor/Switch, 1m Long	
<input type="checkbox"/> I105932-99-0	3.5mm Ext. Cable for Sensor/Switch, 2m Long	
<input type="checkbox"/> I105933-99-0	3.5mm Ext. Cable for Sensor/Switch, 3m Long	
<input type="checkbox"/> I105934-99-0	3.5mm Y-splitter Cable for Sensor/Switch <i>Allows two Proximity Sensors/Mechanical Switches to connect to the same direction jack.</i>	

Additional Order Instructions:

Chair Order Policy:

You are highly encouraged to contact Permobil Sales with ANY questions regarding proper completion of our order forms. The order cannot be processed without the following:

- A client name or code, the client's measurements, the client's weight
- The order form must also be signed by a person authorized by your company to acknowledge the items selected. Permobil is not responsible for configuration or size discrepancies resulting from customer errors on the order form.
- Order confirmations will be provided by Permobil summarizing the items selected. Please review this carefully.
- Once the order is shipped, any changes are subject to the returns policy and may be prohibited.

Order Acknowledgement:

I, _____, am an agent of the medical equipment provider named on this order form and I have the authority to contract for the purchase of powered wheelchairs and related parts on behalf of said provider. I acknowledge that I have reviewed this order and that it is complete and accurate to the best of my knowledge. I further acknowledge that any changes to this order after submission of this order form are subject to the returns policy and may carry additional charges.

Signed: _____

Printed Name: _____

Title: _____

Date: _____

NOTE: All specifications and prices are subject to change without notice. Please note that prices displayed are only valid if ordered with the wheelchair. Any parts ordered for service or future alterations will carry different pricing and warranty.