



*** Indicates a Required Field**

Dealer Information

Order ID: _____
 *ATP: _____
 Dealer Code: _____
 *Dealer Name: _____
 Address: _____
 *City: _____
 * State/Zip: _____
 * Phone#: _____
 PO#: _____
 Email Address: _____


Client Information

Permobil recommends that the client is evaluated by a certified rehab specialist.

* First Name: _____
 * Last Name: _____
 Diagnosis: _____
 Funding Source: _____
 Client Age: _____

Client Measurements

* Weight: _____
 *Height: _____
 A) Top of Shoulders: _____
 B) Chest Depth: _____
 C) Chest Width: _____
 D) Seat Depth: _____
 E) Top of Head: _____
 F) Elbow to Hand: _____
 G) Seat Pan to Elbow: _____
 H) Hip Width: _____
 *I) Knee to Foot: _____



Three column pricing on order/quote:

Display HCPCS Codes on order/quote:

Please send order/quote to fax#: **(800) 231-3256**

Email to: sales@permobil.com

Permobil Inc.
 300 Duke Dr.
 Lebanon, TN 37090
 Tel: (800) 736-0925
 Fax: (800) 231-3256
www.permobil.com

Prices effective **April 01, 2018**

K300 PS Jr

Chair Model and Colors

CHOOSE BASE:

| Part Number | Description | Price | HCPCS Code |
|---------------------------------------|---|----------|------------|
| <input type="checkbox"/> I103573-99-0 | K300 Base PS - C300 PS0 <i>Price Includes: FWD C300 Power Base Set Up For No Power Seat Functions, Shock Absorbing Suspension System, R-net Controller, Tie-Down Hardware for Strap Systems, Support Wheels, Charger, Flat Free Tires, and Standard Positioning Belt. Must Choose Flat Free or Pneumatic Tires Below. 5.0mph, Weight Capacity = 165lbs, STFH = 16.5" (17.5" with Power Tilt).</i> | 7,889.00 | K0848 |
| <input type="checkbox"/> I10431 | K300 Base PS - C300 PS1 <i>Price Includes: FWD C300 Power Base Set Up For One Power Seat Function, Shock Absorbing Suspension System, R-net Controller, Tie Down Hardware for Strap Systems, Support Wheels, Charger, Flat Free Tires, and Standard Positioning Belt. Must Choose Flat Free or Pneumatic Tires Below. 5.0mph, Weight Capacity = 165lbs, STFH = 16.5" (17.5" with Power Tilt).</i> | 7,889.00 | K0856 |
| <input type="checkbox"/> I10520 | K300 Base PS - C300 PS2 <i>Price Includes: FWD C300 Power Base Set Up For Two Power Seat Functions, Shock Absorbing Suspension System, R-net Controller, Tie Down Hardware for Strap Systems, Support Wheels, Charger, Flat Free Tires, and Standard Positioning Belt. Must Choose Flat Free or Pneumatic Tires Below. 5.0mph, Weight Capacity = 165lbs, STFH = 16.5" (17.5" with Power Tilt).</i> | 7,889.00 | K0861 |
| <input type="checkbox"/> I10637 | 6.5mph High Speed Upgrade Package - C300/K300 <i>Price Includes: ESP Module, High Speed 6.5mph Motor Package, and 120A Power Module.</i> | 3,100.00 | |

CHOOSE DRIVE TIRES:

| Part Number | Description | Price | HCPCS Code |
|---------------------------------------|---|-----------|------------|
| <input type="checkbox"/> I103499-99-0 | Tri-Spoke Split Rim w/ Flat Free Tire 3x8" Gray | No Charge | |
| <input type="checkbox"/> I103498-99-0 | Tri-Spoke Split Rim w/ Pneumatic Tire 3x8" Gray | No Charge | |

OPTIONAL DRIVE TIRES:

| Part Number | Description | Price | HCPCS Code |
|---------------------------------|---|--------|------------|
| <input type="checkbox"/> I10782 | Tri-Spoke Split Rim w/ Snow Tire Pneumatic 3x8" <i>Snow tires are designed for use in cold weather conditions. The nylon compound in these tires remains flexible at low temperatures, but wears more easily. Therefore, the snow tires (sold as wheel assemblies) are not recommended for year-long use and must be ordered along with standard drive wheels. Not available with black tires or black rims.</i> | 470.00 | |

CHOOSE SHROUD COLOR:

| Part Number | Description | Price | HCPCS Code |
|---------------------------------------|---------------------------------|-----------|------------|
| <input type="checkbox"/> I105701-99-0 | C300/K300 Shroud - Cobalt Blue | No Charge | |
| <input type="checkbox"/> I105702-99-0 | C300/K300 Shroud - Purple Rain | No Charge | |
| <input type="checkbox"/> I105703-99-0 | C300/K300 Shroud - Onyx Black | No Charge | |
| <input type="checkbox"/> I105704-99-0 | C300/K300 Shroud - Phantom Gray | No Charge | |
| <input type="checkbox"/> I105705-99-0 | C300/K300 Shroud - Popstar Pink | No Charge | |
| <input type="checkbox"/> I105706-99-0 | C300/K300 Shroud - Volcano Red | No Charge | |

CHOOSE HUBCAP COLOR:

| Part Number | Description | Price | HCPCS Code |
|---------------------------------|---|-----------|------------|
| <input type="checkbox"/> I10871 | Match Hub Cap Color to Shroud: _____ <i>Note: Hub caps will be matched to shroud color if no selection is made.</i> Solid Color Hub Caps - Pair <i>Platinum Silver: _____ Onyx Black: _____ Cobalt Blue: _____ Phantom Grey: _____ Pop Star Pink: _____ Volcano Red: _____ Purple Rain: _____ Titan Blue: _____ Sublime Green: _____ Crystal Blue: _____ Vintage Wine: _____ Extreme Yellow: _____ Chili Pepper Red: _____</i> | No Charge | |

Base Options

MUST CHOOSE EITHER SEAT ELEVATOR OR FIXED SEAT TUBE:

| Part Number | Description | Price | HCPCS Code |
|---------------------------------------|--|-----------|------------|
| <input type="checkbox"/> I101995-99-0 | Fixed Seat Tube | No Charge | |
| <input type="checkbox"/> I106749-99-0 | Power Adjustable Seat Height - 8" Travel | 3,139.00 | E2300 |

K300 PS Jr

Base Options

MUST CHOOSE BATTERIES:

| Part Number | Description | Price | HCPCS Code |
|--|--|--------|------------|
| <input checked="" type="checkbox"/> I107032-99-0 | Batteries, Grp 34 (60Ah) Sealed Gel, x2, Installed | 960.00 | E2359 |

OPTIONAL ITEMS:

| Part Number | Description | Price | HCPCS Code |
|---------------------------------------|--|----------|------------|
| <input type="checkbox"/> I105821-99-0 | Enhanced Steering Perform Unit (ESP), R-net <i>MUST be selected if ordering any ventilator tray (I102769-99-0, I103241-99-0, or I102768-99-0). Recommended item with switched drive inputs for improved tracking.</i> | 1,519.00 | |
| <input type="checkbox"/> I105803-99-0 | LED Light Kit for C/K300 w/ R-net (New Shroud) | 1,231.00 | |
| <input type="checkbox"/> I102768-99-0 | Fixed Vent Shelf for C300 <i>If selecting this part, you MUST also pick: PN I105821-99-0 - ESP Unit. Fits an LP10 Sized Vent & a U1 Battery. Must be ordered with an expandable controller (ex. R-net Electronics).</i> | 479.00 | E1029 |
| <input type="checkbox"/> ID51113 | Programming Key for R-net - Permobil Version | 526.00 | |

Joystick Options

CHOOSE JOYSTICK TYPE:

| Part Number | Description | Price | HCPCS Code |
|---------------------------------------|---|-----------|------------|
| | PLEASE NOTE: <i>For additional input devices, please refer to the Alternative Drive Systems Section.</i> | | |
| <input type="checkbox"/> I108319-99-0 | Permobil Joystick Module <i>Designed exclusively for Permobil, this R-net color joystick features: control of up to 4 Bluetooth devices, Infrared control, programmable buttons and a programmable switch jack; all in a compact, durable Aluminum housing.</i> | 1,082.00 | |
| <input type="checkbox"/> I105246-99-0 | R-net LED Remote Joystick <i>This Joystick does NOT have a Display Screen and does NOT contain mono jacks. Limited to five profiles with this Joystick Module.</i> | No Charge | |
| <input type="checkbox"/> I105700-99-0 | R-net Remote Color JS w/ Lights, Heavy Duty Gimbal <i>Heavy-Duty modification enables gimbal to hold up to extreme operating forces. Reinforced housing is attached to a standard Color LCD R-net Joystick. 3" long joystick knob has 3/4" diameter. Not compatible with alternative joystick handles.</i> | 1,232.00 | |

CHOOSE JOYSTICK HANDLE TYPE:

| Part Number | Description | Price | HCPCS Code |
|---------------------------------------|--|-----------|------------|
| <input type="checkbox"/> I10045 | Cone Shaped Joystick Knob - Standard | No Charge | |
| <input type="checkbox"/> I100960-99-0 | Chin Cup For Joystick Knob | 127.00 | E2324 |
| <input type="checkbox"/> I100961-99-0 | "T" Handle For Joystick Knob | 179.00 | E2323 |
| <input type="checkbox"/> I100962-99-0 | "Large Ball" For Joystick Knob | 127.00 | E2323 |
| <input type="checkbox"/> I100963-99-0 | "Softball" For Joystick Knob | 167.00 | E2323 |
| <input type="checkbox"/> I100964-99-0 | "Mushroom" For Joystick Knob | 167.00 | E2323 |
| <input type="checkbox"/> I100965-99-0 | "Stick" For Joystick Knob | 167.00 | E2323 |
| <input type="checkbox"/> I105442-99-0 | 1.75" Black Foam Ball For Joystick Knob | 127.00 | E2323 |
| <input type="checkbox"/> I105846-99-0 | Bodypoint J/S Handle 3" U-Shaped w/ Flex-Shaft | 167.00 | E2323 |
| <input type="checkbox"/> I105847-99-0 | Bodypoint J/S Handle 4" U-Shaped w/ Flex-Shaft | 167.00 | E2323 |
| <input type="checkbox"/> I105848-99-0 | Bodypoint J/S Handle Dome Shaped, Rubber | 167.00 | E2323 |

CHOOSE JOYSTICK MOUNT:

| Part Number | Description | Price | HCPCS Code |
|---------------------------------------|---|-----------|------------|
| <input type="checkbox"/> I10412 | Joystick Mount Right - Fixed, R-net Only | No Charge | |
| <input type="checkbox"/> I10413 | Joystick Mount Left - Fixed, R-net Only | No Charge | |
| <input type="checkbox"/> I102537-99-0 | Joystick Mount Right - Swing Away | 374.00 | E1028 |
| <input type="checkbox"/> I102538-99-0 | Joystick Mount Left - Swing Away | 374.00 | E1028 |
| <input type="checkbox"/> I108037-99-0 | Slimline Retractable Joystick Mount - Right | 503.00 | E1028 |
| <input type="checkbox"/> I108038-99-0 | Slimline Retractable Joystick Mount - Left | 503.00 | E1028 |
| <input type="checkbox"/> I11232 | Adjustable Height Panel Bracket <i>Provides 1.77" of vertical adjustment below the armrest and 25° of angle adjustment. Compatible with all Permobil joystick mounts except for the M1 Fixed Joystick Mount.</i> | 273.00 | E1028 |

K300 PS Jr

Joystick Options

CHOOSE JOYSTICK MOUNT:

| Part Number | Description | Price | HCPCS Code |
|---------------------------------------|---|--------|------------|
| <input type="checkbox"/> I103472-99-0 | Mount Standard Joystick On Backrest, PS, PSJR, KO <i>Standard R-net Joystick is mounted to backrest. Can be chosen instead of PN ID50882. Mount on Left: _____ Mount on Right: _____</i> | 311.00 | |
| <input type="checkbox"/> I108268-99-0 | Bodypoint Tri-Lock Midline Mount Assembly <i>Includes L-shaped joystick arm and height adjustable mounting clamp. Joystick mount sold separately (select from the mounts listed below). Please indicate if you want midline mounted on left or right armrest: Mount Midline on Left Armrest: _____ Mount Midline on Right Armrest: _____</i> | 642.00 | |
| <input type="checkbox"/> I108280-99-0 | Bodypoint Tri-Lock Mount for R-net/Compact Adv JS | 267.00 | |
| <input type="checkbox"/> I108279-99-0 | Bodypoint Tri-Lock Mount for Compact Joystick | 267.00 | |
| <input type="checkbox"/> I108310-99-0 | Bodypoint Tri-Lock Mount for Omni | 302.00 | E1028 |
| <input type="checkbox"/> I10842 | PERMOFix Mount for Mini/Compact Joystick <i>Includes a shim washer that allows the PERMOFix joint to be mounted on the joystick arm. Provides an adjustable mount for a Mini Joystick or Compact Joystick.</i> | 267.00 | |

OPTIONAL ITEMS:

| Part Number | Description | Price | HCPCS Code |
|---------------------------------------|---|----------|------------|
| <input type="checkbox"/> I10411 | Expandable Controller - R-net <i>This item is required when ordering power tilt, power recline, and power elevating legrests together OR when ordering an Alternative Drive Control.</i> | 699.00 | E2377 |
| <input type="checkbox"/> I10513 | Harness for Expandable Controller <i>Must be ordered with any Expandable Controller.</i> | 620.00 | E2313 |
| <input type="checkbox"/> ID50882 | Attendant Control for R-net Control System Only | 1,372.00 | E2331 |
| <input type="checkbox"/> I107350-99-0 | R-net Bluetooth iDevice Module <i>Allows input device commands (such as joystick nudges or switch clicks) to be used to operate an iPhone, iPad or iPod. By assigning each input to an iOS function, users can make phone calls, send texts and emails and access applications. Kit includes a Bluetooth iDevice Module, an R-net Bus Cable, and two mono jack cables for direct switch inputs (switches not included). Must be ordered with a Color JSM or Omni for mode display.</i> | 962.00 | K0108 |
| <input type="checkbox"/> I107367-99-0 | R-net Bluetooth Mouse Module <i>Enables mouse control through any R-net input device. It can operate the mouse on a computer and on select Android devices (with R-net for Android app). Kit includes a Bluetooth Mouse Module, an R-net Bus Cable, two mono jack cables for L/R click (switches not included) and a USB BT Adapter for user's computer. Note: BT Adapter may not work with all computers. Must be ordered with a Color JSM or Omni for mode display.</i> | 962.00 | K0108 |
| <input type="checkbox"/> I107368-99-0 | Add a second Bluetooth Mouse Module <i>Allows R-net to simultaneously pair to a second compatible computer or Android device. Each module is accessible through a separate R-net mode. Kit includes a unique Bluetooth Mouse Module, an R-net Bus Cable, two mono jack cables for L/R click and a USB BT Adapter.</i> | 962.00 | K0108 |
| <input type="checkbox"/> I10582 | Output Module (IOM) for R-net Control System <i>Enables a user to control an external device (ex: Communication Device, ECU, Door Opener, etc.) through the drive input (JSM, Head Array, etc.). Each module can provide up to 4 switch closures. May require specialty cable or wiring. Up to three output modules can be installed on one wheelchair.</i> | 1,292.00 | K0108 |

Seating System

CHOOSE ONE OF THE FOLLOWING:

| Part Number | Description | Price | HCPCS Code |
|---------------------------------------|--|-----------|------------|
| <input type="checkbox"/> I103712-99-0 | PS Junior Seat Frame - Only <i>Weight limit 165lbs.</i> | No Charge | |
| <input type="checkbox"/> I103715-99-0 | PS Jr Seat w/ 45° Manual Tilt <i>Must be ordered with Stroller Handles (I103811-99-0). Weight limit 165lbs. Note: Cannot be ordered with an 11" seat width.</i> | 746.00 | K0108 |
| <input type="checkbox"/> I103716-99-0 | PS Junior Seat Frame with 45° Power Tilt <i>Weight limit 165lbs.</i> | 6,682.00 | E1002 |

OPTIONAL ITEMS:

| Part Number | Description | Price | HCPCS Code |
|-------------|-------------|-------|------------|
|-------------|-------------|-------|------------|

K300 PS Jr

Seating System

OPTIONAL ITEMS:

| Part Number | Description | Price | HCPCS Code |
|---------------------------------------|---|--------|------------|
| <input type="checkbox"/> I103811-99-0 | Stroller Handle for PS Jr Seating System <i>Angle adjustable. This item must be selected when ordering manual tilt (I103715-99-0).</i> | 241.00 | |

Seating System Options

CHOOSE SEAT FUNCTION CONTROL TYPE:

| Part Number | Description | Price | HCPCS Code |
|---------------------------------|--|-----------|------------|
| <input type="checkbox"/> I10418 | Single Seat Function Control Switch for R-net <i>Please choose either Push Buttons or Toggle Switches below.</i> | No Charge | E2310 |
| <input type="checkbox"/> I10419 | Single Seat Function Control Kit for R-net <i>For operation of a single seat function through the wheelchair electronics via the drive control and/or ICS Switchbox. Order with R-net Controller or remote joystick options. Please indicate below if you want push buttons or toggles on the ICS Switchbox.</i> | 1,932.00 | E2310 |
| <input type="checkbox"/> I10585 | Multiple Seat Function Control Kit for R-net <i>MUST be selected if ordering two or more seat functions. This device is used to operate two or more seat functions through the wheelchair electronics via the drive control and/or ICS switch box. Order with R-net electronics or remote joystick options. Please indicate below if you want push buttons or toggles on the ICS Switchbox.</i> | 3,657.00 | E2311 |

CHOOSE PUSH BUTTONS OR TOGGLES:

| Part Number | Description | Price | HCPCS Code |
|---------------------------------------|--|-----------|------------|
| <input type="checkbox"/> I102798-99-0 | Push Buttons on ICS Switchbox (1.4m) | No Charge | |
| <input type="checkbox"/> I102799-99-0 | Toggle Switches on ICS Switchbox (1.4m) | 327.00 | |
| <input type="checkbox"/> I105409-99-0 | Do Not Install ICS Switchbox, Send as Part | No Charge | |

CHOOSE MOUNTING POSITION OF ICS SWITCHBOX:

| Part Number | Description | Price | HCPCS Code |
|---------------------------------------|---|-----------|------------|
| <input type="checkbox"/> I10626 | Mount ICS Switchbox onto R-net Joystick Bracket <i>ICS Switchbox will be mounted to the same bracket as the R-net Joystick Module. If this item is selected, the ICS Switchbox will swing away with the Joystick Module.</i> | No Charge | |
| <input type="checkbox"/> I104651-99-0 | Mount ICS Switchbox onto Armrest <i>Mount on Left: _____ Mount on Right: _____</i> | No Charge | |
| <input type="checkbox"/> I106256-99-0 | Mount ICS Switchbox on Backrest <i>Mount on Left: _____ Mount on Right: _____</i> | 157.00 | |

OPTIONAL ITEMS:

| Part Number | Description | Price | HCPCS Code |
|---------------------------------------|---|----------|------------|
| <input type="checkbox"/> I104196-99-0 | ICS Alternative Switchbox for R-net Electronics <i>Seat function interface box that accommodates up to 8 momentary switches. Provides an additional method for accessing the seat functions. 1/8" mini-mono switches not included.</i> | 1,278.00 | K0108 |
| <input type="checkbox"/> I106743-99-0 | Remote Stop System - All Electronics <i>Consists of a Radio Frequency (RF) handheld transmitter and an RF receiver mounted on the wheelchair. Custom charges may be required when ordering with a PS/PS Jr seat, a low backrest height or a pediatric wheelchair. If ordering on a K450, push handles must also be selected.</i> | 1,264.00 | E2321 |

Backrest Options

CHOOSE ONE OF THE FOLLOWING:

| Part Number | Description | Price | HCPCS Code |
|--|---|-----------|------------|
| <input checked="" type="checkbox"/> I103717-99-0 | PS Jr Back Canes, Height & Angle Adjustable <i>Must also select PS Jr Style Armrests if ordering this backrest. Choose PS Jr Style Backrest Type (Metal Back or Omit Back) below. Backrest width should match seat width. Back height is adjustable from 16.5"T - 22.5"T, measured from the seat pan to the top of the back.</i> | No Charge | |

CHOOSE ONE OF THE FOLLOWING:

| Part Number | Description | Price | HCPCS Code |
|---------------------------------|---|-----------|------------|
| <input type="checkbox"/> I10448 | PS Jr Curved Metal Backrest <i>Backrest width will match seat width. Back height is adjustable from 16.5"T - 22.5"T, measured from the seat pan to the top of the backrest. Please note: As back height is increased, a proportional gap is created between the seat pan and the bottom of the backrest.</i> | No Charge | |
| <input type="checkbox"/> I10449 | PS Jr - Omit Backrest | No Charge | |

K300 PS Jr

Backrest Options

CHOOSE ONE OF THE FOLLOWING:

| Part Number | Description | Price | HCPCS Code |
|---------------------------------------|--|--------|------------|
| <input type="checkbox"/> I104616-99-0 | Backrest Adjustment Quick Lock, PS / PS Jr | 198.00 | |

Seat Sizing

CHOOSE ONE OF THE FOLLOWING:

| Part Number | Description | Price | HCPCS Code |
|---------------------------------|----------------------|-----------|------------|
| <input type="checkbox"/> I10432 | PS Jr Seat Width 11" | No Charge | |
| <input type="checkbox"/> I10433 | PS Jr Seat Width 12" | No Charge | |
| <input type="checkbox"/> I10434 | PS Jr Seat Width 13" | No Charge | |
| <input type="checkbox"/> I10435 | PS Jr Seat Width 14" | No Charge | |
| <input type="checkbox"/> I10436 | PS Jr Seat Width 15" | No Charge | |
| <input type="checkbox"/> I10437 | PS Jr Seat Width 16" | No Charge | |

CHOOSE ONE OF THE FOLLOWING:

| Part Number | Description | Price | HCPCS Code |
|---------------------------------|--|-----------|------------|
| <input type="checkbox"/> I10985 | SET SEAT DEPTH AT: _____ <i>Seat depth is adjustable from 5" to 18". Please round desired seat depth to nearest 0.5".</i> | No Charge | |

CHOOSE ONE OF THE FOLLOWING:

| Part Number | Description | Price | HCPCS Code |
|---------------------------------|---|--------|------------|
| <input type="checkbox"/> I10276 | Miniflex Seat Cushion Leatherette <i>Approximately 2" thick. Please select cushion size: 10"W x 10"D: _____ 12"W x 12"D: _____ 14"W x 14"D: _____ Custom sizes are available upon request.</i> | 595.00 | |
| <input type="checkbox"/> I10999 | Ergonomic Seat Cushion Leatherette <i>Please select cushion size: 16.5"W x 15"D: _____ 17.5"W x 16.5"D: _____ 17.5"W x 19.5"D: _____ Custom sizes available upon request.</i> | 578.00 | E2601 |

Armrest Options

CHOOSE ARMREST STYLE:

| Part Number | Description | Price | HCPCS Code |
|---------------------------------|--|-----------|------------|
| <input type="checkbox"/> I10824 | Flip Back, Height Adjustable Armrests <i>Mounted off the backrest canes. Includes: 45° of angle adjustability, rotational adjustment, 3" offset adjustment both medially and laterally, and height adjustment of 6" - 14" from the seat pan. In one inch increments, please indicate where you would like the armrests set. SET ARMRESTS AT (6" - 14"): _____</i> | 389.00 | E0973 |
| <input type="checkbox"/> I10825 | Flip Back, Height Adjustable Armrests - NC <i>No charge when ordering power tilt. In one inch increments, please indicate where you would like the armrests set. SET ARMRESTS AT (6" - 14"): _____</i> | No Charge | |

CHOOSE ARMREST PAD RIGHT:

| Part Number | Description | Price | HCPCS Code |
|---------------------------------------|--|-----------|------------|
| <input type="checkbox"/> I105100-99-0 | 4"x10" Corpus Armrest Pad Right - Leatherette | No Charge | |
| <input type="checkbox"/> I105104-99-0 | 4"x13" Corpus Armrest Pad Right - Leatherette | No Charge | |
| <input type="checkbox"/> I107212-99-0 | 4"x10" Corpus Armrest Pad Right - Black Fabric | No Charge | |
| <input type="checkbox"/> I106047-99-0 | 4"x13" Corpus Armrest Pad Right - Black Fabric | No Charge | |
| <input type="checkbox"/> I106486-99-0 | 4"x10" Corpus Arm w/ Stealth Gel Pad - Right | 175.00 | |
| <input type="checkbox"/> I106488-99-0 | 4"x13" Corpus Arm w/ Stealth Gel Pad - Right | 180.00 | |
| <input type="checkbox"/> I104951-99-0 | 3.5"x8" Unitrack Armrest Pad Right (C) | No Charge | |
| <input type="checkbox"/> I104950-99-0 | 3.5"x12.5" Unitrack Armrest Pad Right (G) | No Charge | |

CHOOSE ARMREST PAD LEFT:

| Part Number | Description | Price | HCPCS Code |
|---------------------------------------|---|-----------|------------|
| <input type="checkbox"/> I105098-99-0 | 4"x10" Corpus Armrest Pad Left - Leatherette | No Charge | |
| <input type="checkbox"/> I105102-99-0 | 4"x13" Corpus Armrest Pad Left - Leatherette | No Charge | |
| <input type="checkbox"/> I107213-99-0 | 4"x10" Corpus Armrest Pad Left - Black Fabric | No Charge | |
| <input type="checkbox"/> I106046-99-0 | 4"x13" Corpus Armrest Pad Left - Black Fabric | No Charge | |

K300 PS Jr

Armrest Options

CHOOSE ARMREST PAD LEFT:

| Part Number | Description | Price | HCPCS Code |
|---------------------------------------|---|-----------|------------|
| <input type="checkbox"/> I106487-99-0 | 4"x10" Corpus Arm w/ Stealth Gel Pad - Left | 175.00 | |
| <input type="checkbox"/> I106489-99-0 | 4"x13" Corpus Arm w/ Stealth Gel Pad - Left | 180.00 | |
| <input type="checkbox"/> I104952-99-0 | 3.5"x8" Unitrack Armrest Pad Left (C) | No Charge | |
| <input type="checkbox"/> I104949-99-0 | 3.5"x12.5" Unitrack Armrest Pad Left (G) | No Charge | |

OPTIONAL ITEMS:

| Part Number | Description | Price | HCPCS Code |
|---------------------------------------|--|--------|------------|
| <input type="checkbox"/> I107339-99-0 | Armrest Pouch - Right | 244.00 | |
| <input type="checkbox"/> I107338-99-0 | Armrest Pouch - Left | 244.00 | |
| <input type="checkbox"/> I107840-99-0 | UT Armrest Side Support - Right Outer <i>Includes 2"T x 6.5"W tapered pad.</i> | 352.00 | |
| <input type="checkbox"/> I107863-99-0 | UT Armrest Side Support - Left Outer <i>Includes 2"T x 6.5"W tapered pad.</i> | 352.00 | |
| <input type="checkbox"/> I107864-99-0 | UT Armrest Side Support - Right Inner <i>Includes 2"T x 6.5"W tapered pad. May require an increase to the width between the armrests to ensure armrests will flip-up.</i> | 352.00 | |
| <input type="checkbox"/> I107865-99-0 | UT Armrest Side Support - Left Inner <i>Includes 2"T x 6.5"W tapered pad. May require an increase to the width between the armrests to ensure armrests will flip-up.</i> | 352.00 | |
| <input type="checkbox"/> I107839-99-0 | UT Elbow Support - Right <i>Must also select PN I10881.</i> | 174.00 | K0108 |
| <input type="checkbox"/> I107838-99-0 | UT Elbow Support - Left <i>Must also select PN I10882.</i> | 174.00 | K0108 |
| <input type="checkbox"/> I10881 | Adjustable Elbow Support Hardware - Right | 256.00 | |
| <input type="checkbox"/> I10882 | Adjustable Elbow Support Hardware - Left | 256.00 | |

Legrest Options

CHOOSE ONE OF THE FOLLOWING:

| Part Number | Description | Price | HCPCS Code |
|---------------------------------|---|-----------|------------|
| <input type="checkbox"/> I10586 | 70° Swing Away Legrests w Multi Axis Adj Footplate <i>Includes 5"W x 8"D Multi Axis Footplate. Range of adjustment from the seat pan to the top of each footplate is: 6 - 17".</i> | No Charge | |
| <input type="checkbox"/> I10587 | 90° Swing Away Legrests w Multi Axis Adj Footplate <i>Includes 5"W x 8"D Multi Axis Footplate. Range of adjustment from the seat pan to the top of each footplate is: 6 - 17".</i> | No Charge | |
| <input type="checkbox"/> I10823 | Aeron Legrest Assembly - PS Jr - Narrow 2pc Footpl <i>90° Legrest Hanger with Manual Adjustment Legrest Elevation. Price includes 6"W x 9"D Corpus Style, Angle Adjustable, Two Piece Footplates. Range of adjustment from the seat pan to the top of the footplate is: 12 - 17" (with calf pads) and 6.5 - 17" (without calf pads).</i> | 368.00 | |
| <input type="checkbox"/> I11052 | Aeron Legrest Assembly - PS Jr - Narrow 1pc Footpl <i>90° Legrest Hanger with Manual Adjustment Legrest Elevation. Price includes 12.5"W x 10"D, Angle Adjustable, One Piece Footplate. Range of adjustment from the seat pan to the top of the footplate is: 12 - 17" (with calf pads) and 6.5 - 17" (without calf pads).</i> | 368.00 | |

OPTIONAL ITEMS:

| Part Number | Description | Price | HCPCS Code |
|---------------------------------------|---|-----------|------------|
| <input type="checkbox"/> I104684-99-0 | UT Calf Support Kit (B) w Mesh Cover & Adj Hardware <i>3.5"W x 4.5"T Pad. Recommended pad to be used with lower leg lengths of 12" or less. Bracket designed for use with Aeron Center Mount Legrest Assembly. Cannot be ordered with Swing Away Legs (I10586 & I10587).</i> | No Charge | |
| <input type="checkbox"/> I104686-99-0 | UT Calf Support Kit (J) w Mesh Cover & Adj Hardware <i>5"W x 6.5"T Pad. Bracket designed for use with Aeron Center Mount Legrest Assembly. Cannot be ordered with Swing Away Legs (I10586 & I10587).</i> | No Charge | |
| <input type="checkbox"/> I107503-99-0 | BP Ankle Huggers, Med Side Release, Adj (Pair) <i>This item will be sent as a part with the wheelchair. It will not be installed.</i> | 288.00 | |

Accessories

K300 PS Jr

Accessories

MUST CHOOSE POSITIONING BELT:

| Part Number | Description | Price | HCPCS Code |
|---------------------------------------|---|-----------|------------|
| <input type="checkbox"/> I104518-99-0 | Bodypoint Non-Padded Hip Belt - Medium | No Charge | |
| <input type="checkbox"/> I105849-99-0 | Bodypoint Padded Hip Belt - Extra Small | 153.00 | |
| <input type="checkbox"/> I104520-99-0 | Bodypoint Padded Hip Belt - Medium | 153.00 | |
| <input type="checkbox"/> I100519-99-0 | Positioning Belt Retractable | 304.00 | |

HEADREST OPTIONS:

| Part Number | Description | Price | HCPCS Code |
|---------------------------------------|--|--------|------------|
| <input type="checkbox"/> I100749-99-0 | MF Headrest for Koala/PS Jr - Leatherette <i>Must also select PN I10552.</i> | 342.00 | E0955 |
| <input type="checkbox"/> I10552 | Adj/Removable Headrest Hardware | 280.00 | E1028 |
| <input type="checkbox"/> I107499-99-0 | Universal Headrest Adapter <i>Please Note: Recline is limited to 160° or less (depending on backrest height) when ordering a Universal Headrest Adapter on a 23" or 24" tall Corpus backrest.</i> | 276.00 | |

OPTIONAL ITEMS:

| Part Number | Description | Price | HCPCS Code |
|---------------------------------------|--|--------|------------|
| <input type="checkbox"/> I103077-99-0 | UT Hip Suppt (B) w/ Mesh Cover - Low (Pair) <i>4.5"W x 3.5"T Pad w/ 2" Tall Mounting Bracket. Must also select PN I10550.</i> | 288.00 | E0956 |
| <input type="checkbox"/> I103079-99-0 | UT Hip Suppt (A) w/ Mesh Cover - Low (Pair) <i>6"W x 5"T Pad w/ 2" Tall Mounting Bracket. Must also select PN I10550.</i> | 288.00 | E0956 |
| <input type="checkbox"/> I103095-99-0 | UT Thigh/Hip Suppt (C) w/ Mesh Cover - High (Pair) <i>8"W x 3.5"T Pad w/ 4" Tall Mounting Bracket. Must also select PN I10550.</i> | 288.00 | E0956 |
| <input type="checkbox"/> I10550 | Adj Removable - Thigh Support Hardware (Pair) | 551.00 | E1028 |
| <input type="checkbox"/> I106210-99-0 | RAM Self-Leveling Cup Holder | 85.00 | |
| <input type="checkbox"/> I106806-99-0 | RAM X-Grip Phone Holder | 85.00 | |
| <input type="checkbox"/> I107100-99-0 | RAM X-Grip 10" Tablet Holder w/Swing Away Mount <i>Compatible with most 10" screen tablets with or without a protective case. Includes adjustable swing away mount (I103476-99-0). Customization may be required if ordered with an Omni.</i> | 427.00 | |
| <input type="checkbox"/> I107222-99-0 | USB Charger <i>Provides up to 1.5A (at 5V) to charge most smartphones and tablets. The charger features an isolated circuit design for improved safety; and by pulling power from the ICS, prevents battery drainage when the chair is powered off. Includes mounting hardware and ICS connector. Please select armrest for mounting, keeping in mind that charger may need to be relocated. Mount on left armrest: _____ Mount on right armrest: _____</i> | 212.00 | |
| <input type="checkbox"/> I105805-99-0 | Medical Necessities Bag Hooks <i>Please Note: Recline is limited to 165° or less (depending on backrest height) when ordered with 23" or 24" tall Ergo backrest. Custom required if ordered with a ROHO AGILITY backrest.</i> | 101.00 | K0108 |
| <input type="checkbox"/> I102059-99-0 | Medical Necessities Bag <i>Includes hardware to attach bag to backrest. Please Note: Recline is limited to 165° or less (depending on backrest height) when ordering with a 23" or 24" tall Ergo backrest. Custom required if ordered with a ROHO AGILITY backrest.</i> | 293.00 | |
| <input type="checkbox"/> ISH212 | BP Shoulder Harness Medium, Front Pull <i>This item will be sent as a part with the wheelchair. It will not be installed.</i> | 251.00 | E0960 |
| <input type="checkbox"/> I100578-99-0 | Small Angle Adjustable Tray <i>Must also select PN I10556.</i> | 310.00 | E0950 |
| <input type="checkbox"/> I105221-99-0 | Angle Adjustable Tray - 14" W <i>Must also select PN I10556.</i> | 347.00 | E0950 |
| <input type="checkbox"/> I105222-99-0 | Angle Adjustable Tray - 16" W <i>Must also select PN I10556.</i> | 364.00 | E0950 |
| <input type="checkbox"/> I106863-99-0 | Angle Adjustable UniTrack Tray <i>Side mounted tray that attaches to the UniTrack Accessory Rail. Dimensions are: 19.5"W x 13"D. Can only be ordered with chairs that offer the UniTrack System. Must also select PN I10556.</i> | 376.00 | E0950 |
| <input type="checkbox"/> I10556 | Adj Removable - Tray Hardware (1 Each) | 280.00 | |
| <input type="checkbox"/> I103798-99-0 | After Market Seating System Adapter <i>Used to connect after market seating systems to the UniTrack Accessory Rail on the PS Jr Seat Frame. Includes Mounting Hardware.</i> | 99.00 | K0108 |

K300 PS Jr

Stealth Accessories

OPTIONAL ITEMS:

| Part Number | Description | Price | HCPCS Code |
|---------------------------------------|--|--------|------------|
| <input type="checkbox"/> I106222-99-0 | Stealth 3" x 5" Lateral Supports - PS Jr <i>Must also select I10723. Includes THIN Stealth hardware and specialized mounting bracket for Flip Back Armrest assembly. Only compatible with Flip Back Armrest hardware.</i> | 427.00 | E0956 |
| <input type="checkbox"/> I10723 | Adjustable Stealth Lateral Hardware (Pair) | 579.00 | E1028 |

Alternative Drive Systems

ALTERNATIVE PROPORTIONAL JOYSTICK SYSTEMS:

| Part Number | Description | Price | HCPCS Code |
|---------------------------------------|---|----------|------------|
| <input type="checkbox"/> I10751 | Mini-Joystick for Omni - (9-pin connector) <i>Operating Force = 0.5 oz (13 g) Neutral to Full Deflection = 0.2" (5mm)</i> | 3,314.00 | E2312 |
| <input type="checkbox"/> I103893-99-0 | Compact Joystick for Omni - (9-pin connector) <i>Operating Force = 7.8 oz (220 g) Neutral to Full Deflection = 1.1" (28mm) Joystick includes two 1/8" mono jacks for Mode and On/Off. Includes Standard Joystick Knob and 1.75" (42mm) Foam Ball Knob.</i> | 1,621.00 | E2373 |
| <input type="checkbox"/> I103894-99-0 | Compact Joystick LITE for Omni - (9-pin connector) <i>Operating Force = 1.8 oz (50 g) Neutral to Full Deflection = 0.7" (17mm) Joystick includes two 1/8" mono jacks for Mode and On/Off. Includes Standard Joystick Knob and 1.75" (42mm) Foam Ball Knob.</i> | 1,834.00 | E2373 |

ALT PROPORTIONAL JOYSTICK MOUNTING HARDWARE:

| Part Number | Description | Price | HCPCS Code |
|--|---|--------|------------|
| NOTE: <i>Please see the PERMOfix Mounting Solutions, within the Alt Drive Accessories section, for additional mounting options.</i> | | | |
| <input type="checkbox"/> I104690-99-0 | Alt JS Swing Away Mount Kit - Armrest <i>Mount in front of Left Arm: _____ Mount in front of Right Arm: _____</i> | 284.00 | E1028 |
| <input type="checkbox"/> I107135-99-0 | Alt JS / Omni Combo Swing Away Mount Kit - Armrest <i>Used to mount Omni and Alternative Proportional Joystick to the same armrest. Includes Swing Away Omni/Tablet Mount, with a longer mounting rod, two PERMOfix modules, one 3" rod and one 4" rod. Mount in front of Left Arm: _____ Mount in front of Right Arm: _____</i> | 537.00 | E1028 |
| <input type="checkbox"/> I104695-99-0 | Alt JS Adjustable Mount Kit - Midline <i>Mount on Left Arm: _____ Mount on Right Arm: _____</i> | 833.00 | E1028 |
| <input type="checkbox"/> I104143-99-0 | Alt JS Adjustable Mount Kit - Chin Boom <i>Swing-Away / Removable Chin Boom mounted on right of backrest shell. Custom charge may be applied if ordering lower backrest heights.</i> | 885.00 | E1028 |
| <input type="checkbox"/> I10752 | Haims Harness for Chin Control Mount, Large <i>For neck width up to 7½".</i> | 629.00 | E1028 |
| <input type="checkbox"/> I10753 | Haims Harness for Chin Control Mount, Small <i>For neck width up to 5½".</i> | 629.00 | E1028 |

ALT PROPORTIONAL JOYSTICK ACCESSORIES:

| Part Number | Description | Price | HCPCS Code |
|---------------------------------------|---|--------|------------|
| <input type="checkbox"/> I104904-99-0 | Dual Switches for Compact Joystick - Fixed Mount <i>Includes two pushbutton switches mounted on rigid arms. Switches can be used for Mode & On/Off. Switch operation force = 4.5oz (130 g) Mount at front of Joystick: _____ Mount at rear of Joystick: _____ Mount at left of Joystick: _____ Mount at right of Joystick: _____</i> | 263.00 | |
| <input type="checkbox"/> I104641-99-0 | Dual Switches for Compact Joystick -Gooseneck Mount <i>Includes two pushbutton switches mounted on adjustable gooseneck arms. Switches are positioned on Left and Right sides of Compact Joystick. Switches can be used for Mode & On/Off. Switch operation force = 4.5oz (130 g)</i> | 474.00 | |

SIP & PUFF SYSTEMS:

| Part Number | Description | Price | HCPCS Code |
|---------------------------------|--|----------|-------------|
| <input type="checkbox"/> I10573 | Sip & Puff System for Omni <i>Requires Omni Control Module. Includes breath tube kit.</i> | 1,773.00 | E2325,E2326 |
| <input type="checkbox"/> I10574 | Sip & Puff System with Fiber Optic Opti-Stop <i>Requires Omni Control Module. Includes breath tube kit.</i> | 2,360.00 | E2325,E2326 |

K300 PS Jr

Alternative Drive Systems

SIP & PUFF SYSTEMS:

| Part Number | Description | Price | HCPCS Code |
|---------------------------------------|--|--------|------------|
| <input type="checkbox"/> I106807-99-0 | Sip & Puff Kit for Total Control Head Array <i>For use with Total Control Head Array. Requires Omni Control Module. Includes breath tube kit.</i> | 914.00 | E2325 |
| <input type="checkbox"/> I10762 | Adjustable Mounting Hrdw - Sip & Puff <i>Required when selecting PN I10573, I10574 or I106807-99-0.</i> | 284.00 | |

PROPORTIONAL FOOT CONTROL:

| Part Number | Description | Price | HCPCS Code |
|---------------------------------|--|----------|------------|
| <input type="checkbox"/> I10767 | Proportional Foot Control <i>Mount on Left: _____ Mount on Right: _____</i> | 3,379.00 | |
| <input type="checkbox"/> I10766 | Adjustable Mounting Hrdw - Foot Control <i>Required when selecting PN I10767.</i> | 284.00 | |

Total Control Head Array

CONFIGURE HEAD ARRAY:

| Part Number | Description | Price | HCPCS Code |
|---------------------------------------|--|-----------|------------|
| <input type="checkbox"/> I106805-99-0 | Permobil Head Array - 3-Switch System <i>Features an occipital pad sensor and two stainless steel arms configured with proximity sensors. The arms are mounted to the upper attachment points on the head array. Longer arms are also included to extend the sensor adjustment range. ESP recommended. Choose Omni Control Module (ID51154IR) separately.</i> | 4,893.00 | E2330 |
| <input type="checkbox"/> I11072 | Program Head Array as a 5-Switch System <i>Select additional arms/sensors below and indicate the intended use of each sensor in the Notes & Comments section. Choose Omni Control Module (ID51154IR) separately.</i> | No Charge | |

CHOOSE OCCIPITAL PAD SIZE (REQUIRED):

| | | | |
|---------------------------------------|---|-----------|--|
| <input type="checkbox"/> I105920-99-0 | Small Occipital Pad for Head Array - 4.25"W x 3"H | No Charge | |
| <input type="checkbox"/> I105921-99-0 | Large Occipital Pad for Head Array - 6"W x 4"H | No Charge | |

CHOOSE OCCIPITAL MOUNT (REQUIRED):

| | | | |
|---------------------------------------|---|-----------|-------|
| <input type="checkbox"/> I10847 | Fixed Mount for Head Array Occipital Pad | No Charge | |
| <input type="checkbox"/> I105936-99-0 | Adj. Swivel Mount for Head Array Occipital Pad <i>Enables Occipital Pad to be angled up/down and left/right independent of the lateral arms.</i> | 439.00 | K0108 |

CHOOSE HEAD ARRAY MOUNT HARDWARE (REQUIRED):

| | | | |
|---------------------------------------|---|--------|-------|
| <input type="checkbox"/> I105927-99-0 | Stealth TWB480 Backrest Mounting Hdwr with Links <i>Please select the length of the vertical stem below: Standard length (10"): _____ Short length (6"): _____</i> | 293.00 | E1028 |
| <input type="checkbox"/> I105941-99-0 | Stealth TWB480 Backrest Mnt Hdwr with Flipback | 473.00 | E1028 |
| <input type="checkbox"/> I105925-99-0 | Stealth HMO475 MultiAxis Backrest Mounting Hdwr | 260.00 | E1028 |
| <input type="checkbox"/> I106319-99-0 | Stealth TWB480 Backrest Mnt Hdwr w SubOcptl Mount | 508.00 | E1028 |
| <input type="checkbox"/> I106282-99-0 | Stealth 600 SubOccipital Pad - 8.5" W x 2.75" H <i>Must select when PN I106319-99-0 is ordered.</i> | 104.00 | K0108 |

SELECT ADDITIONAL LEFT ARM PACKAGE (OPTIONAL):

| | | | |
|---------------------------------|--|--------|--|
| <input type="checkbox"/> I10976 | Add Left Lateral Arm w Proximity Sensor <i>Choose Lateral Arm Adjustment Range: Short (5.5" - 10"): _____ Long (5.5" - 14"): _____</i> | 630.00 | |
| <input type="checkbox"/> I10977 | Add Left Lateral Arm w Egg Switch <i>Choose Lateral Arm Adjustment Range: Short (5.5" - 10"): _____ Long (5.5" - 14"): _____</i> | 515.00 | |
| <input type="checkbox"/> I10978 | Add Left Lateral Arm w Pad (No Sensor or Switch) <i>Choose Lateral Arm Adjustment Range: Short (5.5" - 10"): _____ Long (5.5" - 14"): _____</i> | 420.00 | |

SELECT ADDITIONAL RIGHT ARM PACKAGE (OPTIONAL):

K300 PS Jr

Total Control Head Array

CONFIGURE HEAD ARRAY:

| Part Number | Description | Price | HCPCS Code |
|---------------------------------|---|--------|------------|
| <input type="checkbox"/> I10979 | Add Right Lateral Arm w Proximity Sensor <i>Choose Lateral Arm Adjustment Range: Short (5.5" - 10"): _____ Long (5.5" - 14"): _____</i> | 630.00 | |
| <input type="checkbox"/> I10980 | Add Right Lateral Arm w Egg Switch <i>Choose Lateral Arm Adjustment Range: Short (5.5" - 10"): _____ Long (5.5" - 14"): _____</i> | 515.00 | |
| <input type="checkbox"/> I10981 | Add Right Lateral Arm w Pad (No Sensor or Switch) <i>Choose Lateral Arm Adjustment Range: Short (5.5" - 10"): _____ Long (5.5" - 14"): _____</i> | 420.00 | |

Alt Drive Systems Setup

CHOOSE OMNI (REQUIRED):

| Part Number | Description | Price | HCPCS Code |
|---------------------------------------|---|----------|------------|
| <input type="checkbox"/> ID51154IR | Omni Control Module for R-net w/ Infrared <i>Must also select mounting hardware, either I103476-99-0 or I107135-99-0. Please specify which type of alternative drive input you will be using if not ordering it from Permobil. Input Device: _____</i> | 1,843.00 | |
| <input type="checkbox"/> I103476-99-0 | Adjustable Swing-Away Hardware - Omni/Tablet <i>Mount in front of Left Arm: _____ Mount in front of Right Arm: _____</i> | 284.00 | |

CHOOSE MOUNTING POSITION OF SEAT CONTROL SWITCHBOX:

| Part Number | Description | Price | HCPCS Code |
|---------------------------------|--|-----------|------------|
| <input type="checkbox"/> I10080 | Mount Seat Control Switchbox on Backrest | No Charge | |

Alt Drive Accessories

PERMOFIX MOUNTING SOLUTIONS:

| Part Number | Description | Price | HCPCS Code |
|---------------------------------------|---|--------|------------|
| <input type="checkbox"/> I106655-99-0 | PERMOfix UniTrack Fixed Mount Kit <i>Includes two base joints, two 4" rods and one 8" rod.</i> | 286.00 | |
| <input type="checkbox"/> I106656-99-0 | PERMOfix UniTrack Removable Mount Kit <i>Includes quick release mount and thigh support hardware, two base joints and three different length rods.</i> | 502.00 | |
| <input type="checkbox"/> I106657-99-0 | PERMOfix UniTrack Swing Away Mount Kit <i>Includes PERMOfix Swing Away module, two base joints, one 2" rod and one 4" rod.</i> | 719.00 | |
| <input type="checkbox"/> I105907-99-0 | Proximity Sensor + Swivel Assembly <i>Only available with the Total Control Head Array. Swivel assembly allows sensor to be mounted with PERMOfix or directly to lateral arm of Total Control Head Array. Since Proximity Sensors require power, they CANNOT be used as On/Off switches.</i> | 363.00 | E1028 |
| <input type="checkbox"/> I105919-99-0 | Egg Switch + Swivel Mount Assembly <i>Swivel assembly allows switch to be mounted with PERMOfix or directly to lateral arm of Total Control Head Array.</i> | 245.00 | |
| <input type="checkbox"/> I105912-99-0 | Solid Oval Pad + Swivel Assembly <i>Swivel assembly allows pad to be mounted with PERMOfix or directly to lateral arm of Total Control Head Array.</i> | 148.00 | |
| <input type="checkbox"/> I106675-99-0 | Oval Plate + Swivel Assembly - No Switch <i>Blank plate provides mounting option for Micro Lite, Buddy Button or other switch options. Swivel assembly allows pad to be mounted with PERMOfix or directly to lateral arm of Total Control Head Array.</i> | 130.00 | |

OPTIONAL ITEMS:

| Part Number | Description | Price | HCPCS Code |
|---------------------------------------|--|--------|------------|
| <input type="checkbox"/> I101769-99-0 | Micro Lite Switch Qty: _____ | 198.00 | |
| <input type="checkbox"/> I101767-99-0 | Buddy Button Qty: _____ | 167.00 | |
| <input type="checkbox"/> I106293-99-0 | Egg Switch Black Qty: _____ | 167.00 | |
| <input type="checkbox"/> I101771-99-0 | Mini Cup Switch Qty: _____ | 167.00 | |
| <input type="checkbox"/> I107500-99-0 | Proximity Switch Qty: _____ | 662.00 | |
| <input type="checkbox"/> I107495-99-0 | Mini Joystick 5 Switch OC for Omni - (9-pin) | 718.00 | |
| <input type="checkbox"/> I107498-99-0 | Touch Contact with CR2032 Battery Qty: _____ | 408.00 | |

K300 PS Jr

Alt Drive Accessories

OPTIONAL ITEMS:

| Part Number | Description | Price | HCPCS Code |
|---------------------------------------|---|-------|------------|
| <input type="checkbox"/> I105935-99-0 | 3.5mm Ext. Cable for Sensor/Switch, 0.3m Long | 33.00 | |
| <input type="checkbox"/> I105931-99-0 | 3.5mm Ext. Cable for Sensor/Switch, 1m Long | 36.00 | |
| <input type="checkbox"/> I105932-99-0 | 3.5mm Ext. Cable for Sensor/Switch, 2m Long | 42.00 | |
| <input type="checkbox"/> I105933-99-0 | 3.5mm Ext. Cable for Sensor/Switch, 3m Long | 46.00 | |
| <input type="checkbox"/> I105934-99-0 | 3.5mm Y-splitter Cable for Sensor/Switch | 36.00 | |

Allows two Proximity Sensors/Mechanical Switches to connect to the same direction jack.

Additional Order Instructions (for Permobil):

Notes & Comments (not for Permobil):

Chair Order Policy:

You are highly encouraged to contact Permobil Sales with ANY questions regarding proper completion of our order forms at 1-800-736-0925. The order cannot be processed without the following:

- A client name or code, the client's measurements, the client's weight
- The order form must also be signed by a person authorized by your company to acknowledge the items selected. Permobil, Inc. is not responsible for configuration or size discrepancies resulting from customer errors on the order form.
- Order confirmations will be provided by Permobil summarizing the items selected. Please review this carefully.
- Once the order is shipped, any changes are subject to the returns policy and may be prohibited.

Order Acknowledgement:

I, _____, am an agent of the medical equipment provider named on this order form and I have the authority to contract for the purchase of powered wheelchairs and related parts on behalf of said provider. I acknowledge that I have reviewed this order and that it is complete and accurate to the best of my knowledge. I further acknowledge that any changes to this order after submission of this order form are subject to the returns policy and may carry additional charges.

Signed: _____

Printed Name: _____

Title: _____

Date: _____

NOTE: All specifications and prices are subject to change without notice. Please note that prices displayed are only valid if ordered with the wheelchair. Any parts ordered for service or future alterations will carry different pricing and warranty. *The HCPCS codes provided are not intended to be billing or legal advice, rather our interpretation of the code definitions. Use of the codes does not ensure coverage or payment for the item. For coverage information, verify the policy of the appropriate payer.*