

For Internal Use Only	
Customer Order #	
Line #	
CSR	

Cushion:
MSRP Starting At: \$654.00 **HPCPS:** E2609

Back:
MSRP Starting At: \$704.00 **HPCPS:** E2617

Design-A-ROHO® Product Worksheet

Reference Guide

Cell Heights

1 2 3 4

Deleted Cells **Cut-out Cells**
Base Remains **Remove Base**

X #

Compartment Divisions

If more than one air compartment is desired, please follow these steps, using the features in your Adobe Acrobat software:

- Click the "Tools" tab
- Click open underneath the "Comment" Tool
- Click the shape icon that appears in the toolbar and select the line tool
- Use the line tool to mark your desired compartments
- Repeat step 4 for any additional compartments
- Click the 'X' next to the toolbar to close the comment tool and continue filling out the rest of the form.

Valve and Valve Hose

Type "V" in one or more of the corner boxes to indicate desired valve location. Please ensure the number of compartments matches the number of valves marked to the right or below. If more than one valve is desired in one specific corner, type "V V" in the corner box.

(Note: For any valves required outside the corner boxes, please include the location in the special instruction box)

Hook and Loop Fastener

If a hook or loop fastener is desired, please indicate the location in the checklist below.

		Front of Cushion															
		1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16
		A	B	C	D	E	F	G	H	I	J	K	L	M	N	O	P
Length	1	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
	2	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
	3	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
	4	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
	5	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
	6	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
	7	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
	8	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
	9	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
	10	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
	11	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
	12	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
	13	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
	14	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
	15	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
	16	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>

- #### Checklist
- Step 1:** Choose seat cushion or back
 Seat Cushion or Back
- Step 2:** Specify overall dimensions of cushion
 (Length x Width in inches)
- _____
- Step 3:** Specify cell heights (1, 2, 3, or 4) and ensure the proper symbol is utilized for deleted cells (X) and cut out areas (#)
- _____
- Step 4:** Specify number of valves
 1 2 3 4 5 6
- Step 5:** If there is more than one valve, ensure compartments are indicated
- Step 6:** If fastener is required, please indicate:
Specify: Hook or Loop
Placement of Fastener:
 Cushion Only Both
 Cover Only (Special) Standard Placement
- Location: _____
- Step 7:** Cover Required
 Yes or No

Specify Base Model	
HP	_____
MP	_____
LP	_____
QS HP	_____
QS MP	_____
QS LP	_____
ENH	_____
CS	_____
Other	_____

Special Instructions:

Example DESIGN-A-ROHO:
 This chart shows a 2R66 with 3in/7.6cm cells. The cushion has deleted cells in the back leaving the rubber base. The compartments are divided between rows 3 (C) and 4 (D), and the additional valve is placed in the front right corner.

		1	2	3	4	5	6
		A	B	C	D	E	F
V	1	<input type="checkbox"/> 3	<input type="checkbox"/> 3	<input type="checkbox"/> 3	<input type="checkbox"/> 3	<input type="checkbox"/> 3	<input type="checkbox"/> 3
	2	<input type="checkbox"/> 3	<input type="checkbox"/> 3	<input type="checkbox"/> 3	<input type="checkbox"/> 3	<input type="checkbox"/> 3	<input type="checkbox"/> 3
	3	<input type="checkbox"/> 3	<input type="checkbox"/> 3	<input type="checkbox"/> 3	<input type="checkbox"/> 3	<input type="checkbox"/> 3	<input type="checkbox"/> 3
	4	<input type="checkbox"/> 3	<input type="checkbox"/> 3	<input type="checkbox"/> 3	<input type="checkbox"/> 3	<input type="checkbox"/> 3	<input type="checkbox"/> 3
	5	<input type="checkbox"/> 3	<input type="checkbox"/> 3	<input type="checkbox"/> X	<input type="checkbox"/> X	<input type="checkbox"/> 3	<input type="checkbox"/> 3
	6	<input type="checkbox"/> 3	<input type="checkbox"/> 3	<input type="checkbox"/> X	<input type="checkbox"/> X	<input type="checkbox"/> 3	<input type="checkbox"/> 3

Design-A-ROHO® Product Worksheet

Provider/Distributor Name:			
ROHO Account Number:		Choose One: Quote or Order	
Address:			
City:	State:	Zip:	Country:
Phone:		Fax:	
Email:			
Person Ordering:			
Mark For (Reference Name):			
P.O. Number (Required to Place an Order):			

How to use the DESIGN-A-ROHO® Product Worksheet

Specific clients may want a ROHO cushion designed for their unique needs. This worksheet is intended for that purpose.

The grid on the front of this worksheet represents a 16 cell width by a 16 cell length cushion, but cushions needing larger dimensions may also be available. Contact ROHO Customer Support for information on limitations that may exist.

Follow these steps when ordering a ROHO® specialty cushion:

Step 1: Determine whether you need a seat cushion or a back.

Step 2: Determine the overall size of the cushion required (width by depth) e.g. 15in/38cm by 15in/38cm.

Step 3: The rounded squares represent individual cells. Specify the height of each cell by inserting the desired height. The numbers should be either 1, 2, 3, or 4 which represents 1in/2.5cm, 2.25in/5.5cm, 3.25in/8.5cm, 4.25in/10.5cm respectively. If you want certain cells deleted, indicate by placing the letter "X" in the square. When cells are deleted, all that remains is the rubber base. If you want this area removed, please indicate by placing the symbol "#" in the square on worksheet. Please note that ROHO will not make a cushion which is only one row of cells on the outer edge of the cushion or one cell standing alone. ROHO suggests that cell variances of 2 or more inches be placed on separate valves for better immersion.

Step 4: Specify the number of valves. If you would like the valve location somewhere other than in the corners, please specify the location(s) in the special instructions box.

Step 5: If you have specified more than one valve, indicate compartment divisions by using the line tool to mark the different compartment divisions.

Step 6: If you would like hook and loop fastener strips attached to the base of the cushion and/or cover, specify the type and the location on the front of this worksheet. Strips near the edge of the cushion can be both glued and sewn, but strips located toward the center of the cushion can only be glued.

Step 7: Is a cover required?

Standard Cell Area Dimensions

Width or Depth		
Inches	Number of Cells	CM
8.25	5	21
10	6	25.5
11.75	7	29.5
13.25	8	34
15	9	38
16.75	10	42.5
18.5	11	47
20.25	12	51
21.75	13	55.5
23.5	14	60
25.25	15	64

ROHO® Neoprene Cushion Sizing Chart

Fits Wheelchair Size (in. / cm)	Number of Cells
10 / 22 to 26.5	5
11 / 12 to 26.5 to 31	6
13 / 31 to 35	7
14 to 15 / 35 to 39.5	8
16 to 17 / 39.5 to 44	9
18 to 19 / 44 to 48.5	10
20 to 21 / 48.5 to 52.5	11
22 / 52.5 to 57	12
23 to 24 / 57 to 61	13
25 / 61 to 65.5	14
26 to 27 / 65.5 to 70	15