

# ROHO® AirLITE® CUSHION



ROHO has *improved* the AirLITE with you in mind!

- Lower anterior height for improved transfers
- Same great positioning, with even deeper leg troughs
- Built in protection in a non-adjustable positioning cushion
- Softer, more comfortable foam...

*...All at a reduced cost!*

---

## STABILITY WITH SAFETY COMBINED

The AirLITE Cushion improves the ability to transfer by lowering the front and middle of the cushion 3/4 inch (2 cm), while keeping its strong positioning capability by deepening the leg troughs. The lower height also contributes to a greater feeling of security in the chair.

## “BOTTOM-OUT” PROTECTION

The non-adjustable ROHO AIR FLOATATION™ air insert built into the cushion under the ischials has been enhanced to provide better bottom-out protection.

## SOFTER FOAM

Softer foam allows the cushion to conform to the individual's physiology and provide a more comfortable seating surface.

## DURABILITY AND SIMPLICITY

The 24 month warranty shows our confidence in the quality and durability of the product. The AirLITE is an off the shelf, ready to use product. No adjustment is required.

## HEAVY DUTY COVER

The cover design, with full zipper coverage, provides enhanced fluid resistance to protect the foam. At the same time, the cover is 100% PVC free, reinforcing ROHO's commitment to more environmentally friendly products.

*AirLITE's design and technology provides the comfort individuals want and the support they need.*

**ROHO**

# ROHO® AirLITE® CUSHION

## APPLICATIONS

### THE AirLITE CUSHION IS APPROPRIATE FOR INDIVIDUALS WHO MAY:

- be at risk for skin/soft tissue breakdown
- have normal sensation
- require increased stability
- require lower extremity alignment
- need a lightweight, non-adjustable seating system
- require a more stable environment for transfers

## SPECIFICATIONS

### INCLUDES:

AirLITE cushion, cover\*\* and operation manual.

### CONSTRUCTION:

Cushion: Contoured polyurethane foam; polyethylene base

Sealed Air Insert: Non-adjustable, segmented, ROHO AIR FLOATATION insert made of polyurethane

Cover: Fluid-resistant two-way stretch top and sides; non-skid bottom

### SIZES (18 Models):

- 13" x 13"/33 cm x 33 cm
- 13" x 15"/33 cm x 38 cm
- 14" x 14"/35.5 cm x 35.5 cm
- 14" x 16"/35.5 cm x 40.5 cm
- 15" x 15"/38 cm x 38 cm
- 15" x 17"/38 cm x 43 cm
- 16" x 16"/40.5 cm x 40.5 cm
- 16" x 18"/40.5 cm x 45.5 cm
- 17" x 17"/43 cm x 43 cm
- 17" x 18"/43 cm x 45.5 cm
- 17" x 20"/43 cm x 51 cm
- 18" x 16"/45.5 cm x 40.5 cm
- 18" x 18"/45.5 cm x 45.5 cm
- 18" x 20"/45.5 cm x 51 cm
- 19" x 17"/48.5 cm x 43 cm
- 19" x 19"/48.5 cm x 48.5 cm
- 20" x 17"/51 cm x 43 cm
- 20" x 18"/51 cm x 45.5 cm

### WEIGHT:

Varies by size, approximately 2.25 lbs / 1 kg (based on an 18" x 16" size AirLITE Cushion)

### WEIGHT LIMIT:

300 lbs. / 136 kg\*

### WARRANTY:

24-month limited warranty

### U.S. MEDICARE CODE:

E2605



\*Cushion must be properly sized to the individual using the AirLITE.

\*\*The AirLITE is intended to be used with the supplied cover.

## MORE PRODUCTS FROM ROHO



ROHO® MOSAIC® CUSHION



ROHO® PRODIGY® MATTRESS OVERLAY



ROHO® AGILITY® Minimum Contour Back System

100 North Florida Avenue • Belleville, IL 62221-5429 • USA • roho.com  
U.S.: 1 800 851 3449 • 1 618 277 9173 • fax 1 888 551 3449  
Outside the U.S.: 1 618 277 9150 • fax 1 618 277 6518

Products in this brochure may be covered by one or more U.S. or worldwide patents. Visit roho.com for further details. Decisions about how to code a product in the U.S. are the responsibility of the provider/supplier. Since U.S. Medicare coding is subject to change, the provider should always confirm the HCPCS code and coverage criteria as part of the client assessment process. For coverage criteria in other countries, consult your local ROHO distributor.

©2015 ROHO, Inc. The following are registered trademarks of ROHO, Inc.: AGILITY®, AIR FLOATATION™, AirLITE®, MOSAIC®, PRODIGY®, ROHO® and shape fitting technology®. ROHO, Inc. has a policy of continual product improvement and reserves the right to amend this document. The current version of this document is available on our website at roho.com.  
LITG02006, Rev. 12/15 Print 12/15

