

permobil

TILITE MANUAL WHEELCHAIRS

TR

THE BEST FROM EVERY ANGLE

The sweeping curves of the titanium dual-tube frame are the perfect fusion of striking design and driven performance. Lighter, more efficient, and more elegant than any chair anywhere. As our flagship rigid chair, the TR is simply the best of the best.





TIFIT. PERSONAL ERGONOMICS.

The difference between a TiFit® chair and all other chairs is that a TiFit chair is a wheeled prosthetic: Every frame is made to measure. A TiFit chair brings unparalleled performance and function. You will feel the difference in everything you do. TiFit...not standard.

CUSTOM CONFIGURATION.

Custom Configuration allows you to add options and components to your chair to meet your needs and to optimize your riding experience. This includes choices of wheels, handrims, tires, wheel locks, color, tattoos, and more. TiLite offers more options than any other manufacturer!

TECHNICAL SPECIFICATIONS

Frame Features	Titanium Dual-Tube Rigid TiFit Limited Adjustability
Total Chair Width*	22"
Seat Width	12" - 20"
Seat Depth	12" - 20"
Front Seat Height	15" - 21"
Rear Seat Height	12.5" - 21"
Front Frame Angle	60°, 65°, 70°, 75°, 80°, 85°, or 90°
Seat Back Type	<i>Folding Seat Back</i> Adj. Height 8.5" - 21" Adj. Angle 80° - 101° <i>Fixed Seat Back</i> Fixed Height 5" - 20" Fixed Angle 80° - 101°

Camber	0°, 2°, 4°, 6°, 8° or 12°
Center of Gravity Adjustment	5"
Front Seat Width (Seat Taper)	0", 1", 2", 3"
Ergonomic Seat	4" - 10"
Transit Options	Q-Strait Lock Down Approved
Transport Weight*	9.3 lb.
Weight Limit	265 lb.
Warranty	Five Years on Frame One Year on TiLite Components

*Weights and widths based on 16" x 16" chair with 0° camber, minimum rear wheel spacing and all standard features. Rear wheels are not included in transport weight. For reference, a pair of standard 24" rear wheels weighs 7.2 lb.


COLOR OPTIONS

HIGH GLOSS PAINT



TiLite Red Burgundy Neon Pink Super Wet Black

MATTE PAINT



Deep Blue Ice Titanium Grey Black Lava

CANDY PAINT



Ruby Red Azure Blue Electric Plum

NATURAL



Satin Titanium Polished Titanium

METALLIC PAINT



Salsa Red Metallic Tangerine Metallic Acid Green Metallic Ocean Blue Metallic

ANODIZED COLORS



Red Pepper Burnt Orange Sapphire Ultraviolet Fuchsia Black

PEARLESCENT PAINT



Midnight Blue Pearl White Pearl