




Back System

Large Back Cane Hardware




AGCCBA-LG-PR



Back System



Large Back Cane Hardware



AGCCBA-LG-PR

BACK CANE BRACKET COMPATIBILITY		Fits Back Cane Diameters
Circular		(A) >1.00 in. & ≤1.25 in. (>2.5 cm & ≤3.0 cm) (B) >1.00 in. & ≤1.25 in. (>2.5 cm & ≤3.0 cm)
Square Key*		(A) = 1.00 in. (2.5 cm)
		*(C) ≤0.125 in. (≤0.3 cm)** (D) ≤0.5 in. (≤1.3 cm)** **All measurements are approximate.

BACK CANE BRACKET COMPATIBILITY		Fits Back Cane Diameters
Circular		(A) >1.00 in. & ≤1.25 in. (>2.5 cm & ≤3.0 cm) (B) >1.00 in. & ≤1.25 in. (>2.5 cm & ≤3.0 cm)
Square Key*		(A) = 1.00 in. (2.5 cm)
		*(C) ≤0.125 in. (≤0.3 cm)** (D) ≤0.5 in. (≤1.3 cm)** **All measurements are approximate.



ROHO, Inc., 100 North Florida Avenue, Belleville, Illinois 62221-5429 U.S.A.
 U.S.A.: 800-851-3449 • Fax: 888-551-3449 • Customer Care: cc@roho.com
 Outside the U.S.A.: 618-277-9150 • Fax: 618-277-6518 • intl@roho.com • permobilus.com

ROHO, Inc., 100 North Florida Avenue, Belleville, Illinois 62221-5429 U.S.A.
 U.S.A.: 800-851-3449 • Fax: 888-551-3449 • Customer Care: cc@roho.com
 Outside the U.S.A.: 618-277-9150 • Fax: 618-277-6518 • intl@roho.com • permobilus.com