

# TiLite Canadian Seating System Order Form

Effective November 18, 2019 (New Items in Red)

Page 1 of 5

DEALER INFORMATION		RTS INFORMATION	
Dealer:	Account #:	RTS:	Phone:
PO:	Ordered By:	Email:	
Email:		SHIP TO INFORMATION	
Mark For:		Name:	Ship Via:
TiLite Sales Rep:		Address:	
<b>TiLite Information:</b>	Promo Code:	City:	Province: Code:
Ship Date:	Order #:	Phone:	Fax:

**Please Email or Fax to the Permobil Manual Customer Service Team at [tilitesales@permobil.com](mailto:tilitesales@permobil.com) or 866.586.2413 / 509.586.2413**

## INSTRUCTIONS

ROHO® Seating available with new TiLite chair orders only.

### SHIPPING

SSH3  Install ROHO® AGILITY® back (Must select seat back height range and set at measurement or Fixed Height when selected.)..... **NCO**

## BACK SUPPORTS

**ROHO® AGILITY® ACTIVE ALUMINUM MINIMUM CONTOUR** (Not available with Non-Folding Fixed Height Backs or Custom Height Rigidizer Bar.)

Shell Height  8"  10"

BRS14  Foam..... **\$895**

**ROHO® AGILITY® ACTIVE ALUMINUM MID CONTOUR** (Not available with Non-Folding Fixed Height Backs or Custom Height Rigidizer Bar.)

Shell Height  8"  10"

BRS15  Foam..... **\$895**

**ROHO® AGILITY® MINIMUM CONTOUR** (No available with Non-Folding Fixed Height Backs or Custom Height Rigidizer Bar.)

Shell Height  10"  13"  16"  18"

BRS5  Foam..... **\$785**

BRS6  Foam + Lumbar..... **\$940**

BRS7  Foam w/ Air Insert..... **\$1,095**

ACB1  AGILITY® Active Carbon Upgrade (Not available with TWIST, X Series or Bolt-on Push Handles)..... **\$230**

BRS8  Foam w/ Air Insert + Lumbar..... **\$1,095**

Quick Release (Not available with AGILITY® Carbon Upgrade) **OR**  Fixed (Not available with 18" or 20" Shell Height or Folding Chairs.)..... **NCO**

**ROHO® AGILITY® MID CONTOUR** (Not available with Non-Folding Fixed Height Backs or Custom Height Rigidizer Bar.)

BRS9  Foam w/ Air Insert..... **\$1,115**

ACB2  AGILITY® Carbon Upgrade (Not available with TWIST, X Series or Bolt-on Push Handles)..... **\$230**

BRS10  Foam w/ Air Insert + Lumbar..... **\$1,115**

Quick Release (Not available with AGILITY® Carbon Upgrade) **OR**  Fixed (Not available with 18" or 20" Shell Height or Folding Chairs.)..... **NCO**

**ROHO® AGILITY® MAX CONTOUR** (Not available with 10" or 13" Shell Height, Non-Folding Fixed Height Backs or Custom Height Rigidizer Bar.)

BRS11  Foam w/ Air Insert..... **\$1,310**

BRS12  Foam w/ Air Insert + Lumbar..... **\$1,310**

Quick Release **OR**  Fixed (16" Shell Height only. Not available with Folding Chairs.)..... **NCO**

## ROHO® AGILITY™ SYSTEM ACCESSORIES

PSB1  ROHO Privacy Shield for AGILITY® Back..... **\$70**

**COMFORT ACTA-EMBRACE™** (Not available with seat width of 13" or less.)

Shell Height  10"  12"  14"  16"..... **STD**

SBR16  Acta-Embrace Back with Compass Hardware..... **\$775**

CPSB1  Comfort Privacy Shield..... **\$55**

**COMFORT ACTA-RELIEF** (Not available with seat width of 13" or less.)

Shell Height  10"  12"  14"  16"..... **STD**

BRS13  Acta-Relief..... **\$945**

CPSB1  Comfort Privacy Shield..... **\$55**

COMFORT Install  Do Not Install  Install on TiLite Chair..... **NCO**



Permobil.com  
2701 West Court Street | Pasco, Washington 99301  
509.586.6117 Tel | 800.545.2266 Toll-Free  
509.586.2413 Fax | 866.586.2413 Toll-Free Fax

OF0027\_SeatingSystem\_0419RevA

# TiLite Canadian Seating System Order Form

Effective November 18, 2019 (New Items in Red)

Page 2 of 5

## BACK SUPPORTS

### SHORT FOLDING SEAT BACK HEIGHT RANGES AND ROHO® AGILITY™ ORIENTATION

Shell Height	Quick Release	Flipped w/ Quick Release	Fixed	Flipped w/ Fixed	Set at: (Use 1/2" increments.)
10"	<input type="checkbox"/> 16"- 18"	<input type="checkbox"/> 13"- 15"	<input type="checkbox"/> 14"- 19"	<input type="checkbox"/> 13.5"- 15.5"	Set at: _____
13"	<input type="checkbox"/> 17"- 21"	<input type="checkbox"/> 16.5"- 18.5"	<input type="checkbox"/> 17.5"- 20"	<input type="checkbox"/> 17"- 18.5"	Set at: _____
16"	<input type="checkbox"/> 19.5"- 21"	<input type="checkbox"/> 20.5"- 21.5"	<input type="checkbox"/> 19.5"- 21"	<input type="checkbox"/> 19"- 20.5"	Set at: _____
18"	<input type="checkbox"/> 22"- 24"	<input type="checkbox"/> 21.5"	N/A	N/A	Set at: _____
20"	<input type="checkbox"/> 23.5"- 25"	<input type="checkbox"/> 24.5"	N/A	N/A	Set at: _____

### MEDIUM FOLDING SEAT BACK HEIGHT RANGES AND ROHO® AGILITY™ ORIENTATION

Shell Height	Quick Release	Flipped w/ Quick Release	Fixed	Flipped w/ Fixed	Set at: (Use 1/2" increments.)
10"	<input type="checkbox"/> 18"- 22"	<input type="checkbox"/> 14.5"- 18.5"	<input type="checkbox"/> 15.5"- 21"	<input type="checkbox"/> 13"- 18"	Set at: _____
AGILITY Carbon 10"	N/A	N/A	<input type="checkbox"/> 14.5"- 18"	<input type="checkbox"/> 12"- 16"	Set at: _____
13"	<input type="checkbox"/> 18"- 24.5"	<input type="checkbox"/> 16"- 21.5"	<input type="checkbox"/> 16"- 24"	<input type="checkbox"/> 16"- 21"	Set at: _____
AGILITY Carbon 13"	N/A	N/A	<input type="checkbox"/> 17"- 19"	<input type="checkbox"/> 17"- 18"	Set at: _____
16"	<input type="checkbox"/> 19"- 25.5"	<input type="checkbox"/> 19.5"- 24.5"	<input type="checkbox"/> 20"- 25"	<input type="checkbox"/> 20"- 24"	Set at: _____
AGILITY Carbon 16"	N/A	N/A	<input type="checkbox"/> 20"- 21"	<input type="checkbox"/> 19"	Set at: _____
18"	<input type="checkbox"/> 22"- 27"	<input type="checkbox"/> 22"- 25"	N/A	N/A	Set at: _____

### FIXED SEAT BACK HEIGHT AND ROHO® AGILITY™ ORIENTATION (Available only with AGILITY Carbon Upgrade)

Shell Height	Fixed Back with Standard Orientation	Fixed Back with Flipped Orientation
10"	11"- 20"	13"- 20"
13"	14"- 20"	15"- 20"
16"	17"- 20"	19"- 20"

#### Back Orientation

Standard OR  Flipped

#### Fixed Height

10.5"  11"  11.5"  12"  12.5"  13"  13.5"  14"  
 15"  15.5"  16"  16.5"  17"  17.5"  18"  18.5"  
 19"  19.5"  20"

## CUSHIONS

### ROHO® ADJUSTABLE SKIN PROTECTION & POSITIONING (Please see page 3 for size chart.)

CSH4 QUADTRO Select® (Fits Seat Widths 14"-19")  High Profile  Mid Profile  Low Profile.....\$720

### ROHO® ADJUSTABLE SKIN PROTECTION (Please see page 4 for size chart.)

CSH11  Single Compartment (Fits Seat Widths 14"-19")  High Profile  Mid Profile  Low Profile.....\$720

SCD1  Sensor Ready Cushion with Smart Check Device (Not available with Low Profile.).....\$265

### ROHO® AIR / FOAM COMBINATION (Please see page 3 for size chart.)

CSH14  Hybrid Elite® SR single compartment (Fits Seat Widths 14"-19").....\$745

SCD1  Smart Check Device.....\$265

CSH18  AirLITE® (Weight Limit: 275 lbs, Fits chair widths 13"-20").....\$250

### COMFORT CUSHIONS

CSH2  Comfort Curve Cushion (Stretch-Air Cover. Not available with Seat Taper).....\$155

CSH23  Comfort Acta-Embrace™ Comfort-Tek Zero Elevation Cushion (Please see below for size chart.).....\$470

CSH24  Comfort Elements™ 2" Thick Comfort-Tek Cushion (Limited width and depth sizes. Please see below for size chart.).....\$140

# TiLite Seating System Order Form

Effective November 15, 2019 (New Items in Red)

## CUSHIONS

### Acta-Embrace

Rear Seat Width	Seat Depth
12"	10", 11", 12", 13", 14", 15", 16"
13"	10", 11", 12", 13", 14", 15", 16", 17"
14"	10", 11", 12", 13", 14", 15", 16", 17", 18"
15"	14", 15", 16", 17", 18", 19", 20"
16"	14", 15", 16", 17", 18", 19", 20"
17"	15", 16", 17", 18", 19", 20"
18"	15", 16", 17", 18", 19", 20"
19"	16", 17", 18", 19", 20"
20"	16", 17", 18", 19", 20"

### Hybrid Elite®

Rear Seat Width	Seat Depth
14"	14", 15", 16"
15"	15", 16", 17"
16"	15", 16", 17", 18", 20"
17"	16", 17", 18"
18"	16", 17", 18", 20"

### AirLITE®

Rear Seat Width	Seat Depth
13"	13", 15"
14"	14", 16"
15"	15", 17"
16"	16", 18"
17"	17", 18", 20"
18"	16", 18", 20"
19"	17", 19"
20"	17", 18"

### Elements Gel Cushion

Rear Seat Width	Seat Depth
16"	16", 18", 20"
18"	16", 18", 20"
20"	16", 18", 20"

### QUADTRO SELECT®

Rear Seat Width	12" Seat Depth	13" Seat Depth	14" Seat Depth	15" Seat Depth	16" Seat Depth	17" Seat Depth	18" Seat Depth	19" Seat Depth	20" Seat Depth
12"	Low, High	Low, High	N/A	N/A	N/A	N/A	N/A	N/A	N/A
13"	N/A	Low, High	Low, High	Low, High	Low, High	Low, High	N/A	N/A	N/A
14"	N/A	Low, High	Low, High	Low, Mid, High	Low, Mid, High	Low, Mid, High	Low, High	Low, High	N/A
15"	N/A	Low, High	Low, Mid, High	Low, Mid, High	Low, Mid, High	Low, Mid, High	Low, High	Low, High	N/A
16"	N/A	N/A	Low, High	Low, High	Low, Mid, High	Low, Mid, High	Low, Mid, High	Low, Mid, High	Low, High
17"	N/A	N/A	Low, High	Low, High	Low, Mid, High	Low, Mid, High	Low, Mid, High	Low, Mid, High	Low, High
18"	N/A	N/A	Low, High	Low, High	Low, Mid, High	Low, Mid, High	Low, Mid, High	Low, Mid, High	Low, High
19"	N/A	N/A	Low, High	Low, High	Low, Mid, High	Low, Mid, High	Low, Mid, High	Low, Mid, High	Low, High

# TiLite Seating System Order Form

Effective November 15, 2019 (New Items in Red)

## CUSHIONS

### Single Compartment

*Italicized profiles are available without ROHO Smart Check*

**Bolded profiles are available with ROHO Smart Check**

Rear Seat Width	12" Seat Depth	13" Seat Depth	14" Seat Depth	15" Seat Depth	16" Seat Depth	17" Seat Depth	18" Seat Depth	19" Seat Depth	20" Seat Depth
12"	Low, High	Low, High	N/A	N/A	N/A	N/A	N/A	N/A	N/A
13"	N/A	<i>Low, Mid, High</i> <b>High</b>	<i>Low, Mid, High</i> <b>Mid, High</b>	<i>Low, Mid, High</i> <b>Mid, High</b>	<i>Low, High</i> <b>Mid, High</b>	<i>Low, High</i> <b>Mid, High</b>	<i>Low, High</i> <b>High</b>	N/A	N/A
14"	N/A	<i>Low, High</i> <b>High</b>	<i>Low, Mid, High</i> <b>Mid, High</b>	<i>Low, Mid, High</i> <b>Mid, High</b>	<i>Low, Mid, High</i> <b>Mid, High</b>	<i>Low, Mid, High</i> <b>Mid, High</b>	<i>Low, High</i> <b>High</b>	<i>Low, High</i>	N/A
15"	N/A	<i>Low, High</i> <b>High</b>	<i>Low, Mid, High</i> <b>Mid, High</b>	<i>Low, Mid, High</i> <b>Mid, High</b>	<i>Low, Mid, High</i> <b>Mid, High</b>	<i>Low, Mid, High</i> <b>Mid, High</b>	<i>Low, High</i> <b>High</b>	<i>Low, High</i>	N/A
16"	N/A	N/A	<i>Low, High</i> <b>High</b>	<i>Low, High</i> <b>High</b>	<i>Low, Mid, High</i> <b>Mid, High</b>	<i>Low, Mid, High</i> <b>Mid, High</b>	<i>Low, Mid, High</i> <b>Mid, High</b>	<i>Low, Mid, High</i> <b>Mid, High</b>	<i>Low, High</i>
17"	N/A	N/A	<i>Low, High</i> <b>High</b>	<i>Low, High</i> <b>High</b>	<i>Low, Mid, High</i> <b>Mid, High</b>	<i>Low, Mid, High</i> <b>Mid, High</b>	<i>Low, Mid, High</i> <b>Mid, High</b>	<i>Low, Mid, High</i> <b>Mid, High</b>	<i>Low, High</i>
18"	N/A	N/A	<i>Low, High</i> <b>Mid, High</b>	<i>Low, High</i> <b>Mid, High</b>	<i>Low, Mid, High</i> <b>Mid, High</b>	<i>Low, Mid, High</i> <b>Mid, High</b>	<i>Low, Mid, High</i> <b>Mid, High</b>	<i>Low, Mid, High</i> <b>Mid, High</b>	<i>Low, Mid, High</i> <b>Mid, High</b>
19"	N/A	N/A	<i>Low, High</i> <b>Mid, High</b>	<i>Low, High</i> <b>Mid, High</b>	<i>Low, Mid, High</i> <b>Mid, High</b>	<i>Low, Mid, High</i> <b>Mid, High</b>	<i>Low, Mid, High</i> <b>Mid, High</b>	<i>Low, Mid, High</i> <b>Mid, High</b>	<i>Low, Mid, High</i> <b>Mid, High</b>
20"	N/A	N/A	N/A	N/A	<i>Low, High</i> <b>High</b>	<i>Low, High</i> <b>High</b>	<i>Low, Mid, High</i> <b>Mid, High</b>	<i>Low, Mid, High</i> <b>Mid, High</b>	<i>Low, Mid, High</i> <b>Mid, High</b>

## GENERAL

Permobil manufactures TiLite manual wheelchairs in a wide variety of configurations to meet the varied needs of wheelchair users. The final selection of a particular model, its options and accessories, and its measurements rests solely with the wheelchair user and advising healthcare professional. These wheelchairs must be fitted by a qualified healthcare professional.

**Get the latest order forms at permobilus.com.** Prices, standard features, options and accessories are subject to change without notice. All rights reserved. Trademarks are identified by symbols™ and ®. TiLite, TiFit, the TiLite logo, Ti, the Ti logo, and TiShaft are registered trademarks of TiSport, LLC and are licensed to Permobil. Permobil® is a registered trademark of Permobil. Bodypoint® is a registered trademark of Bodypoint Inc. ROHO®, AGILITY®, Low Profile® and Quadro Select® are registered trademarks of ROHO Inc.. The ACTA-BACK®, ACTA-EMBRACE™, COMFORT-TEK™, COMPASS®, CURVE®, ELEMENTS®, and STRETCH AIR™ are registered trademarks of THE COMFORT COMPANY® and are licensed to Permobil. Q'Strain® is a registered trademark of Q'Strain Inc.

## NOTES

---



---



---



---



---