


# M5 Corpus

# permobil



Three column pricing on order/quote:

Please send order/quote to fax#: (800) 231-3256  
Email to: [sales@permobil.com](mailto:sales@permobil.com)

## Dealer Information

\* Contact: \_\_\_\_\_  
Dealer Code: \_\_\_\_\_  
\* Dealer Name: \_\_\_\_\_  
Address: \_\_\_\_\_  
\* City: \_\_\_\_\_  
\* Province/Postal Code: \_\_\_\_\_  
\* Phone#: \_\_\_\_\_  
\* Fax#: \_\_\_\_\_  
PO#: \_\_\_\_\_  
Email Address: \_\_\_\_\_

## Client Information

*Permobil recommends that the client is evaluated by a certified rehab specialist.*

\* First Name: \_\_\_\_\_  
\* Last Name: \_\_\_\_\_  
Diagnosis: \_\_\_\_\_  
Funding Source: \_\_\_\_\_  
Client Age: \_\_\_\_\_

## Client Measurements

\* Weight: \_\_\_\_\_  
\* Height: \_\_\_\_\_  
A) Top of Shoulders: \_\_\_\_\_  
B) Chest Depth: \_\_\_\_\_  
C) Chest Width: \_\_\_\_\_  
D) Seat Depth: \_\_\_\_\_  
E) Top of Head: \_\_\_\_\_  
F) Elbow to Hand: \_\_\_\_\_  
G) Seat Pan to Elbow: \_\_\_\_\_  
H) Hip Width: \_\_\_\_\_  
\*I) Knee to Foot: \_\_\_\_\_

**Permobil Ltd.**  
12 - 75 Mary St.  
Aurora, Ontario  
L4G 1G3  
Tel: (800) 736-0925  
Fax: (800) 231-3256  
[www.permobil.ca](http://www.permobil.ca)

Prices effective **October 21, 2019** and payable in Canadian dollars.

**Chair Model and Colors**

Part Number	Description	MSRP - CA
<input type="checkbox"/> I11210	M5 Base Corpus - M5 SPO <i>Price Includes: MWD M5 Power Base Set Up For One Power Seat Function, FlexLink Pro Suspension, 8" Flat-Free Casters, Tie-Down Hardware for Strap System, Charger and Permobil Connect. Max Speed = 7.5mph, Weight Capacity = 300lbs, STFH = 17.5".</i>	9,353.00
<input type="checkbox"/> I107729-99-0	M5 Base Corpus - M5 MPO <i>Price Includes: MWD M5 Power Base Set Up For Multiple Power Seat Functions, FlexLink Pro Suspension, 8" Flat-Free Casters, Tie-Down Hardware for Strap System, Charger and Permobil Connect. Max Speed = 7.5mph, Weight Capacity = 300lbs, STFH = 17.5".</i>	9,353.00

**CHOOSE DRIVE TIRES:**

Part Number	Description	MSRP - CA
<input type="checkbox"/> I107855-99-0	Flat Free Tire 3"W x 14"T Black - Standard Tread <i>Includes flat free standard tread caster wheels.</i>	No Charge
<input type="checkbox"/> I107854-99-0	Pneumatic Tire 3"W x 14"T Black - Standard Tread <i>Includes flat free standard tread caster wheels.</i>	No Charge
<input type="checkbox"/> I107315-99-0	Flat Free Tire 3"W x 14"T Black - Aggressive Tread <i>Includes flat free aggressive tread caster wheels.</i>	No Charge

**OPTIONAL DRIVE TIRES:**

Part Number	Description	MSRP - CA
<input type="checkbox"/> I107289-99-0	Pneumatic/Snow Tire 3"W x 14"T Black - Aggr Tread <i>Snow tires are designed for use in cold weather conditions. The nylon compound in these tires remains flexible at low temperatures, but wears more easily. Therefore, the snow tires (sold as wheel assemblies) are not recommended for year-long use and must be ordered along with standard drive wheels. Does not include caster wheels.</i>	540.00

**CHOOSE SHROUD/RIM COLOR:**

Part Number	Description	MSRP - CA
<input type="checkbox"/> I107836-99-0	M5 Color Kit - Sunburst Orange	No Charge
<input type="checkbox"/> I107831-99-0	M5 Color Kit - Galactic Green	No Charge
<input type="checkbox"/> I107832-99-0	M5 Color Kit - Midnight Black	No Charge
<input type="checkbox"/> I107833-99-0	M5 Color Kit - Ocean Blue	No Charge
<input type="checkbox"/> I107834-99-0	M5 Color Kit - Power Pink	No Charge
<input type="checkbox"/> I107835-99-0	M5 Color Kit - Radiant Red	No Charge
<input type="checkbox"/> I107850-99-0	M5 Color Kit - Phantom Grey	No Charge
<input type="checkbox"/> I107849-99-0	M5 Color Kit - Titan Blue	No Charge
<input type="checkbox"/> I107899-99-0	M5 Color Kit - Carbon Fiber	No Charge
<input type="checkbox"/> I107900-99-0	M5 Color Kit - Mossy Oak Break-Up Country	No Charge
<input type="checkbox"/> I107983-99-0	M5 Color Kit - CDG Poly	No Charge
<input type="checkbox"/> I108297-99-0	M5 Color Kit - Sugar Skulls	No Charge
<input type="checkbox"/> I108302-99-0	M5 Color Kit - Fireball	No Charge
<input type="checkbox"/> I108307-99-0	M5 Color Kit - Water Lily	No Charge

**Base Options**

**MUST CHOOSE EITHER SEAT ELEVATOR OR FIXED SEAT TUBE:**

Part Number	Description	MSRP - CA
<input type="checkbox"/> I107086-99-0	Active Height Power Adj. Seat Height - 14" Travel <i>Must also select Power Tilt when selecting this option.</i>	3,544.00
<input type="checkbox"/> I11117	Seat Tube Fixed	No Charge

**CHOOSE ONE OF THE FOLLOWING:**

Part Number	Description	MSRP - CA
<input type="checkbox"/> I107435-99-0	Set Seat Height to 17.5"	No Charge

## Base Options

## CHOOSE ONE OF THE FOLLOWING:

Part Number	Description	MSRP - CA
<input type="checkbox"/> I107433-99-0	Set Seat Height at 18.5" <i>If ordering Power Adjustable Seat Height, 1" of elevator travel is used to set seat at 18.5".</i>	No Charge
<input type="checkbox"/> I107434-99-0	Set Seat Height at 19.5" <i>If ordering Power Adjustable Seat Height, 2" of elevator travel is used to set seat at 19.5".</i>	No Charge

## MUST CHOOSE BATTERIES:

Part Number	Description	MSRP - CA
<input checked="" type="checkbox"/> I107031-99-0	BATTERIES OPTIONS: Batteries, Grp 24 (72Ah) Sealed Gel, x2, Installed	965.00

## OPTIONAL ITEMS:

Part Number	Description	MSRP - CA
<input type="checkbox"/> I107998-99-0	Enhanced Steering Perform Unit (ESP), R-net	1,223.00
<input type="checkbox"/> I107459-99-0	Essentials Box with Reflectors <i>Please Note: The total backrest angle will be limited to 150° to avoid contact between backrest components and the Essentials Box.</i>	552.00
<input type="checkbox"/> I107204-99-0	Essentials Box with Lights <i>Please Note: The total backrest angle will be limited to 150° to avoid contact between backrest components and the Essentials Box. Must also select Light Kit.</i>	919.00
<input type="checkbox"/> I107741-99-0	Light Kit for M3 / M5	1,195.00
<input type="checkbox"/> I10715	Compact Laptop Ventilator Tray - Corpus Seat <i>Mounts to rear of Corpus backrest hinge. Accommodates Pulmonetic Systems / CareFusion LTV Series and similar sized Laptop Ventilators. The tray is 330mm Tall, adjustable Width between 200mm - 330mm, and adjustable Depth between 50mm - 165mm. Overall backrest angle limited to 150° when mounting this ventilator tray. Custom Option. Additional drive restrictions may apply.</i>	1,732.00
<input type="checkbox"/> I10787	Battery Tray for Laptop Ventilator Tray <i>Mounts to rear of Laptop Ventilator Tray. Measures 290mm W x 170mm H x 42mm D. Accommodates Pulmonetic Systems Lithium-Ion Battery Pack, CareFusion SprintPack Battery System, or Richardson PowerTech Vent Power Center. Custom Option.</i>	277.00
<input type="checkbox"/> I10788	O2 Holder for Laptop Ventilator Tray <i>Mounts to Left or Right Side of Laptop Vent Tray. O2 Holder is 175mm Tall. Width infinitely adjustable 90mm - 170mm. Depth infinitely adjustable 95mm - 125mm. Custom Option.</i>	391.00
<input type="checkbox"/> ID51113	Programming Key for R-net - Permobil Version	540.00

## Joystick Options

## CHOOSE JOYSTICK TYPE:

Part Number	Description	MSRP - CA
	PLEASE NOTE: <i>For all remote joystick modules shown below, you must also select an expandable controller. For additional input devices, refer to Alternative Drive Systems Section.</i>	
<input type="checkbox"/> I107779-99-0	Color Joystick Module, (CJSM) <i>This R-net color joystick features: A Large High-Resolution 3.5" VGA color display, control of up to 4 Bluetooth devices, Infrared control, programmable buttons and a programmable switch jack; all in a compact, durable Aluminum housing.</i>	1,375.00
<input type="checkbox"/> I107778-99-0	Permobil Joystick Module <i>Designed exclusively for Permobil, this R-net color joystick features: a 2.8" color display, control of up to 4 Bluetooth devices, Infrared control, programmable buttons and a programmable switch jack; all in a compact, durable Aluminum housing.</i>	1,221.00
<input type="checkbox"/> I11186	Permobil Joystick Module, NC <i>Designed exclusively for Permobil, this R-net color joystick features: a 2.8" color display, control of up to 4 Bluetooth devices, Infrared control, programmable buttons and a programmable switch jack; all in a compact, durable Aluminum housing.</i>	No Charge
<input type="checkbox"/> I107829-99-0	3.5mm Y-Splitter Cable for PJSM/CJSM <i>For use with the PJSM or CJSM Mode/Profile jack. Allows 2 mechanical switches to control 4 separate functions. Requires PC Programming to assign the short hold and long hold commands to each switch. Switches sold separately.</i>	43.00

## Joystick Options

## CHOOSE JOYSTICK TYPE:

Part Number	Description	MSRP - CA
<input type="checkbox"/> I105700-99-0	R-net Remote Colour JS w/ Lights, Heavy Duty Gimbal <i>Heavy-Duty modification enables gimbal to hold up to extreme operating forces. Reinforced housing is attached to a standard Colour LCD R-net Joystick. 3" long joystick knob has 3/4" diameter. Not compatible with alternative joystick handles.</i>	2,103.00

## CHOOSE JOYSTICK HANDLE TYPE:

Part Number	Description	MSRP - CA
<input type="checkbox"/> I10045	Cone Shaped Joystick Knob - Standard	No Charge
<input type="checkbox"/> I100960-99-0	Chin Cup For Joystick Knob	123.00
<input type="checkbox"/> I100961-99-0	"T" Handle For Joystick Knob	181.00
<input type="checkbox"/> I100962-99-0	"Large Ball" For Joystick Knob	123.00
<input type="checkbox"/> I100963-99-0	"Softball" For Joystick Knob	165.00
<input type="checkbox"/> I100964-99-0	"Mushroom" For Joystick Knob	165.00
<input type="checkbox"/> I100965-99-0	"Stick" For Joystick Knob	165.00
<input type="checkbox"/> I105442-99-0	1.75" Black Foam Ball For Joystick Knob	123.00
<input type="checkbox"/> I105846-99-0	Bodypoint J/S Handle 3" U-Shaped w/ Flex-Shaft	158.00
<input type="checkbox"/> I105847-99-0	Bodypoint J/S Handle 4" U-Shaped w/ Flex-Shaft	158.00
<input type="checkbox"/> I105848-99-0	Bodypoint J/S Handle Dome Shaped, Rubber	158.00

Part Number	Description	MSRP - CA
<input type="checkbox"/> I10664	Joystick Mount Right - Fixed	No Charge
<input type="checkbox"/> I10665	Joystick Mount Left - Fixed	No Charge
<input type="checkbox"/> I108037-99-0	Slimline Retractable Joystick Mount - Right	466.00
<input type="checkbox"/> I108038-99-0	Slimline Retractable Joystick Mount - Left	466.00
<input type="checkbox"/> I11232	Adjustable Height Panel Bracket <i>Provides 1.77" of vertical adjustment below the armrest and 25° of angle adjustment. Compatible with all Permobil joystick mounts except for the M1 Fixed Joystick Mount.</i>	361.00
<input type="checkbox"/> I105210-99-0	Mount Standard Joystick on Corpus Back Hinge <i>Standard R-net Joystick is mounted to backrest hinge. Can be chosen instead of PN I105408-99-0. Mount on Left: _____ Mount on Right: _____</i>	323.00
<input type="checkbox"/> I108268-99-0	Bodypoint Tri-Lock Midline Mount Assembly <i>Includes L-shaped joystick arm and height adjustable mounting clamp. Joystick mount sold separately (select from the mounts listed below). Please indicate if you want midline mounted on left or right armrest: Mount Midline on Left Armrest: _____ Mount Midline on Right Armrest: _____</i>	711.00
<input type="checkbox"/> I108280-99-0	Bodypoint Tri-Lock Mount for R-net JS	291.00
<input type="checkbox"/> I108279-99-0	Bodypoint Tri-Lock Mount for Compact Joystick	291.00
<input type="checkbox"/> I10842	PERMOFix Mount for Mini/Compact Joystick <i>Includes a shim that allows the PERMOFix joint to be mounted on the joystick arm. Provides an adjustable mount for a Mini Joystick or Compact Joystick.</i>	291.00

Part Number	Description	MSRP - CA
<input checked="" type="checkbox"/> I10411	Expandable Controller - R-net <i>Must be ordered with all expandable/remote joysticks.</i>	No Charge
<input checked="" type="checkbox"/> I10513	Harness for Expandable Controller <i>Must be ordered with any Expandable Controller.</i>	No Charge
<input type="checkbox"/> I105225-99-0	Co-Pilot Attendant Control R-net - Corpus	5,185.00
<input type="checkbox"/> I105408-99-0	Attendant Control R-net	1,346.00

## Joystick Options

Part Number	Description	MSRP - CA
<input type="checkbox"/> I107350-99-0	R-net Bluetooth iDevice Module <i>Allows input device commands (such as joystick nudges or switch clicks) to be used to operate an iPhone, iPad or iPod. By assigning each input to an iOS function, users can make phone calls, send texts and emails and access applications. Kit includes a Bluetooth iDevice Module, an R-net Bus Cable, and two mono jack cables for direct switch inputs (switches not included). Must be ordered with a Color JSM or Omni for mode display.</i>	1,087.00
<input type="checkbox"/> I107367-99-0	R-net Bluetooth Mouse Module <i>Enables mouse control through any R-net input device. It can operate the mouse on a computer and on select Android devices (with R-net for Android app). Kit includes a Bluetooth Mouse Module, an R-net Bus Cable, two mono jack cables for L/R click (switches not included) and a USB BT Adapter for user's computer. Note: BT Adapter may not work with all computers. Must be ordered with a Color JSM or Omni for mode display.</i>	1,087.00
<input type="checkbox"/> I107368-99-0	Add a second Bluetooth Mouse Module <i>Allows R-net to simultaneously pair to a second compatible computer or Android device. Each module is accessible through a separate R-net mode. Kit includes a unique Bluetooth Mouse Module, an R-net Bus Cable, two mono jack cables for L/R click and a USB BT Adapter.</i>	1,087.00
<input type="checkbox"/> I10582	Output Module (IOM) for R-net Control System <i>This allows a user to control an external device (ex: Communication Device, ECU, Door Opener, etc.) through the drive input (JSM, Head Array, etc.). Each module can provide up to 4 switch closures. May require specialty cable or wiring. Up to three output modules can be installed on one wheelchair.</i>	1,801.00

## Seating System

## CHOOSE ONE OF THE FOLLOWING:

Part Number	Description	MSRP - CA
	PLEASE NOTE: <i>Speed limitations are imposed when the total backrest angle is greater than 150° or when the seat is tilted forward more than 5°. Drive is inhibited when the seat is tilted back more than 25°</i>	
<input type="checkbox"/> I11212	Corpus Seat w/ No Tilt & 180° Power Recline	9,436.00
<input type="checkbox"/> I11213	Corpus Seat w/ 50° Tilt	8,461.00
<input type="checkbox"/> I11214	Corpus Seat w/ 50° Power Tilt & 180° Power Recline	13,941.00

## CHOOSE ONE OF THE FOLLOWING:

Part Number	Description	MSRP - CA
	PLEASE NOTE: <i>Active Reach should be recommended by a licensed clinical professional and/or ATP. While operating Active Reach, positioning equipment such as lap belts and chest supports should be considered to meet the user's specific safety needs. Note: A low speed limitation is imposed from 5° to 10° of Active Reach, an extra low speed limitation is imposed from 10° to 20° of Active Reach, and a drive lockout is imposed above 20° of Active Reach.</i>	
<input type="checkbox"/> I11141	5° Active Reach (Anterior Tilt) <i>When ordered with a Power Adjustable Seat Height, tilt range will be -5° to +50°. When ordered without a Power Adjustable Seat Height, tilt range will be -5° to +45°. If this option is selected with a fixed seat, the minimum seat to floor height at 0° of tilt will be approximately 19.5". Power Legrest recommended for full functionality. Minimum legrest angle will be ~ 97° when manual legrest is selected.</i>	No Charge
<input type="checkbox"/> I107090-99-0	10° Active Reach (Anterior Tilt) <i>Please Note: This option is only available on configurations ordered with Power Adjustable Seat Height and Power Tilt. A chest support belt (see Accessories section) is strongly recommended when selecting this option. Power Legrest recommended for full functionality. Minimum legrest angle will be ~ 103° when manual legrest is selected.</i>	No Charge
<input type="checkbox"/> I107091-99-0	20° Active Reach (Anterior Tilt) <i>Please Note: This option is only available on configurations ordered with Power Adjustable Seat Height, Power Tilt, Power Recline and Power Legrest Elevation. Must select Corpus VS Power Elevating Legrest, Corpus VS Footplates and Corpus VS Knee Support from VS Legrest Options section. A chest support belt (see Accessories section) is strongly recommended when selecting this option. Not compatible with retractable positioning belts.</i>	3,047.00

## Seating System Options



## Seating System Options

### CHOOSE SEAT FUNCTION CONTROL TYPE:

Part Number	Description	MSRP - CA
<input type="checkbox"/> I10418	Single Seat Function Control Switch for R-net <i>Please choose either Push Buttons or Toggle Switches below. Does not include thru drive control operation of power seat functions.</i>	No Charge
<input type="checkbox"/> I10419	Single Seat Function Control Kit for R-net <i>For operation of a single seat function through the wheelchair electronics via the drive control and/or ICS Switchbox. Order with R-net Controller or remote joystick options. Please indicate below if you want push buttons or toggles on the ICS Switchbox.</i>	1,711.00
<input type="checkbox"/> I10585	Multiple Seat Function Control Kit for R-net <i>MUST be selected if ordering two or more seat functions. This device is used to operate two or more seat functions through the wheelchair electronics via the drive control and/or ICS switch box. Order with R-net electronics or remote joystick options. Please indicate below if you want push buttons or toggles on the ICS Switchbox.</i>	2,851.00

### CHOOSE PUSH BUTTONS OR TOGGLES:

Part Number	Description	MSRP - CA
<input type="checkbox"/> I102798-99-0	Push Buttons on Seat Function Control Kit	No Charge
<input type="checkbox"/> I102799-99-0	Toggle Switches on Seat Function Control Kit	305.00
<input type="checkbox"/> I105409-99-0	Do Not Install ICS Switchbox, Send as Part	No Charge

### CHOOSE MOUNTING POSITION OF ICS SWITCHBOX:

Part Number	Description	MSRP - CA
<input type="checkbox"/> I10626	Mount ICS Switchbox onto R-net Joystick Bracket <i>ICS Switchbox will be mounted to the same bracket as the R-net Joystick Module. If this item is selected, the ICS Switchbox will swing away with the Joystick Module.</i>	No Charge
<input type="checkbox"/> I104651-99-0	Mount ICS Switchbox onto Armrest <i>Mount on Left: _____ Mount on Right: _____</i>	No Charge
<input type="checkbox"/> I106256-99-0	Mount ICS Switchbox on Backrest <i>Mount on Left: _____ Mount on Right: _____</i>	171.00

### OPTIONAL ITEMS:

Part Number	Description	MSRP - CA
<input type="checkbox"/> I104196-99-0	ICS Alternative Switchbox for R-net Electronics <i>Seat function interface box that accommodates up to 8 momentary switches. Provides an additional method for accessing the seat functions. 1/8" mini-mono switches not included.</i>	1,346.00
<input type="checkbox"/> I107088-99-0	Remote Stop System - F/M Series <i>Consists of a 315MHz Radio Frequency (RF) handheld transmitter and an RF receiver mounted on the wheelchair.</i>	1,482.00
<input type="checkbox"/> I10750	Memory Seat Program, Corpus w/ R-net	125.00
<input type="checkbox"/> I106209-99-0	Legrest Movement Linked to Recline Electronically <i>Power Legrest will raise when Power Recline moves backward and will lower when back moves forward. If Power Legrest is operated, Power Recline will NOT automatically reposition. Feature works through ICS Switchbox and Driver Control.</i>	50.00
<input type="checkbox"/> I10898	Enable Independent Repositioning Mode <i>IRM allows for sequential movement of power seat functions into a posterior position (35° tilt, 115° recline, 135° legrest) and returns to upright position (5° tilt, 100° recline, 110° legrest). This can be accessed using a standard drive control or a single switch. On the Permobil Joystick Module, IRM is accessible through the mode/profile jack (long switch hold). Please see the Alt Drive Accessories section for switch options.</i>	No Charge

## Backrest Options

### CHOOSE ONE OF THE FOLLOWING:

Part Number	Description	MSRP - CA
<input type="checkbox"/> I105067-99-0	Sliding Backrest Kit - 5" Mechanical Shear	No Charge
<input type="checkbox"/> I105066-99-0	Fixed Backrest Kit - No Mechanical Shear	55.00

### CHOOSE ONE OF THE FOLLOWING:

Part Number	Description	MSRP - CA
-------------	-------------	-----------

**Backrest Options**

**CHOOSE ONE OF THE FOLLOWING:**

Part Number	Description	MSRP - CA
<b>BACKREST HEIGHT NOTE:</b>		
<i>Backrest heights are measured from the seat pan to the top of the backrest shell with the backrest set at 90°. Corpus Ergo Backs include lateral wedges and foam lumbar pads. Solid Backs do not include foam lumbar pads.</i>		
<input type="checkbox"/> I106545-99-0	Corpus Ergo Back - 14"W x 20"T (C3GERGBCA) <i>Power Recline will be limited to 150° when ordering this backrest.</i>	1,104.00
<input type="checkbox"/> I106546-99-0	Corpus Ergo Back - 16"W x 20"T (C3GERGBCA) <i>Power Recline will be limited to 150° when ordering this backrest.</i>	1,104.00
<input type="checkbox"/> I106547-99-0	Corpus Ergo Back - 18"W x 20"T (C3GERGBCA) <i>Power Recline will be limited to 150° when ordering this backrest.</i>	1,104.00
<input type="checkbox"/> I106548-99-0	Corpus Ergo Back - 20"W x 20"T (C3GERGBCA) <i>Power Recline will be limited to 150° when ordering this backrest.</i>	1,104.00

**PLEASE NOTE:**

*Taller backrests listed below include: lateral wedges, foam lumbar pads\*, UT rail mounted at the top for lateral adjustment of the headrest, and split backrest feature with height adjustment in 1" increments. Please indicate below where you want the height set at so that a properly sized backrest cushion can be provided. \* Foam lumbar pads not included with Solid Backs.*

<input type="checkbox"/> I106532-99-0	Corpus Ergo Back - 14"W x 23" - 28"T (C3GERGBCA)	1,104.00
<input type="checkbox"/> I106533-99-0	Corpus Ergo Back - 16"W x 23" - 28"T (C3GERGBCA)	1,104.00
<input type="checkbox"/> I106534-99-0	Corpus Ergo Back - 18"W x 23" - 28"T (C3GERGBCA)	1,104.00
<input type="checkbox"/> I106535-99-0	Corpus Ergo Back - 20"W x 23" - 28"T (C3GERGBCA)	1,104.00

**MUST ALSO SELECT BACK HEIGHT:**

- Set Corpus Backrest at 23" \_\_\_\_\_
- Set Corpus Backrest at 24" \_\_\_\_\_
- Set Corpus Backrest at 25" \_\_\_\_\_
- Set Corpus Backrest at 26" \_\_\_\_\_
- Set Corpus Backrest at 27" \_\_\_\_\_
- Set Corpus Backrest at 28" \_\_\_\_\_

**CHOOSE ONE OF THE FOLLOWING:**

Part Number	Description	MSRP - CA
<input type="checkbox"/> I11335	Leatherette Cushion - Corpus Ergo Back <i>Includes Dual Density Backrest Foam.</i>	No Charge
<input type="checkbox"/> I11336	Stretch-Air Cushion - Corpus Ergo Back <i>Includes Dual Density Backrest Foam with Stretch-Air breathable fabric cover.</i>	No Charge

**CHOOSE ONE OF THE FOLLOWING:**

Part Number	Description	MSRP - CA
<b>ROHO AGILITY NOTE:</b>		
<i>The ROHO AGILITY Back Systems are available in widths of 14", 16", 18" and 20" and shell heights of 16", 18" and 20". Direct Mount system includes backrest shell, cushion, air inserts, UT rail for headrest and Direct Mount adapter for Corpus.</i>		
<input type="checkbox"/> I11258	ROHO AGILITY Max Contour Back System - 16" Tall <i>Backrest can be installed at a height (seat pan to top of AGILITY) between 22" and 25". Specify Shell Width (14", 16", 18", 20"W): _____</i>	1,155.00
<input type="checkbox"/> I11259	ROHO AGILITY Max Contour Direct Mount - 18" Tall <i>Backrest can be installed at a height (seat pan to top of AGILITY) between 22" and 26". Specify Shell Width (14", 16", 18", 20"W): _____</i>	1,155.00
<input type="checkbox"/> I11260	ROHO AGILITY Max Contour Direct Mount - 20" Tall <i>Backrest can be installed at a height (seat pan to top of AGILITY) between 24" and 27". Specify Shell Width (14", 16", 18", 20"W): _____</i>	1,155.00

**MUST ALSO SELECT BACK HEIGHT:**

- Set AGILITY Backrest at 22" \_\_\_\_\_
- Set AGILITY Backrest at 23" \_\_\_\_\_

## Backrest Options

### CHOOSE ONE OF THE FOLLOWING:

Part Number	Description	MSRP - CA
	Set AGILITY Backrest at 24" _____	
	Set AGILITY Backrest at 25" _____	
	Set AGILITY Backrest at 26" _____	
	Set AGILITY Backrest at 27" _____	

### CHOOSE ONE OF THE FOLLOWING:

Part Number	Description	MSRP - CA
<input type="checkbox"/> I107488-99-0	Permobil 3" x 5" Lateral Supports - Pair <i>Must also select I107490-99-0. Only compatible with Corpus Ergo backrest shells. Backrest Lateral Wedges will be omitted.</i>	473.00
<input type="checkbox"/> I107489-99-0	Permobil 5" x 7" Lateral Supports - Pair <i>Must also select I107490-99-0. Only compatible with Corpus Ergo backrest shells. Backrest Lateral Wedges will be omitted.</i>	473.00
<input type="checkbox"/> I107490-99-0	Adj/Removable Lateral Link Hardware - Pair	648.00

### OPTIONAL ITEMS:

Part Number	Description	MSRP - CA
<input type="checkbox"/> I105782-99-0	Tubular Backrest Frame - Corpus (1" Canes) <i>Provides 1" diameter, vertical canes for mounting aftermarket backrests. Frame can be adjusted to support backrest widths from 14" to 20". Canes are attached with 2" long UniTrack linkage disks that can rotate to provide space for backrest hardware. The tops of the canes can be adjusted up to a height of 27" above the seat pan. Depending on the aftermarket backrest, several inches of seat depth may be lost when mounting this system.</i>	491.00
<input type="checkbox"/> I106072-99-0	Direct Backrest Frame - Corpus <i>Used to mount ROHO Agility, Jay 3, Matrix Elite, Varilite Icon and Ride Java backrests ranging in width from 14" to 20". Slotted brackets and 1" UniTrack linkage disks provide a close, adjustable mount. Backrest pitch can be tuned by rotating the linkage disks. Some seat depth may be lost due to the position and orientation of the backrest.</i>	491.00
<input type="checkbox"/> I10725	Split Backrest Cushion - Corpus <i>Custom Option.</i>	366.00
<input type="checkbox"/> I10786	Tailored Backrest Cushion - Corpus <i>Sides of cushion are trimmed to provide additional clearance for lateral trunk supports. Must also select a backrest cushion when ordering this custom option.</i>	114.00

## Seat Sizing

### CHOOSE ONE OF THE FOLLOWING:

Part Number	Description	MSRP - CA
<input type="checkbox"/> I107023-99-0	Set Seat Width at 17": _____	No Charge
<input type="checkbox"/> I107024-99-0	Set Seat Width at 19": _____	No Charge
<input type="checkbox"/> I107025-99-0	Set Seat Width at 21": _____	No Charge
<input type="checkbox"/> I107026-99-0	Set Seat Width at 23": _____	No Charge

### CHOOSE ONE OF THE FOLLOWING:

Part Number	Description	MSRP - CA
	PLEASE NOTE: <i>There are three seat depth kits, each providing three seat depth settings: (1) 14" - 16", (2) 17" - 19" and (3) 20" - 22". Please select the required seat depth and the appropriate kit will be provided.</i>	
<input type="checkbox"/> I107013-99-0	Seat Depth 14"	No Charge
<input type="checkbox"/> I107014-99-0	Seat Depth 15"	No Charge
<input type="checkbox"/> I107015-99-0	Seat Depth 16"	No Charge
<input type="checkbox"/> I107016-99-0	Seat Depth 17"	No Charge
<input type="checkbox"/> I107017-99-0	Seat Depth 18"	No Charge
<input type="checkbox"/> I107018-99-0	Seat Depth 19"	No Charge
<input type="checkbox"/> I107019-99-0	Seat Depth 20"	No Charge
<input type="checkbox"/> I107020-99-0	Seat Depth 21"	No Charge



## Seat Sizing

## CHOOSE ONE OF THE FOLLOWING:

Part Number	Description	MSRP - CA
<input type="checkbox"/> I107021-99-0	Seat Depth 22"	No Charge

## CHOOSE ONE OF THE FOLLOWING:

Part Number	Description	MSRP - CA
NOTE: A Corpus II Ergonomic Seat Cushion will be provided for select seat sizes.		
<input type="checkbox"/> I111337	Corpus Ergo Cushion - Leatherette Cover	623.00
<input type="checkbox"/> I111338	Corpus Ergo Cushion - Stretch-Air Cover	623.00
<input type="checkbox"/> I10965	Corpus Ergo Air Cushion (C3GRHOSCA) <i>Features ROHO Dry Floatation technology. The removable ROHO cell pad provides adjustable skin protection while the contoured foam and thigh support wedges help maintain proper positioning and stability. Includes removable stretch fabric cover. Available in widths of 17", 19" and 21" and depths of 17", 19" and 21", as well as 17" x 23". Not Available in 21" x 17"</i>	661.00
<input type="checkbox"/> I111113	Incontinence Cover for Ergo Cushion <i>Made from a stretch Dartex fabric and suitable for all Ergo cushions including the Corpus Ergo Air. The cover is cut to the contours of the Ergo cushion and is designed to be used as an outer layer of protection. Cover size will be matched to cushion size. Sent as a part. Part Number will be provided on the quote / order confirmation based on the Seat Width and Seat Depth selected.</i>	198.00
<input type="checkbox"/> I111166	ROHO QUADTRO SELECT HIGH PROFILE	769.00
<input type="checkbox"/> I111168	ROHO HIGH PROFILE Single Compartment	630.00
<input type="checkbox"/> I111170	ROHO HIGH PROFILE Sensor Ready <i>Single compartment cushion designed to work with Smart Check Device I107505-99-0. Not available in all seat widths and depths.</i>	630.00
<input type="checkbox"/> I107505-99-0	ROHO Smart Check Device <i>Sensing device that allows users to setup and maintain proper inflation pressure in ROHO Sensor Ready cushions. Only available when ordered with a ROHO Sensor Ready cushion.</i>	295.00
<input type="checkbox"/> I111173	Incontinence Cover for ROHO Cushion <i>Cover size will be matched to ROHO cushion size. Sent as a part.</i>	144.00
<input type="checkbox"/> I105118-99-0	Omit Seat Cushion	No Charge

## Armrest Options

## CHOOSE ONE OF THE FOLLOWING:

Part Number	Description	MSRP - CA
PLEASE NOTE: <i>Narrow Armbar Kit can be adjusted from 15"W to 19"W between the armrest pads without any additional parts. Wide Armbar Kit can be adjusted from 19"W to 23"W between the armrest pads without any additional parts. Please select both Armbar Kit and Width Setting below.</i>		
<input type="checkbox"/> I106543-99-0	Narrow Armbar Kit, Inf Adjustable 15"W - 19"W	No Charge
<input type="checkbox"/> I106544-99-0	Wide Armrest Kit, Inf Adjustable 19"W - 23"W <i>Required when selecting a 20" wide backrest.</i>	No Charge
Set the Width Between the Armrests to: _____ <i>Must be wider than backrest width to ensure proper fit. Please round to nearest 1".</i>		

## CHOOSE ONE OF THE FOLLOWING:

Part Number	Description	MSRP - CA
<input type="checkbox"/> I105069-99-0	Height Adjustable Armrests - NC SET ARMREST HEIGHT AT: _____ <i>Corpus Armrests are infinitely adjustable from 8" - 15", measured from the seat pan to the top of the armrest pad. No additional parts are needed to get this full range of adjustment. In one inch increments (ex: 8", 9", etc.), please indicate where you would like the armrests set.</i>	No Charge

## CHOOSE ONE OF THE FOLLOWING:

Part Number	Description	MSRP - CA
<input type="checkbox"/> I105071-99-0	Left Armrest Bracket - Fixed <i>Medial swivel cannot be enabled when selecting this option. Must be selected when ordering arm trough adapters.</i>	No Charge

## Armrest Options

### CHOOSE ONE OF THE FOLLOWING:

Part Number	Description	MSRP - CA
<input type="checkbox"/> I105070-99-0	Right Armrest Bracket - Fixed <i>Medial swivel cannot be enabled when selecting this option. Must be selected when ordering arm trough adapters.</i>	No Charge
<input type="checkbox"/> I105073-99-0	Left Armrest Bracket - Medial Swivel Enabled	No Charge
<input type="checkbox"/> I105072-99-0	Right Armrest Bracket - Medial Swivel Enabled	No Charge

### CHOOSE ONE OF THE FOLLOWING:

Part Number	Description	MSRP - CA
<input type="checkbox"/> I105100-99-0	4"x10" Corpus Armrest Pad Right - Leatherette	No Charge
<input type="checkbox"/> I105104-99-0	4"x13" Corpus Armrest Pad Right - Leatherette	No Charge
<input type="checkbox"/> I105105-99-0	4"x16" Corpus Armrest Pad Right - Leatherette	No Charge
<input type="checkbox"/> I105999-99-0	4"x18" Corpus HD Armrest Pad Right - Leatherette	No Charge
<input type="checkbox"/> I107212-99-0	4"x10" Corpus Armrest Pad Right - Black Fabric	No Charge
<input type="checkbox"/> I106047-99-0	4"x13" Corpus Armrest Pad Right - Black Fabric	No Charge
<input type="checkbox"/> I107214-99-0	4"x16" Corpus Armrest Pad Right - Black Fabric	No Charge
<input type="checkbox"/> I105995-99-0	4"x18" Corpus HD Armrest Pad Right - Black Fabric	No Charge
<input type="checkbox"/> I107591-99-0	4"x10" Corpus Arm w/ Permobil Gel Pad - Right	197.00
<input type="checkbox"/> I107593-99-0	4"x13" Corpus Arm w/ Permobil Gel Pad - Right	203.00
<input type="checkbox"/> I107595-99-0	4"x16" Corpus Arm w/ Permobil Gel Pad - Right	208.00
<input type="checkbox"/> I107597-99-0	4"x18" Corpus HD Arm w/ Permobil Gel Pad - Right	216.00
<input type="checkbox"/> I108861-99-0	Moldable Forearm Support - Right <i>11.3"L x 5"W. Not compatible with armrest pouch, trays (except Small Tray) or the VS Chest Support.</i>	343.00
<input type="checkbox"/> I108863-99-0	Contoured Hand Pad Medium - Right <i>6.4"L x 5.5" W. For use with Moldable Forearm Supports only. Not compatible with Slimline or Fixed joystick mounts or Adjustable Height Panel Bracket.</i>	170.00
<input type="checkbox"/> I108865-99-0	Contoured Hand Pad Large - Right <i>7.7"L x 6.1" W. For use with Moldable Forearm Supports only. Not compatible with Slimline or Fixed joystick mounts or Adjustable Height Panel Bracket.</i>	180.00

### CHOOSE ONE OF THE FOLLOWING:

Part Number	Description	MSRP - CA
<input type="checkbox"/> I105098-99-0	4"x10" Corpus Armrest Pad Left - Leatherette	No Charge
<input type="checkbox"/> I105102-99-0	4"x13" Corpus Armrest Pad Left - Leatherette	No Charge
<input type="checkbox"/> I105107-99-0	4"x16" Corpus Armrest Pad Left - Leatherette	No Charge
<input type="checkbox"/> I105998-99-0	4"x18" Corpus HD Armrest Pad Left - Leatherette	No Charge
<input type="checkbox"/> I107213-99-0	4"x10" Corpus Armrest Pad Left - Black Fabric	No Charge
<input type="checkbox"/> I106046-99-0	4"x13" Corpus Armrest Pad Left - Black Fabric	No Charge
<input type="checkbox"/> I107215-99-0	4"x16" Corpus Armrest Pad Left - Black Fabric	No Charge
<input type="checkbox"/> I105994-99-0	4"x18" Corpus HD Armrest Pad Left - Black Fabric	No Charge
<input type="checkbox"/> I107590-99-0	4"x10" Corpus Arm w/ Permobil Gel Pad - Left	197.00
<input type="checkbox"/> I107592-99-0	4"x13" Corpus Arm w/ Permobil Gel Pad - Left	203.00
<input type="checkbox"/> I107594-99-0	4"x16" Corpus Arm w/ Permobil Gel Pad - Left	208.00
<input type="checkbox"/> I107596-99-0	4"x18" Corpus HD Arm w/ Permobil Gel Pad - Left	216.00
<input type="checkbox"/> I108860-99-0	Moldable Forearm Support - Left <i>11.3"L x 5"W. Not compatible with armrest pouch, trays (except Small Tray) or the VS Chest Support.</i>	343.00
<input type="checkbox"/> I108862-99-0	Contoured Hand Pad Medium - Left <i>6.4"L x 5.5" W. For use with Moldable Forearm Supports only. Not compatible with Slimline or Fixed joystick mounts or Adjustable Height Panel Bracket.</i>	170.00
<input type="checkbox"/> I108864-99-0	Contoured Hand Pad Large - Left <i>7.7"L x 6.1" W. For use with Moldable Forearm Supports only. Not compatible with Slimline or Fixed joystick mounts or Adjustable Height Panel Bracket.</i>	180.00

### Armrest Options

#### OPTIONAL ITEMS:

Part Number	Description	MSRP - CA
<input type="checkbox"/> I105463-99-0	Corpus Arm Trough Adapter 1" - Pair <i>For aftermarket armrest pads. Please indicate below if they should be installed or on chair sent as parts. Please Note: Medial Swivel must be must be locked when selecting this option.</i>	357.00
<input type="checkbox"/> I105462-99-0	Corpus Arm Trough Adpater 7/8" - Pair <i>For aftermarket armrest pads. Please indicate below if they should be installed or on chair sent as parts. Please Note: Medial Swivel must be must be locked when selecting this option.</i> Install Arm Trough Adapters: _____ Send Arm Trough Adapters as Parts: _____	357.00
<input type="checkbox"/> I107339-99-0	Armrest Pouch - Right	260.00
<input type="checkbox"/> I107338-99-0	Armrest Pouch - Left	260.00
<input type="checkbox"/> I107840-99-0	UT Armrest Side Support - Right Outer <i>Includes 2"T x 6.5"W tapered pad.</i>	372.00
<input type="checkbox"/> I107863-99-0	UT Armrest Side Support - Left Outer <i>Includes 2"T x 6.5"W tapered pad.</i>	372.00
<input type="checkbox"/> I107864-99-0	UT Armrest Side Support - Right Inner <i>Includes 2"T x 6.5"W tapered pad. May need to increase width between the armrests to ensure armrests will flip-up.</i>	372.00
<input type="checkbox"/> I107865-99-0	UT Armrest Side Support - Left Inner <i>Includes 2"T x 6.5"W tapered pad. May need to increase width between the armrests to ensure armrests will flip-up.</i>	372.00
<input type="checkbox"/> I107839-99-0	UT Elbow Support - Right <i>Must also select PN I10881.</i>	200.00
<input type="checkbox"/> I107838-99-0	UT Elbow Support - Left <i>Must also select PN I10882.</i>	200.00
<input type="checkbox"/> I10881	Adjustable Elbow Support Hardware - Right	291.00
<input type="checkbox"/> I10882	Adjustable Elbow Support Hardware - Left	291.00

### Legrest Options

#### CHOOSE ONE OF THE FOLLOWING:

Part Number	Description	MSRP - CA
<input type="checkbox"/> I107763-99-0	CHOOSE LEGREST TYPE: Corpus Manual Elevating Legrest	No Charge
<input type="checkbox"/> I107926-99-0	Power Elevating Legrest - 170°	3,655.00
<input type="checkbox"/> I11197	CHOOSE FOOTPLATE STYLE: 2pc Narrow Corpus Footplates 9"D x 6"W	No Charge
<input type="checkbox"/> I11196	2pc Wide Corpus Footplates - 9"D x 7.5"W	No Charge
<input type="checkbox"/> I106760-99-0	1pc Narrow Corpus Footpl - 9.75"D x 12.5"W <i>Footplate is 12.5"W in front, 8"W in back.</i>	390.00
<input type="checkbox"/> I105252-99-0	CALF SUPPORT OPTIONS: UT Calf Support Kit (B) w Mesh Cover & Adj Hardware <i>3.5"W x 4.5"T Pad.</i>	No Charge
<input type="checkbox"/> I105258-99-0	UT Calf Support Kit (J) w Mesh Cover & Adj Hardware <i>5"W x 6.5"T Pad.</i>	No Charge
<input type="checkbox"/> I105250-99-0	UT Calf Support Kit (A) w Mesh Cover & Adj Hardware <i>5.5"W x 6"T Pad.</i>	No Charge
<input type="checkbox"/> I105256-99-0	UT Calf Support Kit (H) w Mesh Cover & Adj Hardware <i>7"W x 7"T Curved Pad.</i>	404.00
<input type="checkbox"/> I105254-99-0	UT Calf Support Kit (D) w Mesh Cover & Adj Hardware <i>8"W x 8"T Pad.</i>	404.00

## Legrest Options

## CHOOSE ONE OF THE FOLLOWING:

Part Number	Description	MSRP - CA
	<b>CORPUS VS LEGREST AND KNEE SUPPORT:</b>	
	<i>NOTE: The Corpus VS Legrest is only available with 20°, 30° and 45° of Active Reach.</i>	
<input type="checkbox"/> I107925-99-0	Corpus VS Power Elevating Legrest - 180°	3,754.00
<input type="checkbox"/> I107038-99-0	Knee Support - Corpus VS	215.00
	<i>Select when ordering Corpus VS Legrest. Must also select I11129.</i>	
<input type="checkbox"/> I11129	Adj Removable Knee Support Hardware (1 Each)	314.00
	<b>CHOOSE FOOTPLATE STYLE:</b>	
<input type="checkbox"/> I11198	2pc Narrow Corpus VS Footplates - 9"D x 6"W	No Charge
	<b>CALF SUPPORT OPTIONS:</b>	
<input type="checkbox"/> I107211-99-0	VS Calf Support Kit (J) w Mesh Cover & Adj Hardware <i>5"W x 6.5"T Pad. Compatible with footplate heights of 10" - 18".</i>	No Charge
<input type="checkbox"/> I106190-99-0	VS Calf Support Kit (A) w Mesh Cover & Adj Hardware <i>5.5"W x 6"T Pad. Compatible with footplate heights of 10" - 18".</i>	No Charge
<input type="checkbox"/> I106188-99-0	VS Calf Support Kit (H) w Mesh Cover & Adj Hardware <i>7"W x 7"T Curved Pad. Compatible with footplate heights of 11" - 18".</i>	433.00
<input type="checkbox"/> I106186-99-0	VS Calf Support Kit (D) w Mesh Cover & Adj Hardware <i>8"W x 8"T Pad. Compatible with footplate heights of 12" - 18".</i>	433.00

## CHOOSE ONE OF THE FOLLOWING:

Part Number	Description	MSRP - CA
	<b>SET FOOTPLATE HEIGHT AT: _____</b>	
	<i>This is measured from the SEAT PAN to the top of the FOOTPLATE, when the legrest assembly is at 90°. DO NOT INCLUDE CUSHION THICKNESS IN YOUR MEASUREMENT. In one inch increments, please indicate desired footplate height. NOTE: Corpus footplate height range is 6"-17" and Corpus VS height range is 9"-17". Height may need to be reduced for shipping. Calf Pads are not available at heights less than 9".</i>	

## CHOOSE ONE OF THE FOLLOWING:

Part Number	Description	MSRP - CA
<input type="checkbox"/> I106530-99-0	Stealth Gel Pads - Narrow Footplates (Pair)	411.00
<input type="checkbox"/> I106225-99-0	Stealth Gel Pads - Wide Footplates (Pair)	423.00

## Accessories

<input type="checkbox"/> I104518-99-0	Bodypoint Non-Padded Hip Belt - Medium	No Charge
<input type="checkbox"/> I104520-99-0	Bodypoint Padded Hip Belt - Medium	148.00
<input type="checkbox"/> I104514-99-0	Bodypoint Non-Padded Hip Belt - Large	119.00
<input type="checkbox"/> I104516-99-0	Bodypoint Padded Hip Belt - Large	203.00
<input type="checkbox"/> I100519-99-0	Positioning Belt Retractable	307.00
<input type="checkbox"/> I100600-99-0	Positioning Belt Retractable Long	333.00
<input type="checkbox"/> I108267-99-0	Bodypoint Padded Hip Belt "Airplane" Latch - Med	260.00
<input type="checkbox"/> I108266-99-0	Bodypoint Padded Hip Belt "Airplane" Latch - Large	275.00

## HEADREST OPTIONS:

Part Number	Description	MSRP - CA
<input type="checkbox"/> I107470-99-0	Permobil 7"W Headrest w/ Link Hardware <i>Must also select I107480-99-0 or I107481-99-0.</i>	452.00
<input type="checkbox"/> I107469-99-0	Permobil 10"W Headrest w/ Link Hardware <i>Must also select I107480-99-0 or I107481-99-0.</i>	452.00
<input type="checkbox"/> I107468-99-0	Permobil 14"W Headrest w/ Link Hardware <i>Must also select I107480-99-0 or I107481-99-0.</i>	452.00
<input type="checkbox"/> I107480-99-0	Adj/Removable Headrest Link Hardware - Standard	314.00
<input type="checkbox"/> I107481-99-0	Adj/Removable Headrest Link Hardware - Extended <i>Includes an additional 3" link to extend headrest pad range.</i>	339.00

## Accessories

## HEADREST OPTIONS:

Part Number	Description	MSRP - CA
<input type="checkbox"/> I108721-99-0	Adj/Removable BodiLink Headrest Hardware Max8 <i>Must be selected when ordering any of the Swing Away Chin Boom assemblies.</i>	440.00
<input type="checkbox"/> I11305	Swing Away Egg Switch Assembly <i>Mount on Upper Left:_____ Mount on Upper Right:_____</i>	614.00
<input type="checkbox"/> I11306	Swing Away Microlight Assembly <i>Mount on Upper Left:_____ Mount on Upper Right:_____</i>	705.00
<input type="checkbox"/> I11307	Swing Away Facial Pad Assembly 4"L x 2"T <i>Mount on Upper Left:_____ Mount on Upper Right:_____</i>	500.00
<input type="checkbox"/> I11308	Swing Away Facial Pad Assembly 6"L x 2"T <i>Mount on Upper Left:_____ Mount on Upper Right:_____</i>	500.00
<input type="checkbox"/> I11276	Provide Headrest Receiver for AGILITY Backrest <i>Select for AGILITY backs and third party backs ordered off chair.</i>	No Charge
<input type="checkbox"/> I106331-99-0	Permobil 10"W Headrest, Leatherette <i>Not available with ROHO AGILITY Back Systems. Must also select PN I10552.</i>	640.00
<input type="checkbox"/> I10552	Adj/Removable Headrest Hardware	137.00
<input type="checkbox"/> I107499-99-0	Universal Headrest Adapter <i>Please Note: Recline is limited to 160° or less (depending on backrest height) when ordering a Universal Headrest Adapter on a 23" or 24" tall Corpus backrest.</i>	258.00

## OPTIONAL ITEMS:

Part Number	Description	MSRP - CA
<input type="checkbox"/> I103070-99-0	UT Amp. Support Pad (H) w/ Mesh Cover - (1 Each) <i>7"W x 7"T Curved Pad. Must also select PN I10553.</i>	384.00
<input type="checkbox"/> I103071-99-0	UT Amp. Support Pad (D) w/ Mesh Cover - (1 Each) <i>8"W x 8"T Flat Pad. Must also select PN I10553.</i>	384.00
<input type="checkbox"/> I10553	Adj Removable - Amp Support Hardware (1 Each) <i>Please indicate if you would like the amputation support mounted on the right or left side of the seat frame: Mount Amp Suppt on Left:_____ Mount Amp Suppt on Right:_____</i>	137.00
<input type="checkbox"/> I106104-99-0	UT Hip Suppt (B) w/ Mesh Cover - High (Pair) <i>4.5"W x 3.5"T Pad w/ 4" Tall Mounting Bracket. Must also select PN I10550.</i>	313.00
<input type="checkbox"/> I106103-99-0	UT Hip Suppt (A) w/ Mesh Cover - High (Pair) <i>6"W x 5.5"T Pad w/ 4" Tall Mounting Bracket. Must also select PN I10550.</i>	313.00
<input type="checkbox"/> I103094-99-0	UT Thigh/Hip Suppt (J) w/ Mesh Cover - High (Pair) <i>6.5"W x 5"T Pad w/ 4" Tall Mounting Bracket. Must also select PN I10550.</i>	508.00
<input type="checkbox"/> I103095-99-0	UT Thigh/Hip Suppt (C) w/ Mesh Cover - High (Pair) <i>8"W x 3.5"T Pad w/ 4" Tall Mounting Bracket. Must also select PN I10550.</i>	508.00
<input type="checkbox"/> I102977-99-0	UT Thigh/Hip Suppt (G) w/ Mesh Cover - High (Pair) <i>12"W x 3.5"T Pad w/ 4" Tall Mounting Bracket. Must also select PN I10550.</i>	508.00
<input type="checkbox"/> I10550	Adj Removable - Thigh Support Hardware (Pair)	137.00
<input type="checkbox"/> I106762-99-0	Abductor Kit Corpus - Leatherette <i>Must also select PN I10555.</i>	380.00
<input type="checkbox"/> I10555	Adj Removable - Abductor Hardware (1 Each)	137.00
<input type="checkbox"/> I105260-99-0	Transfer Handles, Corpus Seat	528.00
<input type="checkbox"/> I105292-99-0	Push Handles, Corpus Seat	No Charge
<input type="checkbox"/> I106210-99-0	RAM Self-Leveling Cup Holder	91.00
<input type="checkbox"/> I106806-99-0	RAM X-Grip Phone Holder	91.00
<input type="checkbox"/> I107100-99-0	RAM X-Grip 10" Tablet Holder w/Swing Away Mount <i>Compatible with most 10" screen tablets with or without a protective case. Includes adjustable swing away mount. Note: Customization may be required if ordered along with an Omni.</i>	504.00



## Accessories

## OPTIONAL ITEMS:

Part Number	Description	MSRP - CA
<input type="checkbox"/> I107222-99-0	USB Charger <i>Provides up to 1.5A (at 5V) to charge most smartphones and tablets. The charger features an isolated circuit design for improved safety; and by pulling power from the ICS, prevents battery drainage when the chair is powered off. Includes mounting hardware and ICS connector. Please select armrest for mounting, keeping in mind that charger may need to be relocated. Mount on left armrest: _____ Mount on right armrest: _____</i>	215.00
<input type="checkbox"/> I105805-99-0	Medical Necessities Bag Hooks <i>Please Note: Recline is limited to 165° or less (depending on backrest height) when ordering with 23" or 24" tall Ergo backrest. Custom required if ordered with a ROHO AGILITY backrest. May ship as a part at shorter backrest heights.</i>	99.00
<input type="checkbox"/> I102059-99-0	Medical Necessities Bag <i>Includes hardware to attach bag to backrest. Please Note: Recline is limited to 165° or less (depending on backrest height) when ordering with a 23" or 24" tall Ergo backrest. Custom required if ordered with a ROHO AGILITY backrest.</i>	280.00
<input type="checkbox"/> I105840-99-0	Bodypoint Monoflex Medium Chest Support Belt Kit <i>18.25"L x 4"W Pad Size. Item will be shipped as a part.</i>	193.00
<input type="checkbox"/> I105839-99-0	Bodypoint Monoflex Large Chest Support Belt Kit <i>21.5"L x 4.5"W Pad Size. Item will be shipped as a part.</i>	193.00
<input type="checkbox"/> I105850-99-0	Bodypoint Monoflex X-Lrge Chest Support Belt Kit <i>24.5"L x 5"W Pad Size. Item will be shipped as a part.</i>	193.00
<input type="checkbox"/> I105262-99-0	Joystick Cutout Tray - 16" W <i>Must also select PN I10556.</i>	469.00
<input type="checkbox"/> I105263-99-0	Joystick Cutout Tray - 18" W <i>Must also select PN I10556.</i>	469.00
<input type="checkbox"/> I10843	Small Angle Adjustable Tray <i>Must also select PN I10556.</i>	469.00
<input type="checkbox"/> I10556	Adj Removable - Tray Hardware (1 Each) <i>Mount on Left: _____ Mount on Right: _____</i>	137.00
<input type="checkbox"/> I108230-99-0	Permobil Adjustable Tray - 13"W	345.00
<input type="checkbox"/> I108228-99-0	Permobil Adjustable Tray - 15"W	345.00
<input type="checkbox"/> I108227-99-0	Permobil Adjustable Tray - 17"W	345.00
<input type="checkbox"/> I108229-99-0	Permobil Adjustable Tray - 19"W	345.00
<input type="checkbox"/> I108232-99-0	Permobil Tray Hardware Right Side Adjustable	250.00
<input type="checkbox"/> I108233-99-0	Permobil Tray Hardware Left Side Adjustable	250.00
<input type="checkbox"/> I105224-99-0	Oxygen holder	341.00
<input type="checkbox"/> I105223-99-0	Crutch holder	354.00

## Alternative Drive Systems

## ALTERNATIVE PROPORTIONAL JOYSTICK SYSTEMS:

Part Number	Description	MSRP - CA
<input type="checkbox"/> I10751	Mini-Joystick for Omni - (9-pin connector) <i>Operating Force = 13 g Neutral to Full Deflection = 5mm</i>	3,122.00
<input type="checkbox"/> I103893-99-0	Compact Joystick for Omni - (9-pin connector) <i>Operating Force = 220 g Neutral to Full Deflection = 28mm Joystick includes two 3.5mm mono jacks for Mode and On/Off. Includes Standard Joystick Knob and 1.75" (42mm) Foam Ball Knob.</i>	1,523.00
<input type="checkbox"/> I103894-99-0	Compact Joystick LITE for Omni - (9-pin connector) <i>Operating Force = 50 g Neutral to Full Deflection = 17mm Joystick includes two 3.5mm mono jacks for Mode and On/Off. Includes Standard Joystick Knob and 1.75" (42mm) Foam Ball Knob.</i>	1,743.00

## ALT PROPORTIONAL JOYSTICK MOUNTING HARDWARE:

Part Number	Description	MSRP - CA
-------------	-------------	-----------

## Alternative Drive Systems

### ALT PROPORTIONAL JOYSTICK MOUNTING HARDWARE:

Part Number	Description	MSRP - CA
	NOTE: Please see the PERMOfix Mounting Solutions, within the Alt Drive Accessories section, for additional mounting options. The Swing Away Chin Boom Assemblies are Compatible with Permobil Headrests I107470-99-0, I107469-99-0 or I107468-99-0 and must be selected with I108721-99-0.	
<input type="checkbox"/> I107946-99-0	Swing Away Chin Boom Assm for Compact Lite - Left	565.00
<input type="checkbox"/> I107948-99-0	Swing Away Chin Boom Assm for Compact Lite - Right	565.00
<input type="checkbox"/> I107857-99-0	Swing Away Chin Boom Assm for Mini Joystick - Left	565.00
<input type="checkbox"/> I107947-99-0	Swing Away Chin Boom Assm for Mini Joystick - Right	565.00
<input type="checkbox"/> I104690-99-0	Alt. JS Swing Away Mount Kit - Armrest Mount in front of Left Arm: _____ Mount in front of Right Arm: _____	296.00
<input type="checkbox"/> I107135-99-0	Alt JS/Omni Combo Swing Away Mount Kit - Armrest Used to mount Omni and Alternative Proportional Joystick to the same armrest. Includes Swing Away Omni/Tablet Mount, with a longer mounting rod, two PERMOfix modules, one 3" rod and one 4" rod. Mount in front of Left Arm: _____ Mount in front of Right Arm: _____	599.00
<input type="checkbox"/> I104143-99-0	Alt. JS Adjustable Mount Kit - Chin Boom Swing-Away / Removable Chin Boom mounted on right of backrest shell. Custom charge may be applied if ordering lower backrest heights.	959.00
<input type="checkbox"/> I10752	Haims Harness for Chin Control Mount, Large For neck width up to 19cm.	685.00
<input type="checkbox"/> I10753	Haims Harness for Chin Control Mount, Small For neck width up to 14cm.	685.00

### ALT PROPORTIONAL JOYSTICK ACCESSORIES:

Part Number	Description	MSRP - CA
<input type="checkbox"/> I104904-99-0	Dual Switches for Compact Joystick - Fixed Mount Includes two pushbutton switches mounted on rigid Switches can be used for Mode & On/Off. Switch operation force = 130 g Mount at front of Joystick: _____ Mount at rear of Joystick: _____ Mount at left of Joystick: _____ Mount at right of Joystick: _____	284.00
<input type="checkbox"/> I104641-99-0	Dual Switches for Compact Joystick - Gooseneck Mnt Includes two pushbutton switches mounted on adjustable gooseneck arms. Switches are positioned on Left and Right sides of Compact Joystick. Switches can be used for Mode & On/Off. Switch operation force = 130 g	517.00

### SIP & PUFF SYSTEMS:

Part Number	Description	MSRP - CA
<input type="checkbox"/> I10573	Sip & Puff System for Omni Requires Omni Control Module. Includes breath tube kit.	1,831.00
<input type="checkbox"/> I10574	Sip & Puff System with Fiber Optic Opti-Stop Requires Omni Control Module. Includes breath tube kit.	2,437.00
<input type="checkbox"/> I106807-99-0	Sip & Puff Kit for Total Control Head Array For use with Total Control Head Array. Requires Omni Control Module. Includes breath tube kit.	1,011.00
<input type="checkbox"/> I10762	Adjustable Mounting Hrdw - Sip & Puff Required when selecting PN I10573, I10574 or I106807-99-0.	296.00

### PROPORTIONAL FOOT CONTROL:

Part Number	Description	MSRP - CA
<input type="checkbox"/> I10767	Proportional Foot Control Not available with Active Reach Packages. Mount on Left: _____ Mount on Right: _____	3,215.00
<input type="checkbox"/> I10766	Adjustable Mounting Hrdw - Foot Control Required when selecting PN I10767.	296.00

## Total Control Head Array

## Total Control Head Array

### CONFIGURE HEAD ARRAY:

Part Number	Description	MSRP - CA
<input type="checkbox"/> I106805-99-0	Permobil Head Array - 3-Switch System <i>Features an occipital pad sensor and two stainless steel arms configured with proximity sensors. The arms are mounted to the upper attachment points on the head array. Longer arms are also included to extend the sensor adjustment range. ESP recommended. Choose Omni Control Module separately.</i>	5,650.00
<input type="checkbox"/> I11072	Program Head Array as a 5-Switch System <i>Select additional arms/sensors below and indicate the intended use of each sensor in the Notes &amp; Comments section. Choose Omni Control Module separately.</i>	No Charge
CHOOSE OCCIPITAL PAD SIZE (REQUIRED)		
<input type="checkbox"/> I105920-99-0	Small Occipital Pad for Head Array - 4.25"W x 3"H	No Charge
<input type="checkbox"/> I105921-99-0	Large Occipital Pad for Head Array - 6"W x 4"H	No Charge
CHOOSE OCCIPITAL PAD MOUNT (REQUIRED):		
<input type="checkbox"/> I10847	Fixed Mount for Occipital Pad - Head Array	No Charge
<input type="checkbox"/> I105936-99-0	Adj. Swivel Mount for Head Array Occipital Pad <i>Enables Occipital Pad to be angled up/down and left/right independent of the lateral arms.</i>	506.00
CHOOSE HEAD ARRAY MOUNT HARDWARE (REQUIRED):		
<input type="checkbox"/> I105927-99-0	Stealth TWB480 Backrest Mounting Hdwr with Links <i>Please select the length of the vertical stem below: Standard length (10"): _____ Short length (6"): _____</i>	314.00
<input type="checkbox"/> I105941-99-0	Stealth TWB480 Backrest Mnt Hdwr with Flipback	519.00
<input type="checkbox"/> I105925-99-0	Stealth HMO475 MultiAxis Backrest Mounting Hdwr	266.00
SELECT LEFT ARM PACKAGE (OPTIONAL)		
<input type="checkbox"/> I10976	Add Left Lateral Arm w Proximity Sensor <i>Choose Lateral Arm Adjustment Range: Short (5.5" - 10"): _____ Long (5.5" - 14"): _____</i>	428.00
<input type="checkbox"/> I10977	Add Left Lateral Arm w Egg Switch <i>Choose Lateral Arm Adjustment Range: Short (5.5" - 10"): _____ Long (5.5" - 14"): _____</i>	285.00
<input type="checkbox"/> I10978	Add Left Lateral Arm w Pad (No Sensor or Switch) <i>Choose Lateral Arm Adjustment Range: Short (5.5" - 10"): _____ Long (5.5" - 14"): _____</i>	174.00
SELECT RIGHT ARM PACKAGE (OPTIONAL):		
<input type="checkbox"/> I10979	Add Right Lateral Arm w Proximity Sensor <i>Choose Lateral Arm Adjustment Range: Short (5.5" - 10"): _____ Long (5.5" - 14"): _____</i>	428.00
<input type="checkbox"/> I10980	Add Right Lateral Arm w Egg Switch <i>Choose Lateral Arm Adjustment Range: Short (5.5" - 10"): _____ Long (5.5" - 14"): _____</i>	285.00
<input type="checkbox"/> I10981	Add Right Lateral Arm w Pad (No Sensor or Switch) <i>Choose Lateral Arm Adjustment Range: Short (5.5" - 10"): _____ Long (5.5" - 14"): _____</i>	174.00

## Alt Drive Systems Setup

### CHOOSE OMNI (REQUIRED):

Part Number	Description	MSRP - CA
<input type="checkbox"/> I108500-99-0	Omni2 Control Module for R-net - Corpus <i>Two-section design features 3.5" VGA high resolution display, USB charger port, built in IR control and integrated Bluetooth technology for control of up to 4 Bluetooth devices. Must also select mounting hardware below. Please specify which type of alternative drive input device you will be using. Input Device: _____</i>	1,995.00
<input type="checkbox"/> I103476-99-0	Adjustable Swing-Away Hardware - Omni/Tablet <i>Mount in front of Left Arm: _____ Mount in front of Right Arm: _____</i>	290.00

## Alt Drive Systems Setup

### CHOOSE MOUNTING POSITION OF SEAT CONTROL SWITCHBOX:

Part Number	Description	MSRP - CA
<input type="checkbox"/> I10080	Mount Seat Control Switchbox on Backrest	No Charge

## Alt Drive Accessories

### PERMOFIX MOUNTING SOLUTIONS:

Part Number	Description	MSRP - CA
<input type="checkbox"/> I106655-99-0	PERMOfix UniTrack Fixed Mount Kit <i>Includes two base joints, two 4" rods and one 8" rod.</i>	318.00
<input type="checkbox"/> I106656-99-0	PERMOfix UniTrack Removable Mount Kit <i>Includes quick release mount and thigh support hardware, two base joints and three different length rods.</i>	559.00
<input type="checkbox"/> I106657-99-0	PERMOfix UniTrack Swing Away Mount Kit <i>Includes PERMOfix Swing Away module, two base joints, one 2" rod and one 4" rod.</i>	806.00
<input type="checkbox"/> I105907-99-0	Proximity Sensor + Swivel Assembly <i>Only available with the Total Control Head Array. Swivel assembly allows sensor to be mounted with PERMOfix or directly to lateral arm of Total Control Head Array. Since Proximity Sensors require power, they CANNOT be used as On/Off switches.</i>	416.00
<input type="checkbox"/> I105919-99-0	Egg Switch + Swivel Mount Assembly <i>Swivel assembly allows switch to be mounted with PERMOfix or directly to lateral arm of Total Control Head Array.</i>	277.00
<input type="checkbox"/> I105912-99-0	Solid Oval Pad + Swivel Assembly <i>Swivel assembly allows pad to be mounted with PERMOfix or directly to lateral arm of Total Control Head Array.</i>	169.00
<input type="checkbox"/> I106675-99-0	Oval Plate + Swivel Assembly - No Switch <i>Blank plate provides mounting option for Micro Lite, Buddy Button or other switch options. Swivel assembly allows pad to be mounted with PERMOfix or directly to lateral arm of Total Control Head Array.</i>	146.00

### OPTIONAL ITEMS:

Part Number	Description	MSRP - CA
<input type="checkbox"/> I101769-99-0	Microlight Switch Qty: _____	183.00
<input type="checkbox"/> I101767-99-0	Buddy Button Qty: _____	138.00
<input type="checkbox"/> I106293-99-0	Egg Switch Black Qty: _____	152.00
<input type="checkbox"/> I101771-99-0	Mini Cup Switch Qty: _____	135.00
<input type="checkbox"/> I107500-99-0	Proximity Switch Qty: _____	640.00
<input type="checkbox"/> I107495-99-0	Mini Joystick 5 Switch OC for Omni - (9-pin)	1,087.00
<input type="checkbox"/> I10772	Dome Switch, Backrest Mnt w Gooseneck Qty: _____ <i>Mount on Left: _____ Mount on Right: _____ Mounts at top of backrest shell. Custom charge may be applied if ordering with short backrest heights. Kit Includes: One PERMOfix joint with Flat Plate, two PERMOfix joints with 45° swivel adjustment, one 150mm rod, one 200mm rod and one 225mm gooseneck switch with 1.2m cable and mounting clamp.</i>	570.00
<input type="checkbox"/> I10846	Dome Switch, Backrest Swg Awy Mnt w Gsneck Qty: _____ <i>Mount on Left: _____ Mount on Right: _____ w 90° Swing-Away Function. Mounts at top of back. Custom charge applies if ordered with short back. Kit Includes: One PERMOfix joint with Flat Plate, one PERMOfix joint with 45° swivel adjustment, two 6" rods, one 90° Swing-Away joint, one 9" gooseneck switch w 50" cable and mounting clamp.</i>	855.00
<input type="checkbox"/> I105935-99-0	3.5mm Ext. Cable for Sensor/Switch, 0.3m Long	39.00
<input type="checkbox"/> I105931-99-0	3.5mm Ext. Cable for Sensor/Switch, 1m Long	43.00
<input type="checkbox"/> I105932-99-0	3.5mm Ext. Cable for Sensor/Switch, 2m Long	47.00
<input type="checkbox"/> I105933-99-0	3.5mm Ext. Cable for Sensor/Switch, 3m Long	52.00
<input type="checkbox"/> I105934-99-0	3.5mm Y-splitter Cable for Sensor/Switch <i>Allows two Proximity Sensors/Mechanical Switches to connect to the same direction jack.</i>	43.00

Notes & Additional Instructions:

**Chair Order Policy:**

You are highly encouraged to contact Permobil Sales with ANY questions regarding proper completion of our order forms at 1-800-736-0925. The order cannot be processed without the following:

- A client name or code, the client's measurements, the client's weight
- The order form must also be signed by a person authorized by your company to acknowledge the items selected. Permobil, LTD. is not responsible for configuration or size discrepancies resulting from customer errors on the order form.
- Order confirmations will be provided by Permobil summarizing the items selected. Please review this carefully.
- Once the order is shipped, any changes are subject to the returns policy and may be prohibited.

**Order Acknowledgment:**

I, \_\_\_\_\_, am an agent of the medical equipment provider named on this order form and I have the authority to contract for the purchase of powered wheelchairs and related parts on behalf of said provider. I acknowledge that I have reviewed this order and that it is complete and accurate to the best of my knowledge. I further acknowledge that any changes to this order after submission of this order form are subject to the returns policy and may carry additional charges.

Signed: \_\_\_\_\_

Printed Name: \_\_\_\_\_

Title: \_\_\_\_\_

Date: \_\_\_\_\_

NOTE: All specifications and prices are subject to change without notice. Please note that prices displayed are only valid if ordered with the wheelchair. Any parts ordered for service or future alterations will carry different pricing and warranty.



# BodiLink® Corpus Accessories Order Form - Canada

## BODILINK HEAD SUPPORT

### HEADREST HARDWARE:

PART NUMBER	DESCRIPTION	MSRP - CA
I108721-99-0	Adj/Removable BodiLink Head Support Hardware - MAX 8 <i>Must be selected when ordering &gt;2 Swing-Away headrest accessories or Chin Boom assemblies. NOTE: Must choose Head Support Pad on page 13.</i>	440.00
<b>SELECT HEADREST BRACKET:</b>		
I108722-99-0	BodiLink Headrest Bracket, Corpus <i>For 23" - 28" Corpus Ergo backrest heights.</i>	141.00
I108723-99-0	BodiLink Headrest Bracket, Corpus - Short <i>For 20" Corpus Ergo backrests &amp; AGILITY Max Contour Back Systems.</i>	141.00

## BODILINK LATERAL TRUNK SUPPORT

### PREMIUM LATERAL SUPPORT PADS:

PART NUMBER	DESCRIPTION	MSRP - CA
I109029-99-0	BodiLink Premium 3" x 5" Lateral Pad - Pair <i>Must also select I109031-99-0. Includes foam wedge for additional support. Only compatible with Corpus Ergo backrest shells.</i>	448.00
I109030-99-0	BodiLink Premium 4" x 6" Lateral Pad - Pair <i>Must also select I109031-99-0. Includes foam wedge for additional support. Only compatible with Corpus Ergo backrest shells.</i>	448.00

### LATERAL TRUNK SUPPORT HARDWARE:

PART NUMBER	DESCRIPTION	MSRP - CA
I109031-99-0	Swing-Away Adjustable BodiLink Lateral Hardware, Corpus - Pair	690.00

## BODILINK LATERAL PELVIC/THIGH SUPPORT (LPTS)

### PREMIUM LATERAL PELVIC/THIGH SUPPORT PADS:

PART NUMBER	DESCRIPTION	MSRP - CA
I109040-99-0	BodiLink Premium LPTS Pad (B), Stretch-Air - Pair <i>4.5"W x 3.5"T Pad w/additional foam insert, Stretch-Air cover and mounting hardware. Must also select PN I109044-99-0.</i>	270.00
I109043-99-0	BodiLink Premium LPTS Pad (A), Stretch-Air - Pair <i>6"W x 5.5"T Pad w/additional foam insert, Stretch-Air cover and mounting hardware. Must also select PN I109044-99-0.</i>	270.00
I109041-99-0	BodiLink Premium LPTS Pad (C), Stretch-Air - Pair <i>8"W x 3.5"T Pad w/additional foam insert, Stretch-Air cover and mounting hardware. Must also select PN I109044-99-0.</i>	270.00
I109042-99-0	BodiLink Premium LPTS Pad (G), Stretch-Air - Pair <i>12"W x 3.5"T Pad w/additional foam insert, Stretch-Air cover and mounting hardware. Must also select PN I109044-99-0.</i>	270.00

### LATERAL PELVIC/THIGH SUPPORT HARDWARE:

PART NUMBER	DESCRIPTION	MSRP - CA
I109044-99-0	Adj/Removable BodiLink LPTS Hardware Med, Corpus - Pair <i>Medium TT Power Mount hardware can be extended upward to clear cushions of 3.0-5.5" thick.</i>	667.00